



www.kiranmedicareequipments.com
www.kiransurgicalsystems.com



ARTHOSCOPY PRODUCTS

KSS-A205

LOADED SUTURE ANCHOR SCREW



Dia. : 3 mm & 5 mm Sterile

KSS-A206

DISPOSABLE CANNULA



Size : 8.25 mm x 70 mm
8.25 mm x 90 mm

KSS-A207

ENDO BUTTON WITH LOOP LOCK



Size :
10 mm / 15 mm / 20 mm /
25 mm / 30 mm / 35 mm /
40 mm / 45 mm (Sterile)

KSS-A207

LIGAMENTAL SPIKE STAPPLE



Size : Small / Medium / Large

KSS-A208

TIBIAL POST SCREW



Length : 20 mm to 40 mm


KSS-A209 ROUND HOLE SUTURE DISK



KSS-A210 D' HOLESUTURE DISK



KSS-A211 LOOP BUTTON



KSS-A212 SPIKED WASHER




KSS-A213 LOW-PROFILE CANCELLOUS SCREW 4.0 mm



Size : 25 mm to 45 mm
(each at interval of 5 mm)

KSS-A214 LOW-PROFILE CANCELLOUS SCREW 6.5 mm



Size : 25 mm to 45 mm
(each at interval of 5 mm)

KSS-A216 BIOABSORBABLE SCREW




Diameter : 7,8,9,10 mm
Length : 25 mm, 30 mm

Diameter : 11 mm
Length : 35 mm

Available In : PLA Material & PLA + β TCP Material

PCM101

MINI STRAIGHT PLATE - WITHOUT GAP



CAT No. : PCM101
 Size : 2 Hole to 50 Hole
 1.5 mm, 2.0 mm, 2.5 mm

PCM102


MINI STRAIGHT PLATE - WITH GAP



CAT No. : PCM102
 Size : 2 Hole to 8 Hole
 1.5 mm, 2.0 mm, 2.5 mm

PCM103


MINI 'T' SHAPE PLATE - WITH CAP / WITHOUT CAP



CAT No. : PCM103
 Size : Small, Big, Extra Big
 1.5 mm, 2.0 mm, 2.5 mm

PCM104

MINI 'L' SHAPE PLATE - 90°




Left

Right

1.5 mm, 2.0 mm, 2.5 mm

PCM105

MINI 'L' SHAPE PLATE - 120°




Left

Right

1.5 mm, 2.0 mm, 2.5 mm

PCM106 MINI 'H' SHAPE PLATE




CAT No. : PCM106
Size : 1.5 mm, 2.0 mm

PCM107 MINI 'I' SHAPE PLATE




CAT No. : PCM107
Size : 1.5 mm, 2.0 mm, 2.5 mm

PCM108 MINI 'X' SHAPE PLATE




CAT No. : PCM108
Size : 1.5 mm, 2.0 mm, 2.5 mm

PCM109 MINI 'Y' SHAPE PLATE




CAT No. : PCM109
Size : 1.5 mm, 2.0 mm, 2.5 mm

PCM110 MINI 'Z' SHAPE PLATE



CAT No. : PCM110
Size : 1.5 mm, 2.0 mm, 2.5 mm

PCM111 MINI ORBITAL CAP WITH GAP




CAT No. : PCM111
Size : 4 Hole, 6 Hole, 8 Hole
1.5 mm, 2.0 mm, 2.5 mm

PCM112 MINI ORBITAL CAP W/O GAP




CAT No. : PCM112
Size : 4 Hole, 6 Hole, 8 Hole
1.5 mm, 2.0 mm, 2.5 mm

PCM118 MINI 3D PLATE




CAT No. : PCM116
Size : 4 Hole, 6 Hole, 8 Hole
2.0 mm

PCM124 WIRE MESH (REGULAR)



CAT No. : PCM124
Size : 1x1, 2x2, 2x4, 2x6, 2x8, 4x4 inch
1.5 mm, 2 mm

PCM125 WIRE MESH (PRE BENT)



CAT No. : PCM125
Size : 40x40mm, 60x60mm, 60x80mm, 100x100mm,
1200x120mm, 120x150mm, 140x180mm,

PCM-ASD SCREW DRIVER



PCM-201 MINI SCREW CROSS SLOTTED

PCM-202 PCM-201 FOR 1.5 mm
SIZE : 3 mm to 12 mm







PCM-203 PCM-202 FOR 2 mm
SIZE : 4 mm to 16 mm

PCM-203 FOR 2.5 mm
SIZE : 4 mm to 16 mm



COLOR CODING

The Color Coding of the Screw Diameter

-  1.2 mm Screw Diameter (Emergency Screw)
-  1.5 mm Screw Diameter
-  1.7 mm Screw Diameter (Emergency Screw)
-  2.0 mm Screw Diameter
-  2.3 mm Screw Diameter (Emergency Screw)
-  2.5 mm Screw Diameter
-  2.7 mm Screw Diameter (Emergency Screw)
-  2.7 mm Screw Diameter

The Color Coding of the Plate Thickness

-  0.8 mm Plate Thickness
-  1.0 mm Plate Thickness
-  1.2 mm Plate Thickness
-  1.5 mm Plate Thickness

BH-301 CANNULATED DUAL FUNCTION DRILL

FEATURES

- Dual Combination Of Drilling & Reaming
- Completely Autoclavable Handpiece
- 4.0 mm Hollow Cannulation For All Brand Guide Needle To Pass Through
- Available With Various Attachments Like K-wire, hudson,reamer & 3.5,4.5 Screw Driver
- Dual Locking Base For Battery Locking
- Battery Of 14.4 V



BH-202 OSCILLATING SAW

FEATURES

- Quick Locked Structure Which Gives You Very Precise Cuts During Joint & Thoracic Operations
- Completely Autoclavable Handpiece
- Dual Locking Base For Battery
- Blades Available In Various Sizes
- Battery Of 14.4 V



ATTACHMENTS



S.S. DRILL CHUCK



SCREW DRIVER



FLEXIBLE REAMING ADAPTOR



K-WIRE

ACCESSORIES



BLADES



BATTERY



CHARGER



CASE

BH-302 SMALL BONE DRILL & BH-302C CANNULATED SMALL BONE DRILL

FEATURES

- SUITABLE FOR HAND & FOOT TRAUMA OPERATION
- AVAILABLE IN CANNULATED & NON CANNULATED OPTIONS WITH VARIABLE SPEED CONTROL
- BATTERY OF 7.2 V



BH-203 SMALL BONE SAW

FEATURES

- SUITABLE FOR HAND & FOOT TRAUMA OPERATION
- BATTERY OF 7.2 V



BLADE OPTIONS

LENGTH X WIDTH X THICKNESS (mm)

- 37 x 12 x 0.3
- 40 x 6.6 x 0.5
- 40 x 8 x 0.5
- 46 x 10 x 0.5
- 55 x 15 x 0.5
- 22 x 7.9 x 0.3
- 22 x 11.7 x 0.3
- 22 x 08 x 0.3
- 22 x 15 x 0.3
- 27 x 14.3 x 0.5
- 28 x 7 x 0.5



BH-403 CRANIOTOMY DRILL & MILL

FEATURES

- Multi Attachment System Comes With Drill & Mill Attachment
- Cranial Drill Attachment With Auto Stop Function Speed Of 0-1000 Rpm
- Cranial Mill (bur) Attachment With Rotation Speed Of 0-40000 Rpm
- Dual Locking Base For Battery
- Battery Of 14.4 V



ORDERING INFO

(SYSTEMS) :



- BH-301 : CANNULATED DRILL (1200 RPM)
- BH-302 : MICRO BONE DRILL (0-1100 RPM)
- BH-302 C : MICRO BONE DRILL CANNULATED (0-1200 RPM)

- BH-304 : SLOW CANNULATED DRILL (700 RPM)
- BH-202 : OSCILLATING SAW (0-18,000 OSC/MIN)
- BH-203 : SMALL BONE SAW (0-15,000 OSC/MIN)
- BH-403 : DRILL : 0-900 RPM/MIN & MILL : 0-40000 OSC/MIN

ATTACHMENTS & ACCESSORIES :

- BH-1101 : DRILL CHUCK WITH KEY
- BH-1102 : FLEXIBLE REAMING ADAPTOR
- BH-1103 : HUDSON ATTACHMENT
- BH-1104 : K-WIRE ATTACHMENT
- BH-1105 : 3.5 SCREW DRIVER

- BH-1106 : 4.5 SCREW DRIVER
- BH-1107 : TRANSFER RINGS
- BH-1108 : ALUMINUM BAG
- BH-1109 : BATTERY 14.4 V
- BH-1110 : BATTERY 7.2 V
- BH-1111 : CHARGER



NEXT GEN POWER TOOLS



MONO AXIAL SCREW

Code No.	Total Length
P101.25	MONO 4.5 - 25MM
P101.30	MONO 4.5 - 30MM
P101.35	MONO 4.5 - 35MM
P101.40	MONO 4.5 - 40MM
P101.45	MONO 4.5 - 45MM



POLY AXIAL SCREW

Code No.	Total Length
P104.25	POLY 4.5 - 30MM
P104.30	POLY 4.5 - 35MM
P104.35	POLY 4.5 - 40MM
P104.40	POLY 4.5 - 45MM
P104.45	POLY 4.5 - 50MM



MONO AXIAL SCREW

Code No.	Total Length
P102.30	MONO 5.5 - 30MM
P102.35	MONO 5.5 - 35MM
P102.40	MONO 5.5 - 40MM
P102.45	MONO 5.5 - 45MM
P102.50	MONO 5.5 - 50MM



POLY AXIAL SCREW

Code No.	Total Length
P105.30	POLY 5.5 - 30MM
P105.35	POLY 5.5 - 35MM
P105.40	POLY 5.5 - 40MM
P105.45	POLY 5.5 - 45MM
P105.50	POLY 5.5 - 50MM



MONO AXIAL SCREW

Code No.	Total Length
P103.30	MONO 6.5 - 30MM
P103.35	MONO 6.5 - 35MM
P103.40	MONO 6.5 - 40MM
P103.45	MONO 6.5 - 45MM
P103.50	MONO 6.5 - 50MM



POLY AXIAL SCREW

Code No.	Total Length
P106.30	POLY 6.5 - 30MM
P106.35	POLY 6.5 - 35MM
P106.40	POLY 6.5 - 40MM
P106.45	POLY 6.5 - 45MM
P106.50	POLY 6.5 - 50MM

5.5 REDUCTION POLY

Code No.	Total Length
P107.40	REDUCTION 5.5 POLY - 40MM
P107.45	REDUCTION 5.5 POLY - 45MM
P107.50	REDUCTION 5.5 POLY - 50MM

6.5 REDUCTION POLY

Code No.	Total Length
P107.40	REDUCTION 6.5 POLY - 40MM
P107.45	REDUCTION 6.5 POLY - 45MM
P107.50	REDUCTION 6.5 POLY - 50MM



CONNECTING ROD (5.5MM)

Code No.	Total Length
P109.50	50MM
P109.60	60MM
P109.70	70MM
P109.80	80MM
P109.90	90MM
P109.100	100MM
P109.120	120MM
P109.140	140MM
P109.160	160MM
P109.180	180MM
P109.200	200MM
P109.400	400MM



TLIF CAGE



Code No.	Diameter*Length
P111.8*25	8*25MM
P111.9*25	9*25MM
P111.10*25	10*25MM
P111.11*25	11*25MM
P111.12*25	12*25MM

TLIF CAGE

Code No.	Diameter*Length
P111.8*25	8*30MM
P111.9*25	9*30MM
P111.10*25	10*30MM
P111.11*25	11*30MM
P111.12*25	12*30MM



PLIF CAGE

Code No.	Diameter*Length
P112.8*25	8*26MM
P112.9*25	9*26MM
P112.10*25	10*26MM
P112.11*25	11*26MM
P112.12*25	12*26MM



EXPANDABLE CAGE

Code No.	Diameter*Length
P210.20*25	10*(20MM TO 25MM)
P210.25*35	10*(25MM TO 35MM)
P210.35*55	10*(35MM TO 55MM)
P210.55*75	10*(55MM TO 75MM)
P212.20*25	12*(20MM TO 25MM)
P212.25*35	12*(25MM TO 35MM)
P212.35*55	12*(35MM TO 55MM)
P212.55*75	12*(55MM TO 75MM)
P214.20*25	14*(20MM TO 25MM)
P214.25*35	14*(25MM TO 35MM)
P214.35*55	14*(35MM TO 55MM)
P214.55*75	14*(55MM TO 75MM)
P216.20*25	16*(20MM TO 25MM)
P216.25*35	16*(25MM TO 35MM)
P216.35*55	16*(35MM TO 55MM)
P216.55*75	16*(55MM TO 75MM)
P218.20*25	18*(20MM TO 25MM)
P218.25*35	18*(25MM TO 35MM)
P218.35*55	18*(35MM TO 55MM)
P218.55*75	18*(55MM TO 75MM)



MESH CAGE

Code No.	Diameter*Length
P203.10	10MM *(20MM-100MM)
P203.12	12MM *(20MM-100MM)
P203.14	14MM *(20MM-100MM)
P203.16	16MM *(20MM-100MM)
P203.18	18MM *(20MM-100MM)
P203.20	20MM *(20MM-100MM)

A.C.P (ANTERIOR CERVICAL PLATING)



Code No.	Total Length
----------	--------------

P322.00	22 MM
P324.00	24 MM
P326.00	26 MM
P328.00	28 MM
P330.00	30 MM
P332.50	32.5 MM
P335.00	35 MM
P337.50	37.5 MM
P340.00	40 MM
P342.50	42.5 MM
P345.00	45 MM
P347.50	47.5 MM
P350.00	50 MM
P352.50	52.5 MM
P355.00	55 MM
P357.50	57.5 MM
P360.00	60 MM
P365.00	65 MM

SELF TAPPING SCREW FOR A.C.P 4.0 MM

Code No.	Total Length
----------	--------------

P301.10	10MM
P301.12	12MM
P301.14	14MM
P301.16	16MM
P301.18	18MM



SELF TAPPING SCREW FOR A.C.P 4.3 MM FOR REVISION

Code No.	Total Length
----------	--------------

P302.10	10MM
P302.12	12MM
P302.14	14MM
P303.16	16MM
P303.18	18MM

SPECIALITY PRODUCT

- Cervical Plate with cage System
- 4.0 mm poly axial scoliosis system
- Posterior Cervical Fixation system

**MINIMALLY INVASIVE SPINE SYSTEM - MIS SYSTEM
(LONG FLANGE)**

5.5 MM POLY

Code No.	Total Length
P401.35	35 MM
P401.40	40 MM
P401.45	45 MM
P401.50	50 MM



6.5 MM POLY


Code No.	Total Length
P402.35	35 MM
P402.40	40 MM
P402.45	45 MM
P402.50	50 MM

**PRE-BENT ROD
6.0MM**

Code No.	Total Length
P450.00	50 MM
P460.00	60 MM
P470.00	70 MM
P475.00	75 MM
P480.00	80 MM
P485.00	90 MM
P490.00	95 MM
P410.00	100 MM
P411.00	110 MM



System now come with stainless steel boll bearing for long maintenance free life. Also system body made from solid aluminum bar not from casting for last long anodizing coating and strong body. And new planetary type gear box for high torque and for better performance.



Driving Unit

220 Volts/5 Amp. A.C. supply.
Heavy Duty motor.
Foot Control for Speed Regulate.
Stand For Driving Unit
Suitcase Bag



Cannulated Drill Handpiece

Steel Boll bearing
Steel Drill Chuck (0-1/4")
Light Wight (600 gms. approx)
Maximum Speed 1200 r.p.m
5.5 mm cannulated
Autoclavable
HardAnodized for long life

Reaming Handpiece

Steel Boll bearing
Quick Coupling for Flex Reamers
Light Wight (600 gms. approx)
Maximum Speed 400 r.p.m
5.5 mm cannulated
Autoclavable
HardAnodized for long life



Saggital Saw Handpiece

Steel Boll bearing
Oscillation Max up to 1300 per min
Light Wight (700 gms. approx)
Five types of blades
Autoclavable
HardAnodized for long life

Flexible Shaft

Steel Boll bearing
Length 1.5 mtrs.
Wt. 1000 gms. approx.
Push-Pull type ends.
Autoclavable



Bone Cutting, Drilling & Reaming System Come with Drill Handpiece With Drill Chuck And Chuck Key, Reaming Handpiece with Quick Coupling Chuck, Saggital Saw with five blade and spanners, Flexible Shaft, Foot Control, Driving Unit With Circuit breaker & stand. And all the above items will fit in suitcase bag.

KSS-A217

SHAVER BLADES

We have full range of compatible disposable Shaver Blades & Bur of Various Companies.



KSS-A218

ACL / PCL INSTRUMENT SET



www.kiranmedicareequipments.com

www.kiransurgicalsystems.com



KIRAN TECHNO SERVICES
PRIVATE LIMITED



Delhi : A - 39, Street No.5, Chander Industrial Estate, I.P Extension, Delhi - 110 092

Gujarat : 308/4, GIDC Estate, Makarpura Industrial Estate, Baroda – 390 010,

Andhra Pradesh : No. 6-5-351/5, Astal Heights, Road No. 1, Banjara Hills, Hyderabad – 500 034

Tamil Nadu : No.33-B4, Sidco Industrial Estate Phase - II, Sundarapuram, Coimbatore - 641 021

Kerala : 333/9, Para Road, Near NIDA, Kanjikode, Palakkad - 678 621

Helpline : +91 9751725000, +91 9751625000, +91 9751425000, +91 9751215000, +91 9751245000

Email : kirangroup108@gmail.com, kiranindia108@gmail.com, kiranbiz108@gmail.com

Website : www.kirangroupofcompanies.com, www.kiran.biz, www.kiranshop.in, www.kiranindia.in,
www.kiransurgicalsystems.com, www.kiranmedicareequipments.com, www.kiranpump.com,
www.kiransolarsystems.com, www.kiranled.com, www.kiranincinerationsystems.com
www.kiranequipments.com, www.kirancompressors.com, www.kiranappliances.com
www.kiranairventilators.com, www.kiranfrpproducts.com, www.kiranpebsystems.com,
www.kiranpanels.com, www.kiranfabssystems.com