

HF2051 EMERGENCY TROLLEYS.



Standard Features:

- > Trolley designed to manage patient in emergency.
- > Size: 650× 475× 970mm
- > ABS molded, 304 grade S.S surface (Corrosion resistance)
- > 2 drawer and cupboard for storage.
- > With Defibrillator Shelf, Extension Table, IV Pole with hooks, a waste bin, CPR Board, External Power Holder and Plug Board, Oxygen tank Holder, PC Transparent Divider, Wire Basket, Non Slip Mat, Seal Lock or Centre Lock and Soft closing Drawer.
- > Self Locking ABS handle.
- > Casters 10 cm Dia with brake

HF2051a : Same as HF2051 but with Size: 750x475x970mm

HF2051b : Same as HF2051 but has 3 drawers.

HF2051c : Same as HF2051b but with Size: 750x475x970mm



HF2052 EMERGENCY TROLLEYS.

Standard Features:

- > Size: 650× 475× 970mm.
- > 5 Full Aluminum-modular Colorful drawer .
- > ABS molded, 304 grade S.S surface (Corrosion resistance)
- > With Defibrillator Shelf, Extension Table, IV Pole with hooks, a waste bin, CPR Board, External Power Holder and Plug Board, Oxygen tank Holder, PC Transparent Divider, Wire Basket, Non Slip Mat, Seal Lock or Centre Lock and Soft closing Drawer.
- > Self Locking ABS handle.
- > Casters of 10 cm Dia with brake

HF2052a : Same as HF2052 but with size 750x475x970mm

HF2052b : Same as HF2052 but in Single color

HF2052c : Same as HF2052b but with size: 750x475x970mm



HF2053 ANESTHESIA TROLLEY



Standard Features:

- > Size: 650×475×970mm
- > ABS molded, 304 grade S.S surface (Corrosion resistance)
- > 3 Modular colorful drawer and cupboard for storage, Self-locking.
- > With Transparent Dispenser /container Box, Extension Table, PC Transparent Divider, Wire Basket, Non Slip Mat, Central Lock and soft closing drawer.
- > Casters of 10 cm Dia with brake

HF2053a : Same as HF2053 but different Size: 750x475x970mm

HF2053b : Same as HF2053 but with 2 drawers and single color.

HF2053c : Same as HF2053b but different size 750x475x970mm



HF2053b

HF2055 ANESTHESIA TROLLEY

Standard Features:

- > Size: 650×475×970mm
- > ABS molded, 304 grade S.S surface (Corrosion resistance)
- > 4 Aluminum-molded drawer.
- > With Transparent Dispenser /container Box, Extension Table, PC Transparent Divider, Wire Basket, Non Slip Mat, Central Lock and soft closing drawer.
- > Casters of 10 cm Dia with brake

HF2055b : Same as HF2055 but with 5 drawers.

HF2055a : Same as HF2055 but in different size 750x475x970mm

HF2055c : Same AS HF2055b but in different size 750x475x970mm



HF2055b



HF2056 MEDICINE TROLLEY

Standard Features:

- > Used for supplying medicine
- > Size: 650×475×900mm
- > ABS molded, 304 grade S.S surface (Corrosion resistance)
- > 5 Aluminum-moded drawer.
- > With Extension Table, Waste Bin, PC Transparent Divider, Wire Basket, Non Slip Mat, Central Lock and soft closing drawer.
- > Casters of 10 cm Dia with brake.

HF2056a : Same as HF2056 but with size 750×475×970mm

HF2056b : Same as HF2056 but with 4 colorful drawers

HF2056c : Same as HF2056b but with size 750×475×970mm



HF2057 TREATMENT TROLLEY

Standard Features:

- > Size : L560×W475×H900mm
- > ABS molded, 304 grade S.S surface (Corrosion resistance)
- > One Drawer and one shelf.
- > With Extension Table, Waste Bin, PC Transparent Divider
- > Casters of 10 cm Dia with brake

HF2057a : Same as HF2057 but with two drawers, Wire Basket, Non Slip Mat, Central Lock and soft closing drawer.



HF2063 DRESSING CHANGE TROLLEY



Standard Features:

- > Size : L560×W475×H900mm
- > ABS molded, 304 grade S.S surface (Corrosion resistance)
- > One drawer design.
- > Two waste bins are attached inwards under the drawer.
- > Castor of 10 cm with brakes.

HF2063a : Same as HF2063 but with size 750x475x900mm

HF2058 INFUSION TROLLEY

Standard Features:

- > Size : L560×W475×H900mm
- > ABS molded, 304 grade S.S surface (Corrosion resistance)
- > One drawer design with shelf, Self-locking., colorful
- > 10 hooks hanging device mechanism on the upper part of trolley
- > With Needle Box, PC Transparent Divider, Wire Basket, Non Slip Mat, Extension Table and soft closing drawer
- > Castors of 10 cm dia with brake.

HF2058a : Same as HF2058 but with different size 750x475x970

HF2058b : Same as HF2058 but with two drawers, Central Lock

HF2058c : Same as HF2058b but with different

Size 750x475x970mm



HF2058B

□ AIRWAY MANAGEMENT/HOUSE KEEPING TROLLEY ULTRA SERIES



HF2059 AIRWAY MANAGEMENT TROLLEY

Standard Features:

- > Size : 800x475x970/1065mm
- > ABS molded, 304 grade S.S surface (Corrosion resistance)
- > 5 drawer design
- > Cart side install fiber lens sideboard, volume 3, top design ASB spring fixture, fixed plastic pipe below, facilitate fiber lens security deposit.
- > Castors of 10 cm dia with brake.

HF2059a : Same as HF2059 with Size 445x110x960/1280mm

HF2059b : Same as HF2059 but with 6 drawers and different color

HF2059c : Same as HF2059b with size 445x110x960/1280mm



HF2059b



HF2061 HOUSE KEEPING TROLLEY



Standard Features:

- > Size: 1050x475x900mm
- > ABS molded, 304 grade S.S surface (Corrosion resistance)
- > Table surface and the below shelf is installed with stainless steel three-side fence
- > Two drawer design
- > Bed linen collection bag is on the other side.
- > With Non Slip Mat and Soft Closing Door
- > Castors of 10cm dia with brakes.

HF2062 HOUSE KEEPING TROLLEY

Standard Features:

- > Size : 1050x475x900mm
- > ABS molded, 304 grade S.S surface (Corrosion resistance)
- > Table surface and the below two shelves are installed with stainless steel three-side fence
- > Bed linen collection bag is on the other side.
- > Castors with 10cm dia with brakes.



HF2200 MEDICINE TROLLEY



Standard Features:

- > Two side Trolley designed to supply medicines. (Front & back same)
- > Size: 850 × 520 × 1045 mm.
- > 2 drawers with dividers suitably designed for keeping medicines.
- > 3 Multi purpose drawers.
- > Equipped with Wastebin, needle disposable container, file cassette & guard rails.
- > Equipped with lock key system.
- > 4 swivel noiseless castor with brakes.



Needle Disposal Container




Drawer & ABS Basket




Guard Rail



Aluminium Column



File Cassette



waste Bin

HF2203 MEDICINE TROLLEY



Standard Features:

- > Trolley made of Alloy, ABS plastic and stainless steel
- > Size : 850mm x 520mm x 1010mm
- > 6 drawers, can be pulled out from both sides, self - locked when closed. 50 medicine boxes in all drawers,
- > History card on each medicine box.
- > With waste bins-2 pcs, Transparent file box-1 pcs, sliding side board, and five hooks hanging device mechanism.
- > 4 swivel noiseless castor of 5" (2 with brakes).



□ EMERGENCY / ANESTHESIA TROLLEYS



HF2204 EMERGENCY TROLLEY



Standard Features:

- > Trolley designed to manage patient in emergency.
- > Size: 850 × 520 × 1045 mm.
- > 5 drawers made of ABS with dividers.
- > Equipped with Wastebin, needle disposable container, multibin container, file cassette, guard rails, Sliding shelf & defibrillator shelf.
- > Equipped with Oxygen Cylinder Holder, I.V. Pole & CPR Board.
- > Equipped with lock key system. (Central locking mechanism)
- > 4 swivel noiseless castor of 5" (2 with brakes).




I.V. Pole



Needle Disposal Container




Sliding Side Shelf




Removable Drawer totes allow for replenishment system



Defibrillator Shelf



File Cassette



Waste Bin

ANESTHESIA TROLLEY



HF2208

Standard Features :

- > Made of mild steel epoxy coated finish.
- > Size: 650 × 420 × 1020 mm.
- > 8 drawers for keeping medicines & goods.
- > One sliding shelf.
- > Equipped with one wastebin, Douche can (SS) & oxygen holder.
- > ABS push handle.
- > 4 Swivel castor all with brakes.
- > Equipped with lock & key.

ANESTHESIA TROLLEY



HF2212

Standard Features :

- > Size: 800 × 490 × 1490 mm.
- > 4 transparent boxes at upper part of trolley.
- > 2 concealed working tops under table top.
- > 5 drawers with locks.
- > With medicine containers/drawers.
- > With 2 multi bin container.
- > 4 Swivel castor all with brakes.

ANESTHESIA TROLLEY



HF2216

Standard Features :

- > Made of mild steel epoxy coated finish.
- > Size: 650 × 420 × 1020 mm.
- > 8 drawers for keeping medicines & goods.
- > Wastebin on both side of trolley.
- > One sliding shelf
- > Douche Can, S.S. & push handle.
- > 4 Swivel castor all with brakes..

MEDICINE TROLLEY



HF2220 Standard Features :

- > For Sickroom Nursing.
- > Size: 650 × 420 × 1020 mm.
- > 6 lockable multipurpose drawers.
- > Adjustable medicine boxes in each drawer.
- > With Wastebin on each side.
- > 4 Swivel castor all with brakes.

Anesthesia Trolley, S.S.



HF2224 Standard Features :

- > Trolley made of stainless steel.
- > Size: 650 × 420 × 1020 mm.
- > 6 drawers suitably designed for keeping medicines.
- > With Wastebin on one side.
- > 4 Swivel castor all with brakes.

MEDICINE TROLLEY



HF2228 Standard Features :

- > For Sickroom Nursing.
- > Size: 610 × 475 × 920 mm.
- > Lockable multipurpose drawers.
- > Adjustable medicine boxes in each drawer.
- > 4 Swivel castor (2 with brakes).

TREATMENT TROLLEY



HF2234 Standard Features :

- > 90% made of ABS.
- > Size: 818 × 525 × 985 mm.
- > 2 concealed working tops under table top.
- > 2 drawers & 2 cupboards.
- > With medicine containers/drawers.
- > With 2 wastebins & I.V. Pole.
- > 4 Swivel castor (2 with brakes).
- > ABS table top, Sz. 700 x490mm

TREATMENT TROLLEY



HF2237 Standard Features :

- > Suitable for emergency area.
- > Size: 650 × 420 × 1020 mm.
- > ABS table top
- > 4 lockable drawers.
- > Supplied with 3 file cassettes.
- > 4 Swivel castor all with brakes.
- > Finish with bumper design.

HF2237a : Same as HF2237, but with S.S. Top.

MEDICINE TROLLEY



HF2240 Standard Features :

- > For Sickroom Nursing.
- > Size: 800 × 490 × 910 mm.
- > Lockable multipurpose drawers.
- > Adjustable medicine boxes in each drawer.
- > 4 Swivel castor (2 with brakes).

EMERGENCY TROLLEY



HF2245 Standard Features :

- > Two side Trolley designed to supply medicines. (Front & back same)
- > Size: 850 × 515 × 985 mm.
- > 4 drawers with dividers suitably designed for keeping medicines.
- > 1 Multi purpose drawers.
- > Equipped with IV pole, Wastebin, needle disposable container, file cassette & guard rails.
- > Equipped with lock key system.
- > 4 Swivel castor (2 with brakes).

EMERGENCY TROLLEY



HF2249 Standard Features :

- > 90% made of ABS.
- > ABS Table Top Sz. 520x610mm.
- > Size: 850 × 528 × 995mm.
- > ABS drawers for keeping medicine & instruments.
- > With Wastebin & I.V Pole.
- > 4 Swivel castor (2 with brakes).

EMERGENCY TROLLEY



HF2253 Standard Features :

- > For first aid patient in emergency.
- > 90% made of ABS.
- > ABS Table Top Sz. 700x490mm.
- > Size: 818 × 525 × 985 mm.
- > 5 drawers of ABS with dividers.
- > 2 concealed working tops under table top.
- > With Wastebin, file cassette & I.V Pole.
- > With lock key system.
- > 4 Swivel castor (2 with brakes).

EMERGENCY TROLLEY



HF2256 Standard Features :

- > For Emergency Ward & Sickroom Nursing.
- > Size: 800 × 490 × 960mm.
- > Table Top can be pulled horizontally towards 2 sides.
- > With Wastebin, file cassette & I.V Pole.
- > 4 Swivel castor (2 with brakes).

EMERGENCY TROLLEY, S.S.



HF2019 Standard Features :

- > Trolley made of stainless steel.
- > Size: 675 × 410 × 895 mm.
- > 2 SS Drawers & Openable Top for ease.
- > 2 Cupboards below the drawers.
- > 4 Swivel castor (2 with brakes).

EMERGENCY TROLLEY, S.S.



HF2013 Standard Features :

- > Trolley made of stainless steel.
- > Size: 650 × 400 × 790 mm.
- > 1 big lockable drawer & 1 small drawer.
- > 2 Cupboards below the drawers.
- > Wastebin on one side of trolley.
- > Supplied with I.V. Rod.
- > 4 Swivel castor all with brakes.

MEDICINE TROLLEY, SUPERIOR, S.S.



- HF2025** Standard Features :
- > Trolley made of stainless steel.
 - > Size: 640 × 450 × 850 mm.
 - > Top surface with 3 side guard railings.
 - > 2 S.S. lockable drawers.
 - > 1 S.S. Shelf below the drawers.
 - > 4 Swivel castor (2 with brakes).

DRESSING TROLLEY, ABS



- HF2268** Standard Features :
- > Size: 820 × 510 × 970 mm.
 - > 2 ABS shelves with guard rails on 3 sides.
 - > Main frame work made of Stainless Steel.
 - > With 1 wastebin & 1 S.S. Bowl.
 - > 1 drawer below the first shelf.
 - > 4 Swivel castor (2 with brakes).

HF2268a : Same as HF2268, but with dressing tray at the top & epoxy coated framework.

DRESSING TROLLEY, S.S.



- HF2072a** Standard Features :
- > Trolley made of stainless steel.
 - > Size: 660 × 440 × 925 mm.
 - > 2 S.S. shelves with guard rails on 3 sides.
 - > With 1 wastebin & 1 S.S. Bowl.
 - > S.S. Dressing tray attached at the top.
 - > 1 S.S. drawer below the first shelf.
 - > 4 Swivel castor (2 with brakes).

INSTRUMENT TROLLEY, ABS



- HF2272** Standard Features :
- > Trolley made of ABS.
 - > Size: 820 × 510 × 970 mm.
 - > Main framework made of S.S.
 - > 2 ABS shelves with guard rails on 3 sides.
 - > 4 Swivel castor (2 with brakes).

INSTRUMENT TROLLEY, S.S.



- HF2096b** Standard Features :
- > Trolley made of stainless steel.
 - > Size: 750 × 470 × 925 mm.
 - > 2 S.S. shelves with guard rails on 3 sides.
 - > 4 Swivel castor (2 with brakes).

INSTRUMENT TROLLEY, S.S.



- HF2105** Standard Features :
- > Trolley made of stainless steel.
 - > Size: 660 × 440 × 870 mm.
 - > 3 SS shelves with guardrails on 3 sides.
 - > 4 Swivel castor

□ INSTRUMENT TROLLEY / OVER BED TABLES

INSTRUMENT TROLLEY, S.S.



HF2091 Standard Features :

- > Trolley made of stainless steel.
- > Size: 640 × 450 × 990 mm.
- > 2 S.S. shelves with guard rails on 3 sides.
- > 1 S.S. drawer below the second shelf.
- > 4 Swivel castor all with brakes.

TROLLEY FOR RECORDS, Epoxy Steel



HF2276a Standard Features :

- > Trolley framework made up of Epoxy coated steel with 4 Aluminium column.
- > SS Working surface with ABS push handle.
- > S.S Guardrail on 3 side.
- > Suitable for 40 clips.
- > Size 850x520x1120mm.
- > 1 Drawer at the bottom.
- > Equipped with Lock & Key.
- > 4 Swivel castor all with brakes.

TROLLEY FOR RECORDS, ABS



HF2159

HF2280a Standard Features :

- > Trolley framework made up of Epoxy coated steel with 4 Aluminium column.
- > ABS Working surface & base.
- > Suitable for 40 clips.
- > Size 650x400x1025mm.
- > Equipped with Lock & Key.
- > 4 Swivel castor all with brakes.

HF2285 : Same as HF2280a, but with 25 clips and size 350x400x1025mm.

HF2159 : Same as HF2280a, but made of stainless steel.

OVER BED TABLE, Height Adjustable



HF2139 Standard Features :

- > Over bed table with Laminated Wooden Top.
- > Size: 795 × 400 × 710 / 970 mm.
- > Framework is made of Epoxy steel.
- > Adjustable height from 710 - 970 mm by gas spring.
- > 4 Swivel castor.

OVER BED TABLE, Height Adjustable



HF2284 Standard Features :

- > Over bed table with ABS top.
- > Size: 770 × 380 × 750 / 1010 mm.
- > Framework is made of Epoxy steel.
- > Adjustable height from 750 - 1010 mm by gas spring.
- > 4 Swivel castor all with brakes.

OVER BED TABLE, Height Adjustable



HF2292 Standard Features :

- > Over bed table with ABS Top.
- > Tilt of Top upto 15°.
- > Size: 770 × 380 × 750 / 1010 mm.
- > Adjustable height from 750 - 1010 mm by gas spring.
- > 4 Swivel castor.

MAYO TYPE TROLLEY, Height Adjustable



HF2126 Standard Features :

- > Trolley with Stainless Steel framework.
- > Used in the O.T. rooms for keeping instruments during surgery.
- > Size: 600 × 400 × 850 / 1000 mm.
- > Adjustable height from 850 - 1000 mm.
- > 2 Swivel castor.

MAYO TYPE TROLLEY, Height Adjustable



HF2122 Standard Features :

- > Trolley with Stainless Steel framework.
- > Used in the O.T. rooms for keeping instruments during surgery.
- > Size: 700 × 430 × 850 / 1000 mm.
- > Adjustable height from 850mm - 1000mm.
- > 4 Swivel castor.

OVER BED TABLE, Height Adjustable



HF2288 Standard Features :

- > Over bed table with Laminated Wooden Top.
- > Size: 750 × 350 × 820 / 1050 mm.
- > Framework made of Stainless steel.
- > Adjustable height from 820 - 1050 mm by gas spring.
- > 4 Swivel castor.

TROLLEY FOR DIRTY LINEN & WASTE ... S.S.



HF2042 Standard Features :

- > Trolley made of stainless steel.
- > Size: 950 × 500 × 900 mm.
- > 3 S.S. shelves with guard rails.
- > Supplied with 1 bag for dirty linen.
- > 4 Swivel castor.

TROLLEY FOR DIRTY LINEN & WASTE, S.S.



HF2048 Standard Features :

- > Trolley made of stainless steel.
- > Size: 800 × 530 × 900 mm.
- > Supplied with one bag for dirty linen fitted on stainless steel framework.
- > 4 Swivel castor (2 with brakes).

TROLLEY FOR DIRTY LINEN & WASTE, S.S.



HF2296 Standard Features :

- > Trolley made of stainless steel.
- > Size: 590 × 480 × 900 mm.
- > Supplied with one bag for dirty linen fitted on stainless steel framework.
- > 4 Swivel castor (2 with brakes).