

TROLLEY FOR WASTE



HF2049 Standard Features :

- > Frame work made of Epoxy powder coated mild steel.
- > Bag made from Fibreflex Fabric.
- > Top cover of Bag made of Plastic.
- > Double bin Sz. 920 x 550 x 865mm.

HF2049a : Single bin Sz. 500 x 550 x 865mm.

TROLLEY FOR CONTAMINANT TUB



HF2300 Standard Features :

- > Framework made of Stainless Steel.
- > Size: 420 x 290 mm.
- > Supplied with plastic contaminant tub.
- > 4 Swivel castor.

MOP WRINGER BUCKET PP



HF2310 Standard Features

- > Made of Polypropylene
- > Size : 460 x 270 x 290mm

NURSE STOOL, S.S.



HF2167 Standard Features :

- > Framework made of Stainless Steel.
- > Cushion seat & back
- > Gas spring height adjustment

REVOLVING STOOL, S.S.



HF2167a Standard Features :

- > Framework made of Stainless Steel.
- > Cushion Seat Top.
- > Gas spring height adjustment

STOOL, S.S.



HF2169 Standard Features :

- > Framework made of Stainless Steel.
- > Top surface is made of Stainless Steel.
- > Size: 320 mm Dia Top x 370 mm Dia Base x 540 mm Height.

HF2168 : Same as HF2169, but with rexine top.

SALINE STAND



- HF567** Standard Features :
- > Saline Stand / IV Stand, Economy, S.S. tubular pipe.
 - > Five legged heavy polypropylene base mounted on 5 cm castor.
 - > 4 hook adjustable S.S. I.V. Rod with a bakelite knob.
 - > Freight saving knock down construction.

SALINE STAND / IV STAND, Economy



- HF568** Standard Features :
- > Saline Stand / IV Stand, Economy, S.S. tubular pipe.
 - > Five legged heavy polypropylene base mounted on 5 cm castor.
 - > 4 hook adjustable S.S. I.V. Rod with a bakelite knob.
 - > Freight saving knock down construction.

HF568a : Same as HF568, but with SS Base.

BEDSIDE SCREEN



- Standard Features:
- > Size : 168H x 244W cms.
 - > Supplied with Rod & Curtains.
 - > Special Knockdown construction.

Frame work made of Epoxy Coated Mild Steel with castors.

- HF626** : BEDSIDE SCREEN - 4 Folds.
- HF628** : BEDSIDE SCREEN - 3 Folds.

Frame work made of Stainless Steel with castors.

- HF629** : BEDSIDE SCREEN - 4 Folds.
- HF623** : BEDSIDE SCREEN - 3 Folds.

STEP STOOL / FOOT STOOL, SS



- Standard Features :
- > Size : Single Step 410x275x210 mm.
 - > Double Step 590x410x200/400 mm.
 - > ABS Platform

Framework	Single Step	Double Step
Stainless Steel	HF583a	HF584a

BOWL STAND, Single



- HF718** Standard Features :
- > Bowl Stand, Single, SS tubular pipe construction mounted on heavy polypropylene base with castor.
 - > Supplied complete with S.S. basin of 35 cm. dia.
 - > Knockdown construction.

BOWL STAND, Double



- HF723** Standard Features :
- > Bowl Stand, Double, SS tubular pipe construction mounted on heavy polypropylene base with castor.
 - > Supplied complete with S.S. basin of 35 cm. dia.
 - > Knockdown construction.