


#### HF5000 STRETCHER SINGLE FOLD WITH TWO WHEELS



**Standard Features:**

- > Light Weight, easy clean, anodized aluminium alloy framework with 2 wheels.
- > Washable, water proof & anti-fungus fabric top with safety belts.
- > Unfolded Size (LxWxH): 185x52x23cm.,  
Folded Size (LxWxH): 93x52x10cm.
- > Self-weight: 5 kg, Load Bearing: 159 kg.



#### HF5007 STRETCHER SINGLE FOLD WITH TWO WHEELS



**Standard Features:**

- > Light Weight, easy clean, anodized aluminium alloy framework with 2 wheels.
- > Washable, water proof & anti-fungus fabric top with safety belts.
- > Unfolded Size (LxWxH): 186x52x19cm,  
Folded Size (LxWxH): 93x52x10cm.
- > Self-weight: 5 kg, Load Bearing: 159 kg.



#### HF5014 STRETCHER SINGLE FOLD WITH ADJUSTABLE BACK REST



**Standard Features:**

- > Light Weight, easy clean, anodized aluminium alloy framework with adjustable back rest & 2 wheels.
- > Washable, water proof & anti-fungus fabric top with safety belts.
- > Unfolded Size (LxWxH): 185x54x24cm,  
Folded Size (LxWxH): 92x54x14cm.
- > Self-weight: 7 kg, Load Bearing: 159 kg.



## □ EMERGENCY SERIES ..... FOLDABLE STRETCHERS




### HF5021 STRETCHER SINGLE FOLD



#### Standard Features:

- > Light Weight, easy clean, anodized aluminium alloy framework.
- > Washable, water proof & anti-fungus fabric top with safety belts.
- > Unfolded Size (LxWxH): 185x52x5cm,  
Folded Size (LxWxH): 93x52x6cm.
- > Self-weight: 3 kg, Load Bearing: 159 kg.

### HF5028 STRETCHER SINGLE FOLD WITH ALUMINIUM TOP



#### Standard Features:

- > Light Weight, easy clean, anodized aluminium alloy framework.
- > Washable, water proof & anti-fungus aluminium top with safety belts.
- > Unfolded Size (LxWxH): 185x52x5cm,  
Folded Size (LxWxH): 94x50x7cm.
- > Self-weight: 7 kg, Load Bearing: 159 kg.

### HF5036 STRETCHER ARMY 3 FOLD WITH TWO WHEELS



#### Standard Features:

- > Light Weight, easy clean, anodized aluminium alloy framework with 2 wheels.
- > Washable, water proof & anti-fungus fabric top with safety belts.
- > Unfolded Size (LxWxH): 188x45x22cm,  
Folded Size (LxWxH): 46x45x15cm.
- > Self-weight: 8 kg, Load Bearing: 159 kg.
- > Supplied with shoulder bag.

### HF5043 STRETCHER ARMY 3 FOLD




**Standard Features:**

- > Light Weight, easy clean, anodized aluminium alloy framework.
- > Washable, water proof & anti-fungus fabric top with safety belts.
- > Unfolded Size (LxWxH): 188x45x4cm,  
Folded Size (LxWxH): 46x45x15cm.
- > Self-weight: 6 kg, Load Bearing: 159 kg.
- > Supplied with carry bag.



### HF5049 STRETCHER DOUBLE FOLD




**Standard Features:**

- > Light Weight, easy clean, anodized aluminium alloy framework.
- > Washable, water proof & anti-fungus fabric top with safety belts.
- > Unfolded Size (LxWxH): 215x54x15cm,  
Folded Size (LxWxH): 108x17x9cm.
- > Self-weight: 5 kg, Load Bearing: 159 kg.
- > Supplied with carry bag.



### HF5056 STRETCHER 4 FOLD



**Standard Features:**

- > Light Weight, easy clean, anodized aluminium alloy framework.
- > Washable, water proof & anti-fungus fabric top with safety belts.
- > Unfolded Size (LxWxH): 215x55x14cm,  
Folded Size (LxWxH): 56x22x20cm.
- > Self-weight: 8 kg, Load Bearing: 159 kg.
- > Supplied with carry bag.



# EMERGENCY SERIES ..... FOLDABLE STRETCHERS



## HF5064 STRETCHER DOUBLE FOLD



### Standard Features:

- > Light Weight, easy clean, anodized aluminium alloy framework.
- > Washable, water proof & anti-fungus fabric top.
- > Unfolded Size (LxWxH): 221x53x15cm,  
Folded Size (LxWxH): 108x18x11cm.
- > Self-weight: 5 kg, Load Bearing: 159 kg.
- > Supplied with carry bag.

**HF5065 :** Same as HF5064 but with Nylon Mesh water proof, easy washing cover.



## HF5071 STRETCHER 4 FOLD



### Standard Features:

- > Light Weight, easy clean, anodized aluminium alloy framework.
- > Washable, water proof & anti-fungus fabric top.
- > Unfolded Size (LxWxH): 210x53x12cm,  
Folded Size (LxWxH): 61x30x12cm.
- > Self-weight: 6 kg, Load Bearing: 159 kg.
- > Supplied with shoulder bag.



## HF5169 STRETCHER VACUUM MATTRESS



### Standard Features:

- > Made of high quality & high resistant PVC material.
- > Adjust perfectly to the shape of body.
- > It allows radiological X-Ray.
- > Unfolded Size (LxWxH): 200x98x7cm, Folded Size (LxWxH): 30x30x50cm.
- > First aid personnel can rapidly and safely adjust the mattress to be soft or hard by using double action pump.
- > Self-weight: 5.5 kg, Load Bearing: 159 kg.



**HF5084 STRETCHER CUM WHEEL CHAIR SINGLE FOLD, with 2 Wheels**




**Standard Features:**

- > Designed for easily carrying patients in narrow and hard passing places.
- > Light Weight, easy clean, anodized aluminium alloy framework.
- > Washable, water proof & anti-fungus fabric top with safety belts.
- > Unfolded Size (LxWxH): 210x55x12cm, Folded Size (LxWxH): 102x55x16cm.
- > Folded in Wheelchair Back-H:138cm Seat-H:53cm
- > Self-weight: 9 kg, Load Bearing: 159 kg.

**HF5085** : Same as HF5084 but with PVC material and Unfolded Size (LxWxH): 190x58x20cm, Folded Size (LxWxH): 105x55x20cm.

**HF5091 STRETCHER ARMY SINGLE FOLD, with Telescopic Lifting Handles**



**Standard Features:**

- > Light Weight, easy clean, aluminium alloy framework with telescopic lifting handles.
- > Washable, water proof & anti-fungus fabric top with safety belts.
- > Unfolded Size (LxWxH): 220x54x15cm, Folded Size (LxWxH): 180x12x15cm.
- > Self-weight: 8 kg, Load Bearing: 159 kg.

**HF5098 STRETCHER ARMY DOUBLE FOLD, with Telescopic Lifting Handles**



**Standard Features:**

- > Light Weight, easy clean, aluminium alloy framework.
- > Washable, water proof & anti-fungus fabric top with safety belts.
- > Unfolded Size (LxWxH): 220x54x15cm, Folded Size (LxWxH): 95x20x10cm.
- > Self-weight: 8 kg, Load Bearing: 159 kg.
- > Supplied with carry bag.

**HF5105 STRETCHER ARMY 4 FOLD, with Telescopic Lifting Handles**



**Standard Features:**

- > Light Weight, easy clean, aluminium alloy framework.
- > Washable, water proof & anti-fungus fabric top with safety belts.
- > Unfolded Size (LxWxH): 220x54x15cm,  
Folded Size (LxWxH): 49x24x19cm.
- > Self-weight: 11 kg, Load Bearing: 159 kg.
- > Supplied with carry bag.

**Hf5113 STRETCHER FOR AMBULANCE CAR**



**Standard Features:**

- > The Bed Surface of the stretcher is made of foamed cushion with adjustable backrest which makes the patient feel more comfortable.
- > It has low frame structure, so that it can move on the ground and can also be lifted as a stretcher.
- > The both ends of the stretcher have opening circles which are used for transporting the patients. Salvage stretcher can be put inside the stretcher.
- > Unfolded Size (LxWxH): 190x55x42cm, Self-weight: 25 kg, Load Bearing: 159 kg.
- > Framework made of light weight aluminium. Maximum backrest adjustable: 65°.
- > Suitable for the ambulances where the wheel height is less than 29cm.



**HF5300** : Fixing/Locking Device for Ambulance

**Imp. Note :** HF5113 & HF5115 must be supplied with HF5300 (Fixing/Locking Device for Ambulance).

**HF5115 :** Same as HF5113 but with Unfolded Size (LxWxH): 189x55x42cm, Self-weight: 25 kg, Load Bearing: 159 kg.


**HF5231 Camping Bed / Battlefield Stretcher**



**Standard Features:**

- > Light Weight, easy clean, aluminium alloy framework.
- > It is mainly used for battle field and outdoor carrying patients and wounded person. When unfolded, it can be used as a simple bed in the tent.
- > Unfolded Size : 190L x 65 W x 42H cm
- > Folded size : 93L x 20W x 20H cm
- > Self Weight : 8 Kg.
- > Load bearing : 159 Kg.
- > Self Weight : 8 Kg

### HF5148 SALVAGE STRETCHER YELLOW



HF5149

**Standard Features:**

- > Designed for fractured and immobile serious patients.
- > Light Weight, easy clean, anodized aluminium alloy framework with safety belts.
- > Size (LxWxH): Max. unfolded Size (LxWxH): 214x42x7cm, Min. unfolded Size (LxWxH): 167x42x7cm.
- > Detachable and folding. This Salvage stretcher can be used with head immobilizer.
- > Length of stretcher adjustable according to patient's body.
- > The foot-end is designed as a narrow-frame structure .
- > Self-weight: 8 kg, Load Bearing: 159 kg.

**HF5149** : Same as HF5148 but in silver color.

**HF5150** : Same as HF5148 but in blue color.


### HF5147 ELECTRIC MOBILE PATIENT LIFT



**Standard Features:**

- > Made up of powder Coated Steel.
- > Electric Pump with battery.
- > Base can be adjusted from 65cm to 120cm
- > Widely used for transporting patients in hospital and home.
- > Max loading bearing: 159Kg
- > Four Swivel castor, two with brakes.


### SPINE BOARD



HF5204

HF5205

HF5206



HF5207

HF5208


HF5209

**Standard Features:**

- > Made of PE material.
- > 100% X-Ray translucent, washable & disinfectable.
- > Compatible with all type of head immobilizer.
- > Self-weight: 7 kg, Load Bearing: 159 kg.

**HF5207** : Same as HF5208 but with Self Wt.: 4 kg and Load bearing capacity: 79 Kg.

### HF5290 Head Immobilizer



**Standard Features:**

- > Made of high-density plastic material, large ear holes for monitoring the patient's ear canal.
- > It is used with Salvage stretcher & Spine Stretchers.
- > Minimal interference with X-ray, MRI or CT scanning procedures
- > Waterproof plastic coating makes it easy to clean & prevents bacterial growth on it.

### HF5293 Belt for Spine Board



# □ AUTOMATIC STRETCHER TROLLEYS ....

## HF5141 STRETCHER CUM CHAIR AUTOMATIC LOADING FOR AMBULANCE



### Standard Features:

- > When loading the stretcher in ambulance, the minimum height is 62cm.
  - > Bed surface of the stretcher is made of foamed cushion, the material is waterproof, the backrest is adjustable.
  - > Mechanical folding legs are controlled by the handles on both sides.
  - > Stretcher can be locked up after being aboard through the chute.
  - > Stretcher can be arbitrarily adjusted by pneumatic spring, it can converted into a chair.
- Framework made of alluminium.
- > It is mainly used in limited space transporting patients, such as lift in hospital, ambulance, city road etc.
  - > Size of castor:Φ125mm
- High position : 191×53×82cm, Low position : 191×53×25cm, Maximum Angle of the Back : 65°, Self-weight : 39 Kg, Load Bearing : 181 Kg.

**HF5142 :** Same as HF5141 but with different design.



**Imp. Note :** HF5141 & HF5142 must be supplied with HF5303 - Fixing/Locking Device for Ambulance



## HF5155 STRETCHER CHAIR / STAIRCASE STRETCHER (FOUR WHEELS)



### Standard Features:

- > Designed for easily carrying patients in places such as stairs, corridors & lifts.
- > Light Weight, easy clean, anodized aluminium alloy framework.
- > Under the stretcher, two flexible uplift rods.
- > With safety belts for patient's safety during transport.
- > The back of the stretcher designed with 4 folded handles.
- > Unfolded Size (LxWxH): 70x52x91cm, Folded Size (LxWxH): 90x54x17cm.
- > Wheel size :Φ100mm.
- > Self Weight: 8 Kg., Load bearing Capacity: 159 Kg.

**HF5156 :** Same as HF5155 but with 2 wheels.

**HF5157 :** Same as HF5155 but the back of the stretcher is supplied with top plastic board.

**HF5158 :** Same as HF5155 but with Black Canvas, Unfolded Size (LxWxH): 50x41x85cm, Folded Size (LxWxH): 19x45x85cm, Self-weight: 11 kg.

### HF5120 STRETCHER AUTOMATIC LOADING FOR AMBULANCE



HF5303-- : Fixing/Locking Device for Ambulance

**Imp. Note :** HF5120 must be supplied with HF5303 - Fixing/Locking Device for Ambulance

#### Standard Features:

- > Light Weight, easy clean, anodized aluminium alloy framework with adjustable Back Angle 65°
- > It has 3 step and 2 wheels of big and small according to any height of ambulance carriage bottom board.
- > Bed surface of stretcher made of foamed cushion with safety belts which makes the patient secure and comfortable.
- > The stretcher has safety lock.
- > The mechanical folding legs are controlled by the handles on both sides.
- > Only one person can load the patient to the ambulance (with this trolley) easily.
- > Size (LxWxH): High Position: 190x55x85cm, Low Position: 190x60x26cm.
- > Mechanical folding legs are controlled by the handles.
- > Stretcher can be locked after being aboard in Ambulance to avoid movement.
- > Size (LxWxH): High Position: 195x55x85cm, Low Position: 195x55x25cm.
- > Self-weight: 38 kg, Load Bearing: 181 kg.

### HF5122 STRETCHER AMBULANCE WITH ADJUSTABLE BACK REST & PULL ALONG HANDLES



HF5303 : Fixing/Locking Device for Ambulance

**Imp. Note :** HF5122 must be supplied with HF5303 - Fixing/Locking Device for Ambulance

#### Standard Features:

- > Light Weight, easy clean, anodized aluminium alloy framework with back rest adjustable to maximum 65° angle.
- > This stretcher design is three-step according to the height of ambulance carriage bottom board. When loading the stretcher up to the ambulance, the lowest height is 64cm.
- > Bed surface of the stretcher is made of hollow plastic board, light weighted, the structure is very sturdy & durable.
- > Stretcher with safety lock, it must be locked up while transporting the patients, and opened out when getting on the ambulance.
- > Mechanical folding legs are controlled by the handles on both sides.
- > Only one person can to load the patient (with this trolley) to the ambulance.
- > The stretcher can be locked up after being aboard through the chute.
- > It's mainly used for hospitals, ambulance and city roads carrying patients and wounded person.
- > The size of castor: Ø150mm wide castors.
- > Size (LxWxH): High Position: 195x55x85cm, Low Position: 195x55x25cm.
- > Self-weight: 40 kg, Load Bearing: 181 kg.

**HF5125 :** Same as HF5122 but bed surface comes with 2 section mattress.

[www.kiransurgicalsystems.com](http://www.kiransurgicalsystems.com) [www.kiranmedicareequipments.com](http://www.kiranmedicareequipments.com), [www.kisusy.com](http://www.kisusy.com), [www.d4h.co.in](http://www.d4h.co.in), [www.d4s.co.in](http://www.d4s.co.in)

## □ STRETCHER TROLLEYS ...

### HF5130 STRETCHER TROLLEY ALUMINIUM WITH ADJUSTABLE HEIGHT, Adjustable Back rest



#### Standard Features:

- > Light Weight, easy clean, anodized aluminium alloy framework with adjustable back rest.
- > Washable, water proof & anti-fungus fabric top.
- > Size (LxWxH): 190x55x88cm, Size of Castors : Ø160mm.
- > Height adjustable with crank mechanism.
- > Supplied with I.V. Rod / Saline Rod.
- > Range from bed surface to ground : 50cm.
- > Self-Weight : 36 kg, Load Bearing: 159 kg.

### HF5135 STRETCHER TROLLEY STAINLESS STEEL



#### Standard Features:

- > Framework made of Stainless Steel.
- > Size (LxWxH): 190x54x85cm.
- > Self-weight: 33 kg, Load Bearing: 159 kg.
- > Removable Table Top.
- > Drop down side rails.
- > Two Diagonally oppsite castor with brakes.

### FIXING DEVICES FOR AMBULANCE STRETCHERS



**Imp. Note** : HF5113 & HF5115 must be supplied with HF5300 - Fixing/Locking Device for Ambulance.



**Imp. Note** : HF5120, HF5122, HF5141 & HF5142 must be supplied with HF5303 - Fixing/Locking Device for Ambulance