

WH2000 WHEEL CHAIR, Economy Type



Standard Features:

- > Chromed Steel Frame.
- > Fixed Armrest, Fixed Footrest.
- > Solid Castor, Solid Rear Wheel With Spokes.
- > Patient Controlled Brakes.

WH2006 : Same as WH2000 but with Cushioned Seat.

WH2012 : Same as WH2000 but with Steel Powder Coating.

WH2018 WHEEL CHAIR, Children



Standard Features:

- > Chromed Steel Frame.
- > Fixed Armrest, Fixed Footrest.
- > Solid Castor, Solid Rear Wheel With Spokes.
- > Patient Controlled Brakes.

WH2024 WHEEL CHAIR, Heavy Duty



Standard Features:

- > Chromed Steel Frame, Double Cross bar.
- > Fixed Armrest, Fixed Footrest.
- > Solid Castor, Solid Rear Wheel Plastic.
- > Patient Controlled Brakes.

WH2036 : Same as WH2024 but with solid rear wheel with Spokes and 51cm seat width.










WH2030 WHEEL CHAIR, Heavy Duty



Standard Features:

- > Chromed Steel Frame, Double Cross bar.
- > Fixed Armrest, Fixed Footrest. Drop Back Handle.
- > Solid Castor, Solid Rear Wheel With Spokes.
- > Patient Controlled Brakes.

See quick comparison chart on page no. 39.111

MAIN PARAMETERS	Over all width	Seat width	Rear Wheel Dia	Front Castor Dia.	Seat Ht. from Floor	Over all height	Loading Capacity	Packed Dimension	Net Weight/ Gross weight	Qty. per 20' Container	Qty. per 40' Container
MEASUREMENTS cm MODEL NO.										CNTS 20' FCL	CNTS 40' FCL
WH2000/WH2012	60(65)	41(46)	59	20	50	87	100	94x23x89	19 21	144	300
WH2018	55	35	55	15	47	87	75	86x23x89	17 19	164	332
WH2024	60(65)	41(46)	59	20	47	91	110	94x28x92	20 22	108	220
WH2036	70	51	59	20	50	91	120	94X28X92	19 21	108	220
WH2030	60(65)	41(46)	59	20	47	91	110	94X28X78	20 22	144	300

□ WHEEL CHAIRS Foldable



Wh2042 WHEEL CHAIR Detachable Foot Rest



Standard Features:

- > Steel Powder Coating.
- > Flip Up Armrest, Detachable Footrest.
- > Solid Castor, Solid Rear Wheel Plastic.
- > Patient Controlled Brakes.

WH2048 : Same as WH2042 but with Solid Rear Wheel With Spokes.

WH2054 : Same as WH2042 but with chromed Steel Frame & Pneumatic Rear Wheel and Drop Back Handle.

WH2060 WHEEL CHAIR Detachable Arm & Foot Rest



Standard Features:

- > Steel Powder Coating.
- > Detachable Armrest, Detachable Footrest.
- > Solid Castor, Solid Rear Wheel With Spokes.
- > Patient Controlled Brakes.

WH2066 : Same as WH2060 but with Chromed Steel Frame & with Fixed Armrest.

WH2072 WHEEL CHAIR Detachable Arm Rest



Standard Features:

- > Chromed Steel Frame.
- > Detachable Armrest, Elevating Footrest.
- > PU Castor, Pneumatic Rear Wheel With Spokes.
- > Patient Controlled Brakes.

WH2078 : Same as WH2072 but with Detachable Foot Rest.

WH2084 WHEEL CHAIR, Heavy Duty



Standard Features:

- > Chromed Steel Frame, Double Cross Bar.
- > Fixed Armrest, Fixed Footrest, Drop Back Handle.
- > Solid Castor, Pneumatic Rear Wheel Plastic.
- > Patient & Attendant Controlled Brakes.

See quick comparison chart on page no. 39.111

MAIN PARAMETERS	Over all width	Seat width	Rear Wheel Dia	Front Castor Dia.	Seat Ht. from Floor	Over all height	Loading Capacity	Packed Dimension	Net Weight/ Gross weight	Qty. per 20' container	Qty. per 40' Container
MEASUREMENTS cm MODEL NO.										CNTS 20' FCL	CNTS 40' FCL
WH2042	57(68)	41(46)	60	20	50	93	110	83x29x95	20 22	130	260
WH2048	56(68)	41(46)	60	20	50	90	110	83x28x93	19 21	144	300
WH2054	58(68)	41(46)	60	15	50	93	100	80x28x76	13 15	168	348
WH2060	60(65)	41(46)	59	20	50	90	100	80x28x91	20 22	126	260
WH2066	60(65)	41(46)	59	20	50	87	100	80x24x89	18 20	172	320
WH2072/WH2078	60(65)	41(46)	61	20	50	90	100	80x28x91	20 22	126	260
WH2084	60(65)	41(46)	40	15	45	89	100	82x28x76	11 13	168	342

WH2090 WHEEL CHAIR Detachable Foot Rest



Standard Features:

- > Steel Powder Coating with Super Seat Width & Double Cross Bars
- > Flip Up Armrest, Detachable Footrest.
- > PU Castor, Pneumatic Rear Wheel Plastic.
- > Patient Controlled Brakes.

WH2096 WHEEL CHAIR Detachable Arm & Foot Rest



Standard Features:

- > Chromed Steel Frame with Double Cross Bars.
- > Detachable Armrest, Detachable Footrest.
- > Solid Castor, Pneumatic Rear Wheel With Spokes.
- > Patient Controlled Brakes.

WH2102 WHEEL CHAIR, Reclining High Back



Standard Features:

- > Steel Powder Coating. Reclining High Back.
- > Detachable Armrest, Elevating Footrest.
- > Solid Castor, Solid Rear Wheel With Spokes.
- > Patient & Attendant Controlled Brakes.

WH2114 WHEEL CHAIR, Aluminium












Standard Features:

- > Aluminium Light Weight Frame.
- > Fixed Armrest, Fixed Footrest.
- > Solid Castor, Pneumatic Rear Wheel With Spokes..
- > Patient Controlled Brakes.

WH2120 : Same as WH2114 but with Patient & Attendant Controlled Brake.

WH2126 : Same as WH2114 but with Drop Back Handle and Patient & Attendant Controlled Brake.

See quick comparison chart on page no. 39.111

MAIN PARAMETERS	Over all width	Seat width	Rear Wheel Dia	Front Castor Dia.	Seat Ht. from Floor	Over all height	Loading Capacity	Packed Dimension	Net Weight/ Gross weight	Qty. per 20' Container	Qty. per 40' Container
MEASUREMENTS cm MODEL NO.										CNTS 20' FCL	CNTS 40' FCL
WH2090	71(68)	56(53)	60	20	54	94	125	83x35x93	21 23	112	232
WH2096	70	51	60	20	50	92	120	80x28x94	20 22	124	260
WH2102/WH2108	60(65)	41(46)	60	20	50	120	90	95x34x90	26 28	92	186
WH2114	61(68)	41(46)	60	20	89	89	110	94x29x92	13 15	108	220
WH2120	63(68)	41(46)	56	15	46	88	100	92x29x88	13 15	108	272
WH2126	63(68)	41(46)	56	15	46	88	100	94x28x76	13.6 15.8	164	300

□ WHEEL CHAIRS Foldable

WH2132 WHEEL CHAIR, Electronic



Standard Features:

- > Wheel Chair Electronic With Safety Lamp.
- > Adjustable Armrest Pad, Detachable Sideboard.
- > Adjustable Angle Backrest.
- > Maximum Distance covered on Single Charge: 20 km.
- > Speed Range: 1-6 km/hour. Climbing Power: 12°.
- > Battery: 28AH, 12V x 2., Motor: 168W x 2.

WH2138 COMMUNE WHEEL CHAIR



Standard Features:

- > Chromed Steel Frame.
- > Fixed Armrest, Fixed footrest.
- > Solid castor, Solid Rear Wheel With Spokes.
- > Plastic Commode seat.
- > Patient Controlled Brakes.

WH2144 COMMUNE WHEELCHAIR



Standard Features:

- > Commode Chair Foldable.
- > Chromed Steel Frame mounted on 4 Wheels.
- > Height Adjustable 75-85 cm.
- > Seat Width 46 cm.

WH2150 COMMUNE CHAIR



Standard Features:

- > Commode Chair Foldable.
- > Chromed Steel Frame mounted on feet.
- > Height Adjustable 75-85 cm.
- > Seat Width 43 cm.










WH5090 EXERCISER



Standard Features :

- > Exerciser.
- > Epoxy Coated Steel Framework.
- > Width 35 cm.
- > Weight: 3.6 Kgs.

See quick comparison chart on page no. 39.111

MAIN PARAMETERS	Over all width	Seat width	Rear Wheel Dia	Front Castor Dia.	Seat Ht. from Floor	Over all height	Loading Capacity	Packed Dimension	Net Weight/ Gross weight	Qty. per 20' Container	Qty. per 40' Container
MEASUREMENTS cm MODEL NO.										CNTS 20' FCL	CNTS 40' FCL
WH2132 (Electronic)	65	43	31	20	52	99	110	82x68x89	50 58	50	104
WH2138	64	44	58	20	49	87	100	94x29x89	23 25	120	244
WH2144	49	46	7.5	7.5	45-55	74	90	52x23x76	9 10	335	670
WH2150	43	49	---	---	45-55	75-85	100	52x23x78	8 9	335	670

Quick Reference of Wheel Chair's Features

Item Code	Material			Type		Arm Rest		Foot Rest		Rear Wheel			Front Castors		Brake		Other Features						
	Chromed Steel Frame	Steel Powder Coating	Aluminium	Standard	Heavy Duty	Reclining High Back	Fixed	Detachable	Fixed	Detachable	Elevating	Solid	Pneumatic	Plastic	Solid	PU	Patient Controlled	Attendant Controlled	Seat Width (cm)	Commode	Loading Capacity (Kg.)	Drop Back Handle	
WH2000	Yellow	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White
WH2006	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White
WH2012	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White
WH2018	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White
WH2024	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White
WH2036	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White
WH2030	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White
WH2042	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White
WH2048	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White
WH2054	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White
WH2060	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White
WH2066	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White
WH2072	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White
WH2078	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White
WH2084	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White
WH2090	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White
WH2096	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White
WH2102	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White
WH2114	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White
WH2120	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White
WH2126	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White
WH2138	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White	White

Yellow box indicates that the feature is available in that model. White box indicates that the feature is not available in that model.