

□ COMMODE CHAIRS & SHOWER CHAIRS



COMMODE CHAIR



WH3000 Standard Features :

- > Commode Chair Foldable.
- > Chromed Steel Frame.
- > Cushioned Top.
- > Seat Width 44 cm.

COMMODE CHAIR



WH3004 Standard Features :

- > Commode Chair Fixed.
- > Epoxy coated Steel Frame.
- > Height Adjustable 43-53 cm.
- > Seat Width 52 cm.

COMMODE FOLDABLE SEAT



WH3008 Standard Features :

- > Foldable Commode Stand.
- > Epoxy coated Steel Frame.
- > Seat Width 36 cm.

SHOWER CHAIR, Aluminium



WH3050 Standard Features :

- > Shower Chair Fixed.
- > Aluminium Framework.
- > Plastic Back & Seat.
- > Height Adjustable 74-84 cm.
- > Seat Width 51 cm.

SHOWER CHAIR, Aluminium



WH3054 Standard Features :

- > Shower Chair Fixed.
- > Aluminium Framework.
- > Plastic Back & Seat.
- > Height Adjustable 73-83 cm.
- > Seat Width 40 cm.

SHOWER CHAIR, Aluminium



WH3058 Standard Features :

- > Shower Chair Fixed.
- > Aluminium Framework.
- > Plastic Seat.
- > Height Adjustable 39-49 cm.
- > Seat Width 51 cm.

MAIN PARAMETERS	Overall Depth	Seat width	Seat Height	Overall Height	Seat Depth	Net Weight Per pc	Carton Size	Pcs. per carton	Gross Wt. Per Carton
MEASUREMENTS cm MODEL NO.						 N.WT(kg)	 (cm)	 PCS CN	 G.W. CN (kg)
WH3000	48	44	47	84	43	8	50x19x75	1	11
WH3004	48	52	43-53	59	40	6	60x38x51	2	13
WH3008	40	38	38	---	40	1.5	68x38x25	6	10
WH3050	49	50	35-45	74-84	30	2.8	69x20x40	2	6
WH3054	49	40	34-44	73-83	37	2.6	55x54x26	2	6
WH3058	40	51	34-44	39-49	30	2	16x18x40	2	5

WALKER FOLDABLE, with Two Wheels



WH4000 Standard Features :

- > Walker foldable with two wheels.
- > Aluminium Framework.
- > Height Adjustable 78-96 cm.
- > Seat Width 44 cm.

WALKER ROLLATOR FOLDABLE, Aluminium



WH4003 Standard Features :

- > Walker Rollator foldable with two wheels.
- > Aluminium Framework.
- > Height Adjustable 81-91 cm.
- > Seat Width 47 cm.

WALKER FOLDABLE, Aluminium



WH4006 Standard Features :

- > Walker foldable.
- > Aluminium Framework.
- > Height Adjustable 77-95 cm.
- > Seat Width 47 cm.

WALKER FOLDABLE, Aluminium



WH4009 Standard Features :

- > Walker foldable.
- > Aluminium Framework.
- > Height Adjustable 77-95 cm.
- > Seat Width 45 cm.

WALKER ROLLATOR FOLDABLE, Aluminium



WH4012 Standard Features :

- > Walker Rollator foldable with 4 wheels.
- > Epoxy Coated Steel Framework.
- > Height Adjustable 81-93 cm.
- > Seat Width 64 cm.

WALKER ROLLATOR FOLDABLE, Epoxy Steel



WH4015 Standard Features :

- > Walker Rollator foldable with 4 wheels.
- > Aluminium Framework.
- > Height Adjustable 81-93 cm.
- > Seat Width 58 cm.

MAIN PARAMETERS	Overall Depth	Seat width	Seat Height	Overall Height	Seat Depth	Net Weight Per pc	Carton Size	Pcs. per carton	Gross Wt. Per Carton
MEASUREMENTS cm MODEL NO.						 N.WT(kg)	 (cm)	 PCS CN	 G.W. CN (kg)
WH4000	49	44	---	83-98	---	2.6	58x18x81	2	7
WH4003	53	47	55-66	82-92	---	4	61x36x94	2	9
WH4006	47	46	---	80-98	---	2.4	54x18x80	2	6
WH4009	47	46	---	78-96	---	2.5	54x18x80	2	6
WH4012	71	64	61	81-93	---	8	60x40x78	1	9
WH4015	78	46	53	78-89	---	8	62x24x82	1	9

□ WALKING STICKS



WALKING STICKS ALUMINIUM, Height Adjustable



WALKING STICKS ALUMINIUM Height Adjustable & with Chromed Steel base.



FOLDABLE WALKING STICKS CUM CHAIR Aluminium framework with Plastic Seat



WH5095 UNDER ARM CRUTCHES, Aluminium



- > Made of Aluminium.
- > Height Adjustable.
Large Size: 134-154 cm.
Medium Size: 134-144 cm.
Small Size: 94-114 cm.

WALKING STICKS ALUMINIUM Height Adjustable & with Epoxy Steel base.



MAIN PARAMETERS	Overall Height	Net Weight Per pc	Carton Size	Pcs per Carton	Gross Weight per carton
MEASUREMENTS cm MODEL NO.	HEIGHT (cm)	N.W. (Kg)	(cm)	PCS CN	G.W. CN (kg)
WH5000	72-94	0.4	66x16x29	20	8.3
WH5002	95.5-118.5	0.55	100x30x30	20	12
WH5004	92-122	0.55	94x30x31	20	12
WH5006	66-88	0.3	64x16x29	20	7
WH5008	77.5-100	0.35	97x34x19	20	7
WH5010	75-97	0.3	68x35x18	20	7
WH5012	120	0.3	54x13x27	40	7
WH2014	78-88	0.3	82x17x29	20	7
WH5020	84	0.9	86x31x41	10	10
WH5050	75-97	1.1	60x28x47	10	12
WH5053	75-95	0.9	76x34x39	10	10
WH5056	75-97	1.1	79x30x44	10	12
WH5059	72-95	1.1	79x34x43	10	12
WH5080	79-94	1.1	79x34x43	10	12
WH5083	65-87	0.9	73x31x35	10	9
WH5086	72.5-95	0.9	75x34x39	10	10