

HF1841 ORTHOPAEDIC BED, FOWLER



Standard Features:

- > Orthopaedic bed with mechanically operated back rest & knee rest all by Super Smooth crank mechanism.
- > Size: 2150 × 940 × 500 mm.
- > Epoxy Coated Steel frame work and 4 section perforated top.
- > Stainless Steel Orthopaedic attachment (without weight chain and weight hangers).
- > Head panel and foot panel equipped with safety lock and roller bumpers.
- > Mounted on protective stumps.
- > Provision for I.V. Rod.
- > Freight saving knock-down construction.

HF1865 ORTHOPAEDIC BED, with ABS Panels

Standard Features:

- > Orthopaedic bed with mechanically operated back rest & knee rest by super-smooth crank mechanism.
- > Size: 2150 × 940 × 500 mm.
- > Epoxy coated mild steel frame work and 4 section perforated top.
- > Removable high quality ABS head & foot panels.
- > Mounted on protective stumps.
- > Bumper Protection at four corners.
- > Provision for I.V. Rod.
- > Freight saving knock-down construction.

HF1868 : Same as HF1865, but with castor.



HF1853 ORTHOPAEDIC BED (with Aluminium Balkan Frame & ABS Panels)

Standard Features:

- > Orthopaedic bed with mechanically operated back rest & knee rest by crank mechanism.
- > Size: 2140 × 940 × 500 mm.
- > Epoxy Coated frame work with beautiful ABS panels. 4 section perforated top.
- > Balkan Frame (made of aluminium) with monkey bars and monkey chains.
- > Collapsible side railings.
- > Four swivel castor (all with individual brakes).
- > Provision for I.V. Rod.
- > Freight saving knock-down construction.



0 to 70°



0 to 35°

HF1851 ORTHOPADIC BED (WITH SS FRAME & ABS PANEL)

Standard Features:

- > Orthopedic traction bed with cranks operated back rest tilting 0-70°, Knee rest tilting 0-40°, Height Adjustment (730mm – 530mm). Left and right leg part can be raised up separately, convenient for treatment and use.
- > Size: L 2235mm × W 960mm × H 530-730mm (Adjustable Height).
- > Epoxy coated Mild steel frame work and 4 section perforated top..
- > Removable High Quality ABS engineering Plastic head & foot panel. Size: 1910mm × 830mm
- > Collapsible aluminum side rails and stainless steel frame with pulley.
- > 125mm dia noiseless castors, 2 with brakes.



0 to 70°



0 to 40°



530 - 730 mm