


**Foam Mattress Plain**

**HF905** : Foam Mattress Plain, 4" (10cm.)  
**HF908** : Foam Mattress Plain, 5" (12cm.)



**Foam Mattress 2 Section**

**HF912** : Foam Mattress 2 Section, 4" (10cm.)  
**HF917** : Foam Mattress 2 Section, 5" (12cm.)




**Foam Mattress 4 Section**

**HF925** : Foam Mattress 4 Section, 4" (10cm.)  
**HF932** : Foam Mattress 4 Section, 5" (12cm.)



**Foam Mattress 5 Section**

**HF940** : Foam Mattress 5 Section, 4" (10cm.)  
**HF948** : Foam Mattress 5 Section, 5" (12cm.)



**WH6000 Anti Bedsore Mattress**

- > Anti Bedsore Mattress based on alternative pressure pad & pump system.
- > Power supply: 220V/50Hz, 1A/250V, Air output: 6-8 liter/min. > Pressure range: 30 mmHg – 100 mmHg.
- > Top cover: Nylon / PU with button
- > Size:190cm(L) x 85cm(W) x 11cm (H).
- > Exchangeable air cell, water proof. Low Noisy, Long life time benefit.

**WH6001** : Same as WH6000, but with CPR valve.



**WH6005 Anti Bedsore Mattress**

- > Anti Bedsore Mattress (Air Jet Type).
- > Tiny Holes on the surface for decreasing humidity.
- > Replaces the foul air & reduce growth of microbes.
- > Power supply: 220V/50Hz, 1A/250V, Air output: 3.5-4 liter/min.,
- > Pressure range: 50 mmHg – 90 mmHg.
- > Material: PVC
- > Size:200cm(L) x 90cm(W) x 6.5cm (H)
- > Low Noisy, Long life time benefit.