

□ BLOOD DONOR CHAIRS



HF2193 BLOOD DONOR CHAIR

Standard Features:

- > Framework made of mild steel epoxy coated finish.
- > Chair has two positions, sitting and lying on the back.
- > Size: 1000 × 660 × 1100 mm.



HF2196 BLOOD DONOR CHAIR CUM BED

Standard Features:

- > Framework made of mild steel epoxy coated finish.
- > In folding position, it can be used as chair.
- > When chair is spread, it can be used as a bed and the patient can lie on it.
- > Size: 1960 × 720 × 400 mm.



HF2197 BLOOD DONOR CHAIR CUM BED

Standard Features:

- > Framework made of mild steel epoxy coated finish.
- > In folding position, it can be used as chair.
- > When chair is spread, it can be used as a bed and the patient can lie on it.
- > Size: 1900 × 530 × 240/400 mm.

