



ORTHOTIC SUPPORTS & REHABILITATION AIDS

...**we are** behind your
every movement



inside

- ▶ CERVICAL AIDS
- ▶ BODY BELTS & BRACES
- ▶ WRIST & FORE ARM PRODUCT
- ▶ FINGER SPLINTS
- ▶ KNEE SUPPORT
- ▶ ANKLE SUPPORTS
- ▶ TRACTION KITS
- ▶ GENERAL AIDS
- ▶ SILICONE FOOT CARE PRODUCTS
- ▶ UTILITY SPLINTS

CERVICAL AIDS



RH101 : CERVICAL PILLOW
Size : UNIVERSAL

Feature :

- > Resilient, firm and washable foam
- > Combined with cotton to give comfortable support to the Neck.

Application :

- > Offers support to patients suffering from cervical spondylosis
- > Provides proper extension for faster recovery
- > Enhances full support to the neck allowing shoulder & upper back muscle to relax
- > Gives relief from stress in case of back pain



RH104 : CERVICAL COLLAR HARD ADJUSTABLE
Size : S, M, L

Feature :

- > High-density polythene sheet
- > Proper foam padding for comfort
- > Double layer facilitates height adjustment
- > Holes for ventilation
- > Velcro closure for easy application

Application :

- > For cervical spine immobilization
- > Used during Pre and post surgical cases
- > Provides rigid support and strong immobilization



RH106 : CERVICAL COLLAR SOFT EYELET
Size : S, M, L, XL

Feature :

- > Washable and high-density PU foam
- > Contoured at the below chin level for proper support to lower cervical spine
- > Holes for better ventilation

Application :

- > Semi-Rigid support to the lower cervical spine
- > Used during cervical spondylosis
- > During whiplash injuries



RH108 : PHILADELPHIA COLLAR
Size : S, M, L, XL

Feature :

- > Two Piece design with proper velcro provision for comfortable fitting. Holes for ventilation
- > Provides Preformed Chin support
- > Front opening for tracheotomy
- > Washable foam to maintain hygiene conditions

Application :

- > Pre and Post operative cases
- > Severe Spondylosis
- > Rigid support to the mild and lower cervical spine



RH109 : CERVICAL COLLAR (BONELESS)
Size : Child, S, M, L

Feature :

- > Soft foam collar to provide mild support to the neck
- > Washable foam and fabric to maintain hygiene

Application :

- > Provide proper posture with mid support
- > Prevents any muscle pull or hyper flexion or extension while travelling
- > Preventive care for neck while reading, watching T.V.
- > Used as semi- rigid collar for daily routine jobs
- > During whiplash injuries or hyper extension injuries



RH110 : AMBULANCE COLLAR
Size : S, M, L

Feature :

- > Single piece design with semi-rigid plaster padded
- > Foam to provide better immobilization
- > Large opening allows tracheotomy

Application :

- > Used in emergency to immobilize neck for accident victims.

BODY BELTS & BRACES



RH201 : TAYLOR'S BRACE/DORSO LUMBAR BRACE, Size : S, M, L, XL

Feature :

- > Elastic abdominal panels provide rigid compression
- > Thoracic straps provide comfortable fixation of brace to the body.
- > Adjustable para-spinal and lateral bars are pre-contoured to keep spine in it's right position
- > Clavicle straps with cotton stockinet for right shoulder posture
- > Foam pads for good cushioning

Application :

- > Used for post-operative cases
- > Mild thoraces and Lumber injuries
- > Herniated Disc problems.



RH207 : CONTOURED LUMBO SACRAL SUPPORT
Size : S, M, L, XL, XXL, XXXL

Feature :

- > Highly porous orthopaedic webbing
- > Hook / Loops provided for perfect fitting and strong immobilization
- > Easy application and removal

Application :

- > Degenerative Lumbar disc disease
- > Osteoporotic pain of the lumbar spine
- > Lumbosacral strain



RH209 : LUMBO SACRAL SUPPORT
Size : S, M, L, XL, XXL, XXXL

Feature :

- > High grade elastic for flexibility
- > Double support for proper fitting and durability
- > Skin friendly fabric and elastic.

Application :

- > Low back pain cases
- > Used for general purpose with mild compression and support



RH302: ABDOMINAL SUPPORT 8"
Size : S, M, L, XL, XXL

Feature :

- > Fully elasticated belt ensures better compression
- > Special stretch construction covers fully the area of incision and firmly holds dressings etc.
- > Elastic made of cotton and imported rubber.

Application :

- > Post surgical cases for mild compression
- > Post delivery cases to support abdominal muscles
- > Abdominal hernia support.



RH305: ABDOMINAL SUPPORT 9"
Size : S, M, L, XL, XXL

Feature :

- > Fully elasticated belt ensures better compression
- > Special stretch construction covers fully the area of incision and firmly holds dressings etc.
- > Elastic made of cotton and imported rubber.

Application :

- > Post surgical cases for mild compression
- > Post delivery cases to support abdominal muscles
- > Abdominal hernia support.



RH308: HERNIABELT (DOUBLE PAD)
Size : Universal

Feature :

- > Lightweight made from good quality surgical elastic with elasticated straps
- > Easy application
- > Provides comfortable pad pressure
- > Pads can be adjusted

Application :

- > Compression support over herniated tissue in the pelvic region before surgery
- > Applied where surgery is contra-indicated.



RH810: UNIVERSAL SHOULDER IMMOBILISER
Size : S, M, L, XL

Feature :

- > Durable fabric based material
- > Double Straps fit over shoulders, cross at the back and knot around waist for even weight distribution
- > Conventional type of design to give maximum immobilization

Application :

- > Shoulder dislocation or fractures
- > Post surgery care.



RH815: CLAVICLE BRACE
Size : S, M, L

Feature :

- > Foam with stockinet for better comfort
- > Straps are provided with metal fasteners.
- > Triangular pattern for two way adjustment at back.

Application :

- > Management of clavicular injuries.
- > Acromio-clavicular and scapular separations due to sprain
- > Posture correction in case of stooped shoulders.



RH816: RIB BELT
Size : S, M, L, XL, XXL

Feature :

- > Washable and porous elastic that enhances comfort
- > Velcro support for proper fitting
- > Anterior position cut in form of two notches for right posture

Application :

- > Used for rib fractures or injuries
- > Post surgery cases
- > During recovery period of fractures / injuries to control sudden activities like coughing, sneezing.

WRIST & FORE ARM PRODUCT



RH416: WRIST BRACE WITH THUMB
Size : UNIVERSAL

Feature :

- > Thumb hole present to fit through the thumb to the wrist region
- > Sets of velcro present on elastic for adjustment.

Application :

- > Sprains and strains
- > Post trauma compression support.



RH417: WRIST BRACE WITH DOUBLE LOCK
Size : S, M, L, XL

Feature :

- > Elastic strap with two sets of velcro closures for better compression
- > Easy to apply and peel to open.

Application :

- > Sprain, strains post trauma compression support
- > Preventive support during sports activities.



WRIST AND FOREARM SPLINT
RH418 : LEFT, RH419 : RIGHT
Size : S, M, L, XL

Feature :

- > Pre shaped splint for immobilization in dorsiflexion
- > High quality fabric and foam padding for comfort wear
- > Thumb flap with velcro closure.

Application :

- > Strain and sprain in the wrist
- > Early cast removal
- > Minor fractures of the wrist and distal forearm.

Orthotic Supports & Rehabilitation Aids



RH423 : TENNIS ELBOW SUPPORT
Size : S, M, L, XL

- Feature :**
- > Foam padded synthetic strap with velcro and buckle provision for proper fittings.
- Application :**
- > Tennis elbow syndrome or lateral epicondylitis
 - > Preventive care support for tennis elbow syndrome.



RH424 : ELBOW BRACE
Size : S, M, L, XL

- Feature :**
- > Soft yarn and heat resistant rubber thread for comfortable and prolonged use
 - > Gives graduated compression.
- Application :**
- > Relieves discomfort from tired, aching or swollen elbow



RH812 : POUCH ARM SLING
Size : Child, S, M, L, XL

- Feature :**
- > Cotton fabric and wide straps to come round the neck to provide least single point pressure
 - > Perfect design to hold the arm in desired position.
- Application :**
- > Shoulder dislocation
 - > Bursitis, soft tissue injuries of the shoulder joint
 - > Support to plastered arm.

FINGER SPLINTS



RH401 : FINGER BASEBALL SPLINT
Size : S M L

- Feature :**
- > For protection and support to tip of finger
 - > Its flaps fold to maintain the finger in fixed position.
- Application :**
- > To prevent contractures due to burns
 - > In case of fractures and sprains



RH402 : FINGER COT
Size : S, M, L

- Feature :**
- > V shaped malleable Aluminium splints designed for easy application
 - > Poly foam padding for comfort
 - > Fits any finger, meant to keep straight for given period of time.
- Application :**
- > Fractures of distal inter phalangeal joints
 - > Collateral ligament injury of the distal inter phalangeal joint.
 - > Protection from Nail bed injuries of fingers.



RH403 : FROG SPLINT
Size : S, M, L

- Feature :**
- > Foam padded aluminium splint in unique shape to cover the dorsal surface of finger to remain open
 - > Four extended arms to ensure good fitting.
- Application :**
- > Hyper flexion injuries
 - > Swan neck deformity
 - > To hold distal inter phalangeal joints in correct alignment
 - > Collateral ligament injuries of distal inter phalangeal joint



RH405 : FINGER EXT. SPLINT
Size : UNIVERSAL

- Feature :**
- > Malleable aluminium used to keep interphalangeal joints in right position
 - > Polyfoam padding for comfort
 - > Fits any finger
 - > Can be moulded by bare hands
- Application :**
- > Extensor tendon injuries of the finger
 - > Collateral ligament injuries of the inter phalangeal joint
 - > Boutonniere's deformity
 - > Immobilization to the fractured finger



RH409 : THUMB SPICA SPLINT
Size : UNIVERSAL

- Feature :**
- > High quality fabric with suede lining for comfort and durability
 - > Allows performance of routine tasks in comfort as it limits extreme motion of the joint of the thumb.
- Application :**
- > Designed for patients with painful wrists and thumb tendinitis
 - > Splint provides light compression to assist in decreasing pain and relieving inflammation



RH425 : FINGER MALLET SPLINT
Size : UNIVERSAL

- Feature :**
- > Splint keeps the finger joints in hyper extension
 - > Aluminium made with padding for comfort wear.
- Application :**
- > Hyper flexion injuries like mallet finger

KNEE SUPPORT



RH501 : TUBULAR KNEE SUPPORT

Size : S, M, L, XL

Feature :

- > Soft yarn for comfort
- > Imported heat resistant rubber used for its prolonged uses

Application :

- > Provides mild support and compression
- > Can be used for arthritic conditions
- > Aid for stiff, swollen and painful knee



RH503 : TUBULAR KNEE SUPPORT WITH HINGES

Size : S, M, L, XL

Feature :

- > High quality and long staple cotton yarn
- > Hinges on the lateral sides provide extra support to injured knee muscles

Application :

- > Support and warmth to weak or painful knees
- > Control painful knee movement
- > Post operative care
- > Preventive care in sports activities.



RH504 : KNEE SUPPORT WITH HINGES

Size : S, M, L, XL

Feature :

- > Wide elastic base with velcro / buckle combination for better compression
- > Patellar opening provided to avoid pressure at patella
- > Two lateral hinges provided for rigid compression

Application :

- > Rigid compression support for strains and sprains
- > Post operative support
- > Can be used for arthritis patients 4. Sports injuries



RH505 : TUBULAR KNEE SUPPORT WITH CENTER HOLE

Size : S, M, L, XL

Feature :

- > Soft yarn with elastic combination for mild compression
- > Heat resistant rubber used
- > Patellar opening provided to avoid pressure at patella

Application :

- > Mild support to injured knee in cases of strains and sprains
- > Soft tissue injuries
- > For stiff swollen and painful knee conditions



RH509 : THIGH BAND

Size : UNIVERSAL

Feature :

- > Elastic and velcro combination for better compression
- > Shaped to give tubular structure.

Application :

- > Compression support in case of hamstring pull
- > Quadriceps muscle strain.



RH519 : R.O.M. KNEE BRACE

Size : UNIVERSAL

Feature :

- > Dial screw to give desired angle for range of motion
- > Foam padded splints to give comfortable fitting

Application :

- > Pre and post operative and rehabilitation care
- > Anterior cruciate ligament (ACL) and Posterior cruciate ligament tear (PCL)
- > Used for Total Knee Replacement (TKR) rehabilitation care.



RH522 : KNEE BRACE SHORT TYPE

Size : S, M, L, XL

Feature :

- > Rigid splints provided on lateral and dorsal side
- > Elastic and Velcro combination gives required compression
- > High quality fabric used for comfort

Application :

- > Post-operative knee brace
- > Support after cast removal



RH523 : KNEE BRACE LONG TYPE

Size : S, M, L, XL

Feature :

- > Rigid splints provided on lateral and dorsal side
- > Elastic and Velcro combination gives required compression
- > High quality fabric used for comfort

Application :

- > Post-operative knee brace
- > Support after cast removal



RH524 : KNEE IMMOBILIZER

Size : S, M, L, XL, XXL

Feature :

- > High quality foam, fused fabric with metallic splinting for comfort and rigid support
- > Velcro closures for easy and quick application
- > Anatomically shaped splinting for proper fit.

Application :

- > An excellent product for support and stability to the injured knee
- > Emergency care in case of fracture, muscle / ligament / tendon injury or dislocation of the knee
- > Post operative rehabilitative care

Orthotic Supports & Rehabilitation Aids

ANKLE SUPPORTS



RH510 : ANKLET (PAIR)
Size : S, M, L, XL

Feature :

- > Soft yarn with elastic combination for comfort and mild compression
- > Heat resistant rubber used.

Application :

- > Mild support to ankle in case of strains or sprains
- > Soft tissue injuries 3. For stiff swollen and painful ankle



RH512 : ANKLE BRACE
Size : S, M, L, XL

Feature :

- > Shoe shaped design
- > Velcro attachment to fit accordingly
- > Detachable plastic splint to give required rigidity

Application :

- > Ligament and soft tissue injuries
- > Minor fracture bracing
- > Post operative cases
- > Early cast removal.



RH513 : ANKLE BINDER
Size : S, M, L, XL

Feature :

- > Soft yarn with elastic combination for comfort and mild compression
- > Heat resistant rubber used
- > Long elastic strap attached with velcro to give figure of 8 support.

Application :

- > Post surgical rehabilitation care
- > For strains and sprains at ankle joint
- > Mild fractures and ligament injuries.



RH514 : ANKLE SUPPORT
Size : UNIVERSAL

Feature :

- > Two sets of moulded plastics padded with foam
- > Adjustable velcro straps
- > Ankle size adjustment for day-to-day activity and comfort

Application :

- > Post surgical rehabilitation care
- > Strains and sprains at ankle joint
- > Mild fractures and ligament injuries.



RH807 : CAST SHOE
Size : S, M, L, XL

Feature :

- > Durable coated fabric construction in shoe shaped design
- > Foam padded rubber sole for durability
- > Velcro straps for proper fitting

Application :

- > Prevent plaster from chipping, cracking and soiling
- > Facilitates stability while walking.



RH1022 : F.DROP SPLINT RIGHT/LEFT
Size : S, M, L

Feature :

- > Well-padded unbreakable, washable and reusable
- > Light weight splint
- > Flat base preventing lateral dropping of foot.

Application :

- > Used in pre-cast conditions
- > To prevent lateral dropping of foot during sleep
- > De-rotation splint to prevent hip- rotation following surgical procedures.

TRACTION KITS




RH604 : CERVICAL TRACTION KIT (SLEEPING)
Size : UNIVERSAL

Feature :

- > Padding on head halter to give comfortable wear
- > Spreader bar and pulley system provided
- > Heavy-duty weight bag with separate graduations for water and sand
- > Heavy-duty strings provided in the kit.

Application :

- > Traction upon medical practitioners recommendation.



RH605 : CERVICAL TRACTION KIT (SITTING)
Size : UNIVERSAL

Feature :

- > Padding on head halter to give comfortable wear
- > Spreader bar and pulley system provision
- > Heavy-duty weight bag with separate graduations for water and sand
- > Heavy-duty strings provided in the kit.

Application :

- > Cervical spine traction in sitting position
- > For short time cervical traction in muscle spasm
- > Cervical spine dislocation for continuous



RH701 : SKIN TRACTION KIT
Size : Child, Adult

Feature :

- > Foam padded heavy-duty fabric and crepe bandage provided to wrap around for below knee skin traction
- > Heavy-duty strings provided
- > Availability of Pulley system
- > Heavy-duty weight bag with separate graduations for water and sand.

Application :

- > Below knee traction to relieve pain in the back and knee region.

...Rehabilitation Aids



RH709 : PELVIC TRACTION KIT

Size : S, M, L, XL

Feature :

- > Foam padded pelvic traction belt with velcro designed to provide traction in the lower spine
- > Spreader bar, pulley system provided
- > Heavy-duty strings provided 4. Heavy duty weight bag with separate graduations for water and sand.

Application :

- > Pelvic traction to relieve low back pain
- > Lumbar disc dislocations.



RH710 : ANKLE TRACTION KIT

Size : UNIVERSAL

Feature :

- > Cotton fabric designed specially to hold the ankle during traction
- > Buckles and velcro system for proper fitting
- > Spreader bar, pulley system provided
- > Heavy-duty strings provided
- > Heavy-duty weight bag with separate graduations for water and sand.

Application :

- > For ankle traction



RH711 : LEG TRACTION KIT

Size : Child, Adult

Feature :

- > Foam padded heavy duty fabric with velcro straps for easy application
- > Comfortable to wear 3. Heavy-duty strings provided
- > Availability of Pulley system
- > Heavy-duty weight bag with separate graduations for water and sand.

Application :

- > Leg traction to treat lower back pain

GENERAL AIDS



RH204 : BACK REST

Size : UNIVERSAL

Feature :

- > Lightweight and portable
- > Spinal support when placed on any chair, car seat, soft, wheel chair etc.
- > Easy to use during long hours chair sitting

Application :

- > Can be used for slip disc patients
- > Removes discomfort and pain after surgery



RH515 : COMPRESSION STOCKING MID THIGH

Size : S, M, L, XL, XXL

Feature :

- > Two ways stretchable seamless stocking from thigh to ankle
- > Soft yarn and elastic combination provides comfortable and uniform compression

Application :

- > Varicose vein conditions
- > Swollen legs and joints pain



RH517 : COMPRESSION STOCKING BELOW KNEE

Size : S, M, L, XL, XXL

Feature :

- > Two ways stretchable seamless stocking from below knee to ankle
- > Soft yarn and elastic combination provides comfortable and uniform compression

Application :

- > Varicose vein conditions
- > Swollen legs and joints pain



RH907 : ORTHOPAEDIC STOCKINETTE

Size : 5 cm , 7.5 cm, 10 cm, 15 cm

Feature :

- > Made up of finest, soft, non-absorbent tubular rib weaved cotton
- > Available in disposable roll that is easy to preserve and keeps contents free from dust.

Application :

- > Used inside plaster.



RH910 : WEIGHT CUFF

Size : 0.5 Kg, 1 Kg, 2Kg

Feature :

- > Vertically weight pockets to bend freely at wrist or ankle
- > Leather cover for easy application

Application :

- > .During graded physiotherapy exercises



RH818 : EXERCISING BALL

Size : UNIVERSAL

Feature :

- > Soft version of polyurethane foam for proper resistance
- > Support for orthopaedic patients

Application :

- > Wrist and forearm injuries
- > Minor wrist fractures and sprain

SILICONE FOOT CARE PRODUCTS



RH525 : HEEL PAD
Size : S, M, L

Feature :

- > Can be fitted to any footwear instantaneously
- > Easy to use for both feet.

Application :

- > Prevents spreading and flaring away of the heel fat pad, confines it to its natural position and shape
- > Provides natural and effective cushioned interface between the calcaneum and the ground.



RH526 : HEEL CUSHION
Size : S, M, L

Feature :

- > Silicone based material to give proper cushioning
- > Cup shaped design for better compression at heel joint

Application :

- > For spur heel
- > Plantar fasciitis
- > Bursitis of heel



RH527 : ARCH PAD
Size : UNIVERSAL

Feature :

- > Cushion pad made from high grade silicone to support the plantar arch of foot
- > Designed to maintain the plantar arch comfortable

Application :

- > Pain at arch point
- > Repetitive strains and injury at foot arch



RH528 : INSOLE FULL LENGTH
Size : S, M, L & XL

Feature :

- > Anatomically shaped, made from a skin-like medical grade silicone that resists bacterial growth the does not dry out, extending the life of the product for longer lasting comfort.
- > Designed to reduce stress to the joints and spine even under constant pressure.

Application :

- > For spur heel, Plantar fasciitis, Bursitis of heel



RH529 : TOE SEPARATOR
Size : UNIVERSAL

Feature :

- > A corrective shield placed in between toes, to prevent overlapping of the toes.
- > Relieves pain in Halux Valgus.

Application :

- > Overlapping toes, Hammer toes, Bunions,
- > Toes that cross over.



RH904 : META TARSAL SUPPORT
Size : UNIVERSAL


Feature :

- > Padding attached with elastic to give correct contour. Cotton fabric to give comfortable wear.


Application :

- > Flat foot conditions.


UTILITY SPLINTS




RH819 : FRACTURE BRACING JOINT
Size : Universal



RH906 : CRAMER'S WIRE SPLINT
Size : Available in various sizes.



U' SPLINT
RH1012 : Short Arm, RH1013 : Long Arm
RH1014 : Short Leg, RH1015 : Long Leg



RH905 : THOMAS SPLINT
Size : Child, S, M, L, XL



RH1018 : CTEV SPLINT, Size : 1, 2, 3, 4, 5
(Left & Right)



RH1020 : CTEV D.B. SPLINT
Size : 1, 2, 3, 4, 5 (Left, Right)