

GYNECOLOGICAL PRODUCTS.....



VE001a

- > VE001a : Vacuum Extractor Set, Malmstorm type, manual operated. Complete with 3 S.S. Cups, (40, 50, 60mm) PVC pressure tube, traction handle, suction pump, glass jar in between top and bottom metal plates fitted with vacuum gauge and control valve. Packed in plastic case.



VE019

- > VE019 : Vacuum Extractor Set, manual operated. Complete with 3 Silicone Cups (without valves), (40, 50, 60mm) Autoclavable Silicone pressure tube, suction pump, glass jar fitted between M.S. Basket with vacuum gauge and control valve. Packed in plastic case.



VE026

- > VE026 : Vacuum Extractor Set, Bird Type, manual operated. Complete with 3 stainless steel cups (40, 50, 60mm). Autoclavable Silicone pressure tube. Traction handle. Suction pump, glass jar fitted between M.S. Basket with vacuum gauge and control valve.



VE027

- > VE027 : Vacuum Extractor Set, Manual operated. Complete with 3 stainless steel cups (40, 50, 60mm) & 2 Silicone Cups (without valves). Autoclavable Silicone pressure tube. Traction handle. Suction pump, glass jar fitted between M.S. Basket with vacuum gauge and control valve.




VE028

VE028 : VACCUM EXTRACTOR ELECTRIC WITH S.S. CUPS

- > Housing: ABS Plastic Moulded
- > Capacity: -575 mm Hg \pm 10 , 45 Ltrs/ min
- > Pump Type: Diaphragm Pump
- > Jars: 1 x 2 Ltrs Polycarbonate Jar
- > Filter: Bacterial Filter Autoclavable/Reusable
- > Tubing: 8 mm ID x 2 mtr. (Silicone)
- > Vacuum Gauge: 6.25 cm dia 0-760 mm Hg calibration.
- > Power: 220 V AC, 50 / 60 Hz, 90Watt (110V on request)
- > Noise Level: 55 dB A \pm 3
- > Dimension & Weight: 43 x 22 x 37 cms, 8 Kg.
- > Overflow Safety: Mechanical / Electronic Acoustic Alarm
- > Supplied with 3 S.S. cups (40, 50, & 60 mm dia.)

□ GYNECOLOGICAL PRODUCTS.....




VE045

VE045 : EURO-DEL : Ideal Delivery Aid for Obstetricians : Safe for Mother/Baby


Specifications/Details :

- > Housing: ABS Plastic Moulded
- > Capacity: -600 mm Hg \pm 10 , 12 Ltrs/ min
- > Pump Type: Piston Pump
- > Jars: 1 x 1 Ltr Polycarbonate Jar
- > Filter: Bacterial Filter Autoclavable/Reusable
- > Tubing: 8 mm ID x 2 mtr. (PVC)
- > Vacuum Gauge: 5 cm dia 0-760 mm Hg calibration.
- > Power: 220 V AC, 50 / 60 Hz, 90Watt (110V on request)
- > Noise Level: 55 dB A \pm 3
- > Dimension & Weight: 30 x 14 x 33 cms, 4.5 Kg
- > Overflow Safety: Mechanical / Electronic Acoustic Alarm
- > Venteuse cup : Silicone - 50, 60 & 70mm : minimizes risk of trauma and caphalohematome of fetal cranium.



VE090

- > VE070 : Set of 3 S.S. Vacuum Cups with Rubber tubing
- > VE071 : Set of 3 S.S. Vacuum Cups with PVC tubing
- > VE072 : Set of 3 S.S. Vacuum Cups with Silicone tubing
- > VE090 : Manual Vacuum Pump (Suction Pump)
- > VE091 : Glass Bottle for VE001 & VE001a
- > VE092 : Glass Bottle for VE026
- > VE093 : Vacuum Gauge




VE060

VE062

VE061

SILICONE CUP FOR VACUUM EXTRACTOR :

- > VE058 : 40mm opening dia. With release valve
- > VE059 : 50mm opening dia. With release valve
- > VE064 : 60mm opening dia. With release valve.
- > VE065 : 70mm opening dia. With release valve.
- > VE060 : 50mm opening dia. Without release valve
- > VE061 : 60mm opening dia. Without release valve.
- > VE062 : 70mm opening dia. Without release valve




VE120

VE120 : AMNIHOOK

The amnihook amniotic membrane perforator provides a safe, comfortable method to rupture the amniotic membrane the first time

- > 100% laser vision inspected.
- > Rounded tip and protected point designed to rupture the amniotic membranes with out injury to mother and foetus..
- > 10" long, with molded sides and a non slip grip for precise control



VE121


VE121 : CYTOBRUSH /ENDOCERVICAL BRUSH

- > Cyto brush /Endocervical brush is used for Cervical Papanicolaou (Pap) smear screening
- > Consist of brush and the wooden spoon strip.



VE122


- > VE122 : LATEX GYNECOLOGICAL, LONG CUFF GLOVES
- > VE122-S : Small
- > VE122-M : Medium
- > VE122-L : Large



DS190-B

DS190-A


- UMBILICAL CORD CLAMP
- DS190A : Superior, openable
 - DS190B : Regular, non-openable



RU008

VE086

- TAPE, UMBILICAL, Non-Sterile
- > VE084 : 3mm x 20mts.
 - > VE085 : 3mm x 50mts.
 - > VE086 : 3mm x 100mts.
 - > RU008 : Breast Pump (Glass barrel with rubber bulb).



- DS065 : ABSORBABLE MEDICAL UNDERPAD
- Size : 60 x 90cms



DS060


DS061

- MATERNITY PADS :
- Our maternity pads come with excellent absorbent capacity and unique design which enables them to serve its utility to the fullest. The core pad is made of fluff pulp which gives the pad super absorbent capacity. It comes with extra long tabs which are adjusted with elastic belt to tie around the waist.
- DS060 Maternity Pads with Extra Long Tabs, Size:22.5cm X 7.5cm X 0.4cm
 - DS061 Maternity Pads with Extra Long Tabs, Size:25cm X 8.0cm X 0.9cm



VE080

- VE080 : MIDWIFERY CASE
- Case, Aluminium, for midwife kit, approx 305 x 216 x 150 mm (12" x 8-1/2" x 6") with removable tight fitting lid, collapsible carrying handle on cover, and one sturdy spring loaded clamp at each small end as a locking device. All corners fully rounded and all edges smooth finished.




VS002

VS003


- VAGINAL SPECULUM - CUSCO, Stainless steel
- > VS001S : Small, VS001M : Medium
 - > VS001L : Large
- Vaginal Speculum - Duck Bill (Sims), Stainless steel
- > VS002S : Small, VS002M : Medium
 - > VS002L : Large, VS002XL : Extra Large
- Vaginal Speculum - Graves, Stainless steel
- > VS003S : Small, VS003M : Medium
 - > VS003L : Large
- Vaginal Speculum Graves, Brass chrome plated
- > VS004S : Small, VS004M : Medium
 - > VS004L : Large

Manual Vacuum Aspiration Kit (MVA Kit).....




MVA01 : Manual Vacuum Aspiration Kit (MVA Kit), Double Valve, Autoclavable Syringe (Aspirator), Non Sterile, 60 cc double valve, Autoclavable, Karman Cannulas, ETO Sterilized 4, 5, 6, 7, 8, 9, 10, 12 mm (one each), Adaptors 5 Pcs. O-Ring, Collar Stop, 2 cc Silicone Oil.

MVA02: Manual Vacuum Aspiration Kit (MVA Kit), Double Valve, Autoclavable Syringe (Aspirator), Non Sterile 60 cc double valve, Autoclavable Karman Cannulas, ETO Sterilized 6, 7, 8, 9, 10, 12 mm (one each), Adaptors 5 Pcs. O-Ring, Collar Stop, 2 cc Silicone Oil.



MVA03 : Manual Vacuum Aspiration Kit, Double Valve, Non Autoclavable Syringe (Aspirator), 60 cc double valve, Non Sterile Karman Cannulas, ETO Sterilized 4, 5, 6, 7, 8, 9, 10, 12 mm (one each), Adaptors 5 Pcs O-Ring, Collar Stop, 2 cc Silicone Oil.



MVA04 : Manual Vacuum Aspiration Kit, Double Valve, Non Autoclavable Syringe (Aspirator), 60 cc double valve, Non Sterile Karman Cannulas, ETO Sterilized 6, 7, 8, 9, 10, 12 mm (one each) Adaptors 5 Pcs O-Ring, Collar Stop, 2 cc Silicone Oil.




MVA05 : Manual Vacuum Aspiration Kit, Single Valve, Non Autoclavable Syringe (Aspirator), 60 cc single valve, Non Sterile Karman Cannulas, ETO Sterilized 4, 5, 6 mm (one each) O-Ring, Collar Stop, 2 cc Silicone Oil.



MVA06 : Dilator Set, Non Sterile Reusable - Steam Autoclavable or Chemical Sterilization. Set of five (Sizes 5/6 to 13/14 mm). Double ended. Total Length: 29 cm (11.5 inch). Latex Free

* Dilator Set, Sterile (MVA06s)



Karman Cannula without Adapter
ETO Sterilized
Individually packed in peel pouch
Compatible with all MVA syringes.

Code	MVA07-4	MVA07-5	MVA07-6	MVA07-7
Size	4	5	6	7
Code	MVA07-8	MVA07-9	MVA07-10	MVA07-12
Size	8	9	10	12