

SPHYGMO BULB




PVC001 : Sphygmo Bulb
Size 82x43mm, Matt Surface,
Capacity : 40ml

SPHYGMO BULB WITH PRESSURE CONTROL VALVE




PVC005 : Sphygmo Bulb with pressure control valve
Size : 82x43mm
Capacity : 40ml

BREAST PUMP



BREAST PUMP




Code No.	Bulb Size	Bulb Capacity in cm ² (ml)	Raw material milk receiver
PVC1000	Breast pump bulb size 3	41	Polycarbonate
PVC1003	Breast pump bulb size 3	41	Polypropylene
PVC1009	Breast pump bulb size 3	41	Glass - Hand blown
PVC1012	Breast pump bulb size 3	41	Glass - machine made
PVC1015	Breast pump bulb size 4	57	Polycarbonate
PVC1018	Breast pump bulb size 4	57	Polypropylene
PVC1021	Breast pump bulb size 4	57	Glass - hand blown
PVC1024	Breast pump bulb size 4	57	Glass - machine made

PVC EAR AND ULCER SYRINGES WITH RINGS



Code No.	Size	Bulb Capacity in cm ² (ml)	Surface
PVC1050	1 oz	25	with rings
PVC1055	2 oz	55	with rings
PVC1060	3 oz	70	with rings

PVC EAR AND ULCER SYRINGES WITH MATT SURFACE



Code No.	Size	Capacity in cm ² (ml)
PVC1070	3	27
PVC1073	5	41
PVC1076	7	75
PVC1079	9	135
PVC1081	11	206
PVC1084	13	317

PVC RECTAL SYRINGES, MATT SURFACE WITH PE TIP



Code No.	Size	Capacity in cm ² (ml)	Weight in gr.	Length Of PE TIP
PVC1080	2	25	19	60mm
PVC1083	2	25	18	46mm
PVC1086	3	40	23	60mm
PVC1089	3	40	22	46mm
PVC1091	5	89	42	60mm
PVC1093	7	143	65	60mm
PVC1096	9	224	100	60mm
PVC1099	11	347	152	70mm
PVC1102	13	483	182	70mm