

MICROPIPETTES

- » Individual finely polished stainless steel piston ensure Accuracy.
- » Light weight and ergonomic design for user comfort.
- » Quick click volume setting, separate tip ejector and clear digital read out.
- » Each pipette is individually tested and calibrated.
- » Autoclavable tip cone
- » One Year warranty.

FIXED VOLUME PIPETTE

Code No.	Volume Capacity	Accuracy %	Precision (CV%)
MPT680	1 µl	0.4	5
MPT700	2.5 µl	2.4	2.5
MPT720	5 µl	1.5	1.4
MPT740	10 µl	0.9	0.8
MPT760	20 µl	0.6	0.5
MPT780	25 µl	0.6	0.5
MPT800	50 µl	0.6	0.4
MPT820	100 µl	0.5	0.3
MPT840	200 µl	0.8	0.3
MPT860	250 µl	0.8	0.3
MPT880	500 µl	0.8	0.3
MPT900	1000 µl	0.8	0.3
MPT920	2000 µl	1.0	0.2
MPT940	3000 µl	1.0	0.2
MPT960	5000 µl	1.0	0.2
MPT980	10000 µl	1.0	0.2



VARIABLE VOLUME PIPETTE

Code No.	Volume Range	Increment	Accuracy %	Precision (CV%)
MPT440	0.5 - 10 µl	0.1 µl	3.5 - 1.0	3.0 - 0.8
MPT460	2.0 - 20 µl	0.1 µl	3.0 - 1.0	1.5 - 0.4
MPT480	5.0 - 40 µl	1.0 µl	2.0 - 0.6	2.0 - 0.3
MPT500	5.0 - 50 µl	1.0 µl	2.0 - 0.6	2.0 - 0.3
MPT520	10 - 100 µl	1.0 µl	3.0 - 0.8	1.0 - 0.2
MPT540	20 - 200 µl	1.0 µl	1.0 - 0.6	0.7 - 0.2
MPT560	40 - 200 µl	1.0 µl	1.8 - 0.6	0.7 - 0.2
MPT580	50 - 200 µl	1.0 µl	1.8 - 0.7	0.7 - 0.3
MPT600	200 - 1000 µl	10 µl	0.9 - 0.8	0.5 - 0.2
MPT620	100 - 1000 µl	10 µl	1.0 - 0.8	0.5 - 0.2
MPT640	1.0 - 5.0 ML	100 µl	1.8 - 0.5	0.3 - 0.2
MPT660	2.0 - 10 ML	100 µl	1.8 - 0.5	0.3 - 0.2