




□ HEIGHT & WEIGHT SCALES

We offer a variety of equipment and tools used for Anthropometric measurement and growth monitoring. The measurements are reliable and give a correct assessment of nutritional status of individuals. Beware of under-or-over-nutrition.

scaling
new
heights...



inside

- ▶ Personal Weighing Scales
- ▶ Baby Weighing Scales
- ▶ Stature Meter
- ▶ Weighing Trousers
- ▶ Infant Length Measuring Board
- ▶ Child / Adult Height Scales
- ▶ Stature Gauge
- ▶ Digital Scales
- ▶ X-ray Accessories

WEIGHING SCALES



WS542




WS536



WS551



WS539



WS545



WS541

Personal Weighing Scale Mechanical with shock absorbing mechanism.

AIWA GREEN : WS542
> Capacity : 130 kg. x 0.5 kg.

AIWA GREY : WS536
> Capacity : 150 kg. x 0.5 kg.


AIWA PINK : WS551
> Capacity : 150 kg. x 0.5 kg.

AIWA SPORTY : WS539
> Capacity : 150 kg. x 0.5 Kg.


AIWA CLASSIC : WS545
> Capacity : 150 kg. x 0.5 Kg.

Glass Electronic Personal Scale


SOLAR POWER : WS541
> Capacity: 150 kg. X 0.1 kg.



WS700



WS703



WS706



WS708

HEIGHT MEASURING SCALES

- > WS700 : Height Measuring Scale, Floor model without weighing scale
- > WS703 : Height Measuring Scale, Floor model with DIGITAL weighing scale
- > WS706 : Height Measuring Scale, Floor model with MECHANICAL weighing scale
- > WS708 : Height Measuring scale, Wall mounted. Measuring Range : 60-200 cm

*For WS700, WS703, WS706 Height Measuring range : 20-205cm

WEIGHING SCALES...



WS552 : DIGITAL BABY SCALE

- > Large LCD Display
- > Max. Capacity : 20 Kg., Graduation : 5 gm.
- > Tare weight function : Hold function for recording weight.
- > Automatic switch off, Battery included

BABY WEIGHING SCALES (PAN TYPE) WITH PLASTIC PAN

- > WS590 : Capacity: 10 Kg. x 50 Gms.
- > WS591 : Capacity: 16 Kg. x 50 Gms.
- > WS592 : Capacity: 20 Kg. x 50 Gms.




WS2010

WS2010 : WEIGHING SCALE, MECHANICAL COLUMN TYPE FOR ADULTS

Eye level sliding weights
Capacity: 0-160kgs with 100g graduations
Height measuring rod, range 75 to 200cm / division 5mm Large non-slip platform, Robust manufacture, corrosion protected.

WS2012 : Same as above but the capacity 200 Kg. x 100 gms. graduation.



DS660A : Mid-upper arm Circumference measuring tape, Adult
DS660C : Mid-upper arm Circumference measuring tape, Child



WS552

WS590



WS023

WS024

MEASURING TAPE

- > WS023 : Vinyl coated fibre, black marking. Length : 1.5 metre / 60 inches.
- > WS024 : Steel, graduated in cm & inches. Length : 3 metre / 118 inches.



WS045

WS045 : STATURE METER
Measuring Range : 0-200cm
Division Value : 0.1cm
Self Weight : 0.133 kg.
Minimum Order Qty. 200 pcs.

WEIGHING SCALES




WS2002 : INFANT WEIGHING SCALE - BEAM TYPE

Beam type (mechanical) scale for infants, 16 kg x 10 gm.

Technical Specifications:

- > Measuring range up to approx 16 kg
- > Minimum graduation, 10 gm.
- > Readout in kg only.
- > Easy readable in low light working situations
- > Reading time max 15 seconds
- > Includes removable tray for infant
- > With easy zero adjustment possibility
- > Design allows rough handling
- > Smooth surface/finishing allows for easy cleaning/disinfection.
- > All vital parts made of rust proof materials
- > Easy maintenance and repair in low tech settings
- > Splash proof and shock resistant light-weight body
- > Materials; body, metal (white baked enamel) or plastic
- > Sturdy design allows heavy duty use.

WS2004 : Same as WS2002 but capacity 20kg x 10 gm.



WS2004



WS2008

WS2008 : Weighing Scale, Nylon body, color coded with Hammock
Capacity: 5kg x 100 gms

WS2009 : Hammock only



WS032 / WS037

Popular

Baby Weighing Scale, Salter Type (Dial)

- > With two stainless steel suspension hooks & adjustment screw.

Capacity/ Material	25 kg x 100 gm	50 kg x 200 gm
ABS Plastic Body	WS032	
Metal Body	WS037	WS039

Optional : Weighing trouser - WS035 (Order separately).



WS035

WS040

Weighing Trouser for Baby Weighing Scales

With strap for suspension from infant scale.

- > WS035 : One trouser only (without bag).
- > WS038 : Set of 5 weighing trousers with carrying bag.
- > WS040 : Carry bag for weighing trousers.



WS030 : Baby Weighing Balance / Scale (Spring type)

- > Capacity : 5kg x 100gms.
- > Color coded with zero adjustment device.
- > Compact to fit in midwifery kits.
- > Light weight - Nylon.
- > With hook, suspension ring & rustproof firm grip handle.

Optional : Weighing trouser - WS035 (Order separately).
or WS2009 (Order separately).

Weighing Scale, Spring Type (Multi purpose)

- > Plastic, transparent.
- > With zero adjustment for both force & mass measurement.
- > Calibration in newton & grams.

Code	Capacity
WS060	100Gm / 1 Newton
WS063	250Gm / 2.5 Newton
WS066	500Gm / 5 Newton
WS069	1000Gm / 10 Newton
WS072	2000Gm / 20 Newton

Weighing Scale, Spring Type (Multi purpose)

- > Aluminium (multi purpose).
- > With zero adjustment for both force & mass measurement.
- > Calibration in newton & grams.

Code	Capacity
WS075	100Gm / 1 Newton
WS078	250Gm / 2.5 Newton
WS081	500Gm / 5 Newton
WS084	1000Gm / 10 Newton
WS087	2000Gm / 20 Newton
WS090	3000Gm / 30 Newton
WS093	5000Gm / 50 Newton



Height Measuring Scale (Infant) / Infantometer

Sleek broad acrylic base with one sliding side adjustable as per length of baby, dual scale for direct reading in cms / inches upto 90cm, folding sides for easy storage, baby can comfortably lay on the chest of the scale.

- > WS021 : Size 18" x 5" (Length x Width).
- > WS025 : Size 18" x 7" (Length x Width).



- > WS026 : Height Measuring Rod for babies
Measuring Range : 0-96cm

DIGITAL SCALES

Note : Minimum order quantity : 120 pcs. assorted sizes/types

WS2041 - Bathroom Scale electronics flat with tare parent model

- > Maximum Capacity: 300 kg / 660 lbs.
- > Graduation: 0.1 kg / 0.2 lb.
- > Accuracy: $\pm 150g$
- > Resolution: 1:3000
- > Tare Range: Full Range
- > Repeatability(STD dev): $\pm 1d$
- > Linearity: $\pm 3d$
- > Unit of Measure: Kg/lb
- > Key Function: Hold/Unit/Weight, On/Zero/Off, Tare
- > Stabilisation Time: 1-2Sec
- > Operating Temperature: 5°C-35°C or 13° F -95°F
- > Power Supply: 9V battery or AC adaptor
- > Calibration: Automatic External
- > Calibration Mass: 300KG
- > Indicator Display: 1.0" LCD digital display
- > Dimension(w x d x h):
15.5 x 12.5 x 3.3 inch
39.5 x 32 x 8.5 cm
- > Gross Weight : 7.7lb/3.5 Kg




WS2042 - Weighing Scale, Standing, Adult, Digital with Height Rod

- > Capacity: 300 Kg X 0.1Kg
- > Accuracy: $\pm 150g$
- > Weight Unit: Kg/lb
- > LCD Display: 1.0" LCD digital Display
- > Dimension: Platform Size: 310 x 310 x 54mm,
Column Height: 900mm
- > Key Function: Hold/BMI, Unit/Weight, On/Zero/Off, Tare
- > Power Supply: 1.5V 300mA adapter
1.5V AAA battery *6
- > Operation Temperature: 5°C / 35°C
- > Height Rod: 60-205cm / 6ft 9inch



WS2034 : Digital Baby Scale

- > 20kg/grad.: 5g,
- > Use 4 x AA Batteries(Not Included)
- > Use 4 Load Cells.
- > Unit: kg & g
- > Product Dim.: 60.3x38.5x10cm



WS2035 : Stature Gauges

- > Measuring Range : 20-99cm
 - > Value per division : 0.5cm
- Minimum Order Qty. 200 pcs.



WS2019 : Digital Personal/Bathroom Scale

- > 150 Kg. x 100 gms.
- > Overload & Low Voltage Indicator
- > Auto Zero/Auto Off
- > Power: 2x3V CR2032 Lithium Battery



WS2031 : Wheelchair Scale, Electronic


- > Max weight: 200kg
- > Division: 100g
- > Display: LCD/LED
- > G.W.: 35kg N.W.: 25kg
- > Measurement: 77×64×114cm
- > 53pcs/20'Container

Note : Minimum order quantity : 120 pcs. assorted sizes/types



WS2028 : Wheelchair Scale, Mechanical

- > Max weight: 160/200kg
- > Division: 100g
- > Package size: 77×64×114cm
- > G.W.: 35kg N.W.: 25kg
- > 53pcs/20'Container
- > 107pcs/40' Container



WS2040 : Mechanical Personal/Bathroom Scale

- > 120kg/grad.:500g, 260lb/1 lb
- > Unit: kg or kg&lb
- > Product Dim.:27x26.5x5.5cm




WS2043 : Mechanical Personal/Bathroom Scale

- > 120kg/grad.:500g, 260lb/1 lb
- > The Cartoon Can Be Changed
- > Unit: kg or kg&lb
- > Product Dim.:26.3x25.9x4.9cm



WS2055 : Mechanical Personal/Bathroom Scale

- > 120kg/grad.:1kg,260lb/grad:2lb
- > The Color Can Be Changed
- > Unit: kg or kg&lb
- > Product Dim.:29x29.8x5.7cm



WS2061 : Mechanical Personal/Bathroom Scale

- > 120kg/grad.:1kg,260lb/grad:2lb
- > The Color Can Be Changed
- > Unit: kg or kg&lb
- > Product Dim.:28.8x28.8x5.5cm