


X-RAY ACCESSORIES



X-Ray Viewer, florescent tubes, Plastic Body

- > DP301 : X-Ray Viewer Single Film
- > DP311 : X-Ray Viewer Double Film
- > DP313 : X-Ray Viewer Triple Film
- > DP321 : X-Ray Viewer Four Film



X-Ray Viewer, **LED**

- > DP302 : X-Ray Viewer Single Film
- > DP312 : X-Ray Viewer Double Film
- > DP318 : X-Ray Viewer Triple Film
- > DP320 : X-Ray Viewer Four Film



DP326 : LEAD GLOVES

Lead Gloves have seamless surface that does not allow any radiation to leak in. Convenient to wear and discard.



DP322 : LEAD APRON - WRAP AROUND

- > Lead Apron is specially designed for surgeons.
- > It gives protection to the whole body.
- > Provides very less fatigue and freedom of movement.


DP323 : LEAD APRON-VELCRO TYPE

- > It is designed to provide efficacious shield against radiation during X-ray process.
- > Provide very less fatigue and freedom of movement



DP325 : THYROID COLLARS


- > Thyroid collars provide protection for the neck and thyroid area.
- > It helps to reduce patient exposure during intraoral radiography.



DP330 : GONAD SHIELD

DP332 : OVARY SHIELD


- > Ideal to protect the Male/Female reproductive organ
- > Available in one-size-fits-all and child sizes
- > Designed for ease of wear and removal



DP327 : X RAY CASSETTES

Made of specified aluminum alloy frame, specified engineering plastics and high-stamina aluminum plate.

DP327a : 6.5"X8.5" DP327d : 14"X17"
 DP327b : 8"X10" DP327e : 12"X15"
 DP327c : 10"X12" DP327f : 14"X14"



DP328 : X RAY HANGERS

Made of Stainless Steel

DP328a : 6.5"X8.5" DP328d : 14"X17"
 DP328b : 8"X10" DP328e : 12"X15"
 DP328c : 10"X12" DP328f : 14"X14"



DP340 : APRON STAND FOR 5 APRONS

Made of Stainless Steel, Floor Type with wheels

Note: Apron Stand also available for 3, 7 Aprons.
 Wall mounted type Apron Stand also available



DP341 : CHEST STAND FLOOR MODEL

Made of Steel

Note: Wall mounted Chest Stand also available.



DP324

DP324 : X-RAY PROTECTIVE GLASSES

> The glasses are available for protecting the eyes of the practitioner against X-ray during the Radio Graphic Diagnosis and Interventional Treatment.



DP342 : SAFE LIGHT FOR DARK ROOM

*other X-ray accessories also available