


☐ MISCELLANEOUS PRODUCTS.....



Needle & Syringe Destroyer

NS030 : Instant needle burner & syringe destroyer. This model has an electrical instant needle burner and replacable cutter to cut syringe nozzles after use. Size 165 x 115 x 120mm (Handy desktop model) weight 1.7 kg power 230V AC +- 10% 50 hz., single phase (alternative voltages available) protection 2 amperes 3/4

> MZ011: Dropping Bottle, 10ml, Amber Glass.

NS031 : Spare Cutter for Syringe destroyer.



Kerosene Pressure Stove, Iron Made, Side Tank

- > KS020 : Single burner, capacity : 1.5 ltr.
- > KS029 : Single burner, capacity : 6 ltr.

> KS080 : Hurricane Lantern.

> KS094 : Petromax Pressure type, Full Brass Tank.

Kerosene Pressure Stove, Brass Tank

- > KS031 : Single burner, capacity : 1.5 ltr.
- > KS039 : Single burner, capacity : 2 ltr.



MZ015 : Timer, Mechanical, Spring Driven, with clear and sustained ring, dial easy to read, graduated 0-60 minutes in one minute interval, sturdy construction, shock resistant.

> MZ017 : Stop Watch, Mechanical

> MZ018 : Stop Watch, Digital

MZ016 : Timer, Digital,
Timing Function : 99 minutes and 59 seconds

MISCELLANEOUS PRODUCTS.....




MZ031 : Examination Lamp/Light Floor Stand :

- > Approx. 100 cm high, floor stand type, plastic base with castors.
- > Arm: approx. 75 cm articulated, spring loaded arm.
- > Voltage: 110/220v, bulb: max. 100w incandescent bulb, Pin Type.
- > Shade : Approx. 15cm dia, light opening,
- > Color: Off-white. Free cord: approx. 2 yard.

MZ030 : Specifications same as Mz031 but, with cast iron base.

MZ033 : Halogen Examination Light Adjustable height, with step down Adopter, Halogen Bulb 12V, 50W, 10° angle of emission with UV block front glass.

MZ034 : Solar Examination Light Bulb LED 3W, Goose Neck Battery backup : 4 hours Available with Solar Panel Multifacet Reflector Technology




> **MZ100 :** Safety Pin, Medium 40mm for securing wraps and bandages Pack of 12




> **MZ050 :** Double Edge Razor Blade, Pack of 5 blades

NOTE: Reusable Metal Razors/Disposable Plastic Razors available on request



Blanket Survival (Warming Foil) Both Side Silvered

- > MZ003c : Size 1200mm x 900mm (Child)
- > MZ003a : Size 2200mm x 1400mm (Adult)



- > **MZ051 :** Razor Single Edge, Plastic, Single Use.
- > **MZ052 :** Razor Double Edge, Plastic, Single Use.

□ MISCELLANEOUS PRODUCTS.....




Universal pH Indicator Paper


- > MZ105a : pH 0-6
- > MZ105b : pH 0-14




- > MZ090 : Pen Skin Marker with Flexible Ruler, Sterile Single use



- > MZ110 : Plastic Cotton Tip Buds Disposable, Length 7.5mm, Tip Size 2-4mm, Non Sterile



- > MZ115 : Penlight Torch with Pupil Gauge



- > MZ010 : Nail Brushes, Autoclavable.

MZ012 : Sterile Dry Surgical Hand Brush

- > Made of medical grade polyethylene and sponge.
- > Ergonomic brush design (48*78*35 mm).
- > Easy open package.
- > Soft bristles for no skin irritation.
- > Smooth and soft sponge absorbs up to 90 ml liquid.
- > For external use only.
- > Dry Sterile ETO
- > Shelf life: 3 yrs



Minimum Order Qty.: 5000 pcs.