



## □ HOSPITAL HOLLOWWARE

A complete range of durable, corrosion resistant Stainless Steel utensils necessary for all hospitals and health structures.

# ... Utilitarian Utensils



## inside

- ▶ Kidney Trays
- ▶ Instrument Trays
- ▶ Bed Pans
- ▶ Urinals
- ▶ Douche Cans
- ▶ Lotion Bowls
- ▶ Wash Basins
- ▶ Forceps Jars
- ▶ Gallipots & Sputum Mugs
- ▶ Dressing Jars
- ▶ Bucket / Pail
- ▶ Waste Receptacle
- ▶ Surgical Boxes
- ▶ Plastic Holloware

## KIDNEY TRAYS



Kidney Trays ( Stainless steel) without cover					
Approx. Size		Quality / Code			
In mm	In inch	Export Quality Thickness 0.35mm	Deluxe Quality Thickness 0.5mm	Superior Quality Grade 304, 24g Thickness 0.5mm	Super - Deluxe Qly. Grade 304, 24g Thickness 0.63mm
150	6	HH001	HH005	HH010	HH015
200	8	HH030	HH035	HH040	HH045
250	10	HH060	HH065	HH070	HH075
300	12	HH090	HH095	HH100	HH105
250 - Shallow type (25x15x4cm)	10	-	HH142	HH145	HH140

*Note:* Kidney Tray with cover also available upon request. For code : Add 'D' to above codes.

## INSTRUMENT TRAYS



Instrument Tray with Cover (Stainless steel) - (Inner dia.)					
Approx. Size		Quality / Code			
In mm	In inch	Regular 0.35mm	Superior 0.4mm Tray	Deluxe - 0.5mm 0.5mm Tray	18/8 Gauge 0.5mm Tray 304 Grade
200x79x40	8x3x1 <sup>1/2</sup>	<sup>1</sup> HH200	HH202	HH204	HH206
201x150x48	8x6x2	HH210	HH212	HH214	HH216
227x152x48	9x6x2	HH220	HH222	HH224	HH226
254x200x48	10x8x2	HH230	<sup>2</sup> HH232	HH234	HH236
271x177x48	11x7x2	HH240	HH242	HH244	HH246
305x203x48	12x8x2	HH250	HH252	HH254	HH256
305x250x48	12x10x2	HH260	HH262	HH264	HH266
352x251x48	14x10x2	HH270	HH272	HH274	HH276
376x300x48	15x12x2	HH280	HH282	HH284	HH286
450x300x50	18x12x2	HH290	HH292	HH294	HH296
225x125x50	9x5x2	-	-	-	HH298*
310x195x63	12 <sup>1/2</sup> x7 <sup>3/4</sup> x2 <sup>1/2</sup>	-	-	-	HH299*


*Note:* Instrument Tray without cover also available upon request. For code : Add 'C' to above codes. \* These are 0.6mm

CATHETER TRAY




Catheter Tray (Stainless Steel) with cover				
Approx. Size		Quality / Code		
In mm	In inch	Export Quality, Thickness 0.35mm	Deluxe Quality, Thickness 0.5mm	Superior Quality, Thickness 0.5mm, 304 Grade
205 x 80 x 36	8 x 3	HH690	-	-
205 x 85 x 45	8 x 3.25	-	HH692	HH695
222 x 82 x 41	8.75 x 3.25	-	HH693	HH696
405 x 110 x 48	16 x 4	HH691	HH694	HH697

SHALLOW TRAY



Shallow Tray (Stainless Steel)			
Approx. Size		Quality / Code	
In Inch	In mm	Deluxe Quality, Thickness 0.5mm	Superior Quality, Thickness 0.5mm, 304 Grade
8 x 6	200 x 150 x 25	HH710	HH711
9 x 6	230 x 155 x 25	HH712	HH713
10 x 8	255 x 200 x 25	HH714	HH715
11 x 7	280 x 180 x 25	HH716	HH717
12 x 8	300 x 205 x 25	HH718	HH719
12 x 10	305 x 255 x 25	HH720	HH721
14 x 10	355 x 255 x 25	HH722	HH723
15 x 12	380 x 305 x 25	HH724	HH725
18 x 12	455 x 305 x 25	HH726	HH727


MAYO TABLE TRAY




Mayo Table Tray (Stainless Steel)		
Approx. Size		Code
In mm	In inch	Thickness 0.5mm, 304 Grade
480 x 330 x 22	18 x 12	HH776
610 x 422 x 30	23 x 16	HH775

## BED PAN WITH LID (Stainless Steel)

www.narang.tv/2225



HH305



HH355



HH351



HH307

Bed Pan with lid, stainless steel		
Size in mm	Female (Perfection type)	Male (Slipper type, welded)
355 x 290 x 70	HH302 (Jointed)	HH353
375 x 310 x 70	HH307 (Jointed)	HH351
350 x 285 x 97	HH305 (Seamless)	HH355 (Round Type)

## URINAL / URINE POTS



HH405



HH424

Urinal / Urine Pots (Stainless steel)			
Size ml	Male	Size ml	Female
800	HH405	800	HH420
1030	HH407	920	HH422
-	-	1010	HH424

## Medicine Measure graduated



HH689

## DOUCHE CAN (Irrigators)




HH555

Douche Cans (Irrigators), Stainless Steel	
Code	Details
HH535	Capacity: 1 Litre (Approx.)
HH545	Capacity: 1.5 Litre (Approx.)
HH555	Capacity: 2 Litre (Approx.)

Douche Cans with Douche Tube & Douche Fitting Set	
Code	Details
DC050	Capacity: 1 Litre (Approx.)
DC054	Capacity: 1.5 Litre (Approx.)
DC058	Capacity: 2 Litre (Approx.)

## Measuring Spoon Set (Stainless Steel)



HH503



LOTION BOWLS - Stainless Steel

Lotion Bowls		Quality/Code	
Approx. Dia In mm	Approx. Dia In inch	Regular	Deluxe
100	4	HH430	HH435
125	5	HH440	HH445
150	6	HH450	HH455
180	7	HH460	HH465
200	8	HH470	HH475
250	10	HH480	HH485



LOTION BOWLS (Unicef Pattern), Half curved, open rim  
Stainless steel

Lotion Bowls (Unicef pattern)		
Approx. Size (Dia. X Depth)	Capacity	Code
80 x 40 mm	180 ml	HH476
125 x 57 mm	600 ml	HH477
160 x 75 mm	1200 ml	HH478



BASINS

Basins (Stainless Steel)		Quality / Code	
Dia. (Approx.)		Regular Thickness 0.35mm	Deluxe - Deep Grade 304 Thickness 0.5mm
250 mm	10 inch	HH488	HH489
300 mm	12 inch	HH490	HH495
350 mm	14 inch	HH500	HH505
400 mm	16 inch	HH510	HH515
Code - HH517 : Basin Solution, Deep (Unicef pattern), 330 Dia., 155mm Height, 6 Litre			

Note : For more models in Bowls, visit our website at [www.narang.com](http://www.narang.com)



JARS & FEEDING CUP - Stainless Steel

Jars & Feeding Cup (Stainless Steel)	Deluxe Quality/Code
Thermometer Jar, Sz. 30 x 90mm (1 1/4 x 3 1/2")	HH620
Forceps Jar, Light Base, Sz. 50 x 185mm (2 x 7 1/4")	HH625
Forceps Jar, Heavy Base, Sz. 50 x 185mm (2 x 7 1/4")	HH630
Forceps Jar, Heavy Base, Sz. 75 x 200mm (3 x 8")	HH640
Forceps Jar, Heavy Base, Sz. 100 x 200mm (4 x 8")	HH650
Feeding Cup	HH577

# □ HOSPITAL HOLLOWWARE ..... STAINLESS STEEL



## GALLIPOTS & SPUTUM MUG



Galipots Cup Shaped (Stainless steel) with cover			
Capacity (Approx.)	Diameter x Height (Approx.)	Capacity	Quality / Code
In oz.	In mm	In ml	Deluxe - 18/8G 0.5mm
2 oz.	50 x 30	60	HH559
4 oz.	60 x 40	112	HH560
6 oz.	75 x 45	175	HH565
8 oz.	80 x 50	225	HH568

Sputum Mug with cover	
Item Code	Quality
HH570	Regular, 450ml, 100x50mm
HH575	Deluxe, 450ml, 100x50mm

## DRESSING JAR



## PAIL (BUCKET)



Dressing Jar		
Item Code	Size (mm.)	Size (inches)
HH596	100 x 75 (DxH)	4 x 3 (DxH)
HH597	150 x 150 (DxH)	6 x 6 (DxH)
HH598	225 x 200 (DxH)	9 x 8 (DxH)

Pail (Bucket)		
Capacity	With Cover	Without Cover
11 L	HH611	HH612
12 L	HH615	HH617
14 L	HH613	HH619
18 L	HH621	HH622

Note : Buckets with graduation available on request.

## SURGICAL / SYRINGE BOXES




## MEASURING JUG



Syringe Boxes, Aluminium	
Item Code	Size (in mm)
HH765	80x25x20
HH769	160x50x25
HH772	120x80x25

Surgical Boxes, Aluminium	
Item Code	Size (in mm)
HH780	135x60x30
HH785	160x70x34
HH790	175x115x30
HH795	185x80x40

Measuring Jug	
Item Code	Capacity
HH680	500 ml.
HH683	1000 ml.
HH687	1500 ml.



**KIDNEY TRAY, Polypropylene, Autoclavable**

- HH907 : Size 6"
- HH900 : Size 8"
- HH903 : Size 10"
- HH902 : Size 10" With Handle
- HH906 : Size 12"



**INSTRUMENT TRAY WITH COVER, Polypropylene**

- P194.221570 : Size 220x150x70mm
- P194.451570 : Size 450x150x70mm



**HH931 : Polypropylene Tray with compartment**  
Size: 216L x 175W x 42H mm

- HH932a : Polypropylene Gallipot, Dia 90mm, Cap. 280ml
- HH932b : Polypropylene Gallipot, Dia 65mm, Cap. 90ml



**URINALS, Graduated with Lid, made of white HDPE**

- HH925 : Male.
- HH926 : Female.



**URINALS, Graduated with Lid, made of clear PP, Autoclavable**

- HH927 : Male.
- HH928 : Dual use for male & female with easy fit female attachment



**BED PAN for Adult, Polypropylene, Autoclavable**  
400x260x90mm

HH920

**BED PAN for Adult, Polypropylene, Sleek**  
400x260x65mm

HH923 :



**THERMOMETER JAR PP (without cover)**  
 HH965 : Made of Polypropylene, Autoclavable.  
 Internal Dia. : 30 mm.  
 Depth: 110 mm. Overall Height: 120 mm.

**FORCEP JAR PP (with cover)**  
 HH966 : Made of Polypropylene, Autoclavable.  
 Internal Dia. : Approx. 54 mm.  
 Depth: 180 mm. Overall Height: 200 mm.



**Douche Can Set**  
 DC060 : Plastic can 1500 ml with side outlet,  
 plastic tube & douche fitting set.


HH982 : Bowl, Round Plastic, 6 litre



**JUG, Measuring, Polypropylene, Graduated, Autoclavable**  
 P172.0250 : 250 ml.    P172.0500 : 500 ml.  
 P172.1000 : 1000 ml.    P172.2000 : 2000 ml.  
 P172.1015 : 1500 ml.

**PLASTIC BOTTLE, Plastic with screw cap**  
 HH970 : Round form, narrow neck, capacity 1 litre.

**WASHING BOTTLE, Plastic, squeeze type**  
 HH971 : With bent delivery tube and vented screw cap,  
 Capacity: 250ml.



**SPITTING MUG, PLASTIC**  
 HH969 : Colour Available: Blue & White

**FEEDING CUP**  
 HH967 : Adult Polypropylene Cup, Graduated,  
 200 ml., Lid with wide mouthed spout. Autoclavable