



Kiran TECHNO SERVICES
PRIVATE LIMITED



www.kiransurgicalsystems.com
www.kiranmedicareequipments.com

Electric Bed [D4S 111]

Overall Size : 2300mm X 980 mm X 475mm

Height adjustment : 475 mm to 745 mm

- Backrest movement 80 DEGREE
- Trendelenburg movement 12 DEGREE
- Collapsible railings
- Buffers at four corners
- All four movements operated by wired remote
- 150 mm dia, 4 heavy duty wheels, 2 with stopper
- **OPTIONAL: C.P.R. FACILITY, nurse control panel AND BATTERY BACKUP**



NEW CONCEPT FOR OUR MOST VALUABLE CUSTOMERS

DESIGN YOUR OWN ELECTRIC BED

STEP 1

SELECTION OF BED SURFACE



PERFORATED SHEET TOP



TRAPEZOID TYPE TOP



ABS PLASTIC TOP

STEP 2

SELECTION OF HEAD AND FOOT BOARD



ABS PLASTIC BOARDS



S.S. METAL BOWS



M.S. POWDER COATED BOWS

STEP 3

SELECTION OF CASTORS



DELUX CASTOR



STANDARD CASTOR

STEP 4

SELECTION OF RAILINGS



ABS PLASTIC SWING TYPE RAILING



S.S. SWING TYPE RAILING



S.S. COLLAPSIBLE RAILING



ALLUMINIUM COLLAPSIBLE RAILING



**M. S. POWDER COATED
COLLAPSIBLE RAILING**

**SELECT THE REQUIRED PARTS OF THE BED AND FORWARD US BY EMAIL, FAX
AND WE WILL QUOTE YOU THE PRICES ASAP.**

Electric Homecare Beds

SPECIAL FEATURES :

Complete Wooden Body

- Made from Burma Teak Wood one of the highest quality of wood.
- Look like a part of your designed furniture and don't make your home look like a hospital.
- Controlled with remote control.
- Three different comfortable positions, backrest position, leg rest position and back and leg rest position.
- Size??? To suit your Bedroom.



Backrest Position



Legrest Position



Backrest + Legrest Position

HOSPITAL BEDS MECHANICALLY OPERATED



Plain Hospital Bed (DELUX) [D4S 01]

Overall Size : L 2010 mm x W 865 mm x H 590 mm
Frame : 50 mm x 25 mm, CRCA Rectangular tube

- Perforated CRCA sheet top
- S.S. bows of 32 mm dia.
- 4 flush I.V. rod locations
- Pre-treated and epoxy powder coated



Plain Hospital Bed (STANDARD) [D4S 01A]

Overall Size : L 2010 mm x W 865 mm x H 590 mm
Frame : 50 mm x 25 mm, CRCA Rectangular tube

- Perforated CRCA sheet top
- M.S. bows of 32 mm dia.
- 4 flush I.V. rod locations
- Pre-treated and epoxy powder coated



Manual Back Rest Bed [D4S 02]

Overall Size : L 2010 mm x W 865 mm x H 590 mm
Frame : 50 mm x 25 mm, CRCA Rectangular tube

- Two section perforated CRCA sheet top
- Manually operated, 5 position back - rest on ratchet
- S.S. bows of 32 mm dia.
- 4 flush I.V. rod locations
- Pre-treated and epoxy powder coated



Semi-Fowler Bed [D4S 03]

Overall Size : L 2010 mm x W 865 mm x H 590 mm
Frame : 50 mm x 25 mm, CRCA Rectangular tube

- Two section perforated CRCA sheet top
- Back rest operated by screw with folding handles.
- Detachable S.S. bows of 32 mm dia
- 4 flush I.V. rod locations
- Pre-treated and epoxy powder coated



Semi-Fowler Bed [D4S 03A]

Overall Size : L 2010 mm x W 865 mm x H 590 mm
Frame : 50 mm x 25 mm, CRCA Rectangular tube

- Two section perforated CRCA sheet top
- Back rest operated by screw with folding handles
- Detachable S.S. bows of 32 mm dia. With metal panel (Different colour options available)
- 4 flush I.V. rod locations
- Pre-treated and epoxy powder coated

GAT - Global Access System
 I.V. - In-line Venous Connections

Fowler Bed [D4S 04]

Overall Size : L 2060mm x W 865 mm x H 590 mm
Frame : 50 mm x 25 mm, CRCA Rectangular tube

- Four section perforated CRCA sheet top
- Back rest and knee rest operated by attached folding handles
- Detachable s.s. bows with metal panels (Different colour options available)
- S.S. swing type railings: Optional
- 125 mm dia heavy duty wheels; 2 with brakes
- 4 flush I.V. rod locations
- Pre-treated and epoxy powder coated



GAR - Gold Finish Epoxy
EAC - Cold Zinc Coating

Recovery Bed [D4S 05]

Overall Size : L 2010 mm x W 865 mm x H 590 mm

Frame : 50 mm x 25 mm, CRCA Rectangular tube

- Four section perforated CRCA sheet top.
- 3 Separate screws for back rest, knee rest, trendelenburg & reverse trendelenburg position by attached folding handles.
- Detachable s.s. bows with metal panel (different colour options available).
- S.S. side railings on both sides (Swing type or collapsible type) and buffer on all corners.
- 125 mm dia heavy duty wheels; 2 with brakes.
- 4 flush LV. rod locations.
- Pre-treated and epoxy powder coated.



GAR - Global Access System
EAS - Endless Access System

ICU Bed [D4S 06]

Overall Size : L 2060mm x W 905 mm x H 590 mm

Frame : 50 mm x 25 mm, CRCA Rectangular tube

- Four section perforated CRCA sheet top
- 4 Separate screws for back rest, knee rest trendelenburg & reverse trendelenburg and height adjustment position by attached folding handles. (ALSO AVAILABLE WITH FIXED HEIGHT OPTION MODEL D4S 006-A)
- S.S. side railings on both sides (Swing type or collapsible type) and buffer on all corners
- Detachable s.a. bows with metal panel (different colour options available) . .
- 125 mm dia heavy duty wheels; 2 with brakes
- 4 flush I.V. rod locations
- Pre-treated and epoxy powder coated





**PLAIN HOSPITAL BED WITH CRADLE
[D4S 001B]**

Overall Size : L 2010 mm x W 865 mm x H 590 mm
Frame : 50 mm x 25 mm, CRCA Rectangular tube

- Perforated CRCA sheet top
- S.S. bows of 32 mm dia.
- 4 flush I.V. rod locations
- Pre-treated and epoxy powder coated

CRIB SPECIFICATIONS

SIZE (36" x 15")

M.S Tubular pipe constructi

Pre-treated and Epoxy Powder coated



PAEDIATRIC BED [D4S 000]

Size : 137cm x 76cm x 60cm.

- Complete with full-length dropside safety railings
- Head and leg bows of equal height
- Pre-treated and epoxy powder coated

Operation Table (Height Adjustable) [D4S 20]

Overall Size : L 1820mm x W 560 mm x Ht. Adjustable from 750mm to 920 mm

- Stainless steel top with kidney bridge
- Detachable and folding foot end
- Height Adjustment and head low/up position by screw and lever mechanism
- Pre-treated and epoxy powder coated
- Available in stainless steel top & epoxy powder coated base / S.S. base



Labour Table (Foot End Folding) [D4S 22]

Overall Size : L 1830mm x W 870mm x Ht. 800mm

- Two section Stainless steel top with head low/up position with screw mechanism
- Detachable and folding foot end
- 'U' notch with drainage tray
- Available in full stainless steel / stainless steel top & epoxy powder coated base

Labour Table (Telescopic) [D4S 23]

Overall Size : L 1830mm x W 870mm x H 800mm

Frame : 50 mm x 25 mm Rectangular Stainless steel tube

- Three section Stainless steel top
- Back rest on ratchet; head low/up by screw mechanism
- Middle Section with 'U' notch with bucket / bowl
- Sliding foot end
- Railings on three sides
- Available in full stainless steel / stainless steel top & epoxy powder coated base
- Lithotomy Bars Gutter type / Leather Strip type



Examination cum Gynec Table [D4S 31]

Overall Size : L 1820 mm x W 560 mm x H 815 mm

- CRCA sheet / S.S. top
- 'U' notch with foot end folding and detachable
- Lithotomy bars
- Pre-treated and epoxy powder coated / Full S.S.

Patient Examination Couch [D4S 32]

Overall Size : L 1830mm x W 575 mm x H 850 mm
Frame : CRCA Sheets and Rectangular tubes

- 4 Cabinets; 1 with lock, drawer with lock
- 100MM THICK FOAM Mattress
- Sliding Stepping Stool and Writing Pad
- Pre-treated and epoxy powder coated



Gynec Examination Couch [D4S 33]

Overall Size : L 1830mm x W 575 mm x H 850 mm
Frame : CRCA Sheets and Rectangular tubes

- 3 Cabinets; 1 with lock, drawer with lock.
- Foam Mattress with 'U' notch 100 mm thick
- Examination Trolley with three trays
- Sliding Stepping Stool and Writing Pad
- Pre-treated and epoxy powder coated

TROLLEYS

Stretcher Trolley (Folding) [D4S 16]

Overall Size : L 1830mm - 990mm x W 570 mm x H 860 mm

Frame : 50 mm x 25 mm, CRCA Rectangular tube

- Folding Stainless steel top fixed on CRCA tube frame mounted on 150 mm dia. heavy duty wheels
- Provision for fixing I.V. rod
- Pre-treated and epoxy powder coated



Stretcher-On-Trolley [D4S 17]

Overall Size : L 1830 mm x W 530 mm x H 870 mm

Frame : 50 mm x 25 mm, CRCA tube

- Detachable Stainless steel top fixed on CRCA tube frame mounted on 150 mm dia. heavy duty wheels
- Provision for fixing I.V. rod
- Pre-treated and epoxy powder coated



Casualty Trolley [D4S 18]

Overall Size : L 1830 mm x W 630 mm x Ht. adjustable from 640 mm to 810mm
Frame : 50 mm x 25 mm, CRCA tube

- Height adjustable trolley with X-Ray permeable detachable top and full length X-Ray cassette holder
- Height adjustment, trendelenburg & reverse trendelenburg positions by screw and lever mechanism
- 150 mm dia. heavy duty wheels, 2 with brakes
- Provision for fixing I.V. rod
- Swing away type S.S. railings on both side, rubber buffers on all corners
- Pre-treated and epoxy powder coated



Emergency Trolley [D4S 26]

Overall Size : W 720mm x D 650mm x H 1600mm

- 6 Drawer Modular system upper two drawers with medicine containers; pre-treated and epoxy powder coated
- 6 Detachable plastic bins
- 2 cages for oxygen and nitrous cylinders
- Top tray for monitor and pulse oxymeter
- System mounted on 100 mm dia. wheels, with pair of front brakes
- Available in Full S. S. also



Baby Trolley [D4S 36]

Overall Size : W 390 mm x D 800 mm x H 1080 mm

- Stainless steel tray in head low position
- 3 drawer modular system
- Top with medicine containers
- Adjustable I.V. rod & lamp arrangement
- Bottom tray for mounting suction apparatus
- Trolley mounted on 100 mm. dia. wheels
- Pre-treated and epoxy powder coated



Oxygen Trolley [D4S 34]

A) Pull Type with three wheels of 50 mm dia.

B) Push Type with two wheels of 100 mm dia

G.A.P. - Dental Division
U.S.A. - United States

Food Trolley [D4S 14]

Overall Size L 760 mm x W 360 mm x Ht. adjustable from 750 mm to 1050 mm

- Stainless steel top
- Height adjustable by lever gear mechanism
- Mounted on 50 mm dia wheels
- Pre-treated and epoxy powder coated



Monitor Trolley (for scopic surgery) [D4S 29]

Overall size : W 660 mm x D 560 mm x H 1498 mm

- Available in full Stainless steel & epoxy powder coated
- Adjustable Shelves
- 125 mm dia. heavy duty castors; 2 with brakes





Kick Bucket and Bowl Stand [D4S 27]

- Stainless steel bowl dia - 345 mm.
- 2 Stainless steel bowls of dia - 345mm each

Stand available in stainless steel / epoxy powder coated



Instrument Trolley [D4S 24]

Overall Size : L 710mm x W 460mm x Ht. 830mm
Frame Stainless steel square tubes

- Stainless steel bottom and top shelves with railing on top shelf
- Trolley Mounted on 100 mm dia. wheels
- Available in full stainless steel / Powder coated frame with stainless steel shelves and railings



Dressing Trolley with Bowl and Bucket [D4S 35]

Overall Size : L 1080mm x B 460 mm x H 830 mm
Frame Stainless steel square tubes

- 2 stainless steel shelves with railing on top
- Stainless steel bowl dia 345 mm.
- Bucket dia. 300 mm.
- Available in full stainless steel/ Powder coated frames with stainless steel shelves and railings



Mayo's Trolley [D4S 25]

Tray size : L 555mm x W 430mm x Ht. adjustable from 860mm to 1300 mm
Frame Stainless steel rectangular and square tubes

- Trolley Mounted on 50 mm dia. wheels
- Available in full stainless steel / Powder coated frame with stainless steel tray



Multipurpose Bed Side Lockers / Cabinets [D4S 214]

Overall Size : W 410mm X D 370mm X H 900mm

- Complete ABS (Acrylonitrile Butadiene Styrene) plastic
- Different colors to select (only for bulk quantity)
- OEM supplies available by incorporating your logo (only for bulk quantity)
- Can be used in Hospitals / Hotels / Offices / Laboratories etc.



Bed Side Locker (Delux) [D4S 11]

Overall Size : W 395mm x D 305 mm x H 800 mm

- Crca sheet main cabinet and drawer
- Mounted on 50 mm dia twin roller wheels front wheels with break
- S.S. top with three sides raised and front turned down
- Pre-treated and epoxy powder coated



Multipurpose Stool [D4S 28-A]

Overall Size : Top 280 mm x 280 mm x H 445 mm

- Stainless steel top / Full S.S.
- Pre-treated and epoxy powder coated



Bed Side Locker (standard) [D4S 10]

Overall Size : W 370mm x D 305 mm x H 800 mm

- Two Stainless steel tops with 3 sides raised and front turned down
- Railings on 3 sides above main cabinet
- Handle and magnetic latch on door
- CRCA Steel box mounted on 20 mm x 20 mm CRCA sq. tube frames
- Pre-treated and epoxy powder coated



Revolving Stool [D4S 28]

Overall Size : Top 300 mm dia. with Ht. adjustable from 510 mm to 710 mm

- Stainless steel top / Full S. S.
- Pre-treated and epoxy powder coated



Instrument Cabinet [D4S 30]

Overall Size : W 610 mm x D 355 mm x H 1525 mm
Frame : CRCA Rectangular tube, CRCA sheet panels


- Pre-treated and epoxy powder coated /full stainless steel
- 5 Shelves, doors & 2 sides of 5 mm thick fine glass



Bed side screen [D4S 13]

Overall Size : L Open 2360 mm /
close 1200 mm x W 550 mm x H 1730 mm

- CRCA tube structure mounted on 50 mm dia wheels
- 2 Sides folding



Saline Stand [D4S 37]

Ht. adjustable from Min.:
1370 mm to Max. :2435 mm
• Available in full stainless steel /
epoxy powder coated stand & stainless
steel saline rod



KIRAN TECHNO SERVICES
PRIVATE LIMITED



Delhi : A - 39, Street No.5, Chander Industrial Estate, I.P Extension, Delhi - 110 092
Gujarat : 308/4, GIDC Estate, Makarpura Industrial Estate, Baroda – 390 010,
Andhra Pradesh : No. 6-5-351/5, Astal Heights, Road No. 1, Banjara Hills, Hyderabad – 500 034
Tamil Nadu : No.33-B4, Sidco Industrial Estate Phase - II, Sundarapuram, Coimbatore - 641 021
Kerala : 333/9, Para Road, Near NIDA, Kanjikode, Palakkad - 678 621
Helpline : +91 9751725000, +91 9751625000, + 91 9751425000, + 91 9751215000, + 91 9751245000
Email : kirangroup108@gmail.com, kiranindia108@gmail.com, kiranbiz108@gmail.com
Website : www.kirangroupofcompanies.com, www.kiran.biz, www.kiranshop.in, www.kiranindia.in,
www.kiransurgicalsystems.com, www.kiranmedicareequipments.com, www.kiranpump.com,
www.kiransolarsystems.com, www.kiranled.com, www.kiranincinerationsystems.com
www.kiranequipments.com, www.kirancompressors.com, www.kiranappliances.com
www.kiranairventilators.com, www.kiranfrpproducts.com, www.kiranpebsystems.com,
www.kiranpanels.com, www.kiranfabssystems.com