



KIRAN TECHNO SERVICES
PRIVATE LIMITED



KIRAN REHAB PRODUCTS

www.kiransurgicalsystems.com www.kiranmedicareequipments.com



KIRAN TECHNO SERVICES
PRIVATE LIMITED



Delhi : A - 39, Street No.5, Chander Industrial Estate, I.P Extension, Delhi - 110 092
 Gujarat : 308/4, GIDC Estate, Makarpura Industrial Estate, Baroda – 390 010,
 Andhra Pradesh : No. 6-5-351/5, Astal Heights, Road No. 1, Banjara Hills, Hyderabad – 500 034
 Tamil Nadu : No.33-B4, Sidco Industrial Estate Phase - II, Sundarapuram, Coimbatore - 641 021
 Kerala : 333/9, Para Road, Near NIDA, Kanjikode, Palakkad - 678 621
 Helpline : +91 9751725000, +91 9751625000, +91 9751425000, +91 9751215000, +91 9751245000
 Email : kirangroup108@gmail.com, kiranindia108@gmail.com, kiranbiz108@gmail.com
 Website : www.kirangroupofcompanies.com, www.kiran.biz, www.kiranshop.in, www.kiranindia.in,
www.kiransurgicalsystems.com, www.kiranmedicareequipments.com, www.kiranpump.com,
www.kiransolarsystems.com, www.kiranled.com, www.kiranincinerationsystems.com
www.kiranequipments.com, www.kirancompressors.com, www.kiranappliances.com
www.kiranairventilators.com, www.kiranfrpproducts.com, www.kiranpebsystems.com,
www.kiranpanels.com, www.kiranfabssystems.com

Company Profile

About Us

Kiran is a well known company in Indian as well as International Market for its innovative, durable, quality products since 2003. The company has achieved a very good reputation in the field of Patient Serving Industry by manufacturing products which are not only economical and patient friendly but are also the first choice of our distributors. The main reason behind the success is a personal attention and involvement of it's directors in the process of development of the products, manufacturing process, supply of products and last but very important, the feedback from actual users about the products. All our product ranges are designed on the basic principles like requirements of actual users the patients, value oriented prices that gives complete satisfaction to our customers and distributors and committed after sales support.

Kiran Rehabilitation products our main motto is

- To provide the complete range of Rehabilitation products which will help the patients to improve their functionality, which will results into maximum possible independence and bring the happiness in their life.
- To become a market leader in the field of rehabilitation by getting blessing of the patients
- To help our distributors to grow in Patient Serving Industry with business satisfaction
- To decrease our customer base by helping in improving their health, this ultimately will result in reference of new customers.

Index

Cervical Splints	1.....14	Knee, Calf & Ankle Splints	62.....102
Back Splint	15.....30	Fracture Aids	103.....120
Abdominal Splints	31.....37	General Aids	121.....127
Finger Splints	38.....61	Splints	128.....147



Cervical Pillow Regular -D4R-101

Sizes Available | Universal

Application | It Is Used To Support Neck Muscles And Ligaments During Sleeping Hours. It Maintains Proper Spinal Alignment In Back And Side Sleeping Positions.

Features and Benefits |

- Made from high quality of foam of optimum density.
- Special foam provides soft and cushion effect.
- Soft fabric covers for added comfort to neck.
- Provides support for head and neck.
- Light weight and easy to carry while travelling.
- It can be used for relief in different slipping discomforts)



Cervical Pillow Round -D4R-102

Sizes Available | Universal

Application | It Is Used To Support Neck Muscles And Ligaments During Sleeping Hours.it Maintains Proper Spinal Alignment In Back And Side Sleeping Positions

Features and Benefits |

- Round Tube Shaped Pillow To Support The Neck.
- Made From High Quality Of Foam Of Optimum Density.
- Special Foam Provides Soft And Cushion Effect.
- Special Foaming Covers For Added Comfort To Neck.
- Provides Support For Head And Neck.
- Light Weight And Easy To Carry While Travelling.



Cervical Pillow Contoured -D4R-103

Sizes Available | Universal

Application | It Is Ideal For Cervical Spondylosis.it Is Used To Support Neck Muscles And Ligaments During Sleeping Hours And For Cervical Spinal Disorders.

Features and Benefits|

- Design Of Different Heights On Two Sides With Chain Cover.
- Special Foaming Cover For Added Comfort To Neck.
- Light Weight And Easy To Carry While Travelling.



Cervical Pillow Contoured Latex -D4R-103A

Sizes Available | Universal

Application | It Is Ideal For Cervical Spondylosis.it Is Used To Support Neck Muscles And Ligaments During Sleeping Hours And For Cervical Spinal Disorders.

Features And Benefits |

- Special Latex Mateial For Comfort Cushion.
- Design Of Different Heights On Two Sides With Chain Cover.
- Light Weight And Easy To Carry While Travelling.



Cervical Pillow Contoured Latex -D4R-103B

Sizes Available | Universal

Application | It Is Ideal For Cervical Spondylolysis. It Is Used To Support Neck Muscles And Ligaments During Sleeping Hours And For Cervical Spinal Disorders.

Features And Benefits:-

- Special Latex Material For Comfort Cushion.
- Design Of Different Heights On Two Sides With Chain Cover.
- Light Weight And Easy To Carry While Travelling.



Cervical Collar Hard Adjustable -D4R-104

Sizes Available | S / M / L / XL

Application | It Is Used To Support Neck In Different Conditions Like Sprains, Strains, Whiplash Etc. It Is Ideal For Cervical Spondylolysis. Specially Recommended In Rigid Immobilization Of Cervical Region.

Features And Benefits |

- Made Of Polythene Sheet With Proper Foam Padding For Comfortable Wear.
- Double Layer For Easily Adjustable Height.
- Easy Velcro Closure.
- Provides Holes For Air Ventilation. Light In Weight.
- It Can Be Used During Work, Travelling, Reading, And Watching Television.etc.



Cervical Collar Hard -D4R-105

Sizes Available | S / M / L / XL

Application | It Is Used To Support Neck In Different Conditions Like Sprains, Strains, Whiplash Etc. It Is Ideal For Cervical Spondylolysis. Specially Recommended In Rigid Immobilization Of Cervical Region.

Features And Benefits |

- Made Of Polythene Sheet With Proper Foam Padding For Comfortable Wear.
- Easy Velcro Closure.
- Provides Holes For Air Ventilation.
- It Can Be Used During Work, Travelling, Reading, and Watching Television.etc.



Cervical Collar Soft Eyelet-D4R-106

Sizes Available | S / M / L / XL

Application | It Is Used To Support Neck In Different Conditions Like Sprains, Strains, Whiplash Etc; It Provides Soft And Rigid Support To Limit The Neck Movement Or To Prevent Harmful Neck Movements.

Features And Benefits |

- Anatomical Design Ensures Uniform Support To Neck.
- High Quality Foam Ensures Maximum Comfort To The Patient.
- Polypropylene Nonrust Eyelet Ensures Excellent Air Ventilation Which Avoids Skin Irritation.
- Special Soft Fabric Assures Non Allergic, Rashproof Use Even On Sensitive Skin.
- It Can Be Used During Work, Travelling, Reading And Watching Television.etc.



Cervical Immobilizer / Philadelphia Collar -D4R-107

Sizes Available | S / M / L / XI

Application | It Is Recommended For Cervical Spondylolysis. It Ensures Rigid Immobilization And Total Neck Comfort.

Features And Benefits |

- Premolded Chin Assures Rest To Chin And Balloon Design For Rim Contact Ensures Minimal Contact With Neck Surface
- High Quality Foam Ensures Maximum Comfort To The Patient
- Front Opening Ensures Compatibility To Endotracheal Tube And Excellent Air Ventilation Which Avoids Skin Irritation.
- Light Weight, Easy To Clean And Maintain.
- It Is Very Useful In Burn And Trauma Patients.



Philadelphia Collar Economy/ ESL collar -D4R-108

Sizes Available | S / M / L / XI

Application | It Is Recommended For Cervical Spondylolysis. It Ensures Rigid Immobilization And Total Neck Comfort.ideal For Immobilization Of Neck Upto 60%.

Features and Benefits |

- Premolded Chin Assures Rest To Chin And Balloon Design For Rim Contact Ensures Minimal Contact With Neck Surface
- High Quality Foam Ensures Maximum Comfort To The Patient
- Front Opening Ensures Compatibility To Endotracheal Tube And Excellent Air Ventilation Which Avoids Skin Irritation.
- Light Weight, Easy To Clean And Maintain.
- It Is Very Useful In Burn And Trauma Patients.



Cervical Collar Full Foam -D4R-109

Sizes Available | / M / L / XI

Application | It Is Used To Support Neck In Different Conditions Like Sprains, Strains, Whiplash Etc; It Provides Soft And Rigid Support To Limit The Neck Movement Or To Prevent Harmful Neck Movements.

Features And Benefits |

- Special Design Ensures Uniform Support To Neck.
Made Of High-density Foam.
- High Quality Foam Ensures Maximum Comfort To The Patient.
- Special Soft Fabric Assures Non Allergic, Rashproof Use Even On Sensitive Skin.
- It Can Be Used During Work, Travelling, Reading and Watching Television.etc.



Ambulance Collar/ Protector Collar-D4R-110

Sizes Available | Universal

Application | It Is Used As Temporary Support To Neck In Accidental Conditions. It Provides Rigid Support To Limit The Neck Movement Or To Prevent Harmful Neck Movements.

Features And Benefits |

- Single Piece Flat Design.
- Semi_rigid Plastic, Padded With Sturdy Foam Provides Comfortable Immobilization.
- It Is Used For Temporary Cervical Immobilization For Accidental Victims During The Transportation & Diagnostic Studies.



Somi Cervical Brace -D4R-111

Sizes Available | Universal

Application | It Is Recommended For Support And Immobilization To The Fractured Mid Or Lower Cervical Spine During The Pre & Post Operative Conditions.

Features and Benefits |

- Rigid, Two Post, Light Weight Metal Malleable Plates For Chin, Chest & Back For Proper Contouring And Fit On Patient.
- Three Sets Of Straps Provides Proper Immobilization.
- Foam Padding Provides Cushioning Comfort During Extended Use.
- Post-traction, Normal Cervical Spine Alignment Supports In The Case Of Dislocations.



Four Post Cervical Brace-D4R-112

Sizes Available | Universal

Application | It Is Recommended For Support And Immobilization To The Fractured Mid Or Lower Cervical Spine During The Pre & Post Operative Conditions.

Features And Benefits |

- Rigid, Four Post, Light Weight Metal Malleable Plates For Chin, Chest & Back For Proper Contouring And Fit On Patient.
- Three Sets Of Straps Provides Proper Immobilization.
- Eva Padding Provides Cushioning Comfort During Extended Use.
- Post-traction, Normal Cervical Spine Alignment Supports In The Case Of Dislocations.



**Taylor (Spinal) Brace Ordinary/
Dorso Lumbar Brace: -D4R-201**

Sizes Available | S/m/l/xl

Application | It Is Used To Immobilize Complete Spine Without Affecting The Other Body Movements

Features:

- Two Rigid Up-right Para Spinal and Lateral Bars of Light Weight are Pre-contoured and Benefits To Keep Spine In It's Natural Position.
- Easy to Use, Comfortable.
- Shoulder Straps are Adjustable Which Allows Excellent Gripping and Tightening of Brace
- It Is Recommended in Spinal Deformity, Vertebral Fractures or Collapse, osteoporosis of The Spine, Congenital Abnormalities of The Lumbar Sacral Region and Spondylolysis.



**Taylor (Spinal) Brace Deluxe/
Dorso Lumbar Brace-D4R-202**

Sizes Available | S/m/l/xl

Application | It Is Used To Immobilize Complete Spine Without Affecting The Other Body Movements

Features |

- Two Rigid Up-right Para Spinal and Lateral Bars Of Light Weight are Pre-contoured and Benefits To Keep Spine in It's Natural Position.
- Easy To Use, Comfortable.

- Elastic Abdominal Panels Provide The Required Compression.
- Shoulder Straps are Adjustable Which Allows Excellent Gripping and Tightening of Brace
- It Is Ideal For Weakened Dorso Lumbar Spine Due To T.b & Suitable For Post Operative Rehabilitation of Back Ache Patients.



Orthopaedic Back Rest Long Type-D4R-203

Sizes Available | Universal

Application | To Maintain The Natural Curve of Lower Back and Helps In Releasing The Pressure on Vertebrae and Gives Correct Body Posture While Sitting

Features And Benefits |

- Easy To Use
- Special Foam Provides Better Cushioning and Longer Life To The Product.
- It Is Ideal for Slip Disc Patients.reduces Discomfort and Pain After Surgery.
- It Improves Posture When Placed on Chairs.



Orthopaedic Back Rest Short Type-D4R-204

Sizes Available | Universal

Application | To Maintain The Natural Curve Of Lower Back And Helps In Releasing The Pressure On Vertebrae And Gives Correct Body Posture While Sitting

Features and Benefits |

- Easy To Use.
- It Gives Spinal Support When Placed On Chair, Auto Seat & Wheel Chair.
- Special Foam Provides Better Cushioning And Longer Life To The Product.
- It Is Ideal For Slip Disc Patients.reduces Discomfort And Pain After Surgery.
- It Improves Posture When Placed On Chairs.



Orthopaedic Back Rest Regular-D4R-205

Sizes Available | Universal

Application | To Maintain The Natural Curve Of Lower Back And Helps In Releasing The Pressure On Vertibras And Gives Correct Body Posture While Sitting

Features And Benefits |

- Orthopaedically Designed For Sitting Comfort.
- It Gives Spinal Support When Placed On Chair, Auto Seat & Wheel Chair.
- It Is Ideal For Slip Disc Patients.reduces Discomfort And Pain After Surgery.
- It Improves Posture When Placed On Chairs.



Sacro Lumber Brace (Belt) Heavy Duty-D4R-206

Sizes Available | 28" To 44"

Application | To Provide Uniform Support To Lower Back Section With Rigid Immobilization

Features And Benefits |

- Special Elastic And Thin Material Provides Uniform Compression And Rigid Support To Back Section.
- Easy To Apply, And Easy To Adjust Different Sizes As Per User's Requirements.
- No Rash, No Allergy, Easy To Clean And Maintain
- It Is Beneficial For Lower Back Problems, Degenerative Disc Problems, Lumbosacral Spondylolysis, Osteoporotic Pains Of Lumbar Spines.



Sacro Lumber Brace (Belt) Contoured-D4R-207

Sizes Available | S / M / L / XL / Xxl /
Xxxl

Application | To Provide Uniform
Support To Lower Back Section With
Comfortable Immobilization

Features and Benefits |

- Double Elastic Straps And Thin Material Provides Uniform Compression And Comfortable Support To Back Section
- Four Aluminium Maleable Metal Stays Are Removable & Can Be Given Desired Contoured To Support The Specific Vertebrae.
- Easy To Apply, And Easy To Adjust Different Sizes As Per User's Requirements.
- It Is Beneficial For Patients Having Pendulous Abdomen With Lower Back Problems, Degenerative Disc Problems, lumbosacral Spondylolysis, osteoporotic Pains Of Lumber Spines.



Sacro Lumber Brace (Belt) Contoured Economy-D4R-208

Sizes Available | S / M / L / XL / XXI /
XXXI

Application | To Provide Uniform
Support To Lower Back Section With
Comfortable Immobilization

Features And Benefits |

- Double Elastic Straps And Thin Material Provides Uniform Compression And Comfortable Support To Back Section
- Two Aluminium & Two Iron Maleable Metal Stays Are Removable & Gives Effective Support To The Back..
- Easy To Apply, And Easy To Adjust Different Sizes As Per User's Requirements.
- It Is Beneficial For Patients Having Pendulous Abdomen With Lower Back Problems, Degenerative Disc Problems, lumbosacral Spondylolysis, osteoporotic Pains Of Lumber Spines.



**Sacro Lumber Brace (Belt)
Super Fine Double Lock -D4R-209**

Sizes Available |28" To 44"

Application | To Provide Uniform Support To Lower Back Section With Rigid Immobilization.

Features and Benefits |

- It Is Made Of Tampered Steel Strips.additional Straps & Sliding Buckle Closure System Allows Extra Compression For The Abdomen.
- Fully Elasticated To Add Extra Pressure And Support To The Sacro & Lumbar Regions.
- Easy To Apply, And Easy To Adjust Different Sizes As Per User's Requirements.
- It Is Beneficial For Patients With Lower Back Problems, Degenerative Disc Problems, lumbosacral Spondylolysis, osteoporotic Pains of Lumbar Spines

IMAGE COMING SOON

**Sacro Lumber Brace (Belt)
Super Fine-D4R-210**

Sizes Available | 28" To 44"

Application |To Provide Uniform Support To Lower Back Section With Rigid Immobilization.

Features And Benefits |

- It Is Made Of Tampered Steel Strips.
- Fully Elasticated To Add Extra Pressure And Support To The Sacro & Lumbar Regions.
- Easy To Apply, And Easy To Adjust Different Sizes As Per User's Requirements.
- It Is Beneficial For Patients With Lower Back Problems, Degenerative Disc Problems, lumbosacral Spondylolysis, osteoporotic Pains Of Lumbar Spines



Sacro Lumber Brace (Belt) Fine-D4R-211

Sizes Available | 28" To 44"

Application | To Provide Uniform Support To Lower Back Section With Rigid Immobilization.

Features And Benefits |

- It Is Made Of Tempered Steel Strips.
- Fine Quality Elastic To Add Extra Pressure And Support To The Sacro & Lumbar Regions.
- Easy To Apply, And Easy To Adjust Different Sizes As Per User's Requirements.
- It Is Beneficial For Patients With Lower Back Problems, Degenerative Disc Problems, lumbosacral Spondylolysis, osteoporotic Pains of Lumbar Spines



Sacro Lumber Brace (Belt) Deluxe-D4R-212

Sizes Available | 28" To 44"

Application | To Provide Uniform Support To Lower Back Section With Rigid Immobilization.

Features And Benefits |

- Two Rigid Pre Contoured Splints Provide Support To The Weak Back Muscle.
- Velcro Closure For Easy Application & Adjustment To Comfort Fit.
- Belt Is Made Of Breathable, High Quality Fabric & Elastic.
- Easy To Apply, And Easy To Adjust Different Sizes As Per User's Requirements.
- It Is Beneficial For Patients With Lower Back Problems, Degenerative Disc Problems, lumbosacral Spondylolysis, osteoporotic Pains Of Lumbar Spines.



Eco Frame Back Support-D4R-213

Sizes Available | 28" To 44"

Application | To Provide Firm Support To Back And Lumbar Region

Features And Benefits |

- Light Weight Metal Frame With Flexibility To Adjust The Contoured Of The Back Region
- Special Elastic And Thin Material Provides Uniform Compression And Scientific Design Gives Rigid Support To Back Section
- Easy To Apply, And Easy To Remove. Due To Hook And Loop Arrangement
- No Rash, No Allergy,easy To Clean And Maintain
- Complete Round Back Support Of Metal Frame Restricts The Side Movement Of The Back
- It Is Useful In Treatment Of Acute Lumbago,slip Disc Recurrent Lower Back



Ternal Splint (after Heart Surgery) -D4R-214-S

Sizes Available | S,m,l,xl

Application | For Protection Of Sternum After Heart Surgery

Features And Benefits |

- Special Soft Foam Pad Provides Mild Compression Support To Sternum
- Easy To Apply And Remove
- Allows Normal Respiration
- It Is Useful In Protecting Surgical Incision On Chest.



Lumbo Sacral Spinal Support With Cushion-D4R-215

Sizes Available | S/m/l/xl

Application | To Support Lower Back And To Maintain Posture During Prolonged Sitting

Features And Benefits |

- Anatomically Contoured Shape Helps To Maintain Natural Curve Of Spine.
- Additional Strap Provide Extra Compression.
- Release Pressure On Lower Back Region
- Special Cushion Provides Ultimate Comfort
- It Is Very Useful For Lower Back Pains.it Can Be Used While Driving, Working In Office,resting At Home.



Orthopaedic Coccyx Wedge Seat -D4R-216

Sizes available | universal

Application | To release the pressure on lower spinal discs and tail bone(coccyx) due to prolonged sitting

Features and benefits |

- Scientific design release pressure points on coccyx region
- Special foam provides total comfort to complete lower back section
- It is useful to prevent severe back pain in prolonged sitting position.



6" Abdominal Support-D4R-301

Sizes Available | 28" to 44"

APPLICATION | It Is Used To Support Abdomin Or Back Or To Compress The Abdominal Muscles

Features and Benefits |

- Special Elastic And Thin Material Provides Uniform Compression of Abdominal Muscles and Comfortable Support To Back Section
- Easy To Apply, And Easy To Adjust Different Sizes as Per User's Requirements.
- No Rash, No Allergy, easy To Clean And Maintain Beneficial For Patients With Post Operative Supports To Abdominal Region.



8" Abdominal Support Super fine -D4R-302

Application | It Is Used To Support Abdomin Or Back Or To Compress The Abdominal Muscles

Features and Benefits: |

- Superior Quality Elastic And Thin Material Provides Uniform Compression Of Abdominal Muscles And Comfortable Support To Back Section.
- Easy To Apply, And Easy To Adjust Different Sizes As Per User's Requirements.
- No Rash, No Allergy, Easy To Clean And Maintain
- It Is Beneficial For Patients With Back Problems, post Operative Supports To Abdominal Region, for Mothers After Delivery, for Ventral Hernia



8" Abdominal Support fine -D4R-303

Application | It Is Used To Support Abdomin Or Back Or To Compress The Abdominal Muscles.

Features And Benefits:|

- Fine Elastic And Thin Material Provides Uniform Compression Of Abdominal Muscles And Comfortable Support To Back Section
- Easy To Apply, And Easy To Adjust Different Sizes As Per User's Requirements.
- No Rash, No Allergy, Easy To Clean And Maintain
- It Is Beneficial For Patients With Back Problems, post Operative Supports To Abdominal Region, for Mothers After Delivery, for Ventral Hernia.



9" Abdominal support -D4R-304

Application | It Is Used For Post Surgical Compression Support In Abdominal Region.

Features And Benefits |

- Combination Of Breathable Elastic And Cotton Fabric Offers Compression And Warm Feeling.
- Velcro Closure Allows Easy Application And Proper Fit.
- Additional Double Lock Given For Adjustable Support To The Waist.
- No Rash, No Allergy, Easy To Clean And Maintain.
- It Is Beneficial For Patients With Post Operative
- Supports To Abdominal Region, For Mothers After Delivery.



10" Abdominal support -D4R-305

Application | It Is Used To Support Abdomin Or Back Or To Compress The Abdominal Muscles

Features And Benefits |

- Superior Quality, Breathable Elastic Band Of 10" Width Offers Good Compression Support
- Easy To Apply, And Easy To Adjust Different Sizes As Per User's Requirements.
- No Rash, No Allergy, Easy To Clean And Maintain
- It Is Beneficial For Supports To Abdominal Region. Maintaining Cosmetic Waist Line For Females.




Hernia Support Double Pad-D4R-306

Sizes Available | Universal

Application | to Hold The Tissues Affected By Hernia And To Reduce The Pain And Uncomfortable Feeling In The Affected Area.

Features And Benefits |

- Made Of High Quality Elastic Band Of 3" Width.
- Two Soft Pads, Left And Right Side Helps In Reducing Inguinal Hernia By Targeted Compression On Affected Area.
- Simple Velcro Closure Allows Easy Application & Comfortable Fit.
- It Is Useful For Patients Suffering From Painful Inguinal Hernia Or Even After Hernia Surgery.



Hernia Support Single Pad-D4R-307

Sizes Available | Universal

Application | to Hold The Tissues Affected By Hernia And To Reduce The Pain And Uncomfortable Feeling In The Affected Area

Features And Benefits |

- Made Of High Quality Elastic Band Of 3" Width.
- One Soft Pads, Left And Right Side Helps In Reducing Inguinal Hernia By Targeted Compression On Affected Area.
- Simple Velcro Closure Allows Easy Application & Comfortable Fit.
- It Is Useful For Patients Suffering From Painful Inguinal Hernia Or Even After Hernia Surgery.



Finger Baseball Splint-D4R-401

Size | S,M,L

Application | Hyper Flexion Finger Injuries Due To The Direct Hit Of The Balls. In The Game Like Cricket, Baseball Etc.

Features and Benefits |

- Innovative Finger Tip Fold Over Design
- Alluminium Splint Can Be Reshaped As Per Required Finger Position
- Special Foam Padding Provides Comfortable Fitting.



Finger Cot -D4R-402

Sizes Available | S,m,l

Application | To Protect Finger In Finger Tip Injuries,nail Bed Injuries

Features and Benefits |

- Design Allows Application Without Taping
- U Shaped Alluminium Splint Covers Finger Tip And Assures Rigid Immobilization
- Special Foam Padding Provides Comfortable Fitting
- Used For Fingure Tip Injuries And Burn Injuries



Frog Splint D4R-403

Sizes Available | S,m,l

Application | To Protect And Treat Distal Finger Joint Injuries

Features and Benefits |

- Unique Four Leg Construction Allows Self Closing.
- Versatile Design Allows Dorsal Surface Of Finger To Remain Open For Ventilation.
- Keeps Interphalangeal Joint In Correct Position
- Special Foam Padding Provides Comfortable Fitting
- Useful In Hyperextension Injuries Collateral Ligament Injuries Of Distal Interphalangeal Joint/ Useful For Baseball, Cricket Players.



Straight Splint-D4R-404

Sizes Available | S,m,l

Application | Recommended For Support & Protect Dislocated Fractures

Features And Benefits |

- Simple Flat Design To Protect Hole Finger
- Sizes Allows Use For Full Length Of The Finger



Finger Extension Splint-D4R-405

Sizes Available | -s/m/l

Application | Immobilization Of The Fingers Following Fracture.

Features and Benefits |

- Light Weight Splints Provides Effective Immobilization of The Fingers
- Can Be Moulded Repeatedly. Foam Padding Provides Comfort & Ventilation To The Finger
- Easy To Apply By Just Inserting The Finger



Dynamic Finger Flexion Splint -D4R-406

Sizes Available | Universal

Application | Flexes The Pijoint
Without Interfering The Mcp Joint

Features And Benefits |

- Springs Attached To Give Extension.
- For Extra Comfort Proximal And Distal Pads Are Covered With Soft Padding Which Ensure Minimal Finger Slippage.



Finger Protector Splint -D4R-407

Sizes Available | Universal

Application | to Protect Injured
Finger Tip

Features Anf Benefits |

- Can Be Bent At Any Degree
- Easy To Apply And Remove

Use | Useful In Finger Tip Injuries



Thumb Spica Splint -D4R-408

Sizes Available | Universal

Application | To Support The Cmc Joint And Immobilize The Mcp Joint Of Thumb Without Affecting Hand Movement

Features And Benefits |

- Full Hand And Finger Movement Is Possible
- Reduces Pain And Inflammation
- Customised Abduction Of Thumb Is Possible

Use | For Thumb Injuries



Thumb Abduction Splint -D4R-409

Sizes Available | Universal

Application | for Sprains & Strains To The Thumb's Metacarpal Joint.

Features And Benefits |

- Special Design Minimize Thumb Movements
- Easy To Apply And Remove
- Useful For Skier's Thumb, Bowler's Thumb And General Thumb Injuries



Elastic Wrist Splint Long Type-D4R-410

Sizes Available | Universal (right / Left)

Application | For Sprains And Strains In Wrist

Features And Benefits |

- Unique Design Provides Firm Support To Wrist
- Can Be Used In Left And Right Hand
- Easy To Apply And Remove
- Comfortable And Quick Pain Relief
- Useful In Mild To Moderate Sprains And Strains



Elastic Wrist Splint Short Type -D4R-411

Sizes Available | S,m,l,Xl (right / Left)

Application | For Sprains And Strains In Wrist

Features And Benefits |

- Unique Design Provides Firm Support To Wrist
- Can Be Used In Left And Right Hand
- Easy To Apply And Remove
- Comfortable And Quick Pain Relief
- Useful In Mild To Moderate Sprains And Strains



Dynamic Cock Up Splint (Left & Right) -D4R-412

Sizes Available | S.m.l.XI.XXI

Application | Designed For Dynamic Splinting Following Mcp Arthroplasty Or Radial Nerve Injury Or Any Time Extention. To Keep The Forehand In Desired Angular Position

Features And Benefits |

- Unique Spring Mechanism Allows Perfect Angle Of Flexion And Extension
- Preshaped Splint Helps To Adjust Forehand At Desired Position To Achieve Quick Healing.
- Useful In Sprains, Post Fracture, arthritis



Cock Up Splint With Thumb Support (Left & Right) -D4R-413

Sizes Available | S.m.l (right / Left)

Application | To Immobilize Wrist Movement, Finger & Thumb Movement Also.

Features And Benefits |

- Malleable Structure Can Be Moulded To Fit The Contour Of For Hand And Arm
- Dorsiflexion Angle Can Be Fixed By Adjusting Velcro Straps
- Foam Padding Provides Excellent Comfort To The Patient.
- Useful For Immobilization In Post Fracture, Tendon Injuries Of The Wrist.



Elastic Wrist Brace -D4R-414

Sizes Available | Universal

Application | To Immobilize Wrist Movement Without Affecting Finger Movement

Features and Benefits |

- Unique Design And Elasticity Provides Firm Support To Wrist
- Easy To Apply And Remove
- Comfortable And Quick Pain Relief
- Useful In All Types Of Strains In Wrist.also Useful In Bursities Microtrauma And Arthritis



Wrist Binder With Double Lock -D4R-415

Sizes Available | S,m,l,

Application | To Immobilize Wrist Movement Without Affecting Finger Movement

Features and Benefits |

- Unique Design And Elasticity Provides Firm Support To Wrist.double Lock Helps In Firm Gripping
- Easy To Apply And Remove
- Comfortable And Quick Pain Relief
- Useful In All Types Of Strains In Wrist.also Useful In Bursities Microtrauma And Arthritis



Wrist & Forearm Splint long type-D4R-416

Sizes Available | S,M,L (right /left)

Application | It Is Used For Wrist Fracture, Early Cast Remover & Mild Stable Distal Forearm Fracture.

Maleable Light Metal Splint On The Palmar, Dorsal And Ulnar Sides To Encase The Wrist & Forearm & Restrict Movements Completely.

Features and Benefits |

- Three Sets Of Buckles & Velcro Straps Allow Easy Application & Proper Fit.
- Designed Separately For Right & Left Arms.



Carpal Tunnel Splint-D4R-417

Sizes Available | Universal

Application | To Prevention Of Ct Syndrome Occurance.

Features And Benefits |

- Elastic Band With Velcro Closure Keeps Palmar Splint In Position & Gives Required Compression.
- Thumb Loop Allows Self Application Of The Product.
- High Quality Felt Pad With Rigid But Light Metal Splint Act As Palmar Support To Keep The Wrist In Natural Position.



Solar Splint-D4R-418

Sizes Available | Universal

Application | To Prevention Of Ct Syndrome Occurance.

Features and Benefits |

- Elastic Wrap With Buckle & Reserve Lock Velcro Closure Enables Easy Application & Comfortable Fit.
- Unique Adjustable Loop For Middle Finger Control Flexion Movement Of Wrist.



Tennis Elbow Protector -D4R-419

Sizes Available | S,m,l,xl

Application | For Relief From Tennis Elbow, Forearm Pain.

Features and Benefits |

- High Strength Strap Exerts Counter-rotation Force At The Musculo-Tendon Region Of Pronator-supinator.
- Special Elastic Material Ensures Strong Grip And Comfortable Compression
- Also Recommended In Sport Injuries.



Tennis Elbow Protector (Dlx) -D4R-419A

Sizes Available | S,M,I,XI

Application | For Relief From Tennis Elbow, Forearm Pain.

Features and Benefits |

- Extra Foam Padding For Support To The Inflamed Region
- High Strength Strap Exerts Counter-rotation Force At The Musculo- Tendon Region Of Pronator-supinator.
- Special Elastic Material Ensures Strong Grip And Comfortable Compression
- Also Recommended In Sport Injuries.



Elastic Tubular Elbow Support With/ Without Center Hole (Deluxe)-D4R-420

Sizes Available | S,M,I,XI

Application | It Is Used To Reduce Pains, Control Sprains And Strains.provides Support To Elbow Joint And Surrounding Muscles

Features and Benefits |

- Easy To Apply. A Simple Pull Over Application.
- Special Knitting To Cover The Affected Area In All Four Ways, Which Assures Effective Compression, Instant Relief, And Retains The Body Heat.
- Scientific Design Ensures Uniform Gripping, Excellent Aesthetics And Durability.
- Special Fabric Assures Better Comfort Without Rash And Alergy
- Easy To Clean And Maintain
- Easy Movement Of Elbow.
- It Can Be Used During Daily Routines, Excercises, Sports And Post Surgicals Rehabilitation.specially Recommended In Sprains And Strains In Elbow Joints During Sport Injuries.



Elastic Tubular Elbow Support With/ Without Center Hole (Ordinary)-D4R-420A

Sizes Available | S,M,L,XL

Application | It Is Used To Reduce Pains, Control Sprains And Strains.provides Support To Elbow Joint And Surrounding Muscles

Features And Benefits |

- Easy To Apply. A Simple Pull Over Application.
- Special Knitting To Cover The Affected Area In All Four Ways, Which Assures Effective Compression, Instant Relief, And Retains The Body Heat.
- Scientific Design Ensures Uniform Gripping, Excellent Aesthetics And Durability.
- Special Fabric Assures Better Comfort Without Rash And Allergy
- Easy To Clean And Maintain
- Easy Movement Of Elbow.
- It Can Be Used During Daily Routines, Exercises, Sports And Post Surgicals Rehabilitation. specially Recommended In Sprains And Strains In Elbow Joints During Sport Injuries



Finger mallet splint-D4R-421

Sizes Available | Universal

Application | For Correction Of Mallet Finger Deformity

Features and Benefits |

- Easily Malleable To Get Desired Shape
- Provides Strong Immobilization
- Easy To Apply And Remove
- It Is Useful In Treatment Of Mallet Finger Deformity.





Stax Mallet -D4R-422

Sizes Available | Universal

Application | It Is Used In Hyper Flexion Injuries Like Mallet Finger

Features And Benefits |

- Made From Engineering Plastic That Is Durable & Contains Ventilations Holes.
- Splints Keeps The Finger In 10 Degree Of Hyper Extension.

Elastic Tubular Knee Support Deluxe-D4R-501

Sizes Available | S / M / L / XI / Xxl

Application | It Is Used To Support Knee Joint In Different Conditions Like Weak Knee Due To Old Age, Arthritis, Sport Injuries, Support For Heavy Workouts In Sports, Post Orthopedic Surgical Rehabilitation

Features And Benefits |

- Easy To Apply. A Simple Pull Over Application.
- Special Knitting To Cover The Affected Area In All Four Ways, Which Assures Effective Compression, Instant Relief, and Retains The Body Heat.
- Scientific Design Ensures Uniform Gripping, Excellent Aesthetics And Durability.
- Special Fabric Assures Better Comfort Without Rash And Allergy.
- Easy To Clean And Maintain
- It Can Be Used During Daily Routines, Exercises, sports And Post Surgical Rehabilitation.



Elastic Tubular Knee Support With Center Hole -D4R-502

Sizes Available | S / M / L / XL / Xxl

Application | It Is Used To Support Knee Joint In Different Conditions Like Weak Knee Due To Old Age, Arthritis, Sport Injuries, Support For Heavy Workouts In Sports, Post Orthopedic Surgical Rehabilitation

Features And Benefits |

- Hole Provides To Guide The Patella Position.
- Easy To Apply. A Simple Pull Over Application.
- Special Knitting To Cover The Affected Area In All Four Ways, Which Assures Effective Compression, Instant Relief, and Retains The Body Heat.
- Scientific Design Ensures Uniform Gripping, Excellent Aesthetics And Durability.
- Special Fabric Assures Better Comfort Without Rash And Allergy
- Easy To Clean And Maintain
- It Can Be Used During Daily Routines, Exercises, sports And Post Surgicals Rehabilitation.



Elastic tubular knee support with hinges (Single)-D4R-503

Sizes Available | S / M / L / XL / Xxl

Application | It Is Used In Painful Knee Movement Of Arthritis & To Prevent Mild Strain & Sprains That May Occur Due To Activities Of Daily Living.

Features And Benefits |

- Unique Flexible Splints Provide Lateral Support While Controlling Movement Of The Knee
- Product Design Has Velcro Closure To Fit Support Around Knee.
- Special Knitting To Cover The Affected Area In All Four Ways, Which Assures Effective Compression, •Instant Relief, and Retains The Body Heat.
- Scientific Design Ensures Uniform Gripping, Excellent Aesthetics And Durability.
- Special Fabric Assures Better Comfort Without Rash And Allergy.



Elastic knee support with hinges-D4R-504

Sizes Available | S / M / L / XI / Xxl

Application | It Is Used In Painful Knee Movement Of Arthritis & To Prevent Mild Strain & Sprains That May Occur Due To Activities Of Daily Living.

Features And Benefits |

- Excellent Quality Elastic Of 9" Width Forms The Main Body Of The Support.
- Anterior Patella Hole & Posterior Cut Out Are Provided To Avoid Bunching
- Unique Metal Splints Provide Lateral Support While Controlling Movement Of The Knee
- Product Design Has Velcro Closure To Fit Support Around Knee.



Elastic knee support DELUXE -D4R-505

Sizes Available | S / M / L / XI / Xxl

Application | It Is Used In Painful Knee Movement Of Arthritis & To Prevent Mild Strain & Sprains That May Occur Due To Activities Of Daily Living.

Features And Benefits |

- Excellent Quality Elastic Of 9" Width Forms The Main Body Of The Support.
- Product Design Has Velcro Closure To Fit Support Around Knee.



Elastic Knee Support Ordinary-D4R-506

Sizes Available | S / M / L / XL / XXL

Application | It Is Used In Painful Knee Movement Of Arthritis & To Prevent Mild Strain & Sprains That May Occur Due To Activities Of Daily Living.

Features And Benefits |

- Fine Quality Elastic Of 9" Width Forms The Main Body Of The Support.
- Product Design Has Velcro Closure To Fit Support Around Knee.



Functional Knee Support -D4R-507

Sizes Available | S,m,l,xl,xxl

Application | For Controlled Compression Of Knee Joint And Rigid Lateral Support And Immobilization

Features and Benefit |

- Control Compression & Rigidity Of Splint With Normal Flexion & Free Movement.
- Special Fabric Assures Better Comfort Without Rash And Allergy
- Easy To Clean And Maintain
- Allows Free Flexion Movement
- It Can Be Used During Daily Routines. post S u r g i c a l s Rehabilitation.specially Recommended In Sprains And Strains Arthritis In Patients.

IMAGE COMING SOON

Thigh Band-D4R-508

Sizes Available | S,m,l,xl,xxl

Application | Support In Case Of Quadriceps Muscle Strain & Hamstring Pull.

Features and Benefits |

- High Quality 10" Elastic Band Provides Compression & Supports To Thigh Muscles.
- High Tenacity Strap Along The Length Of Product Maintains It Tubular Structure.



Elastic Tubular Anklelet Deluxe-D4R-509

Sizes Available | S / M / L / XI

Application | It Is Used To Support Ankle Joint In Different Conditions Like Weak Ankle Due To Old Age, ,sport Injuries,support For Heavy Workouts In Sports,post Orthopedic Surgical Rehabilitation

Features And Benefits |

- Easy To Apply. A Simple Pull Over Application.
- Special Knitting To Cover The Affected Area In All Four Ways, Which Assures Effective Compression, Instant Relief,and Retains The Body Heat.
- Scientific Design Ensures Uniform Gripping, Excellent Aesthetics And Durability.
- Special Fabric Assures Better Comfort Without Rash And Allergy.
- Easy To Clean And Maintain.
- It Can Be Used During Daily Routines, Exercises, Sports And Post Surgicals Rehabilitation.



Ankle Brace-D4R-510

Sizes Available | S / M / L / XI

Application | It Is Used In Ligament Tears & Hairline Fracture & Early Cast Removal

Features and Benefits |

- Durable Quoted Fabric With Inner Lining Of Form Comfort & Fit Around The Ankle Joint
- High Quality Semi Rigid Plastic Stays On Sides Of Base Control Abduction.
- Restrict Flexion-extension Movements
- Specialvelcro Straps And Plastic Buckles Assures Effective Reverse Locking With Secured Fitting To Ankle Joint.
- Can Be Worn On Either Left Or Right Ankle



Elastic Ankle Binder -D4R-511

Sizes Available | S / M / L / XI / Xxl

Application | It Is Used In A Mild Ankle Sprain During The Activities Of Daily Living.

Features and Benefits |

- Special Elastic Fabric In Figure Eight Wrap Assures Controlled Gripping All Around The Ankle Joint And Pain Relief.
- Convenient To Apply And Can Be Used For Prolonged Period.
- It Is Used In Different Conditions Like Weak Ankle Due To Old Age, Serious Ankle Injuries, Post Operative Orthopedic Surgical Rehabilitation.
- Can Be Worn Inside Shoes.



Ankle Support-D4R-512

Sizes Available | S / M / L / XI

Application | It Is Used To Support Ankle Joint Against Sprains In Sports & During Normal Activities.

Features And Benefits |

- Pre Moulded Plastic Shells With Sturdy FoamPadding Offer Excellent Lateral Support While Allowing Normal Flexion extension Movements Of Ankle.
- Velcro Closure Allows Quick Application & Proper Fit. adjustable (adjustable) Velcro System Accommodates (accommodates) All Ankle Sizes.



Elastic Tubular Vericose Vein Stockings -Thigh Length-D4R-513

Sizes Available | S,m,l,xl

Application | To Releive Discomfort From Tired Swollen Legs

Features and Benefit |

- Provides Support To Thromboembolism
- Seamless Stocking From Mid Thigh To Ankle
- Scientific Design Ensures Uniform Gripping, Excellent Aesthetics And Durability
- Special Fabric Assures Better Comfort Without Rash And Alergy
- Easy To Clean And Maintain.



Elastic Tubular Vericose Vein Stockings -Below Knee -D4R-513A

Sizes Available | S,m,l,xl

Application | To Relieve Discomfort From Tired Swollen Legs

Features and Benefits |

- Provides Support To Thromboembolism
- Seamless Stocking From Below Knee To Ankle
- Scientific Design Ensures Uniform Gripping, Excellent Aesthetics (aesthetics) And Durability
- Special Fabric Assures Better Comfort Without Rash And Allergy
- Easy To Clean And Maintain.




ROM Brace-D4R-514

Sizes Available | Universal

Application | Provides In The Fracture Of Knee Pre Fabricated. To Control And Restrict Flexion-extension Of The Knee Joint

Features And Benefits |

- Scientific Design Allows Controlled Movements In Any Two Angles
- Can Be Adjusted Easily By Changing The Position Of Pin
- Knee Joint Can Be Immobilize At Any Angle
- Transparent Disc With Pin Locking For Added Safety
- It Can Be Used For Maximum orthopedic Problems Related To Knee Joints.



Knee Brace Short Type -D4R-515

Sizes Available | S / M / L / XI

Application | It Is Used In Fractures,ligament Or Muscle Related Injuries And Dislocation Of Knee Joint

Features and Benefits |

- Easy To Apply. Light Weight,and Comfortable
- Covers Maximum Portion, knee And Calf Which Assures Immobilization Of Affected Knee Joint
- Scientific Design Ensures Uniform Gripping, Excellent Asthetics And Durability With Rigid Splinting At Sides .
- Special Fabric Assures Better Comfort Without Rash And Alergy
- Easy To Clean And Maintain
- It Can Be Used During Daily Routines, Post Surgicals Rehabilitation.specially Recommended In Serious Knee Problems



Knee Brace Long Type -D4R-516

Sizes available | s / m / l / xl

Application | It is used in fractures,ligament or muscle related injuries and dislocation of knee joint

Features and benefits |

- Easy to apply. Light weight,and comfortable
- Covers maximum portion ,knee and calf which assures immobilization of affected knee joint
- Scientific design ensures uniform gripping, excellent asthetics (aesthetics) and durability with rigid splinting at sides .
- Special fabric assures better comfort without rash and alergy
- Easy to clean and maintain
- It can be used during daily routines, post surgicals rehabilitation.specially recommended in serious knee problems.



Knee Brace Deluxe(S) -D4R-517

Sizes Available | S / M / L / XI

Application | It Is Used In Fractures,ligament Or Muscle Related Injuries And Dislocation Of Knee Joint

Features and Benefits |

- High Qulaity Breatheable (breathable) Foaming Fa b r i c Provide Excellent Comfort.
- To Lateral And Three Posterior Rigid Yet Light Metal Stays Pre Contoured For Weight Bearing Support.
- The Stays Are Removable And Mouldable (moldable) For Proper Fitting.
- Buckle & Velcro Closures Above & Below Knee Provide Good Fit.



Heel Pad (Silicone) -D4R-518

Sizes Available | S/m/l

Application | To Provide Cushioning Effect To Heel Against The Painful Un-operatable Overgrowth And Development Of Calcaneum Spur (calcaneal Spur)

Features And Benefits |

- Heel Pad Made Of Soft Non-irritant 100% Medical Grade Silicone
- Scientific Design Provides Comfortable Cushioning And Adjusts Automatically As Per The Pressure On Toe Region Which Gives Relief From Pain To The User
- Durable And Long Life
- Easy To Clean And Maintain
- It Is Useful In Common Orthopedic Problems Of Heel Region.



Heel Pad (Gel) -D4R-519

Sizes Available | S/m/l

Application | To Provide Cushioning Effect To Heel Against The Painful Un-operatable Over Growth And Development Of Calcaneum Spur (calcaneal Spur)

Features and Benefits |

- Heel Pad Made Of Gel
- Scientific Design Provides Comfortable Cushioning And Adjusts Automatically As Per The Pressure On Toe Region Which Gives Relief From Pain To The User
- Durable And Long Life
- Easy To Clean And Maintain
- It Is Useful In Common Orthopaedic Problems Of Heel Region.



Rheel Pad Donut-D4R-520

Sizes Available | S/m/l

Application | To Provide Cushioning Effect To Heel Against The Painful Un-operatable Over Growth And Development Of Calcaneum Spur (calcaneal Spur)

Features and Benefit |

- Heel Pad Donut Made Of Gel
- Scientific Design Provides Comfortable Cushioning And Adjusts Automatically As Per The Pressure On Toe Region Which Gives Relief From Pain To The User
- Durable And Long Life
- Easy To Clean And Maintain
- It Is Useful In Common Orthopaedic Problems Of Heel Region.



Heel Pad Female-D4R-521

Sizes Available | S/m/l

Application | To Provide Cushioning Effect To Heel Against The Painful Un-operatable Over Growth And Development Of Calcaneum Spur (calcaneal Spur)

Features And Benefit |

- Heel Pad Female Made Of Gel
- Scientific Design Provides Comfortable Cushioning And Adjusts Automatically As Per The Pressure On Toe Region Which Gives Relief From Pain To The User
- Durable And Long Life
- Easy To Clean And Maintain
- It Is Useful In Common Orthopaedic Problems Of Heel Region.




Heel Cushion (Silicone)-D4R-522

Sizes Available | S.m.l

Application | To Provide Cushioning Effect To Heel Against The Painful Un-operatable Overgrowth And Development Of Calcaneum Spur (calcaneal Spur)

Features and Benefit |

- Heel Pad Made Of Soft Non-irritant 100% Medical Grade Silicone
- Scientific Design Provides Comfortable Cushioning And Adjusts Automatically As Per The Pressure On Toe Region Which Gives Relief
- From Pain To The User
- Durable And Long Life
- Easy To Clean And Maintain
- Use : Useful In Common Orthopedic Problems Of Heel Region



Heel Cushion (Gel)-D4R-523

Sizes Available | S/m/l

Application | To Provide Cushioning Effect To Heel Against The Painful Un-operatable Over Growth And Development Of Calcaneum Spur (calcaneal Spur)

Features and Benefits |

- Heel Pad Made Of Gel.
- Scientific Design Provides Comfortable Cushioning And Adjusts Automatically As Per The Pressure On Toe Region Which Gives Relief From Pain To The User.
- Durable And Long Life.
- Easy To Clean And Maintain.
- It Is Useful In Common Orthopedic Problems Of Heel Region.



Arch Pad (Gel)-D4R-524

Sizes Available | S.m.l.

Application | To Provide Additional Arch Support To Maintain Proper Foot Position.

Features and Benefits: |

- Special Gel Supports Fallen Arches Weak Arches
- High Arches.
- Useful In Common Orthopaedic Problems Of Arch Region.



Full Length Insole (Silicone) -D4R-525

Sizes Available | S.m.l

Application | Made Of Soft Non-irritant 100% Medical Grade Silicone For Complete Heel To Toe Protection And Comfort

Features And Benefits |

- Absorbs Maximum Shock Impact On Heel To Toe Region
- Effective In Foot As Well As Back Pain
- Easy To Use
- Long Life



Full Length Insole (Gel)-D4R-526

Sizes Available | S.m.l

Application | Made Of Gel For Complete Heel To Toe Protection And Comfort

Features and Benefits |

- Absorbs Maximum Shock Impact On Heel To Toe Region
- Effective In Foot As Well As Back Pain
- Easy To Use
- Long Life




¾ th Length Insole -D4R-527

Sizes Available | S.m.l

Application | ¾ Th Lenth Insole Protection And Comfort Of Heel

Features and Benefits |

- Absorbs Maximum Shock Impact On Heel To Toe Region
- Effective In Foot As Well As Back Pain
- Easy To Use
- Long Life



Toe Separator-D4R-528

Sizes Available | Universal

Application | To Relieve Friction Between Overlapping And Pinch Toes.

Features and Benefits |

- Scientific Design Helps In Spreading And Aligning The Toes Along With Absorption Of Pressure.
- Easy To Use Inside Shoes.
- It Is Useful In Pain And Soreness.



Cervical Traction Kit Sleeping Complete Set-D4R-601

Sizes Available | Universal

Application | To Relieve Pressure By Stretching The Musculature (musculature) And Connective Tissues Of Cervical Spine

Features and Benefits |

- Scientific Design Works Wonderfully In Sleeping Position
- Controlled Traction Can Be Achieved For Longer Duration
- Reduces The Pressure On The Anterior Vertebral Bodies By Expanding The Space Between Intervertebral Joints
- It Is Useful For Professional As Well As Personal Use Related To Cervical Problems




Cervical Traction Kit Sitting Complete Set -D4R-602

Sizes Available | Universal

Application | To Relieve Pain Associated With Muscle Contraction Due To Tension Or Surgery

Features and Benefits |

- Logitudinal (longitudinal) Stretching In Sitting Position By Traction Effectively Relieves Pain Associated With Cervical Spondylosis.
- Scientific Design Works Wonderfully In Sitting Position.
- Controlled Traction Can Be Achieved For Longer Duration
- It Is Useful For Professional As Well As Personal Use Related To Cervical Spine Dislocation



Foot/Skin Traction Kit(S)-D4R-701

Sizes Available | S.m.l

Application | To Provide Traction To Limbs After Injury And Surgery.

Features and Benefits |

- Provides Maximum Aeration To Skin
- Creep Bandage And Foam Padding Provides Ultimate Comfort To The User
- Controlled Traction Can Be Achieved For Longer Duration
- It Is Useful In Relieving Pain In Back And Knee Region.



Pelvic Traction Belt--D4R-702

Sizes Available | 28" To 44"

Application | Traction Arrangement To Control Back Pain

Features And Benefits |

- With Heavy Duty Straps And Strong Metal Rings Belt Is Directly Connectable To Traction Kit.
- Easy To Apply And Remove
- Gives Controlled And Uniform Traction
- It Is Useful For Traction In Lumbar Thoracic (Thoracic) And Sacral Area Of Vertebral (vertebral) Column.

IMAGE COMING SOON

Ankle Traction Belt -D4R-703

Sizes Available | S/m/l

Sizes Available | Universal

Application | To Help In Traction Of Ankle Lower Leg Section

Features and Benefits |

- With Heavy Duty Straps And Strong Rings Belt Is Directly Connectable To Traction Kit.
- Easy To Apply And Remove
- Gives Controlled And Uniform Traction.
- It Is Useful For Traction In Ankle Region

IMAGE COMING SOON ...

Leg traction belt -D4R-704

Sizes Available | S.m.l

Application | To Support The Traction Arrangement For Lower Back Pain, Knee Region

Features and Benefits |

- Special Nonslip Foam Lined Stirrup Makes It Useful For Sensitive Skin
- With Heavy Duty Straps And Strong Metal Rings Belt Is Directly Connectable To Traction Kit.
- Easy To Apply And Remove
- Gives Controlled And Uniform Traction
- It Is Useful In Lower Back And Knee Related Pains

IMAGE COMING SOON...

Pelvic Traction Spreader Bar -D4R-705

Sizes Available | Universal

Application | To Connect Pelvic Traction Belt To Traction Weights.

Features and Benefits |

- Strong And Sturdy Construction
- Easy To Apply



IMAGE COMING SOON

Traction SAND BAG-D4R-706

Sizes Available | Universal

Application | To Connect All Traction Belt To Traction Weights.

Features and Benefits |

- Easy To Apply



IMAGE COMING SOON

Pelvic Traction Complete Set - D4R-707

Sizes Available | 28" To 44"

Application | Traction Arrangement
To Control Back Pain

Features and Benefits |

- With Heavy Duty Straps And Strong Metal Rings Belt Is Directly Connectable To Traction Kit. Pelvic Traction Spreader Bar, Traction Pulley & Sand Bag
- Easy To Apply And Remove
- Gives Controlled And Uniform Traction
- It Is Useful For Traction In Lumbar Thoracic (thoracic) And Sacral Area Of Vertibral (vertebral) Column.



Ankle Traction Complete set -D4R-708

Sizes Available | Universal

Application | To Help In Traction Of
Ankle Lower Leg Section

Features and Benefits |

- With Heavy Duty Straps And Strong Rings Belt Is Directly Connectable To Traction Kit, Traction Pulley & Sand Bag
- Easy To Apply And Remove
- Gives Controlled And Uniform Traction.
- It Is Useful For Traction In Ankle Region.



Leg traction Complete set -D4R-709

Sizes Available | S.m.l

Application | To Support The Traction Arrangement For Lower Back Pain, Knee Region

Features and Benefits |

- Special Nonslip Foam Lined Stirrup Makes It Useful For Sensitive Skin
- With Heavy Duty Straps And Strong Metal Rings Belt Is Directly Connectable To Traction Kit, Traction Pulley & Sand Bag
- Easy To Apply And Remove
- Gives Controlled And Uniform Traction
- It Is Useful In Lower Back And Knee Related Pains



Traction Pulley -D4R-710

Sizes Available | Universal

Application | To Connect Traction Belt To Traction Weights.

Features and Benefit |

- Strong And Sturdy Construction
- Easy To Apply

IMAGE COMING SOON

Soft Roll (10 cm) -D4R-801

Sizes Available | 10 Cm & 15 Cm

Use | For Plastering



Cotton Crepe Bandage B.P with Fast Edges (6 cm)-D4R-802

Sizes Available | 6 Cm / 8 Cm / 10 Cm / 15 Cm

Applications | In The Treatment Of Sprains And Dislocated Painful Joints.

Features And Benefits |

- Excellent Elasticity Provides Uniform Gripping And Compression
- Easy To Apply And Remove.
- Multiple Use In Problems Related To Sprains And Strains



Rubber Walking Heel-D4R-803

Size | Universal
Applications | Used With Weight Bearing Foot Cast

Features & Benefits |

- Moulded Out Of Light Weight Durable & Resilient Polymer
- Durable Rubber With Graduations To Provides Resistance Against Slipping
- Arch Provided To Give Proper Stability



Posture Aid-D4R-804

Size | S/m/l/xl

Applications | to Correct Posture For Patients With Stooped Shoulder Holds Shoulder Upright & Keep In Correct Posture With Chest Out Position.

Features And Benefits |

- Comfortable Cotton Based Fabric Provided With Back Support & Velcro Straps.



Hyper Extension Brace/Ash Brace-D4R-805

Sizes Available | Universal

Application | To Perfectly Stabilize The Vertical Column In A Hyper-Lordotic Position

Features And Benefits |

- Mart & Anatomical Design Provides Extension Of The Spine
- Light Weight
- Easy To Apply And Remove
- It Is Used In Tuberculosis Of Dorso Lumbar Spine, Complications In Thoracolumbar Junction, post Operative Spine Immobilization Etc.



Cast Shoe (Dlx)-D4R-806

Sizes Available | S.m.l.xl

Application | Protects The Plaster Cast Of The Foot From Soiling, Wetting & Mechanical Damage.

Features And Benefits |

- Early Ambulation, Better Balance & Perfect Gait
- Light Weight, Easy To Wear
- Protects From Wear And Tear And Damage To Foot Cast
- It Is Useful In Excessive Wear And Tear Of Plaster Cast And To Reduce Impact On Injured Foot



Cuff & Collar-D4R-807

Sizes Available | Universal

Application | To Support Wrist And Elbow Without Obscuring The Limb Under Treatment

Features And Benefits |

- Made Out Of P.U Foam Covered With Cotton Stockinet.
- Multi-utility, Versatile & Comfortable Sling
- It Is Useful To Provide Support To Injured Wrist Or Elbow.



Universal Shoulder Immobilizer Ordinary -D4R-808

Sizes Available | S.m.l.xl

Application | To Immobilize And Support Dislocated Shoulder

Features And Benefits |

- Made Out Of Durable Cotton Fabric
- Adjustable Hook & Loop For Desired Position
- Two Piece Design Helps In Easy To Apply And Remove
- Secures Humerus And Forearm To Torso
- It Is Useful In Fracture Of Clavical Bone , Shoulder Joint Where Immobilization In Highly Needed.



Universal Shoulder-D4R-809

Immobilizer Deluxe

Sizes Available | S.m.l.xl

Application | To Immobilize And Support Dislocated Shoulder

Features And Benefits |

- Made Out Of Durable Three Layered Foaming Fabric & Wrinkle Free
- Adjustable Hook & Loop For Desired Position
- Two Piece Design Helps In Easy To Apply And Remove
- Secures Humerus And Forearm To Torso
- It Is Useful In Fracture Of Clavical Bone , Shoulder Joint Where immobilization In Highly Needed.



Elastic Shoulder Immobilizer Ring Type-D4R-810

Sizes Available | S.m.l.xl

Application | To Immobilize And Support Dislocated Shoulder

Features And Benefits |

- Elastic Belt & Buckled Grips Immobilize The Injured Arm, Smart, Sleek, Well Ventilated & Comfortable
- To Fit Either Left / Right Arm
- It Is Useful In Fracture Of Clavical Bone , Shoulder Joint Where Immobilization In Highly Needed.



Adjustable Pouch Arm Sling Ordinary-D4R-811

Sizes Available | X-s / S / M / L / XI

Application | To Provide Optimum Support To Injured, Painful Arm By Reducing Neck And Shoulder Pressure

Features And Benefits |

- Special Cotton Fabric Offers Comfortable And Soft Feeling To Injured Arm Portion
- Easily Operative Buckle Helps In Easy Application And Closes Pouch Arm Sling
- Thumb Loop Inside For Desired Hand Rest
- Adjustable Special Foam Pad In Shoulder Takes Load Of Complete Arm
- Very Useful In Countries With Hot And Humid Weather
- It Is Useful In Fractures And Dislocations, Sprains And Strains, To Provide Support And Immobilize Post Operative Arm



Adjustable Pouch Arm Sling Deluxe-D4R-812

Sizes Available | X-s / S / M / L / XI

Application v| To Provide Optimum Support To Injured, Painful Arm By Reducing Neck And Shoulder Pressure

Features And Benefits |

- Special Foaming Fabric Offers Comfortable And Soft Feeling To Injured Arm Portion
- Easily Operative Buckle Helps In Easy Application And Closes Pouch Arm Sling.
- Thumb Loop Inside For Desired Hand Rest
- Adjustable Special Foam Pad In Shoulder Takes Load Of Complete Arm
- Very Useful In Countries With Hot And Humid Weather
- It Is Useful In Fractures And Dislocations, Sprains And Strains, To Provide Support And Immobilize Post Operative Arm



Arm Sling Strap-D4R-813

Sizes Available | S.m.l.xl

Application | For Holding Injured Arm, Shoulder And Wrist

Features And Benefits |

- Heavy Duty Straps Assures Sturdy Design
- Tubular Neck Pad Increase Comfort Level
- Slide Buckle Provides Maximum Adjustability
- It Is Used In Immobilization Of Hand, Wrist & Shoulder



Meek Clavicle Brace-D4R-814

Sizes Available | S.m.l.xl

Application | In Fracture Of Clavical Bone

Features And Benefits |

- Magic Figure 8 Shaped Formed With Straps And Lock Applies Pressure On Clavical Bone And Hold It Together soft Padding Of Special Foam Covered With Special Fabric Gives Ultimate Comfort To The User.
- Scientific Design Of Straps Allows Strong Fitting Of Brace Around Clavicular Region.
- It Is Useful In Fracture Of Clavical Bone, Stopping Shoulder And Hump Back



Rib Belt Male-D4R-815

Sizes Available | 28" To 44"

Application | To Limit The Painful Chest Expansion In Case Of Rib Fracture Or Dislocation Of Ribs

Features And Benefits |

- Light Weight And Comfortable
- Easy To Apply And Remove
- Special Elastic Material Allows Equal Compression Overall Rib Region, At The Same Time Allows Normal Respiration
- It Is Useful In Protection Of Surgical Incision Of Chest And Rib Related Problems



Rib Belt Female -D4R-816

Sizes Available | 28" To 44"

Application | To Limit The Painful Chest Expansion In Case Of Rib Fracture Or Dislocation Of Ribs

Features And Benefits |

- Light Weight And Comfortable
- Easy To Apply And Remove
- Special Elastic Material Allows Equal Compression Overall Rib Region, At The Same Time Allows Normal Respiration
- Anterior Cut Down Design Provides Comfortable Fitting For Female Patients
- It Is Useful In Protection Of Surgical Incision Of Chest And Rib Related Problems



Exercise Ball Ortho & Neuro (Gel) -D4R-817

Sizes Available | Universal

Application | Hand Exerciser To Achieve Grip Strength

Features And Benefits |

- Made From Safe & Insert P.U With Resilience And Softness
- Gives Controlled Isotonic Resistance.
- Provides Exercise Plan For Fingers, Hand And Forearm.



FOAM Splint with foot support (fiber)-D4R-818

Sizes Available | Universal

Application | Used As Emergency Splint

Features And Benefits |

- Light Weight And Made Of Fiber.
- Well Padded, Unbreakable, Washable & Reuseable
- Velcro Straps Provided To Give Proper Fitting



IMAGE COMING SOON....

Arm Tourniquets Black -D4R-901

Size | 8" / 10" / 12" / 14"

Application | To Control Or Stop Blood Flow For Temporary Period

Features And Benefits |

- Easy To Apply



Arm Tourniquets COLOURED -D4R-902

Size | 8" / 10" / 12" / 14"

Application | To Control Or Stop Blood Flow For Temporary Period

Features And Benefits |

- Easy To Apply



Thomas Splint -D4R-903

Sizes Available | X-s / S / M / L

Application | Leather Covered Swivel Ring, Tinned, Used For Fractures Of Upper & Lower Thigh.



Cramer's Wire Splint-D4R-904

Sizes Available |
1"x12", 1"x18", 1"x24"/2"x12"
, 2"x18", 2"x24" / 3"x12", 3"x18", 3"x24"
/ 4"x12", 4"x18", 4"x24"

Application | It Is Used As Emergency Splint Before Applying The Cast

Features And Benefits |

- Made Of Moldable Wire And Vinyl Coated For A Soft Pliable Surface Finish
- Light Weight And Malleable



Orthopaedic Stockinet Roll (5.0" x 10mts) (Rib Cloth) -D4R-905

Sizes Available |
2",2.5",3",3.5",4",5",6"

Application | Inside The Plaster

Features And Benefits |

- Made From Very Soft Non Absorbent Tubular Rib Weaved Cotton
- Easy To Apply



The Cuff Weight Belt-D4R-906

Sizes Available | ½ / 1 / 2 / 3 / 4 / 5
Kgs

Application | For Graded
Physiotherapy Exercises

Features And Benefits |

- Velvet Cover Gives Comfort To Patient
- Easy To Apply
- Weight Packets Sewn Vertically Bents Freely And Can Be Used On Wrist And Ankle



Orthopaedic Heating Pad -D4R-907

Sizes Available | Universal

Application | To Provide Relief From Pain For Arthritis Patient, Decrease Joint Stiffness

Features And Benefits |

- Three Different Levels Of Heating Controlled By Sliding Switch To Produce Appropriate Heat As Per Patients Requirement
- Special Outer Cotton Cover Keeps The Belt In Touch With Patient And Transfers The Maximum Heat To The Patient.
- Use: In Back, Thigh, Shoulder Pains



Dynamic Radial Palsy Splint -D4R-1001

Sizes Available | S.m.l.xl - Left & Right

Application | For Exercises Of Finger And Wrist

Features And Benefits |

- Wrist And Finger Tension Is Adjustable By Detaching Or Attaching The Springs
- Detachable Arrangements For Finger
- Very Light Weight And Soft Padded For Complete Comfort
- Waterproof, Unbreakable Design
- It Is Useful In Physiotherapy After Extensor Tendon Repair, Radial Nerve Palsy And Wrist Drop



Dynamic Radial Palsy Splint -D4R-1002

Sizes Available | S.m.l.xl - Left & Right

Application | For Exercises Of Wrist

Features And Benefits |

- Wrist Tension Is Adjustable By Detaching Or Attaching The Springs
- Very Light Weight And Soft Padded For Complete Comfort
- Waterproof, Unbreakable Design
- It Is Useful In Physiotherapy After Extensor Tendon Repair, Radial Nerve Palsy And Wrist Drop



Volkman's Turn Buckle Splint (Left & Right) -D4R-1003

Sizes Available | s.m.l.xl -left & Right

Application | To Correct Contractors Of Soft Tissues

Features And Benefits |

- Reverting On Plastic Material To Protect Against Undue Force While Correcting
- Turn Buckle Attached To Give Proper Handling Of The Splints
- Very Light Weight And Soft Padded For Complete Comfort
- Waterproof, Unbreakable Design
- It Is Useful In Faster Correction Of Contractors Of Soft Tissues



Elbow Turn Buckle Splint -D4R-1004

Sizes Available | S.m.l.xl

Application | To Correct The Contracture Of Elbow
Turn Buckle For Easy Contracture Release

- Plastic Material With Revetting Arrangement Protects Tissues Against Undue Force
- Very Light Weight And Soft Padded For Complete Comfort
- Waterproof, Unbreakable Design
- It Is Useful In Faster Correction Of Contractors Of Elbow



Knuckle Bender Splint-D4R-1005

Sizes Available | S, m, l

Application | To Flex The Stiff Meta Carpophalangeal Hand Joints

Features And Benefits |

- Scientifically Designed Circular Spring Maintains Constant Pressure
- Stainless Steel Springs Maintain The Constant Pressure
- Very Light Weight And Soft Padded For Complete Comfort
- Waterproof, Unbreakable Design



Humerus Fracture Brace-D4R-1006

Sizes Available | S,m,l,xl - Left & Right

Application | To Increase The Range Of Motion Of Fractured Shoulder And Elbow

Features And Benefits |

- Applied Directly On Fracture Like Claw Around Humerus
- Three Independent Levels Of Tension Control
- It Is Useful In Decreasing The Rehabilitation Time With Lessen Muscular Atrophy



Finger Exerciser Table Model -D4R-1007

Sizes Available | Universal

Application | To Provide Strengthening Exercise To Affected Fingers

Features And Benefits |

- Very Light Weight Splint With Base Stand
- Springs Provides Strength To Patient As Required
- It Is Useful To Give Exercise To Affected Fingers Due To Surgery Or Paralysis



Finger Exerciser (Steel)-D4R-1008

Sizes Available | Universal

Application | to Provide Strengthening Exercise To Affected Fingers

Features And Benefits |

- Very Light Weight Splint Made Of Steel
- Springs Provides Strength To Patient As Required
- It Is Useful To Give Exercise To Affected Fingers Due To Surgery Or Paralysis



Finger Exerciser (Plastic) -D4R-1009

Sizes Available | Universal

Application | To Provide Strengthening Exercise To Affected Fingers

Features And Benefits |

- Very Light Weight Splint Made Of Plastic
- Springs Provides Strength To Patient As Required
- It Is Useful To Give Exercise To Affected Fingers Due To Surgery Or Paralysis



'U' Splint Short Arm -D4R-1010

Sizes Available | Universal

Application | Splinting Of
Dislocated Or Fractured Upper Or
Lower Limb In Emergency

Features And Benefits |

- Malleable, Semi-rigid, Light Metal Splint Provides Support in Anatomical Position
- Provide Immobilization To Injured Part



'U' Splint Long Arm-D4R-1011

Sizes Available | Universal

Application | Splinting Of
Dislocated Or Fractured Upper Or
Lower Limb In Emergency

Features And Benefits |

- Malleable, Semi-rigid, Light Metal Splint Provides Support In Anatomical Position
- Provide Immobilization To Injured Part



'U' Splint Short Leg -D4R-1012

Sizes Available | Universal

Application | Splinting Of
Dislocated Or Fractured Upper Or
Lower Limb In Emergency

Features And Benefits |

- Malleable, Semi-rigid, Light Metal Splint Provides Support In Anatomical Position
- Provide Immobilization To Injured Part



'U' Splint Long Leg -D4R-1013

Sizes Available | Universal

Application | Splinting Of
Dislocated Or Fractured Upper Or
Lower Limb In Emergency

Features And Benefits |

- Malleable, Semi-rigid, Light Metal Splint Provides Support In Anatomical Position
- Provide Immobilization To Injured Part



Tibial Functional Brace-D4R-1014

Sizes Available | 1/2/3/4/5- Left & Right

Application | For Fractures Of Tibia

Features And Benefits |

- Four Part Splints When Joined Together By Velcro Straps Forms A Cylinder To Fit Properly To Patients Leg
- An Ideal Brace For Treating Fracture Of Femur.
- More Comfortable Than Plaster Cast



CTEV Splint With Padding -D4R-1015

Sizes Available | 0/1/2/3/4/5 - Left & Right

Application | For Correcting Mild Ctev Deformities

Features And Benefits |

- Washable, Light Weight Design
- Smooth Padding Provides Extra Comfort To Infants
- Use:useful To Avoid Repeated Plastering And Surgery In Infants For Ctev Deformities



CTEV Splint Without Padding-D4R-1016

Sizes Available | 0/1/2/3/4/5- Left & Right

Application | For Correcting Mild Ctev Deformities

Features And Benefits |

- Washable, Light Weight Design
- Use:useful To Avoid Repeated Plastering And Surgery In Infants For Ctev Deformities



CTEV D.B. Splint-D4R-1017

Sizes Available | 1.2.3.4.5

Application | To Prevent Lateral Dropping Of Foot During Sleep

Features And Benefits |

- Unbreakable, Washable And Reusable Design Light Weight
- Flat Base Provides Stability To Foot Preventing Lateral Dropping
- It Useful As An Alternative To Wire Splint Or Pop Slab In Precast Condition.



Foot Ankle Leg Splint-D4R-1018

Sizes Available | S/m/l/xl- Left/
right

Application | It Is Used To Prevent
Lateral Dropping Of Foot During
Sleep.

Features And Benefits |

- Flat Base To Give Stability To The
Foot Preventing Lateral Dropping
Of Foot.
- Well Padded Unbreakable,
Washable, Useable, Economic &
Light Weight Splint

IMAGE COMING SOON

Foot Drop Splint-D4R-1019

Sizes Available | S/m/l/xl- Left/
right

Application | For Stabilization Of
Ankle And Foot In Foot Drop
Condition

Features And Benefits |

- Premolded, Semi-rigid And
Thermoplastic Design
- Light Weight Padded Splint To Give
Comfortable Wear
- Velcro Straps Provided To Give
Proper Fitting
- Used As Night Splint



Foot Drop Splint-Light-D4R-1020

Sizes Available | S/m/l/xl-Left/
right

Application | For Stabilization Of
Ankle And Foot In Foot Drop
Condition

Features And Benefits |

- Premolded, Semi-rigid And Thermoplastic Design
- Light Weight Padded Splint To Give Comfortable Wear
- Velcro Straps Provided To Give Proper Fitting
- It Is Useful While Walking.

