



www.kiransurgicalsystems.com
www.kiranmedicareequipments.com
www.kisusy.com, www.d4h.co.in
www.d4s.co.in



- **O**rthopedic
- **G**ynecology
- **T**ungsten Carbide
- **T**itanium
- **B**ipolar Instruments


Orthopedic



Orthopedic Instruments



Bone Curette with Fibre Handle



Self Retaining Spinal Retractor



Capner's Gouge Straight / Curved / Deep Curved



Durer Elevator




Ring Curette



Spreader



Polyaxial Screw Driver



Rod Holder



Pedical Elevator



Compressor



Bone Tap



Pelvic Reduction Forceps



Bone Awl with Eye



Straight Pelvic Reduction Forceps 190mm




Wire Tightner Cum Twister



Straight Pelvic Reduction Forceps 400mm

Orthopedic Instruments



K Wire Tractor



K Wire Introducer




K Wire Bender



Steinman Pin Introducer



Wire Passer



Wire Cutter 11"



Long Wire For Twisting Wire



Wire Bending Forcep




Coventry Staple Impactor



Fibula Forceps-Straight



Self Centering Bone Holding Forcep



Lane Bone Holding Forceps with Ratchet (200, 250, 300 mm)



Lane Bone Holding Forcep without Ratchet (200, 250, 300 mm)



Patella Bone Holding Forceps, Three Prongs




Burn's Bone Holding Forcep



Lowman Bone Clamp



Bone Holding Forceps Double Action-Straight



Bone Holding Forcep Double Action-Curved



Tuder Edward Bone Cutting Forcep



Tuder Edward Bone Cutting Forcep



Bone Nibbler Forcep (Double Action) Straight




Bone Nibbler Forcep (Double Action) Angular



Plate Holding Forcep



Forgussion Bone Holding Forcep



Bone & Plate Holding Forcep-Mani Type

Orthopedic Instruments



Chisels With Fibre Handle



Osteotomes With Fibre Handle-straight



Osteotomes With Fibre Handle-curved



Gouge With Fibre Handle-straight



Gouge With Fibre Handle-curved



Lagenback Retractor



Volkman Retractor



Right Angle Retractor



Periosteal Elevator-curved



Periosteal Elevator-straight



Self Retaining Retractor



Skin Retractor



Bone Curette-double Ended



Bone Curette-with Fibre Handle



Bone Hook-small/medium/large



Amputation Knife



Mastoid Retractor



Plier



Wire Cutter For 5mm Wire

Orthopedic Instruments



Hudson Brace With 5 Burs



Muller Compression Clamp



Charnley's Compression Clamp



C Type Skull With Two Burs



Gigli Saw Introducer



Skin Grafting Knife-small



Skin Grafting Knife-medium



Skin Grafting Knife-large



Large Distractor for Tibia and Femur

Orthopedic Instruments



Austin Moore Rasp



Rasp With Tommy Bar With Biopolar Prothesis



Thompson Rasp



Murphy Bone Skid



Extractor With Two Hook



Spare Hook For Extractor



Charnley Retractor



T Handle With Quick Coupling



Wrench For Dhs Screw



Dhs Triple Reamer



Dhs Angle Guide-135 Degree



Long Coupling Screw For Dhs Screw Removal

Orthopedic Instruments



Dcs Triple Reramer



Variable Angle Guide



Bone Plate Bender Roller Type



Plate Bender For Small Fragment Plate



Bone Hammer



Bone Hammer-nylon Faced



Large Heaxagonal Screw Driver



Drill Tap With T Handle



Cannulated Drill Bit 2.5mm To 5mm



Battery Operated Drill Hand Piece



Battery Operated Saw Hand Piece



Plaster Cutter

Electric Cutting Reaming & Drilling System



REAMING HAND PIECE

- ◆ Quick Coupling for Flexible Remarks
- ◆ Light Wight (600gms. approx)
- ◆ Maximum Speed 400 RPM
- ◆ 5.5 mm cannulated
- ◆ Autoclavable, Hard Anodized for long Life.



SAGGITAL SAW HAND PIECE

- ◆ Quick Coupling for Flexible Remarks
- ◆ Five types of blades
- ◆ Light Wight (700gms. approx)
- ◆ Autoclavable, Hard Anodized for long Life.



CANNULATED DRILL HAND PIECE

- ◆ Steel Drill Chuck (0-1/4")
- ◆ Light Wight (600gms. approx)
- ◆ Maximum Speed 1200 RPM
- ◆ 5.5 mm cannulated
- ◆ Autoclavable, Hard Anodized for long Life.

DRIVING UNIT

- ◆ 220 Volts/5 Amp. A.C. supply.
- ◆ Heavy Duty motor with completely enclosed
- ◆ Foot Control for On / Off and Speed

STAND

- ◆ Castors for Mobility
- ◆ Can be dismantled & packed into a Suitcase
- ◆ Auto cutoff switch a special feature



Shaver System



Specification :

- Speed range 900 to 4600 RPM
- Online setting allows online increase and decrease of set speed.
- Variable modes FORWARD, REVERSE & OSCILLATING
- Shaver hand piece with cord fully autoclavable
- Attached blades easy blade fitting and removing in the hand piece.
- Foot switch: Foot switch for easy foot control in variable blade modes (FORWARD, REVERSE & OSCILLATING).
- Shaver hand piece with cord fully autoclavable
- Suction control inbuilt in hand piece.
- Power supply 230V (AC) / 50Hz ± 10%.
- In build voltage control rectifier
- Inbuilt safety fuse
- Electric shock-proof handle
- Dimensions: 350mm (l) x 275mm (w) x 140mm(h)
- Weigh: 5 kgs.

Pulse Lavage System



Specification

- Irrigation pressure can be electronically controlled.
- Easy loading of tubes in the machine
- Autoclavable and reusable tubes of medical grade silicon.
- Complete sterility of system is maintained.
- Foot control operated
- Fully Stainless Steel unit for long life.
- Efficient peristaltic pump
- Simple on/off switch.
- Three types of irrigation possible
- Single stream tip For wider stream & general irrigation.
- Shower spray For cancellous bone preparation, trauma, soft tissues debridement & general irrigation.
- Shower For dual purpose irrigation & suction along with special cup type anti-splash shield.
- Side hole cannula for canal wash
- Multiple bottles with constant force
- In built voltage control rectifier
- Inbuilt safety fuse
- Power supply 230V (AC) / 50Mz ± 10%.
- Dimensions: 190 (l) x 230mm (w) x 190mm (h)
- Weight : 6.5 Kg. Approx



Digital Tourniquet



Specification :


Cuff pressure range	: 10 to 450 mmHg
Pressure regulation	: Better than ± 10 mmHg of set point (Maintains set pressure under leak conditions)
Online setting	: Increase and decrease in pressure settings
Timer	: Range from 9 hours to 59 minutes (0 : 05 hours to 9:59 hours)
Alarm	: Audible Buzzer on set time equaling elapsed value, in intervals of 1 min, after alarm time
Backup	: 3 hours battery backup (in case of full charge)
Switch Over	: Automatic switch Over from mains to battery of vice versa
Power	: 230V (AC) / 50Hz \pm 10%
Battery Charger	: Constant Voltage, Current Limited & Auto Cutoff
Digital display	: Digital display of set pressure, actual pressure, elapsed time & set time
Cuffs	: Different sizes of five cuffs -washable & easy fitting (Large, Big, Medium, Small & Pediatric)
Enclosure	: ABS Plastic, with special Handle
Dimensions	: 215 mm (W) X 140mm (H) X 230mm (L)
Weight	: 3.9 Kgs




Gynecology



Gynecology Instruments



Cusco Vaginal Speculum



Cusco Vaginal Speculum with Fiber Optic Light Carrier



Cusco Vaginal Speculum Grave



Uterine Tenaculum Single Tooth Forceps




Uterine Tenaculum Double Teeth Forceps



Uterine Elevating Forceps



Cervical Biopsy Punch Forceps



(Balfour) Abdominal Self-Retaining Retractor



Pozzi Self-retaining Retractor With Three Blade



Midwifery Forceps Simpson



Wriggle Forceps

Tungsten Carbide



Tungsten Carbide Instruments



Metzenbaum Scissors



Mayo Scissors



Iris Scissors



Knapp Scissors



Pott's Angled Scissors



Kilner Needle Holder Fine



Senning Needle Holder



Gillies Needle Holder With Scissor



Needle Holder Fine Tungsten Carbide



Crile Wood Debaquey Needle Holder



Bayonet Packing Forceps



Bab Cock Forceps


Tungsten Carbide Instruments



Dissecting Forceps



Fine Adson Forceps



Allis Tissue Grasping Forceps



Cottle Rasp



Glabella Rasp



Nasal Rasp With
Diamond Cutting-polishing



Top Cutter
Cutting 2.0mm Wire



Bone Cutter
Cutting 2.5mm Wire



Wire Cutter
Cutting 2.0mm Wire

Tungsten Carbide Instruments

Tungsten Carbide



Wire Cutter Fine



Wire Holding Forceps



Meriam Forceps




Wangenstein Forceps



Pin & Wire Puller



Listen Bone Cutter
Straight 8"



Hard-wire Cutter
Cutting 4mm Wire

Titanium



Titanium Instruments



Adson Forceps
Plain - 5"



Adson Forceps
1x2 Tooth - 5"



Dissecting Forceps
Plain



Dissecting Forceps
Tooth



Mosquito Forceps



Micro Spring Needle Holder
Flat With Out Lock Straight



Micro Needle Holder
Bayonet Type Straight

Titanium Instruments




Debaquey Forceps
Straight



Micro Forceps Straight Tip Flat
Handel, Bayonet Shape



Micro Fine
Needle Holder



Mayo Hegars
Needle Holders



Lambert - Kay Aorta
Clamps



Castaneda Clamp
Cooley Jaws

Titanium Instruments



Debaquey Stainsky
Vascular Clamps



Peripheral Clamps
Curved Debaquey Jaws



Debaquey Vascular
Clamps
Curved



Micro Clamp
Angled Debaquey Jaws



Finochietto
Rib Spreader
Child



Haight Finochietto
Rib Spreader
Medium

General Surgery



General Surgery Instrument

)" "8/%\$#". i



Spencer Wells Artery Forceps



Mosquito Forceps




Kocher's Artery Forceps



Dissecting Forceps (Plain)



Dissecting Forceps (Tooth)



Towel Forceps



Towel Clip Cross-action



B.P. Handle No. 3



Bozemann Needle Holders



B.P. Handle No. 4



B.P. Handle No. 7



Mayo Hegars Needle Holders



Kilner Needle Holders Fine



Dressing Scissors Straight



Dressing Scissors Curved

General Surgery Instrument



Iris Scissors



Iris Scissors Super Cut



Knapp Scissors



Knapp Scissors Super Cut



Stitch Ribbon Scissors



Stitch Ribbon Scissors Super Cut



Metzenbaum Scissors



Metzenbaum Scissors Super Cut



Mayo Scissors



Mayo Scissors Super Cut




Heath Suture Forceps



Listre Bandage Forceps



Pott's Angled Forceps



Allis Tissue Forceps



Lanes Tissue Forceps

General Surgery Instrument



Bab Cock Grasping Forceps



Lahey's Forceps



Mixture Ligature Right Angled Forceps



Hysterectomy Maingots Forceps



Hysterectomy Debakey Forceps



Mayo Safety Pin



Mayo Safety Pin



Kilner Skin Retractor



Double Hook Retractor



Czerny Retractor



Langenbeck Retractor



Langenbeck Retractor



Langenbeck Retractor



London Vaginal Retractor



Kelly's Retractor

General Surgery Instrument



Deavers Retractor



Deavers Retractor Adult



Deavers Retractor With Fiber Optic Light Carrier Adult



Volkman Retractor Sharp



Volkman Retractor Blunt



Morris Retractor



Morris Retractor With Fiber Optic Light Carrier




Doyen Retractor



Doyen Retractor With Fiber Optic Light Carriers



Wullstein Self Retaining Retractor



(Balfour) Abdominal Self-Retaining Retractor



Rendalls Renal Forceps Set of 4



Kelly's Rectal Speculum

General Surgery Instrument



Mathieu Rectal Speculum



Haemorrhoidal Ligating Gun



Piles Band Rubber



Doyen Intestinal Clamp



Green Armytage Forceps



Sponge Holding Forceps



Uterine Vulsellum Forceps



Uterine Tenaculum Single Tooth Forceps




Uterine Tenaculum Double Teeth Forceps




Uterine Elevating Forceps



Interior Vaginal Wall Retractor



Sims Speculum



Sims Speculum With Fiber Optic Light Carrier

General Surgery Instrument



Cusco Vaginal Speculum



Cusco Vaginal Speculum With Fiber Optic Light Carrier



Cusco Vaginal Speculum Grave



Kristeller Vaginal Speculum



Hegar Dilators Set of 8



Hegar Dilators Graduated Fine 1/2mm Set of 10

General Surgery Instrument



Uterine Manipulator



Yankauer Tonsil Suction Tube



Uterine Manipulator Ring Handle



Uterine Manipulator Ring Handle Double Action




Female Catheter



Mymo-screw



M.T.P. Cannula Metal Set Of 6




Episiotomy Scissors



Cervical Biopsy Punch Forceps



Midwifery Forceps Simpson



Midwifery Forceps with A.T.



Debaquey Bull-Dog Clamps



Cooley Back Vessel Clamps

General Surgery Instrument (Laryngoscope)



Fiber Optic Flexible Laryngoscope Set of 3



Fiber Optic Laryngoscope Set of 3



Magill Laryngoscope With Blades Set of 4 Blades



Lamp For Laryngoscope



Mcintosh Laryngoscope With Blades Set of 4 Blades

Bipolar Instrument



Bipolar Coagulation Forceps Bayonet Type



Bipolar Coagulation Forceps Bayonet Type With Two Pins



Bipolar Coagulation Forceps Bayonet Type With Suction




Bipolar Coagulation Forceps Bayonet Type With Irrigation



Bipolar Coagulation Forceps (straight)



Bipolar Coagulation Forceps (straight) With Two Pins



Bipolar Cable Two Pin



Bipolar Cable Round Pin



Bipolar Cable Flat Pin



Bipolar Cable Two Pin



Artery Forceps Straight



Artery Forceps Curved



Turner Warwick Scissors Straight



Turner Warwick Scissors Curved



Lateral Retractor Wide Open



Pederson Vaginal Speculum With Smoke Evacuation Tube



Greaves Vaginal Speculum With Smoke Evacuation Tube

Manufacturing Unit



www.kiransurgicalsystems.com www.kiranmedicareequipments.com
www.kisusy.com, www.d4h.co.in, www.d4s.co.in



KIRAN TECHNO SERVICES
PRIVATE LIMITED



Delhi : A - 39, Street No.5, Chander Industrial Estate, I.P Extension, Delhi - 110 092

Gujarat : 308/4, GIDC Estate, Makarpura Industrial Estate, Baroda – 390 010,

Andhra Pradesh : No. 6-5-351/5, Astal Heights, Road No. 1, Banjara Hills, Hyderabad – 500 034

Tamil Nadu : No.33-B4, Sidco Industrial Estate Phase - II, Sundarapuram, Coimbatore - 641 021

Kerala : 333/9, Para Road, Near NIDA, Kanjikode, Palakkad - 678 621

Helpline : +91 9751725000, +91 9751625000, + 91 9751425000, + 91 9751215000, + 91 9751245000

Email : kirangroup108@gmail.com, kiranindia108@gmail.com, kiranbiz108@gmail.com

Website : www.kirangroupofcompanies.com, www.kiran.biz, www.kiranshop.in, www.kiranindia.in, www.kisusy.com
www.kiransurgicalsystems.com, www.kiranmedicareequipments.com, www.kiranpump.com, www.d4h.co.in
www.kiransolarsystems.com, www.kiranled.com, www.kiranincinerationsystems.com, www.d4s.co.in
www.kiranequipments.com, www.kirancompressors.com, www.kiranappliances.com
www.kiranairventilators.com, www.kiranfrpproducts.com, www.kiranpebsystems.com,
www.kiranpanels.com, www.kiranfabssystems.com