

KIRAN ZERO GRAVITY COATS & BEDS

Special Beds for:

- Aged People
- Obese
- Bed Ridden Patients



ABOUT US

Kiran Beds started their venture, with a mission to be unique manufacturer of ultimate wooden electric beds, beds which can be easily operated by a remote. We had a vision of making something new that has never happened before and by following our vision; we created our excellent range of adjustable wooden beds along with memory foam mattresses.

These wooden beds are extremely famous amongst our precious customers for offering complete back and leg rest with just one touch of electric remote. These beds can change their back and leg position for better comfort – an ideal solution for patients and elder people.

We simply believe in offering sound sleep to each of our customer by providing them our comfortable adjustable wooden beds along with regular and memory foam mattresses, which automatically resulted into good health and that leads you to a good life. Hence we say

Sound Sleep = Good Health = Good Life...

www.kiransurgicalsystems.com

www.kiranmedicareequipments.com



- **BACKACHES**

- **BLOOD PRESSURE**

- **STRESS RELIEF**


- **ELDERLY PEOPLE**

- **BEDRIDDEN PATIENTS**

- **SPINE RELAXATION**



LEG REST POSITION



BACK REST POSITION



ZERO GRAVITY POSITION

Motorised Metal Semi Fowler bed



Motorised Semi Fowler bed wooden paneling



Motorised Fowler bed



Motorised Fowler bed wooden paneling



Kiran Zero Gravity Economy



Kiran Zero Gravity Confort



Kiran Zero Gravity Platina



Wooden Double Bed sunmaica salwood





Double Bed Teak wood



Motorised ICU bed wooden panelling



PLATINA



Model : Platina

Raw Material : Teak wood main Frame, plywood sides with imported rexene covering

Motors and Remote : Imported motors and wire remote

Warranty : 2 Years on motor and remote for manufacturing defects

Head Board/Foot Board : Only head board

Warranty for Wooden Parts : 15 Years warranty on teak wood parts only for manufacturing defects

Movements : Back and knee movements by wired remote control

Mattress Optional : Optional standard foam mattress Or memory and standard foam mattress combination at extra cost



COMFORT



Model : COMFORT

Raw Material : Complete Bed In Teak Wood with machine polish, no side covering

Motors and Remote : Imported motors and wire remote

Warranty : 2 Years on motor and remote for manufacturing defects

Head Board/Foot Board : Both head and foot board

Warranty for Wooden Parts : 15 Years warranty on teak wood parts only for manufacturing defects

Movements : Back and knee movements by wired remote control

Mattress Optional : Optional standard foam mattress Or memory and standard foam mattress combination at extra cost

ULTIMA



Model : ULTIMA

Raw Material : Sal wood frame with ply wood top, plywood sides covered with cloth

Motors and Remote : Imported motors and wire remote

Warranty : 2 Years on motor and remote for manufacturing defects

Head Board/Foot Board : Only head board


Warranty for Wooden Parts : 3 Years warranty on sal wood parts only for manufacturing defects

Movements : Back and knee movements by wired remote control

Mattress Optional : Optional standard foam mattress Or memory and standard foam mattress combination at extra cost



ECONOMY



Model : ECONOMY

Raw Material : Sal wood frame with ply wood top, plywood sides covered with Sunmica

Motors and Remote : Imported motors and wire remote

Warranty : 2 Years on motor and remote for manufacturing defects

Head Board/Foot Board : Only head board

Warranty for Wooden Parts : 3 Years warranty on sal wood parts only for manufacturing defects

Movements : Back and knee movements by wired remote control

Mattress Optional : Optional standard foam mattress Or memory and standard foam mattress combination at extra cost

DOUBLE BED



Model : DOUBLE BED

Raw material : Sal wood frame with ply wood top, plywood slides covered with Sunmica/Vineer (at extra cost). Combination of single motorized bed and single plain bed.

Motors and remote : Imported motors and wire remote

Warranty : 2 years on motor and remote for manufacturing defects

Head board/foot board : only head board

Warranty for wooden parts : 3 years warranty on sal wood parts only for manufacturing defects

Movements : back and knee movements by wired remote control

Mattress optional : optional standard foam mattress or memory and standard foam mattress combination at extra cost



- Standard size of all models (outer size without mattress) : l 84" x w 36" x h 18" (customized sizes are available at extra cost) double bed size : l 84" x w 60" x h 18 "
- No guarantee for colours of mattress cover
- Colours available for polish parts : black, brown , natural
- Laminet/sunmica/vinear colours : black , brown (can be customized with laminet company and code number details from customer at extra cost)

- Warning : bed ,not to be operated by child / mentally retarded/ handicapped person
- Warning: use of the bed is at sole risk of buyer. Company do not take any responsibility of any physical damages ,injuries to the person sleeping on bed or any other person , asistantwho comes in contact with moving parts of the bed.
- Warning: beds is made to sleep for single person upto weight of 90 kg. For heavy weight person additional strong design can be made at extra cost.
- Warning: do not sit on back, knee and leg folders of bed in operational movements. This can damage the parts of bed and cause injuries to the person



Delhi : A - 39, Street No.5, Chander Industrial Estate, I.P Extension, Delhi - 110 092
Gujarat : 308/4, GIDC Estate, Makarpura Industrial Estate, Baroda – 390 010,
Andhra Pradesh : No. 6-5-351/5, Astal Heights, Road No. 1, Banjara Hills, Hyderabad – 500 034
Tamil Nadu : No.33-B4, Sidco Industrial Estate Phase - II, Sundarapuram, Coimbatore - 641 021
Kerala : 333/9, Para Road, Near NIDA, Kanjikode, Palakkad - 678 621
Helpline : +91 9751725000, +91 9751625000, + 91 9751425000, + 91 9751215000, + 91 9751245000
Email : kirangroup108@gmail.com, kiranindia108@gmail.com, kiranbiz108@gmail.com
Website : www.kirangroupofcompanies.com, www.kiran.biz, www.kiranshop.in, www.kiranindia.in, www.kiransurgicalsystems.com, www.kiranmedicareequipments.com, www.kiranpump.com, www.kiransolarsystems.com, www.kiranled.com, www.kiranincinerationsystems.com www.kiranequipments.com, www.kirancompressors.com, www.kiranappliances.com www.kiranairventilators.com, www.kiranfrpproducts.com, www.kiranpebsystems.com, www.kiranpanels.com, www.kiranfabssystems.com