




SMALL FRAGMENTS - SAFETY LOCKING IMPLANTS





[YouTube www.TV.ortho.in/27882](https://www.youtube.com/watch?v=www.TV.ortho.in/27882)



PROXIMAL HUMERUS SAFETY LOCK PLATE 3.5

No. of Holes	Length in mm	Code S.S.
5	85	SF-640.005
6	97	SF-640.006
7	109	SF-640.007
8	121	SF-640.008


No. of Holes	Length in mm	Code S.S.
5	85	Ti-SF-640.005
6	97	Ti-SF-640.006
7	109	Ti-SF-640.007
8	121	Ti-SF-640.008

[YouTube www.TV.ortho.in/72](https://www.youtube.com/watch?v=www.TV.ortho.in/72)




FILOS - PROXIMAL HUMERUS SAFETY LOCK PLATE 3.5 - STANDARD

No. of Holes	Length in mm	Code S.S.
3	90	SF-642.003
4	102	SF-642.004
5	114	SF-642.005



FILOS - PROXIMAL HUMERUS SAFETY LOCK PLATE 3.5 - STANDARD


No. of Holes	Length in mm	Code TITANIUM
3	90	Ti-SF-642.003
4	102	Ti-SF-642.004
5	114	Ti-SF-642.005



[YouTube www.TV.ortho.in/48](https://www.youtube.com/watch?v=www.TV.ortho.in/48)

FILOS - PROXIMAL HUMERUS SAFETY LOCK PLATE 3.5 - LONG

No. of Holes	Length in mm	Code S.S.
5	142	SF-642.205
6	160	SF-642.206
8	196	SF-642.208
10	232	SF-642.210
12	268	SF-642.212



FILOS - PROXIMAL HUMERUS SAFETY LOCK PLATE 3.5 - LONG

No. of Holes	Length in mm	Code TITANIUM
5	142	Ti-SF-642.205
6	160	Ti-SF-642.206
8	196	Ti-SF-642.208
10	232	Ti-SF-642.210
12	268	Ti-SF-642.212



Left Plate - Stainless Steel



Right Plate - Stainless Steel

YouTube www.TV.ortho.in/11150



Left Plate - Titanium



Right Plate - Titanium

DISTAL HUMERUS SAFETY LOCK PLATE 2.7/3.5, DORSOLATERAL

No. of Holes	Length in mm	STAINLESS STEEL	
		Left	Right
3	55	SF-601.03L	SF-601.03R
5	90	SF-601.05L	SF-601.05R
7	116	SF-601.07L	SF-601.07R
9	143	SF-601.09L	SF-601.09R
14	208	SF-601.14L	SF-601.14R

No. of Holes	Length in mm	TITANIUM	
		Left	Right
3	55	Ti-SF-601.03L	Ti-SF-601.03R
5	90	Ti-SF-601.05L	Ti-SF-601.05R
7	116	Ti-SF-601.07L	Ti-SF-601.07R
9	143	Ti-SF-601.09L	Ti-SF-601.09R
14	208	Ti-SF-601.14L	Ti-SF-601.14R

SCREWS : Safety Lock - 2.4mm, 2.7mm, 3.5mm
Cortex - 3.5mm
Cancellous - 4.0mm



Left Plate - Stainless Steel



Right Plate - Stainless Steel

YouTube www.TV.ortho.in/11160



Left Plate - Titanium



Right Plate - Titanium

DISTAL HUMERUS SAFETY LOCK PLATE 2.7/3.5 DORSOLATERAL WITH LATERAL SUPPORT

No. of Holes	Length in mm	STAINLESS STEEL	
		Left	Right
3	55	SF-603.03L	SF-603.03R
5	90	SF-603.05L	SF-603.05R
7	116	SF-603.07L	SF-603.07R
9	143	SF-603.09L	SF-603.09R
14	208	SF-603.14L	SF-603.14R

No. of Holes	Length in mm	TITANIUM	
		Left	Right
3	55	Ti-SF-603.03L	Ti-SF-603.03R
5	90	Ti-SF-603.05L	Ti-SF-603.05R
7	116	Ti-SF-603.07L	Ti-SF-603.07R
9	143	Ti-SF-603.09L	Ti-SF-603.09R
14	208	Ti-SF-603.14L	Ti-SF-603.14R



DISTAL HUMERUS SAFETY LOCK PLATE 3.5, EXTRA-ARTICULAR

No. of Holes	Length in mm	STAINLESS STEEL	
		Left	Right
4	122	SF-612.04L	SF-612.04R
6	158	SF-612.06L	SF-612.06R
8	194	SF-612.08L	SF-612.08R
10	230	SF-612.10L	SF-612.10R
12	266	SF-612.12L	SF-612.12R
14	302	SF-612.14L	SF-612.14R



No. of Holes	Length in mm	TITANIUM	
		Left	Right
4	122	Ti-SF-612.04L	Ti-SF-612.04R
6	158	Ti-SF-612.06L	Ti-SF-612.06R
8	194	Ti-SF-612.08L	Ti-SF-612.08R
10	230	Ti-SF-612.10L	Ti-SF-612.10R
12	266	Ti-SF-612.12L	Ti-SF-612.12R
14	302	Ti-SF-612.14L	Ti-SF-612.14R



YouTube www.TV.ortho.in/11181



MEDIAL DISTAL HUMERUS SAFETY LOCK PLATE 2.7/3.5

No. of Holes	Length in mm	STAINLESS STEEL	
		Left	Right
3	58	SF-604.03L	SF-604.03R
5	83	SF-604.05L	SF-604.05R
7	110	SF-604.07L	SF-604.07R
9	149	SF-604.09L	SF-604.09R
14	201	SF-604.14L	SF-604.14R

No. of Holes	Length in mm	TITANIUM	
		Left	Right
3	58	Ti-SF-604.03L	Ti-SF-604.03R
5	83	Ti-SF-604.05L	Ti-SF-604.05R
7	110	Ti-SF-604.07L	Ti-SF-604.07R
9	149	Ti-SF-604.09L	Ti-SF-604.09R
14	201	Ti-SF-604.14L	Ti-SF-604.14R

YouTube www.TV.ortho.in/11347



METAPHYSEAL SAFETY LOCK PLATE 3.5 FOR DISTAL, MEDIAL HUMERUS

No. of Holes	Length in mm	Code S.S.
7	80	SF-611.07
9	106	SF-611.09
11	132	SF-611.11
13	158	SF-611.13
15	184	SF-611.15

No. of Holes	Length in mm	Code TITANIUM
7	80	Ti-SF-611.07
9	106	Ti-SF-611.09
11	132	Ti-SF-611.11
13	158	Ti-SF-611.13
15	184	Ti-SF-611.15



YouTube www.TV.ortho.in/75

SAFETY LOCK 'T' PLATE 3.5 RIGHT ANGLED



No. of Holes	Length in mm	Code S.S.
3	50	SF-141.103
4	57	SF-141.104
5	67	SF-141.105
6	77	SF-141.106

No. of Holes	Length in mm	Code TITANIUM
3	50	Ti-SF-141.103
4	57	Ti-SF-141.104
5	67	Ti-SF-141.105
6	77	Ti-SF-141.106



YouTube www.TV.ortho.in/74

SAFETY LOCK 'T' PLATE 3.5 OBLIQUE ANGLED



No. of Holes	Length in mm	STAINLESS STEEL	
		Left	Right
3	52	SF-142.103L	SF-142.103R
4	63	SF-142.104L	SF-142.104R
5	74	SF-142.105L	SF-142.105R
6	85	SF-142.106L	SF-142.106R

No. of Holes	Length in mm	TITANIUM	
		Left	Right
3	52	Ti-SF-142.103L	Ti-SF-142.103R
4	63	Ti-SF-142.104L	Ti-SF-142.104R
5	74	Ti-SF-142.105L	Ti-SF-142.105R
6	85	Ti-SF-142.106L	Ti-SF-142.106R





YouTube www.TV.ortho.in/9153

LC-DCP SAFETY LOCK PLATE 3.5

No. of Holes	Length in mm	Code S.S.
4	59	SF-147.104
5	72	SF-147.105
6	85	SF-147.106
7	98	SF-147.107
8	111	SF-147.108
9	124	SF-147.109
10	137	SF-147.110
11	150	SF-147.111
12	163	SF-147.112

No. of Holes	Length in mm	Code TITANIUM
4	59	Ti-SF-147.104
5	72	Ti-SF-147.105
6	85	Ti-SF-147.106
7	98	Ti-SF-147.107
8	111	Ti-SF-147.108
9	124	Ti-SF-147.109
10	137	Ti-SF-147.110
11	150	Ti-SF-147.111
12	163	Ti-SF-147.112



YouTube www.TV.ortho.in/70

PROXIMAL TIBIAL SAFETY LOCK PLATE 3.5

No. of Holes	Length in mm	STAINLESS STEEL	
		Left	Right
4	81	SF-607.04L	SF-607.04R
6	107	SF-607.06L	SF-607.06R
8	133	SF-607.08L	SF-607.08R
10	159	SF-607.10L	SF-607.10R
12	185	SF-607.12L	SF-607.12R
14	211	SF-607.14L	SF-607.14R
16	237	SF-607.16L	SF-607.16R

No. of Holes	Length in mm	TITANIUM	
		Left	Right
4	81	Ti-SF-607.04L	Ti-SF-607.04R
6	107	Ti-SF-607.06L	Ti-SF-607.06R
8	133	Ti-SF-607.08L	Ti-SF-607.08R
10	159	Ti-SF-607.10L	Ti-SF-607.10R
12	185	Ti-SF-607.12L	Ti-SF-607.12R
14	211	Ti-SF-607.14L	Ti-SF-607.14R
16	237	Ti-SF-607.16L	Ti-SF-607.16R



DISTAL TIBIA SAFETY LOCK PLATE 3.5 ANTEROLATERAL

No. of Holes	Length in mm	STAINLESS STEEL	
		Left	Right
5	80	SF-650.05L	SF-650.05R
7	106	SF-650.07L	SF-650.07R
9	132	SF-650.09L	SF-650.09R
11	158	SF-650.11L	SF-650.11R
13	184	SF-650.13L	SF-650.13R
15	210	SF-650.15L	SF-650.15R
17	236	SF-650.17L	SF-650.17R
19	262	SF-650.19L	SF-650.19R
21	288	SF-650.21L	SF-650.21R

No. of Holes	Length in mm	TITANIUM	
		Left	Right
5	80	Ti-SF-650.05L	Ti-SF-650.05R
7	106	Ti-SF-650.07L	Ti-SF-650.07R
9	132	Ti-SF-650.09L	Ti-SF-650.09R
11	158	Ti-SF-650.11L	Ti-SF-650.11R
13	184	Ti-SF-650.13L	Ti-SF-650.13R
15	210	Ti-SF-650.15L	Ti-SF-650.15R
17	236	Ti-SF-650.17L	Ti-SF-650.17R
19	262	Ti-SF-650.19L	Ti-SF-650.19R
21	288	Ti-SF-650.21L	Ti-SF-650.21R



MEDIAL DISTAL TIBIA SAFETY LOCK PLATE 3.5 WITH TAB

YouTube www.TV.ortho.in/24630

No. of Holes	Length in mm	STAINLESS STEEL	
		Left	Right
4	116	SF-645.04L	SF-645.04R
6	142	SF-645.06L	SF-645.06R
8	168	SF-645.08L	SF-645.08R
10	194	SF-645.10L	SF-645.10R
12	220	SF-645.12L	SF-645.12R
14	246	SF-645.14L	SF-645.14R

No. of Holes	Length in mm	TITANIUM	
		Left	Right
4	116	Ti-SF-645.04L	Ti-SF-645.04R
6	142	Ti-SF-645.06L	Ti-SF-645.06R
8	168	Ti-SF-645.08L	Ti-SF-645.08R
10	194	Ti-SF-645.10L	Ti-SF-645.10R
12	220	Ti-SF-645.12L	Ti-SF-645.12R
14	246	Ti-SF-645.14L	Ti-SF-645.14R

SCREWS : For Combi holes - 3.5 Locking, 3.5 Cortical, 4.0 Cancellous.
For Round Locking Holes - 3.5 Locking, 2.7 or 3.5 Cortical, 4.0 Cancellous.



Left Plate - Stainless Steel



Right Plate - Stainless Steel



Left Plate - Titanium



Right Plate - Titanium

MEDIAL DISTAL TIBIA SAFETY LOCK PLATE 3.5 WITHOUT TAB

No. of Holes	Length in mm	STAINLESS STEEL	
		Left	Right
4	116	SF-646.04L	SF-646.04R
6	142	SF-646.06L	SF-646.06R
8	168	SF-646.08L	SF-646.08R
10	194	SF-646.10L	SF-646.10R
12	220	SF-646.12L	SF-646.12R
14	246	SF-646.14L	SF-646.14R

No. of Holes	Length in mm	TITANIUM	
		Left	Right
4	116	Ti-SF-646.04L	Ti-SF-646.04R
6	142	Ti-SF-646.06L	Ti-SF-646.06R
8	168	Ti-SF-646.08L	Ti-SF-646.08R
10	194	Ti-SF-646.10L	Ti-SF-646.10R
12	220	Ti-SF-646.12L	Ti-SF-646.12R
14	246	Ti-SF-646.14L	Ti-SF-646.14R

SCREWS : For Combi holes - 3.5 Locking, 3.5 Cortical, 4.0 Cancellous.
For Round Locking Holes - 3.5 Locking, 2.7 or 3.5 Cortical, 4.0 Cancellous.



YouTube www.TV.ortho.in/24631



DISTAL TIBIA SAFETY LOCK PLATE 3.5/5.0

No. of Holes	Length in mm	STAINLESS STEEL	
		Left	Right
4	116	SF-605.04L	SF-605.04R
6	142	SF-605.06L	SF-605.06R
8	168	SF-605.08L	SF-605.08R
10	194	SF-605.10L	SF-605.10R
12	220	SF-605.12L	SF-605.12R
14	246	SF-605.14L	SF-605.14R


No. of Holes	Length in mm	TITANIUM	
		Left	Right
4	116	Ti-SF-605.04L	Ti-SF-605.04R
6	142	Ti-SF-605.06L	Ti-SF-605.06R
8	168	Ti-SF-605.08L	Ti-SF-605.08R
10	194	Ti-SF-605.10L	Ti-SF-605.10R
12	220	Ti-SF-605.12L	Ti-SF-605.12R
14	246	Ti-SF-605.14L	Ti-SF-605.14R



OLECRANON SAFETY LOCK PLATE 3.5

No. of Holes	Length in mm	STAINLESS STEEL	
		Left	Right
2	86	SF-647.02L	SF-647.02R
4	112	SF-647.04L	SF-647.04R
6	138	SF-647.06L	SF-647.06R
8	164	SF-647.08L	SF-647.08R
10	190	SF-647.10L	SF-647.10R
12	216	SF-647.12L	SF-647.12R


No. of Holes	Length in mm	TITANIUM	
		Left	Right
2	86	Ti-SF-647.02L	Ti-SF-647.02R
4	112	Ti-SF-647.04L	Ti-SF-647.04R
6	138	Ti-SF-647.06L	Ti-SF-647.06R
8	164	Ti-SF-647.08L	Ti-SF-647.08R
10	190	Ti-SF-647.10L	Ti-SF-647.10R
12	216	Ti-SF-647.12L	Ti-SF-647.12R




ONE THIRD TUBULAR SAFETY LOCK PLATE 3.5

No. of Holes	Length in mm	Code S.S.
2	28	SF-138.102
3	40	SF-138.103
4	52	SF-138.104
5	64	SF-138.105
6	76	SF-138.106
7	88	SF-138.107
8	100	SF-138.108
9	112	SF-138.109
10	124	SF-138.110
11	136	SF-138.111
12	148	SF-138.112

No. of Holes	Length in mm	Code Titanium
2	28	Ti-SF-138.102
3	40	Ti-SF-138.103
4	52	Ti-SF-138.104
5	64	Ti-SF-138.105
6	76	Ti-SF-138.106
7	88	Ti-SF-138.107
8	100	Ti-SF-138.108
9	112	Ti-SF-138.109
10	124	Ti-SF-138.110
11	136	Ti-SF-138.111
12	148	Ti-SF-138.112



YouTube www.TV.ortho.in/11397



METAPHYSEAL SAFETY LOCK PLATE 3.5/4.5/5.0 FOR DISTAL TIBIA

No. of Holes	Length in mm	STAINLESS STEEL	
		Left	Right
4+4	123	SF-617.04L	SF-617.04R
4+5	141	SF-617.05L	SF-617.05R
4+6	159	SF-617.06L	SF-617.06R
4+7	177	SF-617.07L	SF-617.07R
4+8	195	SF-617.08L	SF-617.08R
4+9	213	SF-617.09L	SF-617.09R
4+10	231	SF-617.10L	SF-617.10R
4+12	267	SF-617.11L	SF-617.11R
4+14	303	SF-617.12L	SF-617.12R
4+16	339	SF-617.13L	SF-617.13R
4+18	375	SF-617.14L	SF-617.14R
4+20	411	SF-617.15L	SF-617.15R

No. of Holes	Length in mm	TITANIUM	
		Left	Right
4+4	123	Ti-SF-617.04L	Ti-SF-617.04R
4+5	141	Ti-SF-617.05L	Ti-SF-617.05R
4+6	159	Ti-SF-617.06L	Ti-SF-617.06R
4+7	177	Ti-SF-617.07L	Ti-SF-617.07R
4+8	195	Ti-SF-617.08L	Ti-SF-617.08R
4+9	213	Ti-SF-617.09L	Ti-SF-617.09R
4+10	231	Ti-SF-617.10L	Ti-SF-617.10R
4+12	267	Ti-SF-617.11L	Ti-SF-617.11R
4+14	303	Ti-SF-617.12L	Ti-SF-617.12R
4+16	339	Ti-SF-617.13L	Ti-SF-617.13R
4+18	375	Ti-SF-617.14L	Ti-SF-617.14R
4+20	411	Ti-SF-617.15L	Ti-SF-617.15R

YouTube www.TV.ortho.in/71



METAPHYSEAL SAFETY LOCK PLATE 3.5/4.5/5.0 - BROAD

No. of Holes	Length in mm	Code S.S.
5+3	118	SF-615.03
5+4	136	SF-615.04
5+5	154	SF-615.05
5+6	172	SF-615.06
5+7	190	SF-615.07
5+8	208	SF-615.08
5+9	226	SF-615.09
5+11	262	SF-615.11
5+13	298	SF-615.13
5+15	334	SF-615.15

No. of Holes	Length in mm	Code TITANIUM
5+3	118	Ti-SF-615.03
5+4	136	Ti-SF-615.04
5+5	154	Ti-SF-615.05
5+6	172	Ti-SF-615.06
5+7	190	Ti-SF-615.07
5+8	208	Ti-SF-615.08
5+9	226	Ti-SF-615.09
5+11	262	Ti-SF-615.11
5+13	298	Ti-SF-615.13
5+15	334	Ti-SF-615.15

YouTube www.TV.ortho.in/24642



METAPHYSEAL SAFETY LOCK PLATE 3.5

No. of Holes	Length in mm	Code S.S.
6	86	SF-613.06
7	99	SF-613.07
8	112	SF-613.08
9	125	SF-613.09
10	138	SF-613.10
11	151	SF-613.11
12	164	SF-613.12
14	190	SF-613.13
16	216	SF-613.14
18	242	SF-613.15

No. of Holes	Length in mm	Code TITANIUM
6	86	Ti-SF-613.06
7	99	Ti-SF-613.07
8	112	Ti-SF-613.08
9	125	Ti-SF-613.09
10	138	Ti-SF-613.10
11	151	Ti-SF-613.11
12	164	Ti-SF-613.12
14	190	Ti-SF-613.13
16	216	Ti-SF-613.14
18	242	Ti-SF-613.15

YouTube www.TV.ortho.in/24640



CLAVICLE HOOK SAFETY LOCK PLATE 3.5

No. of Holes	Hook Depth in mm	STAINLESS STEEL	
		Left	Right
4	12	SF-609.01L	SF-609.01R
4	15	SF-609.03L	SF-609.03R
4	18	SF-609.05L	SF-609.05R
5	12	SF-609.07L	SF-609.07R
5	15	SF-609.09L	SF-609.09R
5	18	SF-609.11L	SF-609.11R
6	15	SF-609.13L	SF-609.13R
6	18	SF-609.15L	SF-609.15R
7	15	SF-609.17L	SF-609.17R
7	18	SF-609.19L	SF-609.19R

No. of Holes	Length in mm	TITANIUM	
		Left	Right
4	12	Ti-SF-609.01L	Ti-SF-609.01R
4	15	Ti-SF-609.03L	Ti-SF-609.03R
4	18	Ti-SF-609.05L	Ti-SF-609.05R
5	12	Ti-SF-609.07L	Ti-SF-609.07R
5	15	Ti-SF-609.09L	Ti-SF-609.09R
5	18	Ti-SF-609.11L	Ti-SF-609.11R
6	15	Ti-SF-609.13L	Ti-SF-609.13R
6	18	Ti-SF-609.15L	Ti-SF-609.15R
7	15	Ti-SF-609.17L	Ti-SF-609.17R
7	18	Ti-SF-609.19L	Ti-SF-609.19R



CLAVICLE SAFETY LOCK PLATE 3.5
ANTERIOR

No. of Holes	Length mm	STAINLESS STEEL	
		Left	Right
6	94	SF-620.06L	SF-620.06R
7	110	SF-620.07L	SF-620.07R
8	120	SF-620.08L	SF-620.08R



No. of Holes	Length mm	TITANIUM	
		Left	Right
6	94	Ti-SF-620.06L	Ti-SF-620.06R
7	110	Ti-SF-620.07L	Ti-SF-620.07R
8	120	Ti-SF-620.08L	Ti-SF-620.08R



CLAVICLE SAFETY LOCK PLATE 3.5
ANTERIOR WITH LATERAL EXTENSION

No. of Holes	Length mm	STAINLESS STEEL	
		Left	Right
3	69	SF-622.03L	SF-622.03R
4	81	SF-622.04L	SF-622.04R
5	94	SF-622.05L	SF-622.05R
6	108	SF-622.06L	SF-622.06R
7	123	SF-622.07L	SF-622.07R
8	135	SF-622.08L	SF-622.08R



No. of Holes	Length mm	TITANIUM	
		Left	Right
3	69	Ti-SF-622.03L	Ti-SF-622.03R
4	81	Ti-SF-622.04L	Ti-SF-622.04R
5	94	Ti-SF-622.05L	Ti-SF-622.05R
6	108	Ti-SF-622.06L	Ti-SF-622.06R
7	123	Ti-SF-622.07L	Ti-SF-622.07R
8	135	Ti-SF-622.08L	Ti-SF-622.08R



YouTube www.TV.ortho.in/24682

CALCANEAL SAFETY LOCK PLATE 3.5


Length in mm	STAINLESS STEEL	
	Left	Right
60	SF-137.160L	SF-137.160R
70	SF-137.170L	SF-137.170R

for Calcaneus Fractures.




Length in mm	TITANIUM	
	Left	Right
60	Ti-SF-137.160L	Ti-SF-137.160R
70	Ti-SF-137.170L	Ti-SF-137.170R

for Calcaneus Fractures.



CALCANEUS SAFETY LOCK PLATE 3.5

Size	Length mm	STAINLESS STEEL	
		Left	Right
Extra Small	64	SF-138.64L	SF-138.64R
Small	69	SF-138.69L	SF-138.69R
Large	76	SF-138.76L	SF-138.76R
Extra Large	81	SF-138.81L	SF-138.81R



Size	Length mm	TITANIUM	
		Left	Right
Extra Small	64	Ti-SF-138.64L	Ti-SF-138.64R
Small	69	Ti-SF-138.69L	Ti-SF-138.69R
Large	76	Ti-SF-138.76L	Ti-SF-138.76R
Extra Large	81	Ti-SF-138.81L	Ti-SF-138.81R

[YouTube www.TV.ortho.in/24650](http://www.TV.ortho.in/24650)



RECONSTRUCTION SAFETY LOCK PLATE 3.5 - STRAIGHT

No. of Holes	Length in mm	Code S.S.
5	70	SF-144.105
6	84	SF-144.106
7	98	SF-144.107
8	112	SF-144.108
9	126	SF-144.109
10	140	SF-144.110
11	154	SF-144.111
12	168	SF-144.112
13	182	SF-144.113
14	196	SF-144.114
15	210	SF-144.115
16	224	SF-144.116
17	238	SF-144.117
18	252	SF-144.118
19	266	SF-144.119
20	280	SF-144.120

No. of Holes	Length in mm	Code TITANIUM
5	70	Ti-SF-144.105
6	84	Ti-SF-144.106
7	98	Ti-SF-144.107
8	112	Ti-SF-144.108
9	126	Ti-SF-144.109
10	140	Ti-SF-144.110
11	154	Ti-SF-144.111
12	168	Ti-SF-144.112
13	182	Ti-SF-144.113
14	196	Ti-SF-144.114
15	210	Ti-SF-144.115
16	224	Ti-SF-144.116
17	238	Ti-SF-144.117
18	252	Ti-SF-144.118
19	266	Ti-SF-144.119
20	280	Ti-SF-144.120

[YouTube www.TV.ortho.in/24654](http://www.TV.ortho.in/24654)



RECONSTRUCTION SAFETY LOCK PLATE 3.5 ROUND HOLES - STRAIGHT

No. of Holes	Length in mm	Code S.S.
5	65	SF-444.105
6	78	SF-444.106
7	91	SF-444.107
8	104	SF-444.108
9	117	SF-444.109
10	130	SF-444.110
11	143	SF-444.111
12	156	SF-444.112
13	169	SF-444.113
14	182	SF-444.114
15	195	SF-444.115
16	208	SF-444.116
17	221	SF-444.117
18	234	SF-444.118
19	247	SF-444.119
20	260	SF-444.120

No. of Holes	Length in mm	Code TITANIUM
5	65	Ti-SF-444.105
6	78	Ti-SF-444.106
7	91	Ti-SF-444.107
8	104	Ti-SF-444.108
9	117	Ti-SF-444.109
10	130	Ti-SF-444.110
11	143	Ti-SF-444.111
12	156	Ti-SF-444.112
13	169	Ti-SF-444.113
14	182	Ti-SF-444.114
15	195	Ti-SF-444.115
16	208	Ti-SF-444.116
17	221	Ti-SF-444.117
18	234	Ti-SF-444.118
19	247	Ti-SF-444.119
20	260	Ti-SF-444.120



LATERAL DISTAL FIBULA SAFETY LOCK PLATE 2.7/3.5

No. of Holes	Length in mm	STAINLESS STEEL	
		Left	Right
3	73	SF-652.03L	SF-652.03R
4	86	SF-652.04L	SF-652.04R
5	99	SF-652.05L	SF-652.05R
6	112	SF-652.06L	SF-652.06R
7	125	SF-652.07L	SF-652.07R

No. of Holes	Length in mm	TITANIUM	
		Left	Right
3	73	Ti-SF-652.03L	Ti-SF-652.03R
4	86	Ti-SF-652.04L	Ti-SF-652.04R
5	99	Ti-SF-652.05L	Ti-SF-652.05R
6	112	Ti-SF-652.06L	Ti-SF-652.06R
7	125	Ti-SF-652.07L	Ti-SF-652.07R

YouTube www.TV.ortho.in/24652



RECONSTRUCTION SAFETY LOCK PLATE 3.5 - CURVED

No. of Holes	Length in mm	Code S.S.
5	70	SF-144.205C
6	82	SF-144.206C
7	94	SF-144.207C
8	106	SF-144.208C
9	118	SF-144.209C
10	129	SF-144.210C
11	139	SF-144.211C
12	149	SF-144.212C
13	156	SF-144.213C
14	166	SF-144.214C
16	180	SF-144.216C
18	190	SF-144.218C

No. of Holes	Length in mm	Code TITANIUM
5	70	Ti-SF-144.205C
6	82	Ti-SF-144.206C
7	94	Ti-SF-144.207C
8	106	Ti-SF-144.208C
9	118	Ti-SF-144.209C
10	129	Ti-SF-144.210C
11	139	Ti-SF-144.211C
12	149	Ti-SF-144.212C
13	156	Ti-SF-144.213C
14	166	Ti-SF-144.214C
16	180	Ti-SF-144.216C
18	190	Ti-SF-144.218C




**RECONSTRUCTION SAFETY LOCK PLATE 3.5
ROUND HOLES - CURVED**

No. of Holes	Length in mm	Code S.S.
5	65	SF-444.205C
6	78	SF-444.206C
7	91	SF-444.207C
8	104	SF-444.208C
9	117	SF-444.209C
10	130	SF-444.210C
11	143	SF-444.211C
12	156	SF-444.212C
13	169	SF-444.213C
14	182	SF-444.214C
16	208	SF-444.216C
18	234	SF-444.218C

No. of Holes	Length in mm	Code TITANIUM
5	65	Ti-SF-444.205C
6	78	Ti-SF-444.206C
7	91	Ti-SF-444.207C
8	104	Ti-SF-444.208C
9	117	Ti-SF-444.209C
10	130	Ti-SF-444.210C
11	143	Ti-SF-444.211C
12	156	Ti-SF-444.212C
13	169	Ti-SF-444.213C
14	182	Ti-SF-444.214C
16	208	Ti-SF-444.216C
18	234	Ti-SF-444.218C

YouTube www.TV.ortho.in/117



CLOVERLEAF SAFETY LOCK PLATE 3.5

No. of Holes	Length in mm	Code S.S.
3	88	SF-146.103
4	104	SF-146.104
5	120	SF-146.105
6	136	SF-146.106

No. of Holes	Length in mm	Code TITANIUM
3	88	Ti-SF-146.103
4	104	Ti-SF-146.104
5	120	Ti-SF-146.105
6	136	Ti-SF-146.106



VOLAR COLUMN DISTAL RADIUS SAFETY LOCK PLATE 2.4 (8 HEAD HOLES)

L/R	No. of Holes	Length mm	Code S.S.
Left	3	51	SF-120.803L
Right	3	51	SF-120.803R
Left	4	60	SF-120.804L
Right	4	60	SF-120.804R
Left	5	69	SF-120.805L
Right	5	69	SF-120.805R

L/R	No. of Holes	Length mm	Code TITANIUM
Left	3	51	Ti-SF-120.803L
Right	3	51	Ti-SF-120.803R
Left	4	60	Ti-SF-120.804L
Right	4	60	Ti-SF-120.804R
Left	5	69	Ti-SF-120.805L
Right	5	69	Ti-SF-120.805R



VOLAR COLUMN DISTAL RADIUS SAFETY LOCK PLATE 2.4 (9 HEAD HOLES)

L/R	No. of Holes	Length mm	Code S.S.
Left	3	51	SF-120.903L
Right	3	51	SF-120.903R
Left	4	60	SF-120.904L
Right	4	60	SF-120.904R
Left	5	69	SF-120.905L
Right	5	69	SF-120.905R

L/R	No. of Holes	Length mm	Code TITANIUM
Left	3	51	Ti-SF-120.903L
Right	3	51	Ti-SF-120.903R
Left	4	60	Ti-SF-120.904L
Right	4	60	Ti-SF-120.904R
Left	5	69	Ti-SF-120.905L
Right	5	69	Ti-SF-120.905R



VOLAR DISTAL RADIUS SAFETY LOCK PLATE 2.4, EXTRA-ARTICULAR, EXTRA-LONG

Head Holes	Shaft Holes	Length mm	Code S.S.
4	8	100	SF-122.108
4	10	120	SF-122.110
4	12	140	SF-122.112

Head Holes	Shaft Holes	Length mm	Code TITANIUM
4	8	100	Ti-SF-122.108
4	10	120	Ti-SF-122.110
4	12	140	Ti-SF-122.112



VOLAR DISTAL RADIUS SAFETY LOCK PLATE 2.4, EXTRA ARTICULAR (5 HEAD HOLES)

L/R	No. of Holes	Length mm	Code S.S.
Left	3	48	SF-121.504L
Right	3	48	SF-121.504R
Left	5	66	SF-121.505L
Right	5	66	SF-121.505R

L/R	No. of Holes	Length mm	Code TITANIUM
Left	3	48	Ti-SF-121.504L
Right	3	48	Ti-SF-121.504R
Left	5	66	Ti-SF-121.505L
Right	5	66	Ti-SF-121.505R



VOLAR DISTAL RADIUS SAFETY LOCK PLATE 2.4, EXTRA ARTICULAR (4 HEAD HOLES)

L/R	No. of Holes	Length mm	Code S.S.
Left	3	47	SF-121.304L
Right	3	47	SF-121.304R
Left	5	65	SF-121.305L
Right	5	65	SF-121.305R

L/R	No. of Holes	Length mm	Code TITANIUM
Left	3	47	Ti-SF-121.304L
Right	3	47	Ti-SF-121.304R
Left	5	65	Ti-SF-121.305L
Right	5	65	Ti-SF-121.305R



VOLAR DISTAL RADIUS SAFETY LOCK PLATE 2.4, JUXTA-ARTICULAR (5 HEAD HOLES)

L/R	No. of Holes	Length mm	Code S.S.
Left	3	43	SF-123.503L
Right	3	43	SF-123.503R
Left	5	61	SF-123.505L
Right	5	61	SF-123.505R

L/R	No. of Holes	Length mm	Code TITANIUM
Left	3	43	Ti-SF-123.503L
Right	3	43	Ti-SF-123.503R
Left	5	61	Ti-SF-123.505L
Right	5	61	Ti-SF-123.505R



VOLAR DISTAL RADIUS BUTTRESS SAFETY LOCK PLATE 2.4, JUXTA-ARTICULAR (5 HEAD HOLES)

L/R	No. of Holes	Length mm	Code S.S.
Left	3	52	SF-124.503L
Right	3	52	SF-124.503R

L/R	No. of Holes	Length mm	Code TITANIUM
Left	3	52	Ti-SF-124.503L
Right	3	52	Ti-SF-124.503R



DORSAL DISTAL RADIUS SAFETY LOCK PLATE 2.4 - STRAIGHT

No. of Holes	Length in mm	Code S.S.
5	48	SF-125.105
6	57	SF-125.106

No. of Holes	Length in mm	Code TITANIUM
5	48	Ti-SF-125.105
6	57	Ti-SF-125.106



DORSAL DISTAL RADIUS SAFETY LOCK L-PLATE 2.4

Angle Type	Head/Shaft Holes	Length mm	Code S.S.
Right	2/3	40	SF-126.203R
Right	2/4	49	SF-126.204R
Left	2/3	40	SF-126.203L
Left	2/4	49	SF-126.204L
Right	3/3	40	SF-126.303R
Right	3/4	49	SF-126.304R
Left	3/3	40	SF-126.303L
Left	3/4	49	SF-126.304L

Angle Type	Head/Shaft Holes	Length mm	Code TITANIUM
Right	2/3	40	Ti-SF-126.203R
Right	2/4	49	Ti-SF-126.204R
Left	2/3	40	Ti-SF-126.203L
Left	2/4	49	Ti-SF-126.204L
Right	3/3	40	Ti-SF-126.303R
Right	3/4	49	Ti-SF-126.304R
Left	3/3	40	Ti-SF-126.303L
Left	3/4	49	Ti-SF-126.304L



DORSAL DISTAL RADIUS SAFETY LOCK L-PLATE 2.4 - OBLIQUE

Angle Type	Head/Shaft Holes	Length mm	Code S.S.
Right	3/3	43	SF-127.303R
Right	3/4	52	SF-127.304R
Left	3/3	43	SF-127.303L
Left	3/4	52	SF-127.304L



Angle Type	Head/Shaft Holes	Length mm	Code TITANIUM
Right	3/3	43	Ti-SF-127.303R
Right	3/4	52	Ti-SF-127.304R
Left	3/3	43	Ti-SF-127.303L
Left	3/4	52	Ti-SF-127.304L



DORSAL DISTAL RADIUS SAFETY LOCK T-PLATE 2.4

Head/Shaft Holes	Length in mm	Code S.S.
3/3	40	SF-128.303
3/4	49	SF-128.304

Head/Shaft Holes	Length in mm	Code TITANIUM
3/3	40	Ti-SF-128.303
3/4	49	Ti-SF-128.304




PILON SAFETY LOCK PLATE 2.7 CRUCIFORM

No. of Holes	Length in mm	Code S.S.
7	147	SF-649.07
9	173	SF-649.09

No. of Holes	Length in mm	Code TITANIUM
7	147	Ti-SF-649.07
9	173	Ti-SF-649.09

YouTube www.TV.ortho.in/24662



**SAFETY LOCK SCREW
2.4MM - SELF TAPPING**

**SAFETY LOCK SCREW
2.7MM - SELF TAPPING**

Length in mm	STAINLESS STEEL
6	SF-100.206
8	SF-100.208
10	SF-100.210
12	SF-100.212
14	SF-100.214
16	SF-100.216
18	SF-100.218
20	SF-100.220
22	SF-100.222
24	SF-100.224
26	SF-100.226
28	SF-100.228
30	SF-100.230

Drill Bit: 1.8mm

Length in mm	TITANIUM
6	Ti-SF-100.206
8	Ti-SF-100.208
10	Ti-SF-100.210
12	Ti-SF-100.212
14	Ti-SF-100.214
16	Ti-SF-100.216
18	Ti-SF-100.218
20	Ti-SF-100.220
22	Ti-SF-100.222
24	Ti-SF-100.224
26	Ti-SF-100.226
28	Ti-SF-100.228
30	Ti-SF-100.230

Drill Bit: 1.8mm

Length in mm	STAINLESS STEEL
20	SF-101.420
22	SF-101.422
24	SF-101.424
26	SF-101.426
28	SF-101.428
30	SF-101.430
32	SF-101.432
34	SF-101.434
36	SF-101.436
38	SF-101.438
40	SF-101.440
42	SF-101.442
44	SF-101.444
46	SF-101.446
48	SF-101.448
50	SF-101.450
52	SF-101.452
54	SF-101.454
56	SF-101.456
58	SF-101.458
60	SF-101.460

Drill Bit: 2.0mm

Length in mm	TITANIUM
20	Ti-SF-101.420
22	Ti-SF-101.422
24	Ti-SF-101.424
26	Ti-SF-101.426
28	Ti-SF-101.428
30	Ti-SF-101.430
32	Ti-SF-101.432
34	Ti-SF-101.434
36	Ti-SF-101.436
38	Ti-SF-101.438
40	Ti-SF-101.440
42	Ti-SF-101.442
44	Ti-SF-101.444
46	Ti-SF-101.446
48	Ti-SF-101.448
50	Ti-SF-101.450
52	Ti-SF-101.452
54	Ti-SF-101.454
56	Ti-SF-101.456
58	Ti-SF-101.458
60	Ti-SF-101.460

Drill Bit: 2.0mm

YouTube www.TV.ortho.in/24658

SAFETY LOCK SCREW 3.5MM - SELF TAPPING



Drill Bit: 2.8mm


Length in mm	STAINLESS STEEL
10	SF-102.210
12	SF-102.212
14	SF-102.214
16	SF-102.216
18	SF-102.218
20	SF-102.220
22	SF-102.222
24	SF-102.224
26	SF-102.226
28	SF-102.228
30	SF-102.230
32	SF-102.232
34	SF-102.234
36	SF-102.236
38	SF-102.238
40	SF-102.240
42	SF-102.242
44	SF-102.244
46	SF-102.246
48	SF-102.248
50	SF-102.250
52	SF-102.252
54	SF-102.254
56	SF-102.256
58	SF-102.258
60	SF-102.260

Length in mm	TITANIUM
10	Ti-SF-102.210
12	Ti-SF-102.212
14	Ti-SF-102.214
16	Ti-SF-102.216
18	Ti-SF-102.218
20	Ti-SF-102.220
22	Ti-SF-102.222
24	Ti-SF-102.224
26	Ti-SF-102.226
28	Ti-SF-102.228
30	Ti-SF-102.230
32	Ti-SF-102.232
34	Ti-SF-102.234
36	Ti-SF-102.236
38	Ti-SF-102.238
40	Ti-SF-102.240
42	Ti-SF-102.242
44	Ti-SF-102.244
46	Ti-SF-102.246
48	Ti-SF-102.248
50	Ti-SF-102.250
52	Ti-SF-102.252
54	Ti-SF-102.254
56	Ti-SF-102.256
58	Ti-SF-102.258
60	Ti-SF-102.260



Drill Bit: 2.8mm

YouTube www.TV.ortho.in/124



SAFETY LOCK SCREW 3.5MM
SELF TAPPING AND SELF DRILLING

CANCELLOUS SAFETY LOCK SCREW 3.5MM
SELF TAPPING

Length in mm	STAINLESS STEEL
10	SF-102.410
12	SF-102.412
14	SF-102.414
16	SF-102.416
18	SF-102.418
20	SF-102.420
22	SF-102.422
24	SF-102.424
26	SF-102.426
28	SF-102.428
30	SF-102.430
32	SF-102.432
34	SF-102.434
36	SF-102.436
38	SF-102.438
40	SF-102.440
42	SF-102.442
44	SF-102.444
46	SF-102.446
48	SF-102.448
50	SF-102.450
52	SF-102.452
54	SF-102.454
56	SF-102.456
58	SF-102.458
60	SF-102.460

Length in mm	TITANIUM
10	Ti-SF-102.410
12	Ti-SF-102.412
14	Ti-SF-102.414
16	Ti-SF-102.416
18	Ti-SF-102.418
20	Ti-SF-102.420
22	Ti-SF-102.422
24	Ti-SF-102.424
26	Ti-SF-102.426
28	Ti-SF-102.428
30	Ti-SF-102.430
32	Ti-SF-102.432
34	Ti-SF-102.434
36	Ti-SF-102.436
38	Ti-SF-102.438
40	Ti-SF-102.440
42	Ti-SF-102.442
44	Ti-SF-102.444
46	Ti-SF-102.446
48	Ti-SF-102.448
50	Ti-SF-102.450
52	Ti-SF-102.452
54	Ti-SF-102.454
56	Ti-SF-102.456
58	Ti-SF-102.458
60	Ti-SF-102.460

Length in mm	Code No. S.S.
10	SF-102.510
12	SF-102.512
14	SF-102.514
16	SF-102.516
18	SF-102.518
20	SF-102.520
22	SF-102.522
24	SF-102.524
26	SF-102.526
28	SF-102.528
30	SF-102.530
32	SF-102.532
34	SF-102.534
36	SF-102.536
38	SF-102.538
40	SF-102.540
42	SF-102.542
44	SF-102.544
46	SF-102.546
48	SF-102.548
50	SF-102.550
52	SF-102.552
54	SF-102.554
56	SF-102.556
58	SF-102.558
60	SF-102.560

Length in mm	Code No. TITANIUM
10	Ti-SF-102.510
12	Ti-SF-102.512
14	Ti-SF-102.514
16	Ti-SF-102.516
18	Ti-SF-102.518
20	Ti-SF-102.520
22	Ti-SF-102.522
24	Ti-SF-102.524
26	Ti-SF-102.526
28	Ti-SF-102.528
30	Ti-SF-102.530
32	Ti-SF-102.532
34	Ti-SF-102.534
36	Ti-SF-102.536
38	Ti-SF-102.538
40	Ti-SF-102.540
42	Ti-SF-102.542
44	Ti-SF-102.544
46	Ti-SF-102.546
48	Ti-SF-102.548
50	Ti-SF-102.550
52	Ti-SF-102.552
54	Ti-SF-102.554
56	Ti-SF-102.556
58	Ti-SF-102.558
60	Ti-SF-102.560

Drill Bit: 2.8mm

Drill Bit: 2.8mm

□ SMALL FRAGMENT SCREWS SET - COMBO



AS138 SMALL FRAGMENT SCREWS SET COMBO WITH SAFETY LOCK SCREWS & STANDARD SCREWS

Set Consisting of:

	Units
Safety Lock Screw, 2.7mm, Self Tapping, SS 3 each of 10,12,14,16,18,20,22,24,26,28,30,35,40,45,50,55,60mm	51
Cancellous Bone Screw, 4.0mm, Short-Thread, SS 2 each of 10,12,14,16,18,20,22,24,26,28,30,35,40,45,50,55,60mm	34
Cancellous Bone Screw, 4.0mm, Full-Thread, SS 2 each of 10,12,14,16,18,20,22,24,26,28,30,35,40,45,50,55,60mm	34
Cortex Screw, 3.5mm, Self Tapping, SS 6 each of 10,12,14mm 12 each of 16,18,20mm 6 each of 22,24,26,28,30,32,34,36,38, 40,42,44,46,48,50,52,54,56,58,60mm	18 36 120
Safety Lock Screw, 3.5mm, Self Tapping, SS 6 each of 10,12,14mm 12 each of 16,18,20mm 6 each of 22,24,26,28,30,32,34,36,38,40,42,44,46,48,50,52,54,56,58,60mm	18 36 120

AS138.000 Empty Container



AS139 SMALL FRAGMENT SCREWS SET COMBO WITH SAFETY LOCK & STANDARD SCREWS - MODIFIED

Set Consisting of:

	Units
Cancellous Bone Screw, 4.0mm, Short-Thread, SS 2 each of 30,35,40,45,50,55,60mm	14
Cancellous Bone Screw, 4.0mm, Full-Thread, SS 2 each of 10,12,14,16,18,20,22,24,26,28,30,35,40,45,50,55,60mm	34
Cortex Screw, 3.5mm, Self Tapping, SS 4 each of 10,12,14,16,18,20,22,24,26,28,30,32,34,36,38mm 2 each of 40,45,50mm	60 6
Safety Lock Screw, 3.5mm, Self Tapping, SS 5 each of 10,12,14,16,18,20,22,24,26,28,30mm 5 each of 35,40,45,50,55,60mm	60 30
Safety Lock Screw, 3.5mm, Self Tapping & Self Drilling SS 2 each of 16,18,20,22,24,26,28,30mm 2 each of 35,40,45,50,55,60mm	16 12

AS139.000 Empty Container



□ SMALL FRAGMENT SAFETY LOCK PLATES & SCREWS SET




AS136 SMALL FRAGMENT SAFETY LOCK PLATES & SCREWS SET

Set Consisting of:

		Units
AS432.005	2.8mm Drill Sleeve for 3.5mm Safety Lock Screws	2
AS432.001	Drill Bit Dia. 2.8mm, Length 165mm, Calibrated, 2-flute, for Quick Coupling	2
306.010	Quick Coupling Handle (T- Type), SS	1
301.060	Screws Driver Shaft Q.C. End, 2.5mm Tip, SS	1
SF-102.210	Safety Lock Screw 3.5mm, SS, Self Tapping, Length 10mm	2
SF-102.212	Safety Lock Screw 3.5mm, SS, Self Tapping, Length 12mm	2
SF-102.214	Safety Lock Screw 3.5mm, SS, Self Tapping, Length 14mm	2
SF-102.216	Safety Lock Screw 3.5mm, SS, Self Tapping, Length 16mm	2
SF-102.218	Safety Lock Screw 3.5mm, SS, Self Tapping, Length 18mm	2
SF-102.220	Safety Lock Screw 3.5mm, SS, Self Tapping, Length 20mm	2
SF-102.222	Safety Lock Screw 3.5mm, SS, Self Tapping, Length 22mm	2
SF-102.224	Safety Lock Screw 3.5mm, SS, Self Tapping, Length 24mm	2
SF-102.226	Safety Lock Screw 3.5mm, SS, Self Tapping, Length 26mm	3
SF-102.228	Safety Lock Screw 3.5mm, SS, Self Tapping, Length 28mm	3
SF-102.230	Safety Lock Screw 3.5mm, SS, Self Tapping, Length 30mm	3
SF-102.232	Safety Lock Screw 3.5mm, SS, Self Tapping, Length 32mm	3
SF-102.234	Safety Lock Screw 3.5mm, SS, Self Tapping, Length 34mm	3
SF-102.236	Safety Lock Screw 3.5mm, SS, Self Tapping, Length 36mm	3
SF-102.238	Safety Lock Screw 3.5mm, SS, Self Tapping, Length 38mm	3
SF-102.240	Safety Lock Screw 3.5mm, SS, Self Tapping, Length 40mm	3
SF-102.242	Safety Lock Screw 3.5mm, SS, Self Tapping, Length 42mm	3
SF-102.244	Safety Lock Screw 3.5mm, SS, Self Tapping, Length 44mm	3
SF-102.246	Safety Lock Screw 3.5mm, SS, Self Tapping, Length 46mm	3
SF-102.248	Safety Lock Screw 3.5mm, SS, Self Tapping, Length 48mm	2
SF-102.250	Safety Lock Screw 3.5mm, SS, Self Tapping, Length 50mm	2
SF-102.252	Safety Lock Screw 3.5mm, SS, Self Tapping, Length 52mm	2
SF-102.256	Safety Lock Screw 3.5mm, SS, Self Tapping, Length 56mm	2
SF-102.260	Safety Lock Screw 3.5mm, SS, Self Tapping, Length 60mm	2
SF-147.104	LC-DCP Safety Lock Plate 3.5, SS, 4 Holes	1
SF-147.105	LC-DCP Safety Lock Plate 3.5, SS, 5 Holes	1
SF-147.106	LC-DCP Safety Lock Plate 3.5, SS, 6 Holes	1
SF-147.107	LC-DCP Safety Lock Plate 3.5, SS, 7 Holes	1
SF-147.108	LC-DCP Safety Lock Plate 3.5, SS, 8 Holes	1
SF-147.109	LC-DCP Safety Lock Plate 3.5, SS, 9 Holes	1
SF-147.110	LC-DCP Safety Lock Plate 3.5, SS, 10 Holes	1
SF-141.103	Safety Lock T' Plate 3.5 - Right Angled, SS, 3 Holes	1
SF-141.104	Safety Lock T' Plate 3.5 - Right Angled, SS, 4 Holes	1
SF-141.105	Safety Lock T' Plate 3.5 - Right Angled, SS, 5 Holes	1
SF-142.103L	Safety Lock T' Plate 3.5 - Oblique Angled, SS, Left, 3 Holes	1
SF-142.104L	Safety Lock T' Plate 3.5 - Oblique Angled, SS, Left, 4 Holes	1
SF-142.105L	Safety Lock T' Plate 3.5 - Oblique Angled, SS, Left, 5 Holes	1
SF-142.103R	Safety Lock T' Plate 3.5 - Oblique Angled, SS, Right, 3 Holes	1
SF-142.104R	Safety Lock T' Plate 3.5 - Oblique Angled, SS, Right, 4 Holes	1
SF-142.105R	Safety Lock T' Plate 3.5 - Oblique Angled, SS, Right, 5 Holes	1
SF-144.106	Reconstruction Safety Lock Plate 3.5 - Straight, SS, 6 Holes	1
SF-144.108	Reconstruction Safety Lock Plate 3.5 - Straight, SS, 8 Holes	1
SF-144.110	Reconstruction Safety Lock Plate 3.5 - Straight, SS, 10 Holes	1
SF-146.103	Cloverleaf Safety Lock Plate 3.5, SS, 3 Holes	1
SF-146.104	Cloverleaf Safety Lock Plate 3.5, SS, 4 Holes	1
SF-146.105	Cloverleaf Safety Lock Plate 3.5, SS, 5 Holes	1
906.001	Empty Container	

□ SMALL FRAGMENT INSTRUMENT SET - COMBO



AS140 SMALL FRAGMENT INSTRUMENT SET COMBO WITH STANDARD INSTRUMENTS & LOCKING INSTRUMENTS

Set Consisting of:

		Units
STANDARD INSTRUMENTS		
260.250	Drill Bit 2.5mm Dia	2
260.350	Drill Bit 3.5mm Dia	2
305.020	Counter Sink for 3.5, 4.0mm Screws (6.0mm Head)	1
264.025	Tap of 3.5mm Cortex Screw	2
264.038	Tap of 4.0mm Cancellous Screw	2
306.010	Quick Coupling Handle (T-Type)	1
307.060	Drill & Tap Sleeve Combined 3.5/2.5mm	1
307.090	Insert Drill Sleeve 2.5mm / 3.5mm	1
301.060	Small Hexagonal Screw Driver Shaft	1
301.040	Small Hexagonal Screw Driver	1
303.010	Depth Gauge Small	1
371.000	Sharp Hook	1
302.010	Small Plate Bender for Plates 2.7 & 3.5mm	1
340.020	Self Centering Bone Holding Forceps- 190mm	1
379.015	Hohman Retractors, Small, 8mm Wide, 160mm	2
381.012	Periosteal Elevator, Straight/Curved 12mm	1
379.020	Hohman Retractors, 15mm Wide, 160mm	2
341.025	Reduction Forceps, Serrated - 140mm	1
342.035	Reduction Forceps, Pointed - 140mm	1
304.010	Neutral and Loaded Drill Guide (Small)	1
LOCKING INSTRUMENTS		
AS432.001	Drill Bit 2.8mm, Length 165mm, Calibrated, 2-flute, for Quick Coupling	2
AS432.005	2.8mm Drill Sleeve for 3.5mm Safety Lock Screws	2
AS432.009	Torque limiting Screw Driver, 2.5mm Tip	1
AS432.021	Extraction Screw (Left Hand Thread), Conical, for Screws 3.5mm	1
AS432.025	Centering Sleeve for Kirschner Wire 1.2mm, Length 85mm	1
AS140.000	Empty Container	

□ SMALL FRAGMENT SAFETY LOCK INSTRUMENTS SET

 www.TV.ortho.in/27599

AS432 SMALL FRAGMENT SAFETY LOCK INSTRUMENTS SET

<i>Set Consisting of:</i>		Units
AS432.001	Drill Bit 2.8mm, Length 165mm, Calibrated, 2-flute, for Quick Coupling	2
AS432.005	2.8mm Drill Sleeve for 3.5mm Safety Lock Screws	2
AS432.009	Torque Limiting Screw Driver, 2.5mm Tip	1
AS432.013	T-Handle for Quick Coupling	1
AS432.017	HSS Drill Bit 2.5mm (Blue Coating)	2
AS432.021	Extraction Screw (Left Hand Thread), Conical, for Screws 3.5mm	1
AS432.025	Centering Sleeve for Kirschner Wire 1.2mm, Length 85mm	1
AS432.029	Empty Container	

**OPTIONAL INSTRUMENTS
(To be Ordered Separately)**

AS432.030	Torque Limiting Screw Driver, 1.5mm Tip	1
AS432.031	Drill Bit 2.0mm, Length 150mm, 2-flute for Quick Coupling	1
AS432.035	2.0mm Drill Sleeve for 2.4mm Safety Locking Screws	1



DRILL BIT Ø 2.8MM, LENGTH 165MM, CALIBRATED, 2-FLUTE, FOR QUICK COUPLING
Code No.: AS432.001



HSS DRILL BIT Ø 2.5MM
Code No.: AS432.017



**EXTRACTION SCREW (LEFT HAND THREAD)
FOR SCREWS Ø 3.5MM**
Code No.: AS432.021



**2.8MM DRILL SLEEVE
FOR 3.5MM SAFETY LOCK SCREWS**
Code No.: AS432.005



**CENTERING SLEEVE FOR KIRSCHNER
WIRE Ø 1.2MM, LENGTH 85MM**
Code No.: AS432.025



**TORQUE LIMITING SCREW DRIVER
2.5MM TIP**
Code No.: AS432.009



**QUICK COUPLING HANDLE
(T-Type)**
Code No.: AS432.013



**TORQUE LIMITING SCREW DRIVER
1.5MM TIP**
Code No.: AS432.030