





LARGE FRAGMENTS - SAFETY LOCKING IMPLANTS



[www.TV.ortho.in/119](http://www.TV.ortho.in/119)



### BROAD LC-DCP SAFETY LOCK PLATE 4.5 /5.0

No. of Holes	Length in mm	Code S.S.
5	98	SF-151.105
6	116	SF-151.106
7	134	SF-151.107
8	152	SF-151.108
9	170	SF-151.109
10	188	SF-151.110
11	206	SF-151.111
12	224	SF-151.112
13	242	SF-151.113
14	260	SF-151.114
15	278	SF-151.115
16	296	SF-151.116
18	332	SF-151.118
20	368	SF-151.120
22	404	SF-151.122
24	440	SF-151.124

No. of Holes	Length in mm	Code TITANIUM
5	98	Ti-SF-151.105
6	116	Ti-SF-151.106
7	134	Ti-SF-151.107
8	152	Ti-SF-151.108
9	170	Ti-SF-151.109
10	188	Ti-SF-151.110
11	206	Ti-SF-151.111
12	224	Ti-SF-151.112
13	242	Ti-SF-151.113
14	260	Ti-SF-151.114
15	278	Ti-SF-151.115
16	296	Ti-SF-151.116
18	332	Ti-SF-151.118
20	368	Ti-SF-151.120
22	404	Ti-SF-151.122
24	440	Ti-SF-151.124

[www.TV.ortho.in/125](http://www.TV.ortho.in/125)



### NARROW LC-DCP SAFETY LOCK PLATE 4.5/5.0

No. of Holes	Length in mm	Code S.S.
5	98	SF-150.105
6	116	SF-150.106
7	134	SF-150.107
8	152	SF-150.108
9	170	SF-150.109
10	188	SF-150.110
11	206	SF-150.111
12	224	SF-150.112
13	242	SF-150.113
14	260	SF-150.114
15	278	SF-150.115
16	296	SF-150.116
18	332	SF-150.118
20	368	SF-150.120
22	404	SF-150.122
24	440	SF-150.124

No. of Holes	Length in mm	Code TITANIUM
5	98	Ti-SF-150.105
6	116	Ti-SF-150.106
7	134	Ti-SF-150.107
8	152	Ti-SF-150.108
9	170	Ti-SF-150.109
10	188	Ti-SF-150.110
11	206	Ti-SF-150.111
12	224	Ti-SF-150.112
13	242	Ti-SF-150.113
14	260	Ti-SF-150.114
15	278	Ti-SF-150.115
16	296	Ti-SF-150.116
18	332	Ti-SF-150.118
20	368	Ti-SF-150.120
22	404	Ti-SF-150.122
24	440	Ti-SF-150.124

[www.TV.ortho.in/27883](http://www.TV.ortho.in/27883)



### PELVIS SAFETY LOCK PLATE

No. of Holes	Code No. S.S.
4	SF-648.004
6	SF-648.006



### PELVIS SAFETY LOCK PLATE

No. of Holes	Code No. TITANIUM
4	Ti-SF-648.004
6	Ti-SF-648.006

YouTube [www.TV.ortho.in/222](http://www.TV.ortho.in/222)



**BROAD CURVED SAFETY LOCK PLATE 4.5/5.0**

No. of Holes	Length in mm	Code S.S.
12	229	SF-619.12
13	247	SF-619.13
14	265	SF-619.14
15	282	SF-619.15
16	300	SF-619.16
17	318	SF-619.17
18	336	SF-619.18

No. of Holes	Length in mm	Code TITANIUM
12	229	Ti-SF-619.12
13	247	Ti-SF-619.13
14	265	Ti-SF-619.14
15	282	Ti-SF-619.15
16	300	Ti-SF-619.16
17	318	Ti-SF-619.17
18	336	Ti-SF-619.18

YouTube [www.TV.ortho.in/118](http://www.TV.ortho.in/118)



**RECONSTRUCTION SAFETY LOCK PLATE 4.5/5.0**

No. of Holes	Length in mm	Code S.S.
3	56	SF-149.103
4	75	SF-149.104
5	94	SF-149.105
6	113	SF-149.106
7	132	SF-149.107
8	151	SF-149.108
9	170	SF-149.109
10	189	SF-149.110
11	208	SF-149.111
12	227	SF-149.112
13	246	SF-149.113
14	265	SF-149.114
15	284	SF-149.115
16	303	SF-149.116

No. of Holes	Length in mm	Code TITANIUM
3	56	Ti-SF-149.103
4	75	Ti-SF-149.104
5	94	Ti-SF-149.105
6	113	Ti-SF-149.106
7	132	Ti-SF-149.107
8	151	Ti-SF-149.108
9	170	Ti-SF-149.109
10	189	Ti-SF-149.110
11	208	Ti-SF-149.111
12	227	Ti-SF-149.112
13	246	Ti-SF-149.113
14	265	Ti-SF-149.114
15	284	Ti-SF-149.115
16	303	Ti-SF-149.116

YouTube [www.TV.ortho.in/114](http://www.TV.ortho.in/114)



**PROXIMAL FEMORAL SAFETY LOCK PLATE 4.5/5.0/7.3**

No. of Holes	Length in mm	STAINLESS STEEL	
		Left	Right
2	139	SF-656.002L	SF-656.002R
4	175	SF-656.004L	SF-656.004R
6	211	SF-656.006L	SF-656.006R
8	247	SF-656.008L	SF-656.008R
10	283	SF-656.010L	SF-656.010R
12	319	SF-656.012L	SF-656.012R
14	355	SF-656.014L	SF-656.014R
16	391	SF-656.016L	SF-656.016R

No. of Holes	Length in mm	TITANIUM	
		Left	Right
2	139	Ti-SF-656.002L	Ti-SF-656.002R
4	175	Ti-SF-656.004L	Ti-SF-656.004R
6	211	Ti-SF-656.006L	Ti-SF-656.006R
8	247	Ti-SF-656.008L	Ti-SF-656.008R
10	283	Ti-SF-656.010L	Ti-SF-656.010R
12	319	Ti-SF-656.012L	Ti-SF-656.012R
14	355	Ti-SF-656.014L	Ti-SF-656.014R
16	391	Ti-SF-656.016L	Ti-SF-656.016R



YouTube [www.TV.ortho.in/11216](http://www.TV.ortho.in/11216)

**DISTAL FEMORAL SAFETY LOCK PLATE 4.5/5.0**

No. of Holes	Length in mm	STAINLESS STEEL	
		Left	Right
5	156	SF-635.005L	SF-635.005R
6	176	SF-635.006L	SF-635.006R
7	196	SF-635.007L	SF-635.007R
8	216	SF-635.008L	SF-635.008R
9	236	SF-635.009L	SF-635.009R
10	256	SF-635.010L	SF-635.010R
11	276	SF-635.011L	SF-635.011R
12	296	SF-635.012L	SF-635.012R
13	316	SF-635.013L	SF-635.013R

No. of Holes	Length in mm	TITANIUM	
		Left	Right
5	156	TI-SF-635.005L	TI-SF-635.005R
6	176	TI-SF-635.006L	TI-SF-635.006R
7	196	TI-SF-635.007L	TI-SF-635.007R
8	216	TI-SF-635.008L	TI-SF-635.008R
9	236	TI-SF-635.009L	TI-SF-635.009R
10	256	TI-SF-635.010L	TI-SF-635.010R
11	276	TI-SF-635.011L	TI-SF-635.011R
12	296	TI-SF-635.012L	TI-SF-635.012R
13	316	TI-SF-635.013L	TI-SF-635.013R



YouTube [www.TV.ortho.in/132](http://www.TV.ortho.in/132)

**CONDYLAR FEMORAL SAFETY LOCK PLATE 4.5/5.0/6.5**

No. of Holes	Length in mm	STAINLESS STEEL	
		Left	Right
6	170	SF-164.006L	SF-164.106R
8	206	SF-164.008L	SF-164.108R
10	242	SF-164.010L	SF-164.110R
12	278	SF-164.012L	SF-164.112R
14	314	SF-164.014L	SF-164.114R
16	350	SF-164.016L	SF-164.116R

No. of Holes	Length in mm	TITANIUM	
		Left	Right
6	170	TI-SF-164.006L	TI-SF-164.106R
8	206	TI-SF-164.008L	TI-SF-164.108R
10	242	TI-SF-164.010L	TI-SF-164.110R
12	278	TI-SF-164.012L	TI-SF-164.112R
14	314	TI-SF-164.014L	TI-SF-164.114R
16	350	TI-SF-164.016L	TI-SF-164.116R



**CONDYLAR SAFETY LOCK PLATE 4.5/5.0/7.3 - CURVED**

Left Right	No. of Holes	Length in mm	Type	Code S.S.	Code Titanium
Right	10	242	Curved	SF-629.10-24R	TI-SF-629.10-24R
Left	10	242	Curved	SF-629.10-24L	TI-SF-629.10-24L
Right	12	278	Curved	SF-629.12-78R	TI-SF-629.12-78R
Left	12	278	Curved	SF-629.12-78L	TI-SF-629.12-78L
Right	14	314	Curved	SF-629.14-31R	TI-SF-629.14-31R
Left	14	314	Curved	SF-629.14-31L	TI-SF-629.14-31L
Right	16	350	Curved	SF-629.16-35R	TI-SF-629.16-35R
Left	16	350	Curved	SF-629.16-35L	TI-SF-629.16-35L
Right	18	386	Curved	SF-629.18-38R	TI-SF-629.18-38R
Left	18	386	Curved	SF-629.18-38L	TI-SF-629.18-38L
Right	20	422	Curved	SF-629.20-42R	TI-SF-629.20-42R
Left	20	422	Curved	SF-629.20-42L	TI-SF-629.20-42L
Right	22	458	Curved	SF-629.22-45R	TI-SF-629.22-45R
Left	22	458	Curved	SF-629.22-45L	TI-SF-629.22-45L



YouTube [www.TV.ortho.in/136](http://www.TV.ortho.in/136)



**SAFETY LOCK 'T' PLATE 4.5/5.0**

No. of Holes	Length in mm	Code S.S.
3	67	SF-159.103
4	83	SF-159.104
5	99	SF-159.105
6	115	SF-159.106
7	131	SF-159.107
8	147	SF-159.108
9	163	SF-159.109

No. of Holes	Length in mm	Code TITANIUM
3	67	Ti-SF-159.103
4	83	Ti-SF-159.104
5	99	Ti-SF-159.105
6	115	Ti-SF-159.106
7	131	Ti-SF-159.107
8	147	Ti-SF-159.108
9	163	Ti-SF-159.109



YouTube [www.TV.ortho.in/11510](http://www.TV.ortho.in/11510)



**PROXIMAL TIBIAL SAFETY LOCK PLATE 4.5/5.0**

No. of Holes	Length in mm	STAINLESS STEEL	
		Left	Right
4	82	SF-654.004L	SF-654.004R
6	118	SF-654.006L	SF-654.006R
8	154	SF-654.008L	SF-654.008R
10	190	SF-654.010L	SF-654.010R
12	226	SF-654.012L	SF-654.012R
14	262	SF-654.014L	SF-654.014R
16	298	SF-654.016L	SF-654.016R

No. of Holes	Length in mm	TITANIUM	
		Left	Right
4	82	Ti-SF-654.004L	Ti-SF-654.004R
6	118	Ti-SF-654.006L	Ti-SF-654.006R
8	154	Ti-SF-654.008L	Ti-SF-654.008R
10	190	Ti-SF-654.010L	Ti-SF-654.010R
12	226	Ti-SF-654.012L	Ti-SF-654.012R
14	262	Ti-SF-654.014L	Ti-SF-654.014R
16	298	Ti-SF-654.016L	Ti-SF-654.016R

YouTube [www.TV.ortho.in/9073](http://www.TV.ortho.in/9073)



**PROXIMAL LATERAL TIBIAL SAFETY LOCK PLATE 4.5/5.0**

No. of Holes	Length in mm	STAINLESS STEEL	
		Left	Right
5	140	SF-165.105L	SF-165.005R
7	180	SF-165.107L	SF-165.007R
9	220	SF-165.109L	SF-165.009R
11	260	SF-165.111L	SF-165.011R
13	300	SF-165.113L	SF-165.013R

No. of Holes	Length in mm	TITANIUM	
		Left	Right
5	140	Ti-SF-165.105L	Ti-SF-165.005R
7	180	Ti-SF-165.107L	Ti-SF-165.007R
9	220	Ti-SF-165.109L	Ti-SF-165.009R
11	260	Ti-SF-165.111L	Ti-SF-165.011R
13	300	Ti-SF-165.113L	Ti-SF-165.013R

YouTube [www.TV.ortho.in/67](http://www.TV.ortho.in/67)



**PROXIMAL LATERAL TIBIAL SAFETY LOCK PLATE LC CUTS 4.5/5.0**

No. of Holes	Length in mm	STAINLESS STEEL	
		Left	Right
5	140	SF-165.405L	SF-165.305R
7	180	SF-165.407L	SF-165.307R
9	220	SF-165.409L	SF-165.309R
11	260	SF-165.411L	SF-165.311R
13	300	SF-165.413L	SF-165.313R

No. of Holes	Length in mm	TITANIUM	
		Left	Right
5	140	Ti-SF-165.405L	Ti-SF-165.305R
7	180	Ti-SF-165.407L	Ti-SF-165.307R
9	220	Ti-SF-165.409L	Ti-SF-165.309R
11	260	Ti-SF-165.411L	Ti-SF-165.311R
13	300	Ti-SF-165.413L	Ti-SF-165.313R

YouTube [www.TV.ortho.in/123](http://www.TV.ortho.in/123)



MEDIAL PROXIMAL TIBIAL SAFETY LOCK PLATE 4.5

No. of Holes	Length in mm	STAINLESS STEEL	
		Left	Right
4	106	SF-653.004L	SF-653.004R
6	142	SF-653.006L	SF-653.006R
8	178	SF-653.008L	SF-653.008R
10	214	SF-653.010L	SF-653.010R
12	250	SF-653.012L	SF-653.012R
14	286	SF-653.014L	SF-653.014R
16	322	SF-653.016L	SF-653.016R

No. of Holes	Length in mm	TITANIUM	
		Left	Right
4	106	Ti-SF-653.004L	Ti-SF-653.004R
6	142	Ti-SF-653.006L	Ti-SF-653.006R
8	178	Ti-SF-653.008L	Ti-SF-653.008R
10	214	Ti-SF-653.010L	Ti-SF-653.010R
12	250	Ti-SF-653.012L	Ti-SF-653.012R
14	286	Ti-SF-653.014L	Ti-SF-653.014R
16	322	Ti-SF-653.016L	Ti-SF-653.016R

YouTube [www.TV.ortho.in/170](http://www.TV.ortho.in/170)



T BUTTRESS SAFETY LOCK PLATE 4.5/5.0

No. of Holes	Length in mm	Code S.S.
4	83	SF-160.104
5	99	SF-160.105
6	115	SF-160.106
7	131	SF-160.107
8	147	SF-160.108
9	163	SF-160.109
10	179	SF-160.110
12	195	SF-160.112

No. of Holes	Length in mm	Code TITANIUM
4	83	Ti-Sf-160.104
5	99	Ti-Sf-160.105
6	115	Ti-Sf-160.106
7	131	Ti-Sf-160.107
8	147	Ti-Sf-160.108
9	163	Ti-Sf-160.109
10	179	Ti-Sf-160.110
12	195	Ti-Sf-160.112




L BUTTRESS SAFETY LOCK PLATE 4.5/5.0

No. of Holes	Length in mm	STAINLESS STEEL	
		Left	Right
3	69	SF-161.133L	SF-161.103R
4	85	SF-161.134L	SF-161.104R
5	101	SF-161.135L	SF-161.105R
6	117	SF-161.136L	SF-161.106R
7	133	SF-161.137L	SF-161.107R
8	149	SF-161.138L	SF-161.108R
9	165	SF-161.139L	SF-161.109R
10	181	SF-161.140L	SF-161.110R

No. of Holes	Length in mm	TITANIUM	
		Left	Right
3	69	SF-161.133L	SF-161.103R
4	85	SF-161.134L	SF-161.104R
5	101	SF-161.135L	SF-161.105R
6	117	SF-161.136L	SF-161.106R
7	133	SF-161.137L	SF-161.107R
8	149	SF-161.138L	SF-161.108R
9	165	SF-161.139L	SF-161.109R
10	181	SF-161.140L	SF-161.110R

YouTube [www.TV.ortho.in/24746](http://www.TV.ortho.in/24746)

YouTube [www.TV.ortho.in/24747](http://www.TV.ortho.in/24747)



**STAINLESS STEEL**

**SAFETY LOCK SCREW 5.0MM**

**TITANIUM**

Length in mm	Self Tapping	Self Tapping & Self Drilling
12	SF-500.012	SF-500.212
14	SF-500.014	SF-500.214
16	SF-500.016	SF-500.216
18	SF-500.018	SF-500.218
20	SF-500.020	SF-500.220
22	SF-500.022	SF-500.222
24	SF-500.024	SF-500.224
26	SF-500.026	SF-500.226
28	SF-500.028	SF-500.228
30	SF-500.030	SF-500.230
32	SF-500.032	SF-500.232
34	SF-500.034	SF-500.234
36	SF-500.036	SF-500.236
38	SF-500.038	SF-500.238
40	SF-500.040	SF-500.240
42	SF-500.042	SF-500.242
44	SF-500.044	SF-500.244
46	SF-500.046	SF-500.246
48	SF-500.048	SF-500.248
50	SF-500.050	SF-500.250
55	SF-500.055	SF-500.255
60	SF-500.060	SF-500.260
65	SF-500.065	SF-500.265
70	SF-500.070	SF-500.270
75	SF-500.075	SF-500.275
80	SF-500.080	SF-500.280
85	SF-500.085	SF-500.285
90	SF-500.090	SF-500.290

Drill Bit: 4.3mm

Length in mm	Self Tapping	Self Tapping & Self Drilling
12	Ti-SF-500.012	Ti-SF-500.212
14	Ti-SF-500.014	Ti-SF-500.214
16	Ti-SF-500.016	Ti-SF-500.216
18	Ti-SF-500.018	Ti-SF-500.218
20	Ti-SF-500.020	Ti-SF-500.220
22	Ti-SF-500.022	Ti-SF-500.222
24	Ti-SF-500.024	Ti-SF-500.224
26	Ti-SF-500.026	Ti-SF-500.226
28	Ti-SF-500.028	Ti-SF-500.228
30	Ti-SF-500.030	Ti-SF-500.230
32	Ti-SF-500.032	Ti-SF-500.232
34	Ti-SF-500.034	Ti-SF-500.234
36	Ti-SF-500.036	Ti-SF-500.236
38	Ti-SF-500.038	Ti-SF-500.238
40	Ti-SF-500.040	Ti-SF-500.240
42	Ti-SF-500.042	Ti-SF-500.242
44	Ti-SF-500.044	Ti-SF-500.244
46	Ti-SF-500.046	Ti-SF-500.246
48	Ti-SF-500.048	Ti-SF-500.248
50	Ti-SF-500.050	Ti-SF-500.250
55	Ti-SF-500.055	Ti-SF-500.255
60	Ti-SF-500.060	Ti-SF-500.260
65	Ti-SF-500.065	Ti-SF-500.265
70	Ti-SF-500.070	Ti-SF-500.270
75	Ti-SF-500.075	Ti-SF-500.275
80	Ti-SF-500.080	Ti-SF-500.280
85	Ti-SF-500.085	Ti-SF-500.285
90	Ti-SF-500.090	Ti-SF-500.290

Drill Bit: 4.3mm

YouTube [www.TV.ortho.in/18213](http://www.TV.ortho.in/18213)


YouTube [www.TV.ortho.in/27307](http://www.TV.ortho.in/27307)

YouTube [www.TV.ortho.in/18281](http://www.TV.ortho.in/18281)



### CANCELLOUS SAFETY LOCK SCREW 5.0MM

Length in mm	Code - Full Thread		Code - 16mm Thread		Code - 32mm Thread	
	S.S.	TITANIUM	S.S.	TITANIUM	S.S.	TITANIUM
25	SF-500.725	Ti-SF-500.725	SF-500.325	Ti-SF-500.325	-	-
30	SF-500.730	Ti-SF-500.730	SF-500.330	Ti-SF-500.330	-	-
35	SF-500.735	Ti-SF-500.735	SF-500.335	Ti-SF-500.335	-	-
40	SF-500.740	Ti-SF-500.740	SF-500.340	Ti-SF-500.340	SF-500.540	Ti-SF-500.540
45	SF-500.745	Ti-SF-500.745	SF-500.345	Ti-SF-500.345	SF-500.545	Ti-SF-500.545
50	SF-500.750	Ti-SF-500.750	SF-500.350	Ti-SF-500.350	SF-500.550	Ti-SF-500.550
55	SF-500.755	Ti-SF-500.755	SF-500.355	Ti-SF-500.355	SF-500.555	Ti-SF-500.555
60	SF-500.760	Ti-SF-500.760	SF-500.360	Ti-SF-500.360	SF-500.560	Ti-SF-500.560
65	SF-500.765	Ti-SF-500.765	SF-500.365	Ti-SF-500.365	SF-500.565	Ti-SF-500.565
70	SF-500.770	Ti-SF-500.770	SF-500.370	Ti-SF-500.370	SF-500.570	Ti-SF-500.570
75	SF-500.775	Ti-SF-500.775	SF-500.375	Ti-SF-500.375	SF-500.575	Ti-SF-500.575
80	SF-500.780	Ti-SF-500.780	SF-500.380	Ti-SF-500.380	SF-500.580	Ti-SF-500.580
85	SF-500.785	Ti-SF-500.785	SF-500.385	Ti-SF-500.385	SF-500.585	Ti-SF-500.585
90	SF-500.790	Ti-SF-500.790	SF-500.390	Ti-SF-500.390	SF-500.590	Ti-SF-500.590
95	SF-500.795	Ti-SF-500.795	SF-500.395	Ti-SF-500.395	SF-500.595	Ti-SF-500.595
100	SF-500.800	Ti-SF-500.800	SF-500.400	Ti-SF-500.400	SF-500.600	Ti-SF-500.600
105	SF-500.805	Ti-SF-500.805	SF-500.405	Ti-SF-500.405	SF-500.605	Ti-SF-500.605
110	SF-500.810	Ti-SF-500.810	SF-500.410	Ti-SF-500.410	SF-500.610	Ti-SF-500.610
115	SF-500.815	Ti-SF-500.815	SF-500.415	Ti-SF-500.415	SF-500.615	Ti-SF-500.615
120	SF-500.820	Ti-SF-500.820	SF-500.420	Ti-SF-500.420	SF-500.620	Ti-SF-500.620



**CANNULATED CANCELLOUS SAFETY LOCK SCREW 5.0MM**

Length in mm	Code - Full Thread		Code - 16mm Thread		Code - 32mm Thread	
	S.S.	TITANIUM	S.S.	TITANIUM	S.S.	TITANIUM
25	SF-501.725	Ti-SF-501.725	SF-503.325	Ti-SF-503.325	-	-
30	SF-501.730	Ti-SF-501.730	SF-503.330	Ti-SF-503.330	-	-
35	SF-501.735	Ti-SF-501.735	SF-503.335	Ti-SF-503.335	-	-
40	SF-501.740	Ti-SF-501.740	SF-503.340	Ti-SF-503.340	SF-503.540	Ti-SF-503.540
45	SF-501.745	Ti-SF-501.745	SF-503.345	Ti-SF-503.345	SF-503.545	Ti-SF-503.545
50	SF-501.750	Ti-SF-501.750	SF-503.350	Ti-SF-503.350	SF-503.550	Ti-SF-503.550
55	SF-501.755	Ti-SF-501.755	SF-503.355	Ti-SF-503.355	SF-503.555	Ti-SF-503.555
60	SF-501.760	Ti-SF-501.760	SF-503.360	Ti-SF-503.360	SF-503.560	Ti-SF-503.560
65	SF-501.765	Ti-SF-501.765	SF-503.365	Ti-SF-503.365	SF-503.565	Ti-SF-503.565
70	SF-501.770	Ti-SF-501.770	SF-503.370	Ti-SF-503.370	SF-503.570	Ti-SF-503.570
75	SF-501.775	Ti-SF-501.775	SF-503.375	Ti-SF-503.375	SF-503.575	Ti-SF-503.575
80	SF-501.780	Ti-SF-501.780	SF-503.380	Ti-SF-503.380	SF-503.580	Ti-SF-503.580
85	SF-501.785	Ti-SF-501.785	SF-503.385	Ti-SF-503.385	SF-503.585	Ti-SF-503.585
90	SF-501.790	Ti-SF-501.790	SF-503.390	Ti-SF-503.390	SF-503.590	Ti-SF-503.590
95	SF-501.795	Ti-SF-501.795	SF-503.395	Ti-SF-503.395	SF-503.595	Ti-SF-503.595
100	SF-501.800	Ti-SF-501.800	SF-503.400	Ti-SF-503.400	SF-503.600	Ti-SF-503.600
105	SF-501.805	Ti-SF-501.805	SF-503.405	Ti-SF-503.405	SF-503.605	Ti-SF-503.605
110	SF-501.810	Ti-SF-501.810	SF-503.410	Ti-SF-503.410	SF-503.610	Ti-SF-503.610
115	SF-501.815	Ti-SF-501.815	SF-503.415	Ti-SF-503.415	SF-503.615	Ti-SF-503.615
120	SF-501.820	Ti-SF-501.820	SF-503.420	Ti-SF-503.420	SF-503.620	Ti-SF-503.620



CANNULATED SAFETY LOCK SCREW 5.0MM

Length in mm	Code - Full Thread		Code - 16mm Thread		Code - 32mm Thread	
	S.S.	TITANIUM	S.S.	TITANIUM	S.S.	TITANIUM
25	SF-506.725	Ti-SF-506.725	SF-506.325	Ti-SF-506.325	-	-
30	SF-506.730	Ti-SF-506.730	SF-506.330	Ti-SF-506.330	-	-
35	SF-506.735	Ti-SF-506.735	SF-506.335	Ti-SF-506.335	-	-
40	SF-506.740	Ti-SF-506.740	SF-506.340	Ti-SF-506.340	SF-506.540	Ti-SF-506.540
45	SF-506.745	Ti-SF-506.745	SF-506.345	Ti-SF-506.345	SF-506.545	Ti-SF-506.545
50	SF-506.750	Ti-SF-506.750	SF-506.350	Ti-SF-506.350	SF-506.550	Ti-SF-506.550
55	SF-506.755	Ti-SF-506.755	SF-506.355	Ti-SF-506.355	SF-506.555	Ti-SF-506.555
60	SF-506.760	Ti-SF-506.760	SF-506.360	Ti-SF-506.360	SF-506.560	Ti-SF-506.560
65	SF-506.765	Ti-SF-506.765	SF-506.365	Ti-SF-506.365	SF-506.565	Ti-SF-506.565
70	SF-506.770	Ti-SF-506.770	SF-506.370	Ti-SF-506.370	SF-506.570	Ti-SF-506.570
75	SF-506.775	Ti-SF-506.775	SF-506.375	Ti-SF-506.375	SF-506.575	Ti-SF-506.575
80	SF-506.780	Ti-SF-506.780	SF-506.380	Ti-SF-506.380	SF-506.580	Ti-SF-506.580
85	SF-506.785	Ti-SF-506.785	SF-506.385	Ti-SF-506.385	SF-506.585	Ti-SF-506.585
90	SF-506.790	Ti-SF-506.790	SF-506.390	Ti-SF-506.390	SF-506.590	Ti-SF-506.590
95	SF-506.795	Ti-SF-506.795	SF-506.395	Ti-SF-506.395	SF-506.595	Ti-SF-506.595
100	SF-506.800	Ti-SF-506.800	SF-506.400	Ti-SF-506.400	SF-506.600	Ti-SF-506.600
105	SF-506.805	Ti-SF-506.805	SF-506.405	Ti-SF-506.405	SF-506.605	Ti-SF-506.605
110	SF-506.810	Ti-SF-506.810	SF-506.410	Ti-SF-506.410	SF-506.610	Ti-SF-506.610
115	SF-506.815	Ti-SF-506.815	SF-506.415	Ti-SF-506.415	SF-506.615	Ti-SF-506.615
120	SF-506.820	Ti-SF-506.820	SF-506.420	Ti-SF-506.420	SF-506.620	Ti-SF-506.620


Drill Bit Cannulated : 4.3mm

YouTube [www.TV.ortho.in/77](http://www.TV.ortho.in/77)




**CANCELLOUS SAFETY LOCK SCREW 6.5MM**

Length in mm	Code - Full Thread		Code - 16mm Thread		Code - 32mm Thread	
	S.S.	TITANIUM	S.S.	TITANIUM	S.S.	TITANIUM
25	SF-504.125	Ti-SF-504.125	SF-501.325	Ti-SF-501.325	-	-
30	SF-504.130	Ti-SF-504.130	SF-501.330	Ti-SF-501.330	-	-
35	SF-504.135	Ti-SF-504.135	SF-501.335	Ti-SF-501.335	-	-
40	SF-504.140	Ti-SF-504.140	SF-501.340	Ti-SF-501.340	SF-501.540	Ti-SF-501.540
45	SF-504.145	Ti-SF-504.145	SF-501.345	Ti-SF-501.345	SF-501.545	Ti-SF-501.545
50	SF-504.150	Ti-SF-504.150	SF-501.350	Ti-SF-501.350	SF-501.550	Ti-SF-501.550
55	SF-504.155	Ti-SF-504.155	SF-501.355	Ti-SF-501.355	SF-501.555	Ti-SF-501.555
60	SF-504.160	Ti-SF-504.160	SF-501.360	Ti-SF-501.360	SF-501.560	Ti-SF-501.560
65	SF-504.165	Ti-SF-504.165	SF-501.365	Ti-SF-501.365	SF-501.565	Ti-SF-501.565
70	SF-504.170	Ti-SF-504.170	SF-501.370	Ti-SF-501.370	SF-501.570	Ti-SF-501.570
75	SF-504.175	Ti-SF-504.175	SF-501.375	Ti-SF-501.375	SF-501.575	Ti-SF-501.575
80	SF-504.180	Ti-SF-504.180	SF-501.380	Ti-SF-501.380	SF-501.580	Ti-SF-501.580
85	SF-504.185	Ti-SF-504.185	SF-501.385	Ti-SF-501.385	SF-501.585	Ti-SF-501.585
90	SF-504.190	Ti-SF-504.190	SF-501.390	Ti-SF-501.390	SF-501.590	Ti-SF-501.590
95	SF-504.195	Ti-SF-504.195	SF-501.395	Ti-SF-501.395	SF-501.595	Ti-SF-501.595
100	SF-504.200	Ti-SF-504.200	SF-501.400	Ti-SF-501.400	SF-501.600	Ti-SF-501.600
105	SF-504.205	Ti-SF-504.205	SF-501.405	Ti-SF-501.405	SF-501.605	Ti-SF-501.605
110	SF-504.210	Ti-SF-504.210	SF-501.410	Ti-SF-501.410	SF-501.610	Ti-SF-501.610
115	SF-504.215	Ti-SF-504.215	SF-501.415	Ti-SF-501.415	SF-501.615	Ti-SF-501.615
120	SF-504.220	Ti-SF-504.220	SF-501.420	Ti-SF-501.420	SF-501.620	Ti-SF-501.620



CANNULATED CANCELLOUS SAFETY LOCK SCREW 6.5MM

Length in mm	Code - Full Thread		Code - 16mm Thread		Code - 32mm Thread	
	S.S.	TITANIUM	S.S.	TITANIUM	S.S.	TITANIUM
25	SF-504.325	Ti-SF-504.325	SF-507.325	Ti-SF-507.325	-	-
30	SF-504.330	Ti-SF-504.330	SF-507.330	Ti-SF-507.330	-	-
35	SF-504.335	Ti-SF-504.335	SF-507.335	Ti-SF-507.335	-	-
40	SF-504.340	Ti-SF-504.340	SF-507.340	Ti-SF-507.340	SF-507.540	Ti-SF-507.540
45	SF-504.345	Ti-SF-504.345	SF-507.345	Ti-SF-507.345	SF-507.545	Ti-SF-507.545
50	SF-504.350	Ti-SF-504.350	SF-507.350	Ti-SF-507.350	SF-507.550	Ti-SF-507.550
55	SF-504.355	Ti-SF-504.355	SF-507.355	Ti-SF-507.355	SF-507.555	Ti-SF-507.555
60	SF-504.360	Ti-SF-504.360	SF-507.360	Ti-SF-507.360	SF-507.560	Ti-SF-507.560
65	SF-504.365	Ti-SF-504.365	SF-507.365	Ti-SF-507.365	SF-507.565	Ti-SF-507.565
70	SF-504.370	Ti-SF-504.370	SF-507.370	Ti-SF-507.370	SF-507.570	Ti-SF-507.570
75	SF-504.375	Ti-SF-504.375	SF-507.375	Ti-SF-507.375	SF-507.575	Ti-SF-507.575
80	SF-504.380	Ti-SF-504.380	SF-507.380	Ti-SF-507.380	SF-507.580	Ti-SF-507.580
85	SF-504.385	Ti-SF-504.385	SF-507.385	Ti-SF-507.385	SF-507.585	Ti-SF-507.585
90	SF-504.390	Ti-SF-504.390	SF-507.390	Ti-SF-507.390	SF-507.590	Ti-SF-507.590
95	SF-504.395	Ti-SF-504.395	SF-507.395	Ti-SF-507.395	SF-507.595	Ti-SF-507.595
100	SF-504.400	Ti-SF-504.400	SF-507.400	Ti-SF-507.400	SF-507.600	Ti-SF-507.600
105	SF-504.405	Ti-SF-504.405	SF-507.405	Ti-SF-507.405	SF-507.605	Ti-SF-507.605
110	SF-504.410	Ti-SF-504.410	SF-507.410	Ti-SF-507.410	SF-507.610	Ti-SF-507.610
115	SF-504.415	Ti-SF-504.415	SF-507.415	Ti-SF-507.415	SF-507.615	Ti-SF-507.615
120	SF-504.420	Ti-SF-504.420	SF-507.420	Ti-SF-507.420	SF-507.620	Ti-SF-507.620



**CANNULATED CANCELLOUS SAFETY LOCK SCREW 7.3MM**

Length in mm	Code - Full Thread		Code - 16mm Thread		Code - 32mm Thread	
	S.S.	TITANIUM	S.S.	TITANIUM	S.S.	TITANIUM
25	SF-505.725	Ti-SF-505.725	SF-505.325	Ti-SF-505.325	-	-
30	SF-505.730	Ti-SF-505.730	SF-505.330	Ti-SF-505.330	-	-
35	SF-505.735	Ti-SF-505.735	SF-505.335	Ti-SF-505.335	-	-
40	SF-505.740	Ti-SF-505.740	SF-505.340	Ti-SF-505.340	SF-505.540	Ti-SF-505.540
45	SF-505.745	Ti-SF-505.745	SF-505.345	Ti-SF-505.345	SF-505.545	Ti-SF-505.545
50	SF-505.750	Ti-SF-505.750	SF-505.350	Ti-SF-505.350	SF-505.550	Ti-SF-505.550
55	SF-505.755	Ti-SF-505.755	SF-505.355	Ti-SF-505.355	SF-505.555	Ti-SF-505.555
60	SF-505.760	Ti-SF-505.760	SF-505.360	Ti-SF-505.360	SF-505.560	Ti-SF-505.560
65	SF-505.765	Ti-SF-505.765	SF-505.365	Ti-SF-505.365	SF-505.565	Ti-SF-505.565
70	SF-505.770	Ti-SF-505.770	SF-505.370	Ti-SF-505.370	SF-505.570	Ti-SF-505.570
75	SF-505.775	Ti-SF-505.775	SF-505.375	Ti-SF-505.375	SF-505.575	Ti-SF-505.575
80	SF-505.780	Ti-SF-505.780	SF-505.380	Ti-SF-505.380	SF-505.580	Ti-SF-505.580
85	SF-505.785	Ti-SF-505.785	SF-505.385	Ti-SF-505.385	SF-505.585	Ti-SF-505.585
90	SF-505.790	Ti-SF-505.790	SF-505.390	Ti-SF-505.390	SF-505.590	Ti-SF-505.590
95	SF-505.795	Ti-SF-505.795	SF-505.395	Ti-SF-505.395	SF-505.595	Ti-SF-505.595
100	SF-505.800	Ti-SF-505.800	SF-505.400	Ti-SF-505.400	SF-505.600	Ti-SF-505.600
105	SF-505.805	Ti-SF-505.805	SF-505.405	Ti-SF-505.405	SF-505.605	Ti-SF-505.605
110	SF-505.810	Ti-SF-505.810	SF-505.410	Ti-SF-505.410	SF-505.610	Ti-SF-505.610
115	SF-505.815	Ti-SF-505.815	SF-505.415	Ti-SF-505.415	SF-505.615	Ti-SF-505.615
120	SF-505.820	Ti-SF-505.820	SF-505.420	Ti-SF-505.420	SF-505.620	Ti-SF-505.620

Drill Bit Cannulated : 4.5mm

□ LARGE FRAGMENT LOCK SCREWS SET



YouTube [www.TV.ortho.in/27595](http://www.TV.ortho.in/27595)

**AS137 LARGE FRAGMENT LOCK SCREWS SET**

Set Consisting of:

			Units
AS435.002	Drill Bit	4.3mm, Length 221mm, Callibrated, 2-flute, for Quick Coupling	2
AS435.006		4.3mm Drill Sleeve for 5.0mm Safety Lock Screws	2
SF-500.014	Safety Lock Screw	5.0mm, Self Tapping, SS, L:14mm	2
SF-500.016	Safety Lock Screw	5.0mm, Self Tapping, SS, L:16mm	2
SF-500.018	Safety Lock Screw	5.0mm, Self Tapping, SS, L:18mm	2
SF-500.020	Safety Lock Screw	5.0mm, Self Tapping, SS, L:20mm	2
SF-500.022	Safety Lock Screw	5.0mm, Self Tapping, SS, L:22mm	2
SF-500.024	Safety Lock Screw	5.0mm, Self Tapping, SS, L:24mm	2
SF-500.026	Safety Lock Screw	5.0mm, Self Tapping, SS, L:26mm	2
SF-500.028	Safety Lock Screw	5.0mm, Self Tapping, SS, L:28mm	2
SF-500.030	Safety Lock Screw	5.0mm, Self Tapping, SS, L:30mm	2
SF-500.032	Safety Lock Screw	5.0mm, Self Tapping, SS, L:32mm	2
SF-500.034	Safety Lock Screw	5.0mm, Self Tapping, SS, L:34mm	4
SF-500.036	Safety Lock Screw	5.0mm, Self Tapping, SS, L:36mm	4
SF-500.038	Safety Lock Screw	5.0mm, Self Tapping, SS, L:38mm	4
SF-500.040	Safety Lock Screw	5.0mm, Self Tapping, SS, L:40mm	4
SF-500.042	Safety Lock Screw	5.0mm, Self Tapping, SS, L:42mm	2
SF-500.044	Safety Lock Screw	5.0mm, Self Tapping, SS, L:44mm	2
SF-500.046	Safety Lock Screw	5.0mm, Self Tapping, SS, L:46mm	2
SF-500.048	Safety Lock Screw	5.0mm, Self Tapping, SS, L:48mm	2
SF-500.050	Safety Lock Screw	5.0mm, Self Tapping, SS, L:50mm	2
SF-500.055	Safety Lock Screw	5.0mm, Self Tapping, SS, L:55mm	2
SF-500.060	Safety Lock Screw	5.0mm, Self Tapping, SS, L:60mm	2
SF-500.065	Safety Lock Screw	5.0mm, Self Tapping, SS, L:65mm	2
SF-500.070	Safety Lock Screw	5.0mm, Self Tapping, SS, L:70mm	2
SF-500.214	Safety Lock Screw	5.0mm, ST/SD, SS, L:14mm	2
SF-500.218	Safety Lock Screw	5.0mm, ST/SD, SS, L:18mm	2
SF-500.222	Safety Lock Screw	5.0mm, ST/SD, SS, L:22mm	2
SF-500.226	Safety Lock Screw	5.0mm, ST/SD, SS, L:26mm	4
SF-500.230	Safety Lock Screw	5.0mm, ST/SD, SS, L:30mm	4
SF-500.236	Safety Lock Screw	5.0mm, ST/SD, SS, L:36mm	4
SF-500.240	Safety Lock Screw	5.0mm, ST/SD, SS, L:40mm	4
SF-500.244	Safety Lock Screw	5.0mm, ST/SD, SS, L:44mm	2
SF-500.250	Safety Lock Screw	5.0mm, ST/SD, SS, L:50mm	2
SF-500.255	Safety Lock Screw	5.0mm, ST/SD, SS, L:55mm	2
925.607	Empty Container		

□ LARGE FRAGMENT LOCK SCREWS SET - COMBO



**AS143 LARGE FRAGMENT SCREWS SET - COMBO WITH SAFETY LOCK SCREWS & STANDARD SCREWS**

Set Consisting of:

	Units
<b>Cortex Screw 4.5mm, Self Tapping, SS</b>	
5 each of 22,24,26,28mm	20
10 each of 30,32,34,36,38,40,42mm	70
5 each of 44,46,48,50,52,54,56,58,60,65,70,75,80,85,90mm	75

<b>Safety Lock Screw 5.0mm, Self Tapping, SS</b>	
4 each of 22,24,26,28	16
8 each of 30,32,34,36,38,40,42	56
4 each of 44,46,48,50,52,54,56,58mm	32
4 each of 60,65,70,75,80,85,90mm	28

<b>Cancellous Bone Screw 6.5mm, 16mm Thread, SS</b>	
2 each of 30,35,40,45,50,55,60,65,70,75,80,85,90,95,100mm	30

<b>Cancellous Bone Screw 6.5mm, 32mm Thread, SS</b>	
2 each of 40,45,50,55,60,65,70,75,80,85,90,95,100,105,110mm	30

<b>Cancellous Bone Screw 6.5mm, Fully Thread, SS</b>	
2 each of 30,35,40,45,50,55,60,65,70,75,80,85,90,95,100mm	30

AS143.000 Empty Container



**AS144 LARGE FRAGMENT SCREWS SET - COMBO WITH SAFETY LOCK & STANDARD SCREWS - MODIFIED**

Set Consisting of:

	Units
<b>Cortex Screw 4.5mm, Self Tapping, SS</b>	
2 each of 18,20,22,24,26,28,30,32,34,36,38,40mm	24
2 each of 42,44,46,48,50,52,54,56,58,60,62,64,66mm	26

<b>Safety Lock Screw 5.0mm, Self Tapping, SS</b>	
6 each of 20,22,24,26,28,30,32,34,36,38,40,42,44,46,48mm	90
4 each of 50,55,60,65,70,75,80,85,90mm	36


<b>Safety Lock Screw 5.0mm, Self Tapping &amp; Self Drilling, SS</b>	
2 each of 20,22,24,26,28,30,32,34,36,38,40,42,44,46,48,50mm	32
2 each of 55,60,65,70,75,80,85,90mm	16

<b>Cancellous Bone Screw 6.5mm, Fully Thread, SS</b>	
2 each of 30,35,40,45,50,55,60,65,70,75,80,85,90,95,100mm	30

AS144.000 Empty Container



□ LARGE FRAGMENT LOCK PLATES SET





**AS127 LARGE FRAGMENT SAFETY LOCK PLATES SET**

Set Consisting of:

		Units
SF-150.106	Narrow LC- DCP 4.5/5.0, Safety Lock Plate, SS, 6 Holes	1
SF-150.107	Narrow LC- DCP 4.5/5.0, Safety Lock Plate, SS, 7 Holes	1
SF-150.108	Narrow LC- DCP 4.5/5.0, Safety Lock Plate, SS, 8 Holes	1
SF-150.109	Narrow LC- DCP 4.5/5.0, Safety Lock Plate, SS, 9 Holes	1
SF-150.110	Narrow LC- DCP 4.5/5.0, Safety Lock Plate, SS, 10 Holes	1
SF-150.111	Narrow LC- DCP 4.5/5.0, Safety Lock Plate, SS, 11 Holes	1
SF-150.112	Narrow LC- DCP 4.5/5.0, Safety Lock Plate, SS, 12 Holes	1
SF-160.104	T Buttress Safety Lock Plate 4.5/5.0, SS, 3 Holes	1
SF-160.105	T Buttress Safety Lock Plate 4.5/5.0, SS, 4 Holes	1
SF-159.104	Safety Lock T Plate 4.5/5.0, SS, 4 Holes	1
SF-159.105	Safety Lock T Plate 4.5/5.0, SS, 5 Holes	1
SF-159.106	Safety Lock T Plate 4.5/5.0, SS, 6 Holes	1
SF-161.134L	L Buttress Safety Lock Plate 4.5/5.0, Left, SS, 4 Holes	1
SF-161.135L	L Buttress Safety Lock Plate 4.5/5.0, Left, SS, 5 Holes	1
SF-161.136L	L Buttress Safety Lock Plate 4.5/5.0, Left, SS, 6 Holes	1
SF-161.104R	L Buttress Safety Lock Plate 4.5/5.0, Right, SS, 4 Holes	1
SF-161.105R	L Buttress Safety Lock Plate 4.5/5.0, Right, SS, 5 Holes	1
SF-161.106R	L Buttress Safety Lock Plate 4.5/5.0, Right, SS, 6 Holes	1
SF-149.106	Reconstruction Safety Lock Plate 4.5/5.0, SS, 6 Holes	1
SF-149.108	Reconstruction Safety Lock Plate 4.5/5.0, SS, 8 Holes	1
SF-149.110	Reconstruction Safety Lock Plate 4.5/5.0, SS, 10 Holes	1
SF-149.112	Reconstruction Safety Lock Plate 4.5/5.0, SS, 12 Holes	1
SF-149.107	Reconstruction Safety Lock Plate 4.5/5.0, SS, 7 Holes	1
SF-149.109	Reconstruction Safety Lock Plate 4.5/5.0, SS, 9 Holes	1
SF-149.111	Reconstruction Safety Lock Plate 4.5/5.0, SS, 11 Holes	1
925.608	Empty Container	

□ LARGE FRAGMENT INSTRUMENT SET - COMBO



**AS145 LARGE FRAGMENT INSTRUMENT SET - COMBO WITH STANDARD INSTRUMENTS & LOCKING INSTRUMENTS**

Set Consisting of:

Units

STANDARD INSTRUMENTS		
260.320	Drill Bit (Quick Coupling), 3.2mm x 130mm	1
260.450	Drill Bit (Quick Coupling), 4.5mm x 130mm	1
305.030	Counter Sink for 4.5 & 6.5mm Screws (8.0mm Head)	1
306.010	Quick Coupling Handle (T-Type)	1
264.040	Bone Tap with Quick Coupling for 4.5mm Cortex Screws, L:130mm	1
264.050	Bone Tap with Quick Coupling for 6.5mm Cancellous Screws	1
307.095	Insert Drill Sleeve 4.5/3.2mm	1
307.070	Drill & Tap Sleeve combined 3.2/4.5mm	1
307.080	Drill & Tap Sleeve combined 3.2/6.5mm	1
301.070	Large Hexagonal Screw Driver Shaft	1
301.050	Large Hexagonal Screw Driver	1
303.020	Depth Gauge for 4.5mm & 6.5mm Screws	1
371.000	Sharp Hook	1
304.020	Large Neutral & Loaded Drill Guide	1
308.000	Pointed Drill Guide	1
302.020	Large Plate Bender (pair)	1
LOCKING INSTRUMENTS		
AS435.002	Drill Bit 4.3mm, Length 221mm, Calibrated, 2-flute, for Quick Coupling	2
AS435.006	4.3mm Drill Sleeve for 5.0mm Safety Lock Screws	2
AS435.010	Torque limiting Screw Driver, 3.5mm Tip	1
AS435.022	Extraction Screw (Left Hand Thread), Conical, for Screws 4.5mm to 6.5mm	1
AS435.026	Centering Sleeve for Kirschner Wire 2.0mm, Length 110mm	1
AS145.000	Empty Container	

□ LARGE FRAGMENT LOCK PLATES INSTRUMENTS SET

YouTube [www.TV.ortho.in/27600](http://www.TV.ortho.in/27600)

**AS435 LARGE FRAGMENT SAFETY LOCK INSTRUMENTS SET**

*Set Consisting of:*

		Units
AS435.002	Drill Bit 4.3mm, Length 221 mm, Calibrated, 2-flute, for Quick Coupling	2
AS435.006	4.3mm Drill Sleeve for 5.0mm Safety Lock Screws	1
AS435.010	Torque Limiting Screw Driver, 3.5mm Tip	1
AS435.014	T-Handle for Quick Coupling	1
AS435.018	HSS Drill Bit 3.5mm (Blue Coating)	2
AS435.022	Extraction Screw (Left Hand Thread), conical, for Screws 4.5mm to 6.5mm	1
AS432.026	Centering Sleeve for Kirschner Wire 2.0mm, Length 125mm	1
AS435.030	Empty Container	

**OPTIONAL INSTRUMENTS**

*(To be Ordered Separately)*

AS435.034	Cannulated Drill Bit 5.0mm, for 7.3 Cannulated Cancellous Lock Screws, SS
AS435.038	Threaded Drill Sleeve for 7.3mm Cannulated Cancellous Lock Screws
AS435.042	Cannulated Drill Bit 4.3mm, for 5.0 Cannulated Lock Screws, SS
AS435.046	Threaded Drill Sleeve for 5.0mm Cannulated Lock Screws
AS435.050	Threaded Drill Sleeve for 6.5mm Cancellous Lock Screws



**DRILL BIT Ø 4.3MM, LENGTH 221MM, CALIBRATED, 2-FLUTE, FOR QUICK COUPLING**  
Code No.: AS435.002



**HSS DRILL BIT Ø 3.5MM (BLUE COATING)**  
Code No.: AS435.018



**EXTRACTION SCREW (LEFT HAND THREAD), CONICAL, FOR SCREWS Ø 4.5MM TO 6.5MM**  
Code No.: AS435.022



**CENTERING SLEEVE FOR KIRSCHNER WIRE Ø 2.0MM, LENGTH 125MM**  
Code No.: AS432.026



**4.3MM DRILL SLEEVE FOR 5.0MM SAFETY LOCK SCREWS**  
Code No.: AS435.006



**TORQUE LIMITING SCREW DRIVER, 3.5MM TIP**  
Code No.: AS435.010



**QUICK COUPLING HANDLE (T-Type)**  
Code No.: AS435.014