



**Kiran TECHNO SERVICES**  
PRIVATE LIMITED



# KIRAN SURGICAL SYSTEMS



[www.kiransurgicalsystems.com](http://www.kiransurgicalsystems.com)

[www.kiranmedicareequipments.com](http://www.kiranmedicareequipments.com)



**Model No. KSS - 101**

### ELECTRIC ICU BED MULTIFUNCTION

#### Specification:


- Frame work made of Rectangular M.S Tube
- MS Mattress platform in four sections
- Back Rest, Knee Rest, TB/RTB and Hi-Lo Positions obtained by electric actuator operator by hand remote system.
- ABS molded head and foot bows.
- ABS track away safety side Railing.
- Telescopic IV Rod with Two Locations.
- Bed Mounted on 125 mm dia Castor, two with brakes.
- Nurse control panel
- CPR system
- Overall Size: 210 L x 90W x 50 -70Hcms
- Finish: Epoxy Powder Coated

#### Optional

- Bed top also available in metal sheet
- Central locking system also available
- Chair position in bed also available



### Positions :





**Model No. KSS - 102**

**ELECTRIC ICU BED (ABS PANELS & RAILING)**

**Specification:**

- Frame work made of Rectangular M.S Tube.
- MS Mattress platform in four sections.
- Back Rest, Knee Rest and Hi-Lo Positions obtained by electric actuator operator by hand remote system.
- ABS molded head and foot bows.
- ABS track away safety side Railing.
- Bed Mounted on 125 mm Dia Castor with two braking wheels.
- Telescopic IV Rod with Two Locations.
- Overall Size: 210 L x 90W x 50 -70Hcms
- Finish: Epoxy Powder Coated

**Optional**

- Central locking system also available
- Nurse control panel
- CPR system



**Model No. KSS - 103**

**ELECTRIC ICU BED (ABS PANELS)**

**Specification:**

- Frame work made of rectangular M.S. tube perforated CRC sheet top in four sections.
- Back rest, knee rest and hi-lo positions obtained by electric actuator operate by hand remote system.
- ABS molded head and foot bows.
- Aluminum collapsible safety side railings.
- Telescopic IV rod with four locations.
- Bed mounted on 125 mm diameter castor-two with brakes.
- Provided with 100 mm. mattress in 40 density.
- Overall size: 210 (L) x 90 (W) x 60-80 (H) cm.
- Finish: Epoxy powder coated.





**Model No. KSS - 104**  
**ICU BED HI-LO HYDRAULIC**

**Specification:**

- Frame work made of Rectangular M.S Tube.
- 4 Section Top made of perforated M.S. Sheet
- BackRest, KneeRest, TB / RTB Positions manoeuvred by Separate screw from Foot End.
- Height Adjustment with hydraulic foot pump.
- ABS moulded head and foot panels.
- Collapsible Aluminium Safety side Railings.
- SS IV Rod with four Locations.
- Bed mounted on 150 mm Dia Castor (2 with brakes).
- Overall size 2150L x 1050W x 460-750H cms.
- Finish: Epoxy Powder Coated



**Model No. KSS - 105**  
**MECHANICAL ICU BED (ABS PANEL)**

**Specification:**

- Frame work made of rectangular M.S. tube 4 section top made of perforated M.S. sheet.
- Backrest, knee rest, TB/RTB and hi-lo positions maneuvered by separate.
- ABS molded Head & foot bows.
- Collapsible safety side railing.
- Telescopic IV rod with four locations.
- Bed mounted on 125 mm diameter castor two with brakes.
- Overall Size: 2100 (L) x 900 (W) x 600-800 (H) mm
- Finish: Epoxy powder coated





**Model No. KSS - 106**  
**MECHANICAL ICU BED**

**Specification:**

- Frame work made of rectangular M.S. tube 4 section top made of perforated M.S. sheet.
- Backrest, knee rest, TB/RTB and hi-lo positions maneuvered by separate screw.
- S.S Head & foot bows.
- Collapsible railing (Aluminium or Iron or Stainless steel).
- Telescopic IV rod with four locations.
- Bed mounted on 125 mm diameter castor two with brakes.
- Overall Size: 2100 (L) x 900 (W) x 60-800 (H) mm
- Finish: Epoxy powder coated

**Model No. KSS - 107**  
**ELECTRIC FOWLER BED**

**Specification:**

- Frame Work made of rectangular M.S. tubes.
- Overall approx size : 2030mm x 900mm x 600mm.
- Perforated CRCA sheet top in four sections.
- Backrest & Knee rest positions obtained by electric actuator operator by hand remote system.
- Telescopic IV Rod with Four Locations
- Bed mounted on 100mm deluxe castors, two with brakes.
- ABS moulded head and foot bows.
- ABS moulded safety side railings or collapsible railing.
- Pre Treated & Epoxy Powder Coated Finish.





**Model No. KSS - 108**  
**FWLER BED DELUXE (ABS PANEL)**

**Specification:**

- Hospital Fowler Beds made of high quality materials, components and accessories.
- Hospital Fowler Beds has four section sheet metal top.
- Hospital Fowler Beds with adjustable back section and knee-rest.
- Manually operated crank system for various positions.
- High quality castors.
- Standard dimension: bed frame 2030mm L x 900 mm W x 600 mm H (approx).
- Pre-treated and powder coated.



**Model No. KSS - 109**  
**FWLER BED (DELUXE)**

**Specification:**

- Fowler Bed made of high quality materials, components and accessories.
- Four section sheet metal top.
- Fowler Bed with adjustable back section and knee-rest.
- Manually operated crank system for various positions.
- Standard: SS panels, sheet metal platform.
- High quality castors.
- Standard dimension: bed frame 2030mm L x 900 mm W x 600 mm H (approx).
- Pre-treated and powder coated Fowler Bed.



**Model No. KSS - 110**  
**FWLER BED**

**Specification:**

- Fowler Beds made of high quality materials, components and accessories.
- Four section sheet metal top.
- Fowler Beds with adjustable back section and knee-rest.
- Manually operated crank system for various positions.
- Standard: P/C Frame, sheet metal platform.
- High quality castors.
- Standard dimension: bed frame 2030mm L x 900 mm W x 600 mm H (approx).
- Pre-treated and powder coated Fowler Beds.



**Model No. KSS - 111**  
**SEMI-FOWLER BED (DLX)**


**Specification:**

- Size: 206 L x 91 W x 61 H Cms.
- Frame Work of specially designed CRC Rectangular heavy Tubes.
- Two Section Top of Uniformly Perforated CRC Sheets.
- Backrest maneuvered by screw, at Legend. , S.S. Bows are of equal height fitted with Laminated Panel.
- Two Location Saline Hole.
- Pretreated and Epoxy Powder Coated.

**Model No. KSS - 112**  
**SEMI FOWLER BED (ABS PANEL)**

**Specification:**

- Frame work made of Rectangular M.S. Tube.
- 2 Section Top made of Perforated M.S. Sheet.
- Back rest Position manoeuvred by Separate screw from foot end.
- ABS Moulded Head & Foot Panels.
- Four Location for I.V Rod.
- Overall Size : 206L x 90 W x 60 H cms.
- Finish : Epoxy Powder Coated.
- Also Available Electrically Operated



**Model No. KSS - 113**  
**SEMI-FOWLER BED (STD)**

**Specification:**

- Size: 203 L x 90 W x 60 H Cm.
- Mainframe of rectangular pipe.
- Two section top of CRCA sheet.
- Backrest section maneuvered by crank mechanism from the foot end.
- Tubular head and foot bows of unequal height mounted on rubber stump.
- Provision for saline location.
- Pre-treated & powder coated

**Model No. KSS - 114**

**ANTI-DECUBITUS MATTRESS (LATTICE STYLE)**



**Features**

- This air mattress is specially made for short-term patients.
- It is simple in structure, easy in operation and is better in prevention and cure of bed sores.

**Model No. KSS - 116**

**MATTRESS FOR HOSPITAL BEDS**

**Specification:**

- 10 cms thick, 32 / 40 density foam mattress covered with rexine closed with zip at one end.



A -Single mattress for plain beds / Attendent beds



B-Two section mattress for semi-fowler beds



C-Four section mattress for fowler/I.C.U. beds



**Model No. KSS - 115**

**X-RAY PERMEABLE BACK RESET SECTION**

**Specification:**

X-Ray permeable back rest with S.S. X-Ray cassette holder on sliding mechanism.



**Model No. KSS - 117**  
**BALKEN BEAM SET**

**Specification:**

Made of MS/SS vertical longitudinal and transverse tubes with clamps and pulleys.



**Model No. KSS - 118**  
**LIFTING POLE**

**Specification:**

M.S./S.S. tubular pole fitted with hand grip. Easy fitted back side of the head bows.



**Model No. KSS - 119**  
**PATIENT CASE SHEET HOLDER**

**Specification:**

Made from Aluminium/S.S. / ABS sheet.

- A) Single with chrome plated/S.S. clip..
- B) Double (file type).



**Model No. KSS - 120**  
**PROTECTIVE RAILING**

**Specification:**

- Made of Aluminium tubes.
- Packing : 1pair 0.022CBM.
- N.W. : 6.4 KGS.
- G.W. : 7.0 KGS.
- Finish : Epoxy Powder Coated.

**Model No. KSS - 121**

**MOSQUITO NET POLES**

**Specification:**

- Four vertical M.S. tubes of 19mm, two longitudinal & two transverse tubes of 15.8mm, knocked down condition, Pre-treated and epoxy powder coated.



**Model No. KSS - 122**

**ABS PANNEL WITH BUFFER (DELUXE)**



**Model No. KSS - 124**

**I.V. ROD FOR HOSPITAL BED**



**Model No. KSS - 123**

**ABS TUCK AWAY RAILING IN 2PCS.**



**Model No. KSS - 125**

**ABS TUCK AWAY RAILING IN 4 PCS**





**Model No. KSS - 126**  
**HOSPITAL BED PLAIN (STD.)**

**Specification:**

- Frame work made of Rectangular square M.S. Tube.
- Top made of Perforated CRCA MS sheet.
- Tubular Head & Foot Bows of Unequal Height.
- Four Location for I.V Rod.
- Overall Size : 206 L x 90 W x 60 H cms.
- Finish : Epoxy Powder Coated.

**Model No. KSS - 127**

**ATTENDANT BED**

**Specification:**

- Overall Size : 180 L X 60 W x 46 H cms.
- Frame work made of rectangular / square MS Tube.
- 5 cm Foam Padded Top Covered with rexine.
- Legs fitted with PVC Stumps.
- Finish : Epoxy Powder Coated.



**Model No. KSS - 128**  
**PEDIATRIC BED**

**Specification:**

- Overall Approx. Dimension : 140 cmL x 76 cmW x 60cmH.
- Frame work made of CRCA rectangular tubes.
- Top made of CRCA perforated Sheet.
- Full length drop side rails on both sides.
- Finish : Epoxy Powder Coated.



**Model No. KSS - 129**

**ATTENDANT COUCH CUM CHAIR  
(WITH ARM REST)**

**Specification:**

**Overall size:**

Chair position : 770 x 680 x 1125 mm

Bed position : 2050 x 680 x 550 mm

- Mild steel tubular frame work mounted
- Easily convertible to chair
- Castors to facilitate couch position
- High quality durable mattress in three sections.
- Compact construction to save space
- High quality form mattress

**Model No. KSS - 130**

**CHOLERA BED**

**Specification:**

- Cholera Bed Frame Work Constructed With Heavy Duty Ms Steel Tube.
- Deck of Cholera bed Made with PVC backing with provision of hole in the middle Of the Deck.
- Carrying Capacity 160kg.
- Fast and Easy Assembly.
- Provided With Sanitary Bucket.
- Provision of iv rod.
- Pre-Treated and Epoxy Powder Coated Finishing.



### Model No. KSS - 131

#### ELECTRIC EXAMINATION TABLE DELUXE

##### Specification:

Electrically controlled, multi function, variable height couch has been designed to meet the ever growing needs of the busy practice. This new design has evolved from extensive user group feedback.

Size: Width : 600 mm, Length: 1830 mm, Height 560-1000 mm

##### Features :

- Variable height
- Electric hand-set control
- Adjustable back rest section through ratchet system
- Swivel castors two with brakes, for ease of manoeuvrability and safety
- Comfortable Upholster waterproof foam surface
- Paper roll holder on foot end
- Fibber covered base

##### Optional Accessories :

- Safety side railing



### Model No. KSS - 132

#### ELECTRIC EXAMINATION TABLE 3 SECTION

##### Specification:

- Height adjustment by electromotive drive.
- The electric motor is controlled by 1 foot switch (possible option: Hand switch).
- Head part step less up - (+50°) and downwards (-45°) adjustable by gas spring.
- Length head part: 450 mm.
- Middle part fixed. Length: 1.290 mm.
- Foot part step less adjustable from downwards vertical position (-90°) to upwards position (+28°) by gas spring.
- Length foot part: 450 mm.
- Upholstery base: Laminated wood.

Size: Width: 600 mm - Length 1830 mm - Height 560-1000 mm





**Model No. KSS - 133**  
**ELECTRIC EXAMINATION TABLE 2 SECTION**

**Specification:**

- Height adjustable from 560-1000mm
- Scissor action provides vertical lift
- Mains powered electric actuator (240 volt electric)
- Hand switch operation
- Head section electrically adjustable from 0-70 deg.
- Support the roll of paper mounted foot end
- 125 mm swivel castors with brake
- 190 kg load maximum

Size: Width: 600 mm - Length 1830 mm - Height 560-1000 mm

**Model No. KSS - 134**  
**ELECTRIC EXAMINATION TABLE**

**Specification:**

- Height adjustable from 560-1000mm
- Scissor action provides vertical lift
- Mains powered electric actuator (240 volt electrical supply)
- Hand or foot switch operation
- 125mm swivel braking castors
- Maximum load 190kg

**Additional accessories:**

- Paper roll holder fitted to foot end

Size: Width: 600 mm - Length 1830 mm - Height 560-1000 mm



**Model No. KSS - 135**  
**HYDRAULIC EXAMINATION TABLE DELUXE**

**Specification:**

- Back rest position can be adjusted by a crank
- Hi-low position can be adjusted by a high quality hydraulic jack
- S.S swing away safety side railing.
- 150mm diameter castors, with directional positions locks
- Capability of attaching IV stand
- Finishing treatment: pre-treated powder coating
- Bumpers at four corners
- Paper holder roller mount on head side
- Polyurethane rebound foam mattress with water proof cover

**Dimensions :**

|                   |   |                    |
|-------------------|---|--------------------|
| Overall           | : | 1900mm L x 830mm W |
| Mattress platform | : | 183mm L x 720mm W  |
| Hi-Low            | : | 590mm - 1020mm     |
| Back rest angle   | : | 75°                |
| Bed weight        | : | 70 kg              |
| Safe working load | : | 180 kg             |
| Castors           | : | 200mm Dim          |





**Model No. KSS - 136**

**ELECTRIC GYNE EXAMINATION TABLE**

**Specification:**

- Overall Size L 1830mm x W 870mm x Ht. 520 - 1010mm
- Two section Upholster cushion top
- Head low/up position manually with the help of ratchet
- Sliding foot end
- Adjustable Height position through hand remote as per convenience
- 'U' notch with drainage tray
- 100mm goods quality castors, two with brakes
- Lithotomy padded cushion rods for leg support
- Epoxy powder coated base

**Positions :**





**Model No. KSS - 137**

**DELIVERY BED HYDRAULIC**


**Specification:**

- Overall Size : 1880Lx900W x 600-800H mm.
- Tubular Frame work mounted on 125 mm dia castor two with brakes.
- Three Section top made of MS pipe with fixed Cushioned mattress.
- Trendlenburg / Reverse Trendlenburg Positions, Back section adjustment through manually.
- Height adjustment through Hydraulic Pump.
- Leg end section can slide under the main section.
- Provided with a pair of upholstery lithotomy position.
- SS telescopic IV Rod
- Detachable SS Tray.

**Optional:**

5 years guarantee replacement of Hydraulic pump

**Positions :**





**Model No. KSS - 138**

**ELECTRIC DELIVERY BED**

**Specification:**

- Overall Size : 1880L x 900W x 600-800H mm.
- Tubular Frame work mounted on 125 mm dia castor two with brakes.
- Three Section top made of MS pipe with fixed Cushioned mattress.
- Height Adjustment, Back section & Trendlenburg / Reverse Trendlenburg Positions adjustment through electric actuator system.
- Leg end section can slide under the main section.
- Provided with a pair of upholstery lithotomic position.
- SS telescopic IV Rod.
- Detachable SS Tray.

**Model No. KSS - 139**

**LABOUR TABLE HYDRAULIC PUMP OPERATED**

**Specification:**

- Approx. Dimensions: 72"L x 27"W x 24"-32"H.
- Tubular frame work mounted on swiveling castors.
- Stainless steel top in three sections.
- Trendlenburg positions adjustable by pneumatic gas spring system.
- Hydraulically height adjustable by foot paddle from side.
- Backrest section adjustable on ratchet.
- Leg end section can slide under the main section
- A pair of knee crutches with height adjustable.
- Finish : Pre treated & Epoxy powder coated.



**Model No. KSS - 140**

**OBSTETRIC LABOUR TABLE TELESCOPIC (FIXED HEIGHT)**

**Specification:**

- Overall Approx. max. Size : 72" x 27" x 30".
- Frame work made of MS Tube.
- Legs fitted with PVC Stumps.
- 3 Section Top made of S.S. Sheet.
- Back rest section adjustable on ratchet.
- Leg end section can slide under the main section.
- The entire table top gives trendlenberg positions by gas entire system.
- Removable safety side railing are provided on 3 sides of the table.
- Height adjustable a pair of foam padded knee crutches.
- Patient bearing down hand grips made of S.S. Tubes.
- Supplied with waste collecting S.S. Bowl on S.S. bracket and one S.S. I.V. Red.
- Finish : Pre treated and Epoxy Powder Coated.
- Also available in S.S. Frame work and S.S. Top.



**Model No. KSS - 141**

**OBSTETRIC DELIVERY TABLE (3 SECTION)**

**Specification:**

- Approx. Dimensions: 72"L x 24"W x 30"H.
- Rectangular & square tube frame work mounted on PVC stumps.
- Three sections cushioned top.
- Backrest, seat, leg section adjustable by hand lever.
- Complete with SS bowls.
- Height adjustable a pair of knee crutches.
- Pre treated & Epoxy Powder Coated.



**Model No. KSS - 142**

**DELIVERY BED (TELESCOPIC)**

**Specification:**

- Bed Labour and Delivery. Made in 2 sections, body section with legs mounted on rubber tips, leg section on swivel castors with brakes.
- 72"L x 27"W x 30" H

**Model No. KSS - 143**

**OBSTETRIC DELIVERY BED IN 2 PARTS (2 SECTION TOP)**

**Specification:**

- Approx. Dimensions : 72"L x 27"W x 30" H.
- Tubular frame work made in two section.
- Body section with legs mounted on PVC stumps.
- Leg section on swivel castors with brakes.
- Two section cushioned top.
- Height adjustable a pair of knee crutches.
- Pre treated & Epoxy powder coated.





**Model No. KSS - 144**

**GYNAECOLOGICAL EXAMINATION TABLE 3 SECTION**

**Specification:**

- Table, Gynecological, 3 sections, adjustable backrest and leg section with mechanism to varying positions, capable of tilting for obstetric and gynecological use.
- 72" L x 24" W x 30" H

**Model No. KSS - 145**

**GYNAE EXAMINATION TABLE INTO**

**Specification:**

- Overall Approx. Size : 1830 mm L x 600 mm W x 760 mm H.
- Frame work made of MS CRC tubes.
- Two section cushioned top.
- Complete with the pair of S.S. lithotomy rod made of 12 mm S.S. Rod.
- Back rest adjustable by hand lever.
- Epoxy powder coated finish.



**Model No. KSS - 146**

**EXAMINATION TABLE 2 SECTION**

**Specification:**

- 183 L x 60 W x 81 H cms.
- Frame work of CRC tubes.
- Legs fitted with PVC stumps.
- 5 cms foam cushioned top covered with rexine.
- Pre-treated and Epoxy Powder Coated.

**Model No. KSS - 147**

**EXAMINATION COUCH WITH CABINETS AND DRAWERS**

**Specification:**

- Upholstered top upper section with three drawers  
Adjustable upholstered back rest.
- Lower Section with lockable slide door cabinet  
Swinging tray for B.P. Appartus.
- Chemically pre treated and Epoxy powder coated /  
oven baked enameled backrest.
- Size 1890 x 510 x 820 (L x W x H mm).
- Platform size - 450 x 510 (L x W mm).



**Model No. KSS - 148**

**EXAMINATION COUCH**

**Specification:**

- Overall Approx. Size : 180L x 50 W x 80 H cms.
- Frame work of CRC tubes.
- 5 cms foam cushioned top covered with rexine.
- 1Pc. cabinet with drawer.
- Legs fitted with PVC stumps.
- Finish : Pre-treated and epoxy powder coated.

**Optional**

- Single foot step in powder coted / S.S.

**Model No. KSS - 149**

**EXMINATION COUCH WITH GYANE CUT**

**Specification:**

- Overall Approx. Size : 183L x 50 W x 82 H cms.
- Top frame work of CRC tubes covered with 5 cms foam  
cushioned and rexine.
- 3Pc. cabinet with 1 drawer.
- Legs fitted with PVC stumps.
- Finish : Pre-treated and epoxy powder coated.

**Optional**

- Single Foot Step in Powder coated.
- Utility trolley in powder coated with S.S. tray.



↑  
**UTILITY TROLLEY**



**Model No. KSS - 150**

**EMERGENCY & RECOVERY TROLLEY  
(HYDRAULIC)**

**Specification:**

- Approx. size: 210L x 65W x 65-95H cms.
- Safety side railing , I.V. Rod, storage tray & Cylinder cage.
- Removable X-Ray translucent top with tray for holding X-ray Cassettes.
- Adjustable Height by hydraulic pump through foot pedal.
- Trendlenburg / Rev. Trendlenburg by gas spring. 15 cms dia castors.
- Corner buffers.
- Pre-treated & Epoxy Powder Coated.

**Model No. KSS - 151**

**PATIENT SHOWER TROLLEY HYDRAULIC**

**Specification:**



- Hydraulically height adjustable allowing a smooth transition for patients and less strenuous activity for carers.
- Able to achieve both Trendelenburg and anti-Trendelenburg positions at the with help of pneumatic pumps.
- Comes with an tilting position from one side with the help of pneumatic pump.
- Safety rails with an easy mechanism for lowering and raising.
- Flexible drain hose to allow water to run straight into drainage system.
- Vinyl mattress and pillow in grey
- Braked castors.
- Dimensions available: Length 2100 mm, Width 800mm, Height 600mm to 800mm.
- Safe working limit: 200kg.
- Pre-treated & epoxy coated



**Model No. KSS - 152**  
**AMBULANCE STRETCHER**

**Specification:**

- Light weight , strong structure, lovely design, good arts and crafts, easy antisepsis, small size convenient to use.
- Utilizes high precision mechanical transmission.
- Can be used automatically folded through the left and right hand control handle.
- Made of high quality aluminium alloy tubes.
- Dimensions:-  
A-High position : 1,900 x 560 x 900 mm  
B-low position : 1,900 x 560 x 280 mm
- Product weight : 40 kg
- Diameter of the Castors : 125 mm
- Applicable for various kind of hospital, emergency centers, ambulances and battle fields.



**Model No. KSS - 153**  
**PATIENT STRETCHER TROLLEY**

**Specification:**

- Approx. size : 195L x 65W x 68-90H cms.
- Frame Work made of MS Tubes.
- Safety side railing I.V. Rod.
- Backrest Adjustable.
- 15 cms dia castors.
- Pre-treated & Epoxy Powder Coated.
- Also available in complete S.S.



**Model No. KSS - 154**  
**STRETCHER TROLLEY**

**Specification:**

- Overall Approx. Size : 213 L x 56 W x 81 H cms.
- Frame work made of 1.25" 16 G vertical and 1" x 18 G horizontal CRC tube.
- Trolley mounted on 15 cm dia castors 2 with brake.
- Removable stretcher top made of MS sheet.
- Pushing handle at both end of the stretcher top.
- Finish : Pre treated and Epoxy Powder Coated.
- Also available in complete S.S.



**Model No. KSS - 155**  
**FOLDING STRETCHER 2 FOLD**

**Specification:**

- Unfolded Size: 221 L x 53 W x 15 H cms.
- Folded Size : 110 L x 18 W x 11 H cms.
- Net Wt : 5 kg.
- Gross Wt : 6 kg.
- Load Bearing Capacity : 159 kg.
- Also Available in Aluminum

**Model No. KSS - 156**  
**MODEL SCOOP STRETCHER (TELESCOPIC/IMMOBILE)**


**Specification:**

A portable first aid stretcher, is commonly used for transporting bone fractured and seriously Injured patients. A clutch zip design at the middle of the model, the stretcher can be divided into left and right parts. Both parts can be drawn out when you need to put the patient on the operating table or the hospital bed quickly without moving the patient. The stretcher length can be directly transferred into the ambulance. The stretcher can fold or be separated into two parts, easy to carry or transport.



Product size (LxWxH) : 168 x 43 x 7cm (min)  
224 x 45 x 10cm (max)  
Net Weight : 9 kg

Gross Weight : 11 kg  
The bearing pressure : < 159 kg  
Dimension of the packing : 127 x 45 x 10 cm



**Model No. KSS - 157**  
**SPINE BOARD**

**Specification:**

Spine Board is mainly used for immobilizing and safely transferring serious patients during rescuing. It's a first-aid transport equipments for medical department. It also could be equipped on ambulances. Besides, it can go into X-ray directly and float on water. This model is made of high strength engineering plastics. It can be enduring used, and has the advantage of antiaging, weight loading and easy stocking specialties.

Size : 183 L x 45 W x 5 H cms  
Load Bearing : 159 kg  
Net Weight : 6.5 kg  
Gross Weight : 7.5 kg

**Model No. KSS - 158**  
**WHEEL CHAIR FOLDING**

**Specification:**

- 16" Seat with Fixed Armrest Front with 4" Solid Castors & Rear 26" Solid Wheels with brakes.
- Size 39.5" L x 28" W x 34.5 H.



**Model No. KSS - 159**  
**WHEEL CHAIR NON FOLDING**

**Specification:**

- 17" cushioned seat & back fixed armrests front with 5" solid castors & rear 26" solid wheels with brakes.
- Size - 41" x 28" x 33" (L x W x H).





**Model No. KSS - 160**  
**MONITOR TROLLEY**

**Specification:**

- Overall Approx. Size : 800 mm.Lx 600 mm W x 1400 mm H.
- Four shelves and one drawers.
- Trolley having two handles.
- The trolley mounted on 125 mm dia castors -2 with brakes.
- Switch and Socket for medical equipment.
- Supplied with S.S. IV Rod and MS Oxygen Gas Cylinder Case.
- Finish : Pre-treated and Epoxy Powder Coated.

**Model No. KSS - 161**  
**CRASH CART**

**Specification:**

- Overall Approx. Size : 960 Lx 500 W x 1545 H mm.
- Frame work made of S.S. Tubes.
- Two light weight polystyrene boxes each with three drawers.
- Upper drawer with medicine container of different sizes.
- Provision to hold oxygen cylinder made of M.S.
- Massage board for cardiac patient are provided.
- Six Number colour hand out small bins at eye level.
- Total three S.S. Shelves with railing to carry Medical Equipments.
- Provided with S.S. I.V. Rod, corner Buffers.
- The cart mounted on 125 mm dia casters 2 with brake.
- Also available in M.S. Frame work and laminated shelves.



**Model No. KSS - 162**  
**OT DRUG TROLLEY**

**Specification:**

- Approx. Overall Size :  
67 cm LX 46 cm W X 153 cm H.
- Complete in S.S. Structure with 24 nos. small hand bins & 12 nos. big hand bins.
- One tray on top.
- 5" dia. Castors two with break.

**Optional**

- Also available in powder coated structure.



**Model No. KSS - 163**

**ELECTRIC BLOOD DONOR COUCH**

**Specification:**

- Motorized Donor Chairs suitable to use for blood banks, infusion treatment centres, and hospital, clinics, and aphaeresis, plasma centers and other similar discipline.
- Remote switch for easy repositioning
- Hard wood frame construction
- It has solid surface seat to ensure easy & hygienic premium rated foam, cushion and spring component for donor comfort.
- With the help of electric power, we can fully adjust the donor position in between head low-feet high to head high-feet lower position and arm rest (left and right) swing out for easy seating as well as up & down of chair such position minimization of donor.
- High quality contract grade vinyl upholstery
- All exposed metal parts are chrome plate for easy cleaning & durability.

**Technical Details**

|                                  |   |                  |
|----------------------------------|---|------------------|
| Seat Height                      | : | 64 cm            |
| Head Height Max                  | : | 114 cm           |
| Foot Height Max                  | : | 89 cm            |
| Chair Length                     | : | 152 cm           |
| Chair Width (W/Armrest Sections) | : | 92 cm            |
| Load Bearing Capacity            | : | 150Kg            |
| Net Weight                       | : | 52 kg            |
| Gross Weight                     | : | 92 kg            |
| Packing Dimensions               | : | 168 x 76 x 92 cm |



**Model No. KSS - 164**  
**BLOOD DONOR CHAIR, 3 SECTION**

**Specification:**

- Overall Size : 150L x 50W x 46H cms
- Framework made of Stainless Steel
- Back rest & leg rest section adjustable on Ratchet / Lever
- Legs fitted with PVC Stumps
- 5 cms foam cushioned top covered with rexine
- Head section & leg section move simultaneously.
- Adjustable Arm Rest

**Model No. KSS - 165**  
**BLOOD DONOR CHAIR**

**Specification:**

- Blood sampling chair with fixed height, Height= 670 mm
- Comfortable upholstery with good density, width 550
- Steel quality, epoxy paint
- Adjusting back rest via 1 gas jack
- Adjusting leg rest via a self locking mechanism system
- Equipped with 2 blood test splints, every position adjustable by serrated blocks system
- Load capacity: Dynamic 130kg



**Model No. KSS - 166**  
**BLOOD TRANSFUSION CHAIR**

**Specification:**

- Blood sampling chair with fixed height, Height= 460 mm
- Comfortable upholstery with good density, width 420
- Steel quality, epoxy paint
- Adjusting back rest via 1 gas jack
- Equipped with 1 blood test splint, every position adjustable by serrated blocks system
- Load capacity: Dynamic 130kg

**Model No. KSS - 167**  
**BLOOD COLLECTION CHAIR**

**Specification:**

- The frame is made of flat oval steel tubing
- The padded seat, back and armrests provide the patient with comfortable seating.
- The seat height is 500 mm.
- The distance between the wedge-shaped padded armrests 500 mm
- Maximum weight load 150 kg





**Model No. KSS - 168**  
**BED SIDE LOCKERS (DELUXE)**

**Specification:**

- Hospital bedside lockers made of high quality sheet metal.
- Die press Stainless steel top.
- Hospital bedside lockers Provided with two storage cabinet and single draw under the top and space for keeping utilities.
- Hospital bedside lockers optionally available without draw under the top
- Standard dimension :  
500mm L x 450mm W x 820 mm H

**Model No. KSS - 169**  
**BED SIDE LOCKER (SUPER DELUXE)**

**Specification:**

- 16" L x 16" W x 32" H
- Machine Pressed CRCA Sheet enclosed on three sides.
- Having one box and one drawer fitted with membrane /laminated pressed top and front side.
- Fitted with 50 mm dia castors.
- Epoxy powder coated.



**Model No. KSS - 170**  
**BED SIDE LOCKER (DELUXE)**

**Specification:**

- CRCA sheet constructions.
- Having one box and one drawer
- S.S. Top with three sides raise edges.
- Epoxy Powder Coated.
- 16" L x 16" W x 32" H.

**Model No. KSS - 171**  
**BED SIDE LOCKER (STD)**

**Specification:**

- Four MS legs fitted with PVC Stumps.
- Box and door made of MS sheet.
- Top made by machine pressed SS Sheet.
- Finish : Epoxy Powder coated.
- 16" L x 16" W x 32" H





**Model No. KSS - 172**  
**REVOLVING STOOL (SS TOP)**

**Specification:**

- Four leg base with PVC Stumps.
- Height adjustment by machine Screw from 18" x 27".
- A foot rest ring support made of MS Tube.
- SS top supported with MS sheet top.
- Finish : Pre-treated and epoxy Powder Coated.
- Also Available in S.S. Frame work and S.S. Top.

**Model No. KSS - 173**  
**LAB TECHNICIAN STOOL**

**Specification:**

- Height adjustable cushioned Seat.
- 5 legs base on castors.
- Also Available with backrest.



**Model No. KSS - 174**  
**REVOLVING STOOL (CUSHIONED TOP)**

**Specification:**

- Four Leg base fitted with PVC Stumps.
- Height adjustment by machine Screw from 18" x 27".
- A foot rest ring support made of MS Tube.
- Cushioned top support with MS Circle.
- Finished : Pre-treated and Epoxy Powder Coated.
- Also available in S.S. Frame work Cushioned top.

**Model No. KSS - 175**  
**OVER BED TABLE PNEUMATIC**

**Specification:**

- Powder coating
- Wooden printing melamine laminate table top (non-tilt).
- 40mmx80mm outer support steel tube.
- 40mm square steel tubing legs.
- 80cmLx38cmW table top.
- Fingertip pressure height adj.: 76cm-106cm.
- Plastic caster x 4 pcs.



**Model No. KSS - 176**  
**OVER BED TABLE WITH ABS TOP**

**Specification:**

- Sturdy yet light weight
- Easy to clean surfaces
- This over bed table provides beauty without having to worry about laminate wood damaging over time with moisture.
- ABS plastic top.
- Easy roll swivel caster facilitate transport

**Model No. KSS - 177**  
**OVER BED TABLE GEAR OPERATED**

**Specification:**

- Powder coating
- Wooden printing melamine laminate table top (non-tilt).
- 40mmx80mm outer support steel tube
- 40mm square steel tubing legs
- 80cmLx38cmW table top
- Gear handle pressure height adj.: 76cm-106cm
- Plastic caster x 4 pcs



**Model No. KSS - 178**  
**OVER BED TABLE MANUAL**

**Specification:**

- Rectangular Frame Work Mounted On 5cms Castors
- 30"×16" Laminated Top Adjustable By Side Knob
- Pre-treated and epoxy powder coated.



**Model No. KSS - 179**  
**OVER BED TABLE FIXED HEIGHT**

**Specification:**

- Overall approx size: 1050 (L) mm x 400 (W) mm x 975 (H) mm
- Rectangular frame work mounted on 50 mm castors.
- 750 mm x 400 mm laminated top.
- Pre treated and Epoxy powder coated.

**Model No. KSS - 180**  
**SALINE STAND**

**Specification:**

- Five Legs base mounted on 5cms dia castors.
- With one approx. 35cms dia S.S. basins.
- Vertical MS tube made of 25mm X 15G.
- Basin holder MS tube made of 16mm X 18G.
- Pre treated and Epoxy Powder Coated.
- Also Available in SS vertical tube and SS basin holder.



**Model No. KSS - 181**  
**TROLLEY FOR OXYGEN CYLINDER**

**Specification:**

- Frame work made of CRC tube mounted on 10 cms/ 15cms castors
- Epoxy powder coated finish.
- A- for B type cylinder.
- B- for D type cylinder.

**Model No. KSS - 182**  
**KICK BOWL / BUCKET**

**Specification:**

- Approx. S.S. Bowl of 350 mm dia/SS Bucket.
- Frame work made of MS tubes.
- Mounted on 75 mm castors.
- Finish : Epoxy Powder Coated.
- Also available in SS frame work and SS Basin/Bucket.





**Model No. KSS - 183**  
**WASH BASIN STAND (DOUBLE)**

**Specification:**

- Five Legs base mounted on 5cms dia castors.
- With two approx. 35cms dia S.S. basins.
- Vertical MS tube made of 25mm X 15G.
- Basin holder MS tube made of 16mm X 18G.
- Pre treated and Epoxy Powder Coated.
- Also Available in SS vertical tube and SS basin holder.



**Model No. KSS - 184**  
**WASH BASIN STAND (SINGLE)**

**Specification:**

- Five Legs base mounted on 5 cms dia castore.
- Pre-teated and Epoxy Powder Coated.
- Also available in S.S. pipe & rings



**Model No. KSS - 185**  
**FOOT STEP DOUBLE**

**Specification:**

- 18"x18"x18".
- M.S. Tubular Frame fitted with PVC stump.
- CRC sheet duly double press bent.
- Pre-treated and Epoxy Powder Coated.
- Also available in Stainless Steel Model.

**Model No. KSS - 186**  
**FOOT STEP SINGLE**

**Specification:**

- 18"x9"x9".
- M.S. Tubular legs fitted with PVC stump.
- Pre-treated and Epoxy Powder Coated.
- M.S. top fixed with rubber matting and aluminium beedings.
- Also available in Stainless Steel Model.



**Model No. KSS - 187**  
**LINEN CHANGE TROLLEY**

**Specification:**

- MS Tubular frame mounted on 10cms Dia castors.
- Three shelves made of MS Sheet.
- A Canvas Bag fitted on front side.
- Epoxy Powder Coated Finish.



**Model No. KSS - 188**

**BABY CRADLE**

**Specification:**

- M.S. tubular pipe. Stand fitted with PVC stumps.
- Wire mesh hanging cradle.
- Pre-treated and Epoxy Powder coated.
- Size - 36"x15"x 39" (L x W x H)

**Model No. KSS - 189**

**BABY BASSINET**

**Specification:**

- Preprex transparent Crib with mattress.
- Trendlenburg / Rev. Trendlenburg 5 cms castors. I.V. Rod.
- Pretreated & Epoxy Powder Coated.
- Also Available S.S.



**Model No. KSS - 190**

**BED SIDE SCREEN (4 PANELS)**

**Specification:**

- Overall Approx. Size : 1680 H x 2450 mm.
- MS tubular construction made of 19 mm x 18G in four section.
- Legs made of 25 mm x 18G fitted with 50 mm dia castors.
- Each span 610 mm wide.
- Supplied with hooks and spring with curtain cloth.
- Finish : pre treated and Epoxy Powder Coated.

**Model No. KSS - 191**

**BED SIDE SCREEN (3 PANELS)**

**Specification:**

- Screen panel Size: (160cm H x 60cm L)
- Curtain Size: (140cm H X 60cm L)
- M.S steel Tubular frame work mounted on 50mm castors with curtains.
- Curtains available in different colors. (white, green, blue)
- Easy assemble & disassemble or storage.
- Pre-treated & epoxy coated





**Model No. KSS - 192**

**Specification:**

- Frame work of MS tubes mounted on 10 cms castors.
- Two SS shelves with three side railing on top shelf side railing on top shelf.
- Two drawer under the each shelf.
- Epoxy Powder Coated Finish.
- Overall Approx. Size: 76L x 46W x 81 H cms.



**Model No. KSS - 193**

**UTILITY TROLLEY**

**Specification:**

- M.S. Frame work.
- S.S. Shelves.
- One M.S. Drawer.
- 10 cms Castors.



**Model No. KSS - 194**  
**INSTRUMENT TROLLEY**  
**(THREE SHELVES)**

**Specification:**

- Frame work made of 1" x 18G vertical MS Tubes.
- Three 22G S.S. Shelves with three sides railing.
- The trolley mounted on 100 mm dia castros.
- Finish : Pre treated and Epoxy Powder Coated.
- A- 24" x 18" x 32" B-30" x 18" x 32" C- 36" x 24" x 32"
- Also Available in S.S. Frame & S.S. Sheet



**Model No. KSS - 195**

**DRESSING TROLLEY**

**Specification:**

- Approx . Size : 30"x18"x32".
- Frame work made of MS tubes.
- Two S.S. Shelves over with S.S. railings on all four sides.
- The trolley fitted S.S. bucket and S.S. bowl.
- Finish : Pre treated and Epoxy Powder Coated.
- Also available in SS frame work and SS shelves.



**Model No. KSS - 196**

**INSTRUMENT TROLLEY**

**Specification:**

- Tubular pipe framework,
- EPC finished Two shelves-Top & underneath of stainless steel sheet.



**Model No. KSS - 197**

**MAYO'S INSTRUMENT TROLLEY WITH S.S. TRAY**

**Specification:**

- Rectangular CRC tube mounted on 5 cms castors, height adjustable by side Knob.
- Pre - treated and Epoxy Powder Coated finish.
- Also available in S.S. frame & tray.
- Tray size S.S. : 56Lx 40 W cms

**Model No. KSS - 198**

**MAYO'S TABLE GEAR OPERATED**

**Specification:**

- Overall Approx. Tray Size : 560 mm Lx 400 mm W.
- Adjustable height of tray Approx. 760 mm to 1270 mm through gear handle.
- Frame work made of 50 mm x 30 mm x 1.22 mm rectangular MS tubes.
- Tray made of 0.9 mm thick S.S. 304 Grade sheet.
- Table mounted on 50 mm dia castors.
- Finish : Pre treated and Epoxy Powder Coated.
- Also available in SS frame work and SS Tray.



**Model No. KSS - 199**

**INSTRUMENT MAYO'S TROLLEY OVER THE O.T TABLE (MECHANICALLY)**

**Specification:**

- Approx Tray S.S. Size : 82 Lx 50 W cms
- Frame Work Of CRC Tubes Mounted on 10cms Castors.
- Height adjustable By Gear Handle.
- Epoxy Powder Coated.
- Also Available in S.S. Frame & tray.

**Model No. KSS - 200**

**MEDICINE TROLLEY 6 DRAWERS**

**Specification:**

- Overall Approx. Size : 750mm L x 380mm W x 930mm H.
- Six Plastic drawers, two drawers with utility container.
- Fixed one S.S. Tray on top & down side.
- M.S. Tubular frame work fitted with 4" dia. castors.
- Pretreated & Powder Coated frame
- Note : Also available in S.S. Framework & S.S. shelves.



**Model No. KSS - 201**

**MEDICINE TROLLEY 12 DRAWERS**

**Specification:**

- Overall Approx. Size : 750mm L X 380mm W X 930mm H.
- Twelve plastic drawers, two drawer with Utility container.
- One S.S. Tray on top.
- M.S. Tubular framework fitted with 4" dia. castors.
- Pretreated & Powder coated Frame.
- Note : Also available in S.S. Framework & S.S. Top.



**Model No. KSS - 202**

**UNIVERSAL HYDRAULIC OPERATING TABLE**



An Universal hydraulic operating table is specially designed to provide complete flexibility for patient positioning required for modern surgeries. This table is constructed of stainless steel and other high quality materials.

**FEATURES**

- Five sectional stainless steel top which gives ease various surgical positions.
- UP & DOWN by Hydraulic System .
- Standard Positioning :- Height UP/DOWN, Trendelenberg / Reverse Trendelenberg, Lateral tilt sides, Flex / Reflex, Chair, Kidney elevation through drive handle.
- Base and column covered with S.S. Top to ensure high durability and Hygiene for the operation room.
- Mobile on Smooth running castors.
- Floor Lock facility.


**STANDARD ACCESSORIES**

- Five section mattress - 1 set
- Anesthetist Screen - 1 Pcs.
- Padded Shoulder support -2 Pcs.
- Padded Lateral Support-2 Pcs.
- Padded Lithotomy Crutches-2 Pcs.
- Padded Arm Rest - 2 Pcs.
- Drive Handle - 1 Pcs.

| Specificaton/Model No.      | ASC- 704<br>(Approx. Dimeisions) | ASC-704 Electric<br>(Approx. Dimeisions) |
|-----------------------------|----------------------------------|--|
| Overall Length of the table | 1950 mm                          | 1950 mm                                  |
| Overall width of the table  | 500 mm                           | 500 mm                                   |
| Minimum height of the table | 780 mm                           | 780 mm                                   |
| Maximum height of the table | 1030 mm                          | 1030 mm                                  |
| Trendlengburg/Reverse       | 30 Degree                        | 30 Degree                                |
| Flex / Reflex               | 80°/220°                         | 80°/220°                                 |
| Side Tilt                   | 22°                              | 22°                                      |
| Kidney Position             | 125 mm manually                  | 125 mm manually                          |
| Leg Section                 | 90° down                         | 90° down                                 |
| Hi-Low Lift                 | Hydraulic                        | Electric                                 |



**Model No. KSS - 203**  
**OPERATING TABLES C-ARM COMPATIBLE (HYDRAULIC)**



**Specification:**

- Four Sections Radiolucent top.
- C-Arm compatible inbuilt kidney bridge.
- Eccentrically positioned top for maximum access for movement of C-Arm image intensifier.
- Heavy & sturdy T shaped base allowing maximum foot room for operating team.
- Mobile on smooth running castors & Floor locking system for stability
- Heavy & Hygienic stainless steel Covered base and lifting column
- All the fittings, side rails and accessories made of stainless steel
- Accessories clamps made of Stainless Steel

**Standard Positioning**

- Hydraulic UP / DOWN
- Precise and efficient crank mechanism to obtain trendelenberg / reverse trendelenberg
- Lateral tilt both sides
- Flex & Reflex, & Chair
- Detachable, adjustable and interchangeable head & leg sections
- Also Available with Kidney Bridge Position with Radiolucent top.

**Trendelenberg**





C-Arm compatible operating table for Orthopaedic with Leg Traction unit attached to the base



Inner thigh support for lateral position



Arm Support Carter Braine with pad (Adjustable) for lateral positioning (Optional)



Knee Rest for tibia and fibula position / arthroscopic position



Lateral Cassette Holder (Optional)



Pelvic rest (adjustable) at extreme right



Steinmann Pin Clamp



Pelvic rest (adjustable) positioned at extremity left



C-arm Compatible Operating Table (Hydraulic)

X-ray translucent,  
C-Arm Compatible.  
Hand operating table  
with mattress  
and telescopic support.

#### Standard Accessories

- Five section mattress - 1 set
- Anesthetist Screen - 1 Pcs.
- Padded Shoulder support - 2 Pcs.
- Padded Lateral Support - 2 Pcs.
- Padded Lithotomy Crutches - 2 Pcs.
- Padded Arm Rest - 2 Pcs.
- Drive Handle - 1 Pcs.



| Specification/Model No.     | ST - 203<br>(Approx. Dimensions) |
|-----------------------------|----------------------------------|
| Overall Length of the table | 1950 mm                          |
| Overall width of the table  | 500 mm                           |
| Minimum height of the table | 780 mm                           |
| Maximum height of the table | 1030 mm                          |
| Trendlengburg/Reverse       | 30 Degree                        |
| Flex / Reflex               | 80°/220°                         |
| Side Tilt                   | 22°                              |
| Kidney Position             | 125 mm manually                  |
| Leg Section                 | 90° down                         |
| Chair Position              | 85°                              |

**Model No. KSS - 204**

**ELECTRIC OPERATION TABLE C. ARM COMPITABLE**



**Specification:**

**With Orthopaedic Attachment**

Electrically driven, finger touch controlled operation tables are specially designed to provide complete flexibility for patient positioning required for modern surgeries.

**SALIENT FEATURES**

- Mains operated electrical system
- In-built rechargeable battery (optional)
- Eccentrically Positioned table top for maximum access for movement of C Arm image intensifier.
- X-ray Translucent four/five sectional top.
- Standard positioning
- Noiseless and smooth operation.
- Mobile base with floor locking facility.
- Sturdy T-type base for easy operations.
- Base and column covered with stainless steel prevent dirt and fluids from getting into base and easy to be cleaned.
- Frame and side rails are made of S.S
- Completely RUST proof design.
- Adaptable to ortho attachment, Neuro Surgery attachment and other special accessories.

**Model:**

Four/five Sectional X-ray translucent top. Finger touch hand set controlled position:-

- (i) Height Up/Down
- (ii) Trendelenburg / Reverse Trendelenburg.
- (iii) Lateral tilt both sides.
- (iv) Flex/Reflex
- (v) 90 degree chair.
- (vi) Available with Kidney Bridge

Position adjustable mechanically. Detachable, adjustable and interchangeable head & leg sections manually.

**Standard Accessories**

- Five section mattress - 1 set
- Anesthetist Screen - 1 Pcs.
- Padded Shoulder support -2 Pcs.
- Padded Lateral Support-2 Pcs.
- Padded Lithotomy Crutches-2 Pcs.
- Padded Arm Rest - 2 Pcs.
- Drive Handle for kidney bridge only- 1 Pcs.

**Optional Accessories**

- Orthopaedic Attachment
  - > Orthopaedic Leg Traction Device.
  - > L-Shaped Knee rest padded.
  - > Steinmen Pin Holder.
  - > Inner thigh support padded.
  - > Raised arm rest padded.
  - > Lateral Cassette holder.
- Hand Surgery Attachment
- Neuro Surgery Attachment
  - > Horse Shoe shaped
  - > Mayfield Three Pin Head Fixator
- Manually operated Kidney Bridge.
- Hand Surgery Table (X-ray translucent) with telescopic support
- Electric Floor Locking System

**Optional Orthopaedics Accessories**



Hand Traction device  
(Complete with X-ray translucent finger plate )



Spinal Frame



Counter Traction post for  
Tibia and Fibula

**Optional Urology / General Surgery Accessories**



Lithotomy position with  
perineal 'U' cut  
and knee krutches.



inbuilt X-ray translucent  
C-Arm Compatible  
Kidney bridge

**Optional Neuro-Surgery Accessories**



Mayfield three Pin Head Fixator for  
Neuro Surgery in SUPine Position  
(wholly made of S.S.)



Horse shoe shaped Neuro Surgery Head  
Rest full adjustable for prone & supine position  
(wholly made of stainless steel)



Mayfield three Pin Head  
Fixator for Neuro Surgery  
in Sitting Position



**Optional Accessories**

**Urology / GE Rnral Surgery**

- Inbuilt kidney bridge X-ray translucent & C-arm Compatible. PCNL attachment (Extra long leg section X-ray translucent with telescopic support.) TUR Drainage Tray.



**Model No. KSS - 205**  
**CEILING OPERATION LIGHT DOUBLE SOME**  
**(4+4) REFLECTORS**

**Specifications**

|                          |   |                                |
|--------------------------|---|--------------------------------|
| Dome                     | : | 4 + 4                          |
| Dome Diameter            | : | 20 " ~ 510 mm                  |
| Reflector                | : | Aluminum / Dichroic            |
| Intensity                | : | 50000 Lux (A) / 100000 Lux (D) |
| Color Temp.              | : | 4300 K (A)                     |
| CRI                      | : | 93(Dichroic)                   |
| Variable Light Intensity | : | 3 stages                       |
| Bulb specification       | : | 24v - 70w x 4                  |
| Power                    | : | AC 230V                        |
| Battery Backup           | : | Optional                       |

**Model No. KSS - 206**  
**CEILING OPERATION LIGHT DOUBLE SOME (7+4) REFLECTORS**

**Specifications**

|                          |   |   |
|--------------------------|---|---|
| Dome                     | : | 7 + 4   |
| Dome Diameter            | : | Aluminum / Dichroic   |
| Reflector                | : | 80000 Lux (A) + 50000 Lux (A)                                 |
| Intensity                | : | 28 " + 20 " ~ 700 mm + 510 mm<br>150000 Lux(D)+100000 Lux (D) |
| Color Temp.              | : | 4300 K (A)  |
| CRI                      | : | 93(Dichroic)  |
| Variable Light Intensity | : | 3 stages  |
| Bulb specification       | : | 24v - 70w x 7 + 70w x 4                                       |
| Power                    | : | AC 230V   |
| Battery Backup           | : | Optional  |



**Model No. KSS - 207**  
**CEILING OPERATION LIGHT DOUBLE DOME**  
**Specification:**

- The Light Head provides the superior performance and is good for all kinds of surgeries. The compact light head with shadow reduction control makes this light suitable for wide variety
- The Shadow less light equipped with the reserve light source, which is activated in case the main light source brakes down.
- Sterilizable handle has been provided on the dome for easy movement.



**Specifications**

|                         |   |                     |
|-------------------------|---|---------------------|
| Dome Head               | : | 521 mm x 2          |
| Halogen Bulbs           | : | 24V x 150 W x 2 nos |
| Light intensity at 1 mt | : | 1,00,000 Lux        |
| Height Adjustment       | : | 1014 mm             |

|                           |   |             |
|---------------------------|---|-------------|
| Color Temperature         | : | 4200 K ±300 |
| Temp rise in field (1 mt) | : | 2° -3° C    |
| Control Unit              | : | Transformer |



**Model No. KSS - 208**  
**CEILING OPERATION LIGHT SINGLE DOME**

**Specifications**

|                          |   |                             |
|--------------------------|---|-----------------------------|
| Dome                     | : | 7                           |
| Dome Diameter            | : | 28 " ~700mm                 |
| Reflector                | : | Aluminum / Dichroic         |
| Intensity                | : | 80000 Lux (A) 150000 Lux(D) |
| Color Temp.              | : | 4300 K (A)                  |
| CRI                      | : | 93(Dichroic)                |
| Variable Light Intensity | : | 3 stages                    |
| Bulb specification       | : | 24v – 70w x 7               |
| Power                    | : | AC 230V                     |
| Battery Backup           | : | Optional                    |



**Model No. KSS - 209**  
**MOBILE OPERATION LIGHT**

**Specifications**

|                    |   |                                 |
|--------------------|---|---------------------------------|
| Dome               | : | 1                               |
| Dome Diameter      | : | 19 " ~ 480mm                    |
| Reflector          | : | Aluminum / Dichroic             |
| Intensity          | : | 50000 Lux (A)<br>100000 Lux (D) |
| Color Temp.        | : | 4300 K (A)                      |
| CRI                | : | 93(Dichroic)                    |
| Bulb specification | : | 24v – 150w                      |
| Power              | : | AC 230V                         |
| UPS System         | : | Optional                        |



**Model No. KSS - 210**  
**MOBILE OPERATION LIGHT 4 REFLECTOR**

**Specifications**

|                    |   |                                |
|--------------------|---|--------------------------------|
| Dome               | : | 4                              |
| Dome Diameter      | : | 20 " ~ 510mm                   |
| Reflector          | : | Aluminum /<br>Dichroic         |
| Intensity          | : | 50000 Lux (A)<br>100000 Lux(D) |
| Color Temp.        | : | 4300 K (A)                     |
| CRI                | : | 93(Dichroic)                   |
| Bulb specification | : | 24v                            |
| Power              | : | AC 230V                        |
| UPS System         | : | Optional                       |

(Battery Back show in picture is at extra cost)



**Model No. KSS - 211**  
**LED OPERATION THEATRE LIGHT (4+4)**

**Specifications :**

|                              |   |                   |
|------------------------------|---|-------------------|
| Intensity                    | : | 1,50,000 Lux +10% |
|                              |   | 1,50,000 Lux +10% |
| Size of light field          | : | 12-13 cm          |
| Colour Temp.                 | : | 4000-5000 K       |
| Colour Reduction index       | : | 90                |
| LED Life                     | : | >45000 Hrs.       |
| Number of LED                | : | 15x4+15x4         |
| Diameter of Light (External) | : | 500mm x2          |
| Brightness Control           | : | Digital 9 Steps   |
| Power Supply                 | : | 220V/50Hz AC      |
| Focusing                     | : | Adjustable        |



**Model No. KSS - 212**  
**LED OPERATION THEATRE LIGHT (4+3)**

**Specifications :**

|                        |   |                   |
|------------------------|---|-------------------|
| Intensity              | : | 1,50,000 Lux +10% |
|                        |   | 1,00,000 Lux +10% |
| Size of light field    | : | 12-30 cm          |
| Colour Temp.           | : | 4000-5000 K       |
| Colour Reduction index | : | 90                |
| LED Life               | : | >45000 Hrs.       |
| Number of LED          | : | 15x4+15x3         |
| Diameter of Light      | : | 500mm x2          |
| Brightness Control     | : | Digital 9 Steps   |
| Power Supply           | : | 220V/50Hz AC      |
| Focusing               | : | Adjustable        |



**Model No. KSS - 213**  
**LED OPERATION THEATRE LIGHT**

**Specifications :**

|                        |   |                         |
|------------------------|---|-------------------------|
| Intensity              | : | 1,20,000 Lux +10%       |
| Size of light field    | : | 12-30 cm                |
| Colour Temp.           | : | 4000-5000 K             |
| Colour Reduction index | : | 90                      |
| LED Life               | : | >45000 Hrs.             |
| Number of LED          | : | 15x4                    |
| Diameter of Light      | : | 500mm x 1 (Entire Dome) |
| Brightness Control     | : | Digital 9 Steps         |
| Power Supply           | : | 220V/50Hz AC            |
| Focusing               | : | Adjustable              |



**Model No. KSS - 214**

**MOBILE LED EXAMINATION LIGHT**

**Specifications :**

|                        |                     |
|------------------------|---------------------|
| Intensity              | : 1,20,000 Lux +10% |
| Size of light field    | : 12-13 cm          |
| Colour Temp.           | : 4000-5000 K       |
| Colour Reduction index | : 90                |
| LED Life               | : >45000 Hrs.       |
| Number of LED          | : 15x4              |
| Diameter of Light      | : 500mm             |
| Brightness Control     | : Digital 9 Steps   |
| Power Supply           | : 220V/50Hz AC      |
| Focusing               | : Adjustable        |

**Model No. KSS - 215**

**LED EXAMINATION LIGHT**

**Specifications :**

|                        |                   |
|------------------------|-------------------|
| Intensity              | : 90,000 Lux +10% |
| Size of light field    | : 12-30 cm        |
| Colour Temp.           | : 4000-5000 K     |
| Colour Reduction index | : 90              |
| LED Life               | : >45000 Hrs.     |
| Number of LED          | : 15x3            |
| Diameter of Light      | : 500mm           |
| Brightness Control     | : Digital 9 Steps |
| Power Supply           | : 220V/50Hz AC    |
| Focusing               | : Adjustable      |



**Model No. KSS - 216**  
**EXAMINATION LIGHT LED**

**Specifications :**

|                        |                   |
|------------------------|-------------------|
| Intensity              | : 40,000 Lux +10% |
| Size of light field    | : 12-20 cm        |
| Colour Temp.           | : 4500 K          |
| Colour Reduction index | : 90              |
| LED Life               | : >45000 Hrs.     |
| Number of LED          | : 12              |
| Diameter of Light      | : 300mm           |
| Brightness Control     | : Digital         |
| Power Supply           | : 220V/50Hz AC    |



**Model No. KSS - 217**

**AEROSOL DISINFECTOR**

**Specification:**

This model, Ideally suited for decontamination of operation theatres, Wards in Hospitals & Nursing Homes, Pharmaceutical Labs & Production units, Animal housed and Poultry sheds, Hatcheries & EGG Storages, Silk worm rearing houses in sericulture industries etc. Instantly converts liquids formalin into true aerosol particles. Due to high aerodynamic stability the formalin droplets remain suspended in the atmosphere for longer duration and because of penetrating properly ensure 100% sterile germ free and Zero-Bacteria Environment. Eliminates the Cumbersome process involved in the conventional method of fumigation. Can be easily operated even by a ward-boy. A single unit can sterilize 250 cu. mtrs. about 9000 cu. ft. in just 30-60 minutes. Suitable to work on 220 V, Single Phase, 50 Hz, Ac Supply. Available in Aluminium and Stainless Steel body in following capacities. (A) 1.5 Litres (B) 5.0 Litres



**Model No. KSS - 218**

**VERTICAL HIGH PRESSURE STEAM STERILIZER AUTOCLAVE**

**Specification:**

High pressure water steam electrically heated vertical sterilizer used for sterilizing of surgical instruments, Dressing materials, Linen, Rubber, Plastic material, Injection liquids by means of steam under pressure of 15 to 20 PSI. (Adjustable) Tripple chamber with inner chamber and steam jacket made of thick stainless steel jacket made of thick stainless steel jacket. All sterilizers are hydraulically tested upto 2.5 times of working pressure. Outer chamber made of stainless steel sheet reinforced with M.S. sheet. Lid is made of thick stainless steel plate single piece and closed by wing nuts arrangement. Jointless neoprene gasket does not allow any leakage. Fitted with double safety valve, Water level indicator, water inlet and drain valves. Supplied complete with cord and plug but without dressing bins, To work on 220 volts, Single Ph, 50 Hz, Ac supply.

**Available in following Sizes:**

| Dia | Depth | Element | Load  | Total Load |
|-----|-------|---------|-------|------------|
| (A) | 300   | 500     | 2Nos. | 2 Kw 4 Kw  |
| (B) | 400   | 600     | 2Nos. | 3 Kw 6 Kw  |

**Optional Accessories:**

- Automatic low water level cut-off device
- Pressure control switch



**Model No. KSS - 219**

**AUTOCLAVE (PORTABLE) SINGLE / DOUBLE DRUM**

**Specification:**

- Seamless construction made of aluminium alloy with turning on lid highly polished.
- It comes with one stop cock continues bleeding valve, safety valve and pressure guage.
- Size : 9"x11", 12"x12", 12"x14", 12"x16", 12"x 22", 14"x16"



**Model No. KSS - 220**

**SCRUB STATION**

**Specification:**

- Scrub station wall mounted Scrub Station, made of heavy duty.
- Stainless Steel 304 grade designed for splash-free usage touching the taps soap dispensers directly with the handle.



**Model No. KSS - 221**  
**SEAMLESS DRESSING**  
**DRUMS (S.S)**

**Specification:**  
304/202 G - S.S.  
Sizes : 6"x4", 6"x6", 9"x5", 9"x9",  
11"x5", 11"x9", 14"x9", 15"x12"



**Model No. KSS - 222**  
**INSTRUMENT STERILIZER**  
**ELECTRIC (S.S.)**

**Specification:**  
304/202 G- S.S.  
Sizes : 12"x6"x5.5", 14"x6"x5.5",  
16"x6"x5.5", 16"x8"x7.5",  
18"x8"x7, 20"x8"x7.5", 24"x6"x7.5",  
24"x12"x7.5", 24"x12"x12"



**Model No. KSS - 223**  
**KIDNEY TRAY (S.S)**

**Specification:**  
304/202 G - S.S.  
Sizes : 6", 8", 10", 12"



**Model No. KSS - 224**  
**BASIN (S.S)**

**Specification:**  
304/202 G - S.S.  
Sizes : 6", 8", 10", 12", 14" 16"



**Model No. KSS - 225**  
**CATHETER TRAY (S.S)**

**Specification:**  
304/202 G - S.S.  
Sizes : 8"x3", 16"x4"



**Model No. KSS - 226**  
**INSTRUMENT TRAY (S.S)**

**Specification:**  
304/202 G - S.S.  
Sizes : 8"x6", 9"x6", 10"x8", 11"x7",  
12"x8", 12"x10", 15"x12", 18"x12"



**Model No. KSS - 227**  
**ENEMA CAN (S.S.)**

**Specification:**  
304/202 G - S.S.



**Model No. KSS - 228**  
**SPUTUM MUG (S.S.)**

**Specification:**  
304/202 G - S.S.



**Model No. KSS - 229**  
**LOTION BOWL (S.S.)**

**Specification:**  
304/202 G - S.S.



**Model No. KSS - 230**  
**BED PAN (S.S.) MALE /**  
**FEMALE WITH COVER**

**Specification:**  
304/202 G - S.S.



**Model No. KSS - 231**  
**URINAL POT (S.S.)**  
**MALE / FEMALE**

**Specification:**  
304/202 G - S.S.



**Model No. KSS - 232**  
**CIDEX TRAY**

**Specification:**  
(A) 10" x 5" x 3" (B) 17" x 5" x 3"  
(C) 26" x 5" x 3", (D) 36" x 5" x 3"



**Model No. KSS - 233**  
SYRINGE & NEEDLE DESTROYER



**Model No. KSS - 234**  
NEEDLE CUTTER STAINLESS STEEL



**Model No. KSS - 235**  
SYRINGE & NEEDLE DESTROYER



**Model No. KSS - 236**  
LED X-RAY VIEWER BOX

**Specification:**

**Body:** Powder Coated MS Cabinet with Aluminium Frame

**View:** Single, Double, Triple View Available

**Sheet:** 2 mm Acrylic Sheet

**Power:** 12 V, 50-60Hz, 2 Amp

**Bulb:** 264 High Glossy Led Bulbs

**Parts:** Heavy Duty Switch & Socket Used

**Safety:** Shock Proof Mechanism

**Holding:** Plastic Clips for X-ray film Holding

**Handling:** Light Weight & Portable

**Total Lux:** 1000

**LED Bulb life:** 10,000 hours

**Single View**

X-ray View size: 14"X17"

Overall size: 17"X20"X2"

**Double View**

X-ray View size: 28"X17"

Overall size: 31"X20"X2"

**Triple View**

X-ray View size: 42"X17"

Overall size: 45"X20"X2"



**Model No. KSS - 237**  
INSTRUMENT CABINET

**Specification:**

- Made form CRC sheet three side closed.
- Front glass door.
- Pre-treated and Epoxy Powder Coated.



**Model No. KSS - 238**  
WAITING AREA CHAIR

**Specification:**

- Perforated Sheet.
- Galvanized arm rest.
- All S.S.

[www.kiranmedicareequipments.com](http://www.kiranmedicareequipments.com)

[www.kiransurgicalsystems.com](http://www.kiransurgicalsystems.com)



**Delhi :** A - 39, Street No.5, Chander Industrial Estate, I.P Extension, Delhi - 110 092

**Gujarat :** 308/4, GIDC Estate, Makarpura Industrial Estate, Baroda – 390 010,

**Andhra Pradesh :** No. 6-5-351/5, Astal Heights, Road No. 1, Banjara Hills, Hyderabad – 500 034

**Tamil Nadu :** No.33-B4, Sidco Industrial Estate Phase - II, Sundarapuram, Coimbatore - 641 021

**Kerala :** 333/9, Para Road, Near NIDA, Kanjikode, Palakkad - 678 621

**Helpline :** +91 9751725000, +91 9751625000, + 91 9751425000, + 91 9751215000, + 91 9751245000

**Email :** kirangroup108@gmail.com, kiranindia108@gmail.com, kiranbiz108@gmail.com

**Website :** www.kirangroupofcompanies.com, www.kiran.biz, www.kiranshop.in, www.kiranindia.in, www.kiransurgicalsystems.com, www.kiranmedicareequipments.com, www.kiranpump.com, www.kiransolarsystems.com, www.kiranled.com, www.kiranincinerationsystems.com, www.kiranequipments.com, www.kirancompressors.com, www.kiranappliances.com, www.kiranairventilators.com, www.kiranfrpproducts.com, www.kiranpebsystems.com, www.kiranpanels.com, www.kiranfabssystems.com